Scintillating bolometers for Double Beta Decay search

L.Gironi^{*,a,b}

^aDipartimento di Fisica - Universit`a di Milano Bicocca, Italy b INFN - Milano Bicocca, Italy

Abstract

In the field of Double Beta Decay (DBD) searches, the use of high resolution detectors in which background can be actively discriminated is very appealing. Scintillating bolometers containing a Double Beta Decay emitter can largely fulfill this very interesting possibility. In this paper we present the latest results obtained with CdWO₄ and CaMoO₄ crystals. Moreover we report, for the first time, a very interesting feature of CaMoO₄ bolometers: the possibility to discriminate $β$ - $γ$ events from those induced by α particles thanks to different thermal pulse shape.

Key words: Double Beta Decay, Bolometers, Scintillating, CdWO ⁴, CaMoO ⁴, Pulse shape discrimination (PSD) *PACS:* 23.40B, 07.57.K, 29.40M

1. Introduction

The experimental evidence for neutrino oscillations clearly showed that the neutrino is a finite-mass particle. Anyway, two big questions concerning the neutrino are still unsolved: its nature (Dirac or Majorana) and the absolute value of its mass. Neutrinoless Double Beta Decay (0 νDBD) is, at present, one of the most sensitive method to study neutrino properties [\[1\]](#page-3-0) and bolometers are - together with germanium diodes - the detectors which have provided the best results so far. Bolometers allows the application of the so-called calorimetric approach, where the detector is composed of the same material candidate to th e decay, and in the mean time allow the study of many isotopes by means of high resolution devices (FWHM around 0.2-0.5 % above 2500 keV), this latter feature being necessary to resolve the searched peak from background.

The purpose of the so-called 0vDBD new generation experiments is to reach a sensitivity on the mass of the neutrino of the order of 50 meV, approaching - but probably not completely covering - the region of the inverse hierarchy [\[2](#page-3-1), [3,](#page-3-2) [4](#page-3-3), [5\]](#page-3-4). Th e high sensitivity required imply excellent energy resolution, a low number of spurious counts within the region of interest (low background) and a high quantity of the isotope on which the study focuses (large mass).

Further improvements cannot rely simply on the mass increase or on a better energy resolution but will require the implementation of innovative techniques for background discrimination.

In the case of a scintillating bolometer the double independent read-out (heat and scintillation) will allow, thanks to the different scintillation Quenching Factor between α and β - γ , the suppression of α background events. These have been identified as the most important background source in bolometers

dedicated to Double Beta Decay searches [\[6,](#page-3-5) [2](#page-3-1)]. Furthermore , the contribution of environmental gammas can be strongly re duced using a scintillating bolometer containing a DBD emitter whose transition energy exceeds the highest gamma line due to natural radioactivity (the 2615 keV line of ²⁰⁸Tl). Examples are ¹¹⁶Cd (Q_{ββ} = 2805 keV), ⁸²Se (Q_{ββ} = 3000 keV) and ¹⁰⁰Mo (Q_{ββ} = 3030 keV). Finally, this technique is also extremely helpfu l for rejecting other unavoidable sources of background such as direct interactions of neutrons.

2. Environmental Background

In a bolometric experiment, the signature of the neutrinoless Double Beta Decay (0vDBD) is a peak at the $Q_{\beta\beta}$ value of the transition, while in ordinary double beta decay (2 νDBD) is a continuous beta-like spectrum. However, thanks to the excellent energy resolutions of bolometers, 2νDBD is seldom a significant source of background.

However, there are various sources that give rise to spuriou s counts in the region of interest such as external γ background, neutrons, surface contaminations, ²³⁸U and ²³²Th internal contaminations and cosmogenic activity. As already discussed, a way to eliminate the problem of γ lines is to study isotopes with transition energy above the 2615 keV line. This is the 208 Tl line that is the highest energy γ -ray line from natural radioactivity. Above this energy there are only extremely rare high energy $γ$ from ²¹⁴Bi: the total Branching Ratio in the energy window from 2615 up to 3270 keV is 0.15 % in the ²³⁸U decay chain. All the Double Beta Decay isotopes with $Q_{\beta\beta} > 2615$ keV are listed in Table [1.](#page-1-0)

Surface contaminations represents the main source of background for bolometric detectors such as CUORICINO [\[6](#page-3-5)]. This source of background plays a role in almost all detectors but turns out to be crucial for fully active detectors, as in the case of bolometers: a radioactive nucleus located within a few μ m of

[∗]Corresponding author: tel: +39 02 64 48 21 07; fax: +39 02 64 48 24 63 *Email address:* luca.gironi@lngs.infn.it (L.Gironi)

Table 1: Double beta decay isotopes with endpoint energies above the ²⁰⁸Tl line.

Isotope	$Q_{\beta\beta}$ [MeV]	natural abundance
116 Cd	2.80	7.5 %
${}^{82}Se$	3.00	9.2 %
100M ₀	3.03	9.6%
96Zr	3.35	2.8%
150 _{Nd}	3.37	5.6 %
	4.27	0.19%

Figure 1: Operating principle of scintillating bolometers. The release of energy inside a scintillating crystal follows two channels: light production and thermal excitation. The heat is read out by a temperature sensor (NTD) glued on the primary crystal while the light is read by a second crystal (light detector) where it is completely converted into heat.

a surface facing the detector can emit an α particle (the energy of which is, in most of the cases, between 4 and 8 MeV). Those particles can lose part, or even all, of their energy in the few microns of this dead layer before reaching the bolometer. The energy spectra read by the bolometer, therefore, will result in a continuum between 0 and 4-8 MeV, covering, unfortunately all the possible $Q_{\beta\beta}$ values. Furthermore, the same mechanism holds in the case of surface contaminations on the bolometer itself.

3. Scintillating Bolometer and Light Detectors

A scintillating bolometer [\[7\]](#page-3-6) is, in principle, a very simple device. It is composed of a calorimetric mass or absorber (a scintillating crystal coupled with a thermometer) and a Light Detector, facing it, and able to measure the emitted photons(see Fig. [1\)](#page-1-1). The driving idea of this hybrid detector is to combine the two information available: the heat (i.e. that fraction of the energy released in the crystal absorber which is converted into phonons) and the emitted scintillation light (i.e. that small fraction of the energy which is converted into photons). Thanks to the different scintillation yield (or scintillation Quenching Factor) of different particles (namely $β - γ$, α and neutrons) they can be very efficiently discriminated.

The idea to use a bolometer as light detector was first developed by Bobin et al. [\[9\]](#page-3-7) and then optimized [\[10,](#page-3-8) [11\]](#page-3-9) for Dark Matter (DM) searches. Starting from that work we developed a thermal light detector to be used for DBD search.

Our light detector is constituted of a "dark" thin bolometer, normally a crystal of Ge or Si, that can absorb scintillation photons and give a measurable thermal signal. This means that it has the characteristic time constant of bolometers (20 - 500 ms) and have an important advantages respect PMTs: they are sensitive over an extremely large band of photon wavelength.

4. Experimental setup and Results

In these years different crystals were tested in dedicated runs to study their characteristics such as their thermal response, light yield and radio purity. The last test was performed at the beginning of 2009 in order to evaluate the efficiency of discrimination of this technique. In this test were measured some new CdWO₄ and CaMoO₄ crystals.

The detectors are mounted in an Oxford 200 3 He/⁴He dilution refrigerator located deep underground in the National Laboratory of Gran Sasso. The cryostat is surrounded by about 20 cm of lead in order to reduce environmental γ radioactivity. Furthermore, the crystal set-up was mounted below about 5.5 cm of Roman lead in order to further decrease the environmental background. The cryostat is also surrounded by about 7 cm of polyethylene $(CH₂)$ to thermalize fast neutrons and about 1 cm of CB_4 that, thanks to the high neutron capture cross section for thermal neutrons of ^{10}B , allows to reduce the neutron flux on the detectors.

The temperature sensors are Neutron Transmutation Doped (NTD) Ge thermistors glued to each crystal. A Si resistor of \sim 300 kΩ is attached to each crystal and acts as a heater to stabilize the gain of the bolometer [\[18\]](#page-3-10). The read-out of the thermistors is performed via a cold (∼110 K) preamplifier stage located inside the cryostat [\[19\]](#page-3-11). The room temperature front-end [\[20\]](#page-3-12) and the second stage of amplification are located on the top of the cryostat. After the second stage there is an antialiasing filter (a 6 pole roll-off active Bessel filter). The ADC is a NI USB device (16 bit 40 differential input channels) located in a small Faraday cage. The signals (software triggered) are sampled in a 512 ms window with a sampling rate of 2 kHz.

The data analysis is completely performed off-line. It uses an Optimal Filter technique [\[2\]](#page-3-1) to evaluate the pulse amplitudes and to compare pulse-shapes with detector response function. Events not caused by interactions in the crystals are recognized and rejected on the basis of this comparison. Pile-up pulses are identified and dealt with. Accepted-pulse amplitudes are then corrected using the variation in the gain measured with the heat pulses from the Si resistors.

Since, for bolometric detectors, surface alpha contaminations represents the main source of background in DBD region [\[6\]](#page-3-5), in this run ²³⁸U sources was faced to the crystals. The purpose of these sources is to have a high number of alpha counts in the 2-3 MeV region in order to evaluate the efficiency of discrimination of this technique. For this reason, after simulations with the GEANT4 [\[12\]](#page-3-13) package and measurements with a Silicon Barrier Detectors, alpha sources were covered with an aluminated

Figure 2: Scatter plot of Heat (new CdWO₄ crystal, h = 50mm, \varnothing = 40mm, 508g) vs Light. Calibration performed on γ peaks.

Mylar foil (thickness \sim 6 µm). In such a way alpha particles lose part of their energy in this dead layer before reaching the bolometer. The energy spectra read by the bolometer, therefore, will result in a continuum between 0 and ∼4 MeV.

*4.1. CdWO*⁴

In 2008 we tested an array of four CdWO₄ Double Beta Decay scintillating bolometers, read by only one light detector [\[13\]](#page-3-14). The array consists of a plane of four $3x3x3$ cm³ crystals and a second plane consisting of a single $3x3x6$ cm³ crystal. This test demonstrate the technical feasibility of this technique through an array of detectors, and perform a long background measurement in the best conditions in order to prove the achievable background in the 0νDBD region.

In 2009 a measurement was done with a new, bigger, cylindrical CdWO₄ crystal (h = 50mm, \varnothing = 40mm, 508g). The main difference between this crystal and the crystals measured previously is that this crystal has perfect optical surfaces. This feature has led to a significant reduction of the width of the alpha band compared to measurements made previously. In figure [2](#page-2-0) it is possible to observe the high number of counts in the region 2-4 MeV due to alpha source facing the detector. Data relates to a measurements made with the new crystal with a live time of ∼ 93 h. This plot shows how it is possible to separate very well the background due to α particles from the β - γ region. In particular it should be noted that α continuum is completely ruled out thanks to the combined measurement of heat and scintillation.

The resolution of the main bolometer is affected by the fluctuation of the fraction of energy going in the heat and light channels. Since the energies of the two channel are anticorrelated, it is possible reconstruct the intrinsic energy resolution of the bolometer. Rotation the axis by θ we obtain a new plot in which $E_{rot} = E_{heat} cos\theta + E_{light} sin\theta$ is the new energy. After rotation, the energy resolution at the 2615 keV line is 12 keV FWHM.

Figure 3: Scatter plot of Heat (CaMoO₄ crystal, h = 40mm, \varnothing = 35mm, 158g) vs Light. Calibration performed on γ peaks.

*4.2. CaMoO*⁴

Recently CaMoO⁴ was intensively studied as possible cryogenic scintillation bolometric detector for experiments to search for DBD and Dark Matter [\[14,](#page-3-15) [15](#page-3-16), [16\]](#page-3-17). Below we discuss possible use of CaMoO₄ scintillators as potential detectors in search for 0ν DBD decay of 100 Mo.

In the case of scintillating bolometer, thanks to the excellent energy resolution achievable respect to scintillators, the 2νDBD of 100 Mo will not contribute to the 0ν DBD peak. Nevertheless the background due to $2v$ DBD of ⁴⁸Ca will limit the reachable sensitivity. In fact contribution from 2ν DBD of ⁴⁸Ca could be very dangerous (measured half-life of this process is $T_{1/2}(2\nu)$ = 4×10^{19} yr [\[17\]](#page-3-18)). ⁴⁸Ca present in natural Ca (a.i.=0.187%) produces a background of about 10−³ counts/keV/kg/y in the $0v$ DBD region of 100 Mo.

Even if sensitivity of $CaMoO₄$ crystal is limited by 2 ν DBD decay of ⁴⁸Ca, this crystals shows, for the first time, the possibility to discriminate β - γ from α thanks to different shape of the thermal pulses.

4.2.1. Bolometric Pulse Shape Analysis

The release of energy inside a scintillating crystal follows two channels: light production and thermal excitation. The heat is read out by a temperature sensor glued on the primary crystal while the light is read by a second crystal (light detector) where it is completely converted into heat. In this run the CaMoO⁴ crystal was measured with an antialiasing filter with a high cutoff frequency of 120 Hz. This feature allow us to study in detail the rise time of the thermal signal in the $CaMoO₄$ crystal.

The light pulses produced by α particles and β - γ in the CaMoO⁴ scintillator present different decay time constants and relative intensities. Moreover it was observed a temperature dependence of averaged decay time of the light pulses [\[21\]](#page-3-19). This difference between the light pulses and the increase of the decay time costant of the light pulses at low temperature modify the shape of the thermal signal measured on the same scintillating crystal.

Figure 4: Scatter plot of Heat vs Rise Time in $CaMoO₄$ crystal.

As can be seen in figure [4](#page-3-20) the difference in the rise time of the thermal pulses measured in the scintillating crystal allows to discriminate β - γ events from those induced by α particles. In this run the efficiency of discriminations of α events in the 2.7 -8 MeV region is 99.7% and can be improved by measurements devoted to study this feature.

The possibility to discriminate α events thanks to the different pulse shape allow to consider an experiment without light detector and reflecting sheets. This will allow not only to construct a 0νDBD experiment in a much easier way but also to develop a new technique in order to investigate with a very high sensitivity surface contaminations. Indeed the surface of a scintillating bolometer that can be faced to a sample is much bigger than the surface of a standard Silicon Barrier Detector. Moreover scintillating bolometers have not dead layer and can reach energy resolution orders of magnitude better than standard Silicon Barrier Detector.

5. Conclusion

We have performed simultaneous detection of heat and light on different Double Beta Decay scintillating crystals. Very good results was obtained with CdWO₄ and CaMoO₄ crystals. CdWO⁴ has excellent characteristics that make it an excellent candidate as a crystal for the study of 0νDBD.

The sensitivity of $CaMoO₄$ is limited by the background due to 2νDBD of ⁴⁸Ca. However, this crystal shows, for the first time, the possibility to discriminate α events through the study of the shape of the thermal signal. This feature can be used to develop detectors for measurements of surface contamination with high sensitivity.

6. Acknowledgments

This work was funded and implemented under the Bolux experiment of INFN.

References

- [1] S.J. Freedman, B. Kayser, [\[arXiv:physics](http://arxiv.org/abs/physics/0411216)/0411216] (2004)
- [2] C. Arnaboldi *et al.*, Nucl. Instr. Meth. A **518** (2004) 775
- [3] S. Scheonert *et al.*, Proc. Neutrino 2006 (Santa Fe, 2006)
- [4] C.E. Aalseth *et al.*, Yad. Fiz. **67**, 11 (2004)
- [5] A. Piepke *et al.*, Proc. Neutrino (Santa Fe, 2006)
- [6] C. Arnaboldi *et al.*, Phys. Rev. Lett. **95** (2005) 142501
- [7] A. Alessandrello *et al.*, Physics Letters B **420** (1998) 109
- [8] S. Pirro *et al.*, AIP Conference Proceedings **785** (2005) 177
- [9] C. Bobin *et al.*, Nucl. Instr. Meth. **A 386** (1994) 453
- [10] G. Angloher *et al.*, Astropart. Phys. **23** (2005) 325
- [11] S. Cebrian *et al.*, Phys. Lett. B **563** (2003) 48
- [12] Geant4 Collaboration, J. Allison et al., IEEE Trans. Nucl. Sci. **53**, (2006) 270.
- [13] L.Gironi *et al.*, Optical Materials, DOI:10.1016/j.optmat.2008.09.014 [arXiv:physics/0809.5126v15]
- [14] A. Senyshyn et al., Phys. Rev. B 73 (2006) 014104.
- [15] V.B. Mikhailik and H. Kraus, J. Phys. D: Appl. Phys. 39 (2006) 1181.
- [16] S. Pirro et al., Phys. Atom. Nucl. 69 (2006) 2109.
- [17] V.I. Tretyak and Yu.G. Zdesenko, At. Data Nucl. Data Tables 80 (2002),
- p. 83.
- [18] A. Alessandrello *et al.*, Nucl. Instr. Meth. **A 412** (1998) 454
- [19] C. Arnaboldi *et al.*, Nucl. Instr. Meth. **A 559** (2006) 826
- [20] C. Arnaboldi *et al.*, IEEE Trans. Nucl. Sci. **49** (2002) 2440
- [21] A.N. Annenkov et al., Nucl. Instr. Meth. A 584 (2008) 334.