A fibered interference scanning optical microscope for living cell imaging

Jean-Baptiste Decombe, Wilfrid Schwartz, Catherine Villard,

Hervé Guillou, Joël Chevrier, Serge Huant, and Jochen Fick*

Institut Néel, CNRS & Université Joseph Fourier, 25 Avenue des Martyrs, BP 166,

38042 Grenoble, France

* Corresponding author: jochen.fick@grenoble.cnrs.fr

We report on the 3-dimensional imaging of biological cells including living neurons by a recently developed fibered interferometric Scanning Optical Microscope. The topography and surface structure of mouse fibroblasts and hippocampal neurons are clearly revealed. This straightforward technique allows fast, high resolution observation of samples in liquids without lengthy alignment procedures or costly components. © 2023 Optical Society of America

OCIS codes: (060.2350) Fiber optics imaging; (170.5810) Scanning microscopy; (180.3170) Interference microscopy

1. Introduction

The fast progress in nanosciences experienced over the last decades was made possible thanks to the development of sub-micron resolution, scanning visualization techniques such as Atomic Force Microscopy (AFM) [1], Scanning Tunneling Microscopy (STM) [2], or Near-field Scanning Optical Microscopy (NSOM) [3,4]. In the field of biology the AFM is nowadays widely used. However routine observation of biological cells in water requires dedicated instruments [5]. Moreover, the force exerted by the AFM on the cell may influence its behaviort. Finally AFM measurements often need lengthy alignment procedures.

Phase contrast microscopy is also widely used in biology. It allows for noncontact and fast observation of cells. Its resolution is, however, limited by the Rayleigh criterion. 3-dimensional quantitative observation is only possible with complex methods such as the Hilbert Phase Microscope [6] or the Digital Holographic Microscope [7,8], requiring well calibrated optical systems or heavy numerical data processing.

In this paper, we present imaging of biological cells by a new kind of high resolution optical microscope, which we call interferometric Scanning Optical Microscope (iSOM) [9]. The topography and surface structure of mouse fibroblasts and living hippocampal neurons are clearly revealed by iSOM. Its working principle is based on the interference of the internally reflected light at an optical fiber tip with the light reflected by the surface facing the tip. iSOM is working in the optical far field, i.e. at tip-sample distances of the order of a few microns, thus avoiding any mechanical contact with the sample. This characteristic is of particular importance for the observation of living cells. Samples with topography in the micron range can thus be imaged at high resolution (typically half the wavelength laterally and better than 10 nm axially with the fiber) without any closed-loop tracking of the fiber tip height as done in AFM or NSOM.

This important feature allows to use less complex electronic control units and much faster recording speeds. In addition, iSOM can be developed at reduced costs and is well suited for integration with other measurement techniques, as only an optical fiber has to be positioned near the sample surface. Therefore, lengthy optical alignments can be avoided. Similarly, its integration ability into a micro-incubator will permit long term living cell observation. Because iSOM can work in various environments, such as liquids, vacuum, or even low temperatures, we expect it to cover a broad range of applications. Here, we focus on cell imaging in liquid.

2. Experimental

2.A. interferometric Scanning Optical Microscope (iSOM)

The heart of iSOM is the optical fiber tip which is scanned relative to the sample at distances in the order of a few microns. In our set-up the fiber is fixed and the sample is moving (Fig. 1). Two different types of linear translation stages are used for positioning and scanning. Translation stages with piezoelectric inertial motors



Fig. 1. Schematic view of the iSOM set-up

(Mechonics MS30) with 15 nm resolution and 18 mm range allow for selecting laterally the region to be imaged. For image scanning a three axis piezoelectric nanopositioning system (Mad City Labs Nano-3D200) is used. Typical images of $200 \times 200 \ \mu m^2$ with 256×256 pixels can be recorded in about 4 minutes using a speed of 0.5 s/line. The Nano-3D200 is also used for fine tuning the fiber tip height. For the coarse approach the fiber holder is mounted on an manual micrometer translation stage.

The optical part of the iSOM is totally fibered. Thus no alignment is necessary and the components can easily be changed. Two different laser diodes were used for imaging: a fiber-coupled laser diode emitting at 670 nm (Schäfter & Kirchhoff) and a pigtailed laser diode with 830 nm wavelength (Thorlabs LPS-830-FC/APC). The output of the two laser diodes are combined by a dielectric mirror mounted inside a fiber coupled beam combining unit from Schäfter & Kirchhoff. Thus the wavelength can be changed without touching the optical set-up. A wide-band polarization insensitive optical circulator (AC photonics PIOC-3-83-P-2-1-1-11) is used to guide the incoming light to the iSOM fiber tip and the reflected intensity to a fiber-coupled amplified photoreceiver (Newfocus 2001-FC). The whole set-up is computer controlled using an available Nanonis SPM controller. A data acquisition computer card with controlling software would, however, be sufficient to run the experiment.

The iSOM probe tip is the key element that determines the quality of the measurements. Their strong divergence prohibits the use of cleaved fibers. Homemade optimized chemically-etched fibers with tip angles of 20° and 60 nm tip aperture size were used [10, 11]. These fibers allow efficient imaging for fiber sample distances of up to 50 μ m.

2.B. Fibroblast

NIH-3T3 fibroblasts are cultivated as described in [12,13]. In short, they are cultivated in DMEM complemented with 10% of inactivated fetal calf serum, streptomycin and penicillin. They are harvested at 80% of confluence with trypsin/EDTA and plated on glass substrates coated with fibronectin fragment type III domain 7-10 [14]. Prior to fibronectin coating, glass substrates are covered with a sputtered 1000 Å titanium layer. The titanium layer is patterned using ion milling to yield a numbered grid structure that eases the imaging of the same cells with various techniques.

After the fibronectin coating, the cells are inoculated onto the substrate and spread for several hours. They are fixed by two different methods. In the first experiment PBS with 3% paraformaldehyde supplemented with 2% sucrose was used for 5 min at room temperature. In the second experiment methanol fixation was used by dipping PBS rinsed substrates into cold methanol at -20°C for 10 minutes.

2.C. Neurons

Neurons are constituted by a cell body (or "soma") and a burgeoning of fine cables (neurites). Neuronal cultures were prepared as described by Dotti et al. [15]. Briefly, hippocampi brain tissues from E18.5 (embryonic day 18.5) mice were removed and digested in 0.25% trypsin in Hepes-buffered Hanks' balanced salt solution at 37°C for 15 min. After manual dissociation, cells were plated at a concentration of 5,000-15,000 cells/cm² in DMEM-10% fetal bovine serum on different substrates. Usual poly-L-lysine (PLL)-coated coverslips lead to uniform cultures with arbitrary cell connections, while patterned surfaces could be used to direct neural growth. In this last case, PLL patterns were transferred on silanized silicon chips [16] using UV classical photolithography steps, including Shipley S1805 photoresist spinning, insulation through a mask, development, PLL deposition (1mg/ml one night) and lift-off in an ultra-sound ethanol bath. One hour after plating, the medium was changed to DMEM containing B27 and N2 supplements. All culture media present a stable pH value of 7.8 under the 5% CO_2 atmosphere provided by commercial incubators regulated at 37°C. Living cells were observed in this work after 21 days of *in vitro* culture (DIV 21) at room temperature and ambient atmosphere. Under these conditions, cell survival

is of a few hours.

Cell development can also be stopped by a process called fixation that kept the initial cellular morphology. Neurons were here treated in fixative solutions containing paraformaldehyde (PFA) 4% and HMDS [17] for an optimal shape preservation.

3. Working principle

iSOM is exploiting the interference between the light intensity reflected at the fiber tip I_f and the sample I_s (Fig. 2.a). In a first approximation, the overall reflected intensity I_r can be calculated by:

$$I_r = I_f + I_s + 2\sqrt{I_f I_s} \cos(\Phi) = I_0 \left(R_f + \tilde{R}_s + 2\sqrt{R_f \tilde{R}_s} \cos(\Phi) \right)$$
(1)

with I_0 the incident intensity, R_f the inner fiber tip reflectance, $\tilde{R}_s = R_s T_{f \Downarrow} T_{f \uparrow}$ the global sample reflectance, $T_{f \Downarrow / \uparrow}$ the transmission for fiber out-/in-coupling, and $\Phi = 4\pi/\lambda z + \Phi_s$ the total phase-shift of I_s with respect to I_f . Φ_s is the phase-shift encountered during the reflection. This term becomes important for samples with more complex structures such as biological cells.

The cosine term in Eq. 1 results in interference fringes in the $I_r(z)$ curves, which we call z-scans. These fringes are used in xy-scans to deduce the actual fiber - sample surface distance from $I_r(x, y)$ plots. Vertical resolutions below 1 nm can be attained for highly reflective surfaces.

The z-scan of a bare glass substrate is shown in Fig. 2.b for two wavelengths. In



Fig. 2. (a) Scheme of the iSOM working principle. (b) z-scan curves on a bare glass substrate for two wavelengths. The inset shows the FFT.

this case, the minimal fiber sample distance was estimated at $z_0 \leq 100$ nm. The interference fringes described by Eq. 1. are clearly visible. The measured periods, calculated using Fast Fourier Transformation (FFT, see inset), are 2.87 and 2.25 μ m⁻¹ for 670 and 830 nm wavelengths, respectively. These values are slightly lower than the nominal values of 2.98 and 2.41 μ m⁻¹.

The mean intensity and the visibility of the interference fringes show a maximum at $(z-z_0) \approx 0.7$ and 1.0 μ m for $\lambda = 670$ and 830 nm, respectively. At large distances the beam divergence leads to lower back-coupled intensity I_s . Consequently, the fringes amplitude decreases and the mean intensity converges to I_f , which is independent on z. The reasons for the intensity decrease at small distance are more complex. In fact, for distances smaller than the wavelength, the optical near-field in the vicinity of the

fiber tip becomes more intense than the radiation field [18, 19]. The description of the reflected intensity I_s has to be modified. The in-depth treatment in this case is, however, not the scope of the present paper.

The evolution of the interference fringes for long distances was also studied on fixed neurons (Fig. 3). These curves were recorded above and beside one soma on a fixed neuron culture made on silicon substrate and covered with water. The reference distance above the substrate $z_0 > 5 \ \mu m$ is equal for both curves, but much larger than in Fig. 2. Hence, the intensity maximum observed in Fig. 2 is not visible. The interference fringes remain visible for distances up to 50 μm .

For the z-scan beside the soma, i.e. on a homogeneous surface, one single interference frequency peak is observed at 2.97 μ m⁻¹ in the FFT spectra. For the case above the soma, supplementary frequency peaks on the low frequency side of the main peak are observed. The mixing of these interference frequencies results in the observed long range oscillations. This feature has the same origin as the fringes visible on xy-scans with $z \ge 2\lambda$, i.e. the interference of light reflected by different, nonhomogeneous parts of the sample surface.

4. Imaging results and discussion

4.A. Fibroblasts

The fixed fibroblasts spread on glass substrates are imaged by iSOM and phase contrast microscopy in liquid environment, and by AFM (Fig. 4). The iSOM image is



Fig. 3. z-scan curves above and beside a neuronal sona in water ($\lambda = 830$ nm). The inset shows the FFT.

built by superposition of about 20 single scans of $40 \times 40 \ \mu m^2$ size. On iSOM and phase contrast images, three cells of different forms can clearly be distinguished. The contained informations are however supplementary. The phase contrast microscopy image clearly shows the cell outlines and some structures in the cytosol in and around the nucleus.

The iSOM image depicts the cells' topography by contour lines. The distance between each dark/ bright line corresponds to half the wavelength. The upper right cell shows 11 lines from its border to the center. The corresponding height of 3.7 μ m is in good agreement with this kind of cells. As the reflection on the cell surface is not total, the iSOM image contains more information than solely the position of the interface. In the figure, this is at best visible for the lower left hand cell. Its inner structure is revealed by dark spots. With models yet to be developed it will be possible to get



Fig. 4. Fixed fibroblast cells observed by (a) phase-contrast microscopy, (b,d) iSOM, and (c) AFM. (c) and (b) show the lower left cell visible in (a) and (b).

more informations on the nature of the material inside the cell volume.

The AFM image shows the lower cell visible in Fig. 4.a and b and zoomed in Fig. 4.d. Similar to the iSOM image the AFM one is built from about 10 single scans. The cell shape with its ramifications is clearly visible and the cell nucleus contours are revealed. Compared to the iSOM image the AFM one shows more details at the border of the cell. Its 3-dimensional perception is, however, less intuitive and the higher resolution of the AFM image is at the cost of longer acquisition time. Although the cell imaged



Fig. 5. Calculation of the cell topography from iSOM measurements made at two wavelengths. The white lines indicate the chosen cross section.

by AFM was dried by blowing clean air, its height as measured by AFM qualitatively match the height measured by iSOM showing no qualitative volume change upon drying.

In order to check the interpretation of the lines as contour lines in the iSOM image, we have made the same measurements on a similar cell, but making use of two different wavelengths. The corresponding topographies were calculated for a cross section as shown on Fig. 5. A very good agreement is obtained. As expected, the comparison of the two iSOM images shows a slightly better resolution in the case of a smaller imaging wavelength.

4.B. Neurons

Neurons grown on two kinds of substrates, silicon chips (Fig. 6) and coverslips (Fig. 7), have been observed. The silicon chips are used to investigate the local electric activity



Fig. 6. Fixed neurons on an electronic substrate observed by (a) SEM and (b) iSOM.

of neurons grown on their surface. In order to locate specific neurite positions right above the silicon sensors, cell growth was guided according to specific geometries [20].

The SEM and iSOM images of the same neuron fixed after 3 days of culture on a silicon chip is shown on Fig. 6. The neuron's soma and neurites can clearly be distinguished on the metallic lines. Some details have been changed between the two images. In the SEM image, the neurite leaving on the lower right part of the soma is suspended over the edge of the metal strip. In the iSOM image this same neurite is broken. Moreover the smaller soma visible on the left side of the SEM image has disappeared on the iSOM image.



Fig. 7. Living neurons observed by (a) phase contrast microscopy and (b-c) iSOM.

Even if the SEM image is somehow more precise, the iSOM image clearly reveals the main features and gives a good three dimensional impression. The topography of the soma is visible as interference fringes.

Compared to SEM, one important advantage of iSOM is the possibility to image living cells. The use of non pH-buffered solutions already allowed to keep neurons alive for about three hours, but an improvement of life conditions could be easily implemented around th fibered microscope. Compared to phase contrast microscopy, iSOM allows to obtain much more detailed images, as visualized on "old" living cells forming mature networks where the different neurons are cross-linked by a great number of neurites, sometimes grouped in thick bundles (Fig. 7).

The iSOM image gives a clear impression of the topography and many details are revealed. However, no fringes are visible, even for nominal cell heights of some microns. In fact, the refractive index difference between the solution and the cell is too small to ensure efficient reflection at the cell surface. The light is thus penetrating the cell and reflected by the substrate. The observed image is due to the phase shift encountered during the propagation inside the cell.

5. Conclusion

We have presented the ability of the recently introduced iSOM [9] for imaging different biological cells including living neurons in liquid. The obtained images reveal many details and render a clear perception of the cells' topography. Furthermore, iSOM minimizes the impact on the cells, as it is a contactless poorly invasive method. In addition iSOM can also operate in aqueous environments. Because the microscope head is very compact and requires only a small free space to be available around the sample, it can easily be integrated into existing systems. For instance, its impact on living cells could be further minimized by the implementation of a mini-incubator for temperature and pH control. Under these circumstances, we expect iSOM to be able to track living cells for hours, allowing for example the study of mechanical deformation associated with either growth or the electrical activity of neurons [21].

Acknowledgments

We thank the Thermometry and Calorimetry team and J.F. Motte from Institut Néel for the fibroblast substrate preparation and fiber tip fabrication, respectively.

References

- G. Binnig, C. F. Quate, and C. Gerber, "Atomic force microscope," Phys. Rev. Lett. 56, 930–933 (1986).
- G. Binnig, H. Rohrer, C. Gerber, and E. Weibel, "Surface studies by scanning tunneling microscopy," Phys. Rev. Lett. 49, 57–61 (1982).
- D. W. Pohl, W. Denk, and M. Lanz, "Optical stethoscopy: Image recording with resolution lambda/20," Appl. Phys. Lett. 44, 651–653 (1984).
- A. Drezet, A. Hohenau, J. Krenn, M. Brun, and S. Huant, "Surface plasmon mediated near-field imaging and optical addressing in nanoscience," Micron 38, 427–437 (2007).
- D. Fotiadis, S. Scheuring, S. A. Müller, A. Engel, and D. J. Müller, "Imaging and manipulation of biological structures with the AFM," Micron 33, 385–397 (2002).
- G. Popescu, Y. Park, R. Ramachandra, K. Badizadegan, and M. S. Feld, "Coherence properties of red blood cell membrabe motions," Phys. Rev. E 76, 031902 (2007).
- F. Charrière, N. Pavillon, T. Colomb, C. Depeursinge, T. J. Heger, E. A. D. Mitchell, P. Marquet, and B. Rappaz, "Living specimen tomography by digital holographic microscopy: morphometry of testate amoeba," Opt. Express 14, 7005–7013 (2006).
- 8. B. Rappaz, F. Charrière, C. Depeursinge, P. J. Magistretti, and P. Marquet,

"Simultaneous cell morphometry and refractive index measurement with dualwavelength digital holographic microscopy and dye-enhanced dispersion of perfusion medium," Opt. Lett. **33**, 744–746 (2008).

- 9. J. Chevrier, S. Huant, W. Schwartz, and A. Siria, French Patent 08/59091, Université Joseph Fourier - CNRS - CEA (2008), http://www.wipo.int/pctdb/en/wo.jsp?WO=2010076540.
- P. Hoffmann, B. Dutoit, and R.-P. Salathé, "Comparison of mechanically drawn and protection layer chemically etched optical fiber tips," Ultramicroscopy 61, 165–170 (1995).
- N. Chevalier, Y. Sonnefraud, J. F. Motte, S. Huant, and K. Karrai, "Aperturesize-controlled optical fiber tips for high-resolution optical microscopy," Rev. Sci. Instr. 77, 063704 (2006).
- H. Guillou, A. Depraz-Depland, E. Planus, B. Vianay, J. Chaussy, A. Grichine, C. Albigès-Rizo, and M. Block, "Lamellipodia nucleation by filopodia depends on integrin occupancy and downstream rac1 signaling," Exp. Cell Res. **314**, 478–488 (2008).
- B. Vianay, J. Käfer, E. Planus, M. Block, F. Graner, and H. Guillou, "Single cell spreading on a protein lattice adopt an energy minimizing shape," Phys. Rev. Lett. 105, 128101 (2010).
- 14. D. Leahy, I. Aukhil, and H. Erickson, "2.0 Å crystal structure of a four-domain

segment of human fibronectin encompassing the rgd loop and synergy region," Cell 84, 155–164 (1996).

- C. Dotti, C. Sullivan, and G. Banker, "The establishment of polarity by hippocampal neurons in culture," J. Neurosci. 8, 1454–1468 (1988).
- Y. Nam, D. W. Branch, and B. C. Wheeler, "Epoxy-silane linking of biomolecules is simple and effective for patterning neuronal cultures," Biosens. Bioelectron. 22, 589–597 (2006).
- D. F. Bray, J. Bagu, and P. Koegler, "Comparison of hexamethyldisilazane (hmds), peldri ii, and critical-point drying methods for scanning electron microscopy of biological specimens," Microsc. Res. Tech. 26, 489–495 (1993).
- K. Karrai, G. Kolb, G. Abstreiter, and A. Schmeller, "Optical near-field induced current microscopy," Ultramicroscopy 61, 299–304 (1995).
- A. Drezet, J. C. Woehl, and S. Huant, "Diffraction by a small aperture in conical geometry: Application to metal-coated tips used in near-field scanning optical microscopy," Phys. Rev. E 65, 046611 (2002).
- L. Abad, M. Petit, G. Bugnicourt, T. Crozes, T. Fournier, and C. Villard, "Neurofets : Field effect nano-transistors fabrication for neural recording," In: Stett A (ed). Proceedings MEA Meeting 2010. Stuttgart: BIOPRO Baden-Württemberg GmbH 2010; 342-343. (2010).
- 21. G. Kim, "A mechanical spike accompanies the action potential in mammalian

nerve terminals," Biophys J $\mathbf{92},$ 3122–3129 (2007).