### A Better Method for Volume Determination of Regularly and Irregularly Shaped Nanoparticles with Enhanced Accuracy

#### Ravi Kiran Attota

Nanoscale Metrology Group (NSMG) Engineering Physics Division (EPD) Physical Measurement Laboratory (PML) National Institute of Standards and Technology (NIST) Gaithersburg, MD 20899, USA <u>Ravikiran.attota@nist.gov</u>; Ph: +1 301 975 5750

Nanoparticles (NPs) are widely used in diverse application areas, such as medicine, engineering, and cosmetics. The size (or volume) of NPs is one of the most important parameters for their successful application. It is relatively straightforward to determine the volume of regular NPs such as spheres and cubes from a one-dimensional or two-dimensional measurement. However, due to the three-dimensional nature of NPs, it is challenging to determine the proper physical size of many types of regularly and irregularly-shaped NPs (IS-NPs) at high-throughput using a single tool. Here, we present a relatively simple method that statistically determines a better volume estimate of many types of NPs by combining measurements from their top-down projection areas and peak-heights using two tools. The proposed method is significantly faster and more economical than the electron tomography method. We demonstrate the improved accuracy of the combined method over scanning electron microscopy (SEM) and atomic force microscopy (AFM) by using both modeling and measurements. This study also shows that SEM provides a more accurate estimate of size than AFM for most IS-NP size measurements. The method provides a much needed, proper high-throughput volumetric measurement method useful for many applications.

#### Introduction

The projected market figure for nanotechnology incorporated in manufactured goods by 2020 is approximately 3,000 billion US dollars worldwide <sup>1,2</sup>. Nanomaterials in the form of nanoparticles (NPs) represent an important part of nanotechnology. NP usage is increasing at a fast pace with the development of new types of nanomaterials. It is reflected by the existence of over 1,300 diverse listed consumer goods <sup>3</sup> (e.g., sunscreens, cosmetics, clothes, foods, drugs, paints, varnishes, and self-cleaning coating for floors, walls and windows) in addition to medical <sup>4</sup> and high technology <sup>5</sup> applications. It has been reported that nanotechnology (through the use of NPs) has the potential to radically change the way cancer is treated <sup>6</sup>. Engineered NPs are a subset of nanomaterials, of which irregularly-shaped NPs (IS-NPs) are a major component. The environmental health and safety (EHS) of engineered NPs have attracted much attention owing to their potential toxicity. This is reflected by the amount of research (over 10,000 publications) that has been conducted on the subject of the environmental and health effects of nanomaterials in the last 15 years <sup>7,8</sup>.

Several publications have highlighted the relationship between the size and toxicity of NPs <sup>1,6-14</sup>. They have also highlighted the lack of proper standards and characterization tools for NPs <sup>1,6-14</sup>. They have also highlighted the lack of proper standards and characterization tools for NPs <sup>1,6-17</sup>. In fact, some have also complained that, "After over a decade of research, answers to the most basic questions are still lacking," and suggest that more coherence in the experimental methods and materials used is needed <sup>7,10</sup>. One of the pressing issues with characterization is the large size discrepancy when the same set of NPs is measured using different tools <sup>1,18-23</sup>. Some publications have described variations in reported size even when NPs were measured using the same types of tool <sup>1,19,20</sup>. This creates confusion, as one of the fundamental starting points for application of NPs is an accurate knowledge of their size.

Electron tomography is a method that has demonstrated reliable 3D-shape reconstructions of nanoparticles<sup>24-26</sup>. But even though significant improvement has been made recently <sup>27</sup>, the technique is impractical and slow for routine and large-scale measurements required to statistically determine an accurate size and distribution of a batch of monodisperse NPs.

Apart from the electron tomography, at present, it is nearly impossible to accurately determine the size of many types of NPs(including types of IS-NPs<sup>1,3</sup>). For this reason, only reference materials (a less rigorous measurement certification process) exist for IS-NPs, such as gold <sup>28,29</sup>. However, for some nearly-spherical nanoparticles such as polystyrene, standard reference materials (a highly rigorous and accurate measurement process) exist <sup>30</sup>. Although several tools are currently available, more tools and methods are being sought <sup>31,32</sup>, mainly due to a lack of consensus in the reported sizes. One of the latest additions to this long list of tools is through-focus scanning optical microscopy (TSOM)<sup>33,34</sup>.

At this juncture, many critical questions remain: Is it possible to accurately measure the true volume of NPs, especially for IS-NPs and at high-throughput? If it is, which tool or method provides the correct spherical-volume-equivalent mean diameter at high throughput? If none, is there any way to determine which tool provides the smallest deviation from the accurate diameter? Is there any way to determine the deviation from the correct mean diameter so we can estimate the errors involved in the measured values? In this paper, we attempt to provide answers to these important fundamental questions.

#### **Background information**

First, we briefly discuss the measurement procedures used to determine the reported diameters of NPs by the most widely used direct imaging tools, such as SEM, and AFM. These methods are usually considered to be the most accurate reference metrology tools, specifically when they are calibrated to an International System (SI) unit <sup>18</sup>. At the sub-100 nm scale, it is challenging to obtain complete 3-D information about IS-NPs. For this reason, certain assumptions are typically made when using these tools.

A two-dimensional, top-down area of projection is used to measure the diameters of IS-NPs using SEM (and also using conventional transmission electron microscopy (TEM)). Information

regarding the height of the NPs is ignored. The calculated projected area equivalent diameter (Dxy) is reported as the SEM measured diameter (it is assumed that the height (or the 3<sup>rd</sup> dimension) is the same as the Dxy). In principle, both SEM and conventional TEM are supposed to produce the same size. However, in practice, SEM and TEM tools produce slightly different diameters due to differences in the edge detection methods.

Accurately calibrated AFM uses a sharp tip to measure the precise peak heights of NPs. However, it is extremely challenging to obtain an accurate contour (or area of projection) of NPs because of limitations in the physical interaction between the AFM tip and NPs. For this reason, only height information is used for the reported diameter and can be referred to as the peak-height equivalent diameter, which is the same as Dz. It is assumed that the AFM-measured height represents the diameter of perfectly spherical NPs and is reported as such.

In the case of SEM (and also for conventional TEM), only the 2-D projection area is exploited, whereas in the case of AFM, only the 1-D height information is exploited. Because none of the three tools use complete 3-D information on the IS-NPs, which are inherently 3-D in nature, they are prone to error. In fact, we should expect an erroneous diameter measurement and systematic difference in the reported diameters when SEM/TEM and AFM measurements are compared to each other. The published literature provides ample support for this observation (Table 1). We can also expect the reported diameters to include a certain deviation from the not yet known correct diameters. However, with the prevailing conventional knowledge, it is challenging to routinely determine the true diameter/volume of many types of NPs.

Material	Mean d	iameters	s, nm	SEM Dia.	SEM Vol.	Ref.
	SEM	TEM	AFM	Devi. (%)	Devi. (%)	
SiO2	39.00	35.10	30.30	28.71	113.24	[1]
SiO2	46.60	42.90	36.20	28.73	113.32	[1]
SiO2	89.80	86.30	80.20	11.97	40.38	[1]
Au	85.48		79.05	8.13	26.44	[2]
Ag	30.97		23.14	33.84	139.74	[11]
Au	9.90		7.20	37.50	159.96	[16,17]
Au	26.90		23.70	13.50	46.22	[16,17]
Au	54.90		53.90	1.86	5.67	[16,17]

**Table 1. Nanoparticle diameters reported in the literature.** A comparison of the reported mean diameters measured from the same batch using SEM, TEM and AFM shows disagreement in their values. SEM-reported diameters are always larger than the AFM-reported diameters. Diameters reported by SEM deviate upto 37 % from the diameters reported by AFM. Volumes calculated by SEM reported diameters deviate upto 150% from the volumes calculated by AFM reported diameters.

Ideally, we would like to have all types of NPs to be perfectly regular in nature such as spherical, cubic, or cylindrical. However, in reality, many types of NPs are irregular as shown in Fig. 1(a). In this study, we model these types of IS-NPs as ellipsoids (Figs. 1(d) to (d)).



**Figure 1. Nanoparticle shapes.** (a) SEM image of nominally 30 nm diameter Au nanoparticles showing their irregular shapes. (b) to (d) Simplified ellipsoid shapes assigned to IS-NPs. (b') to (d'), Schematics of the corresponding areas of projections. Dx, Dy, and Dz are the diameters in the directions of the X, Y and Z axes, respectively.

Conventionally, the IS-NPs are first spread onto a substrate (or a grid) before measuring with SEM, TEM or AFM tools. This process typically aligns the IS-NPs in the most stable orientations on the substrates such that they have a wide base and a low center of gravity (or in this case due to van der Waals'/electrostatic force)<sup>28</sup>. Because of this, it is highly unlikely that the IS-NPs will come to rest in a prolate ellipsoid orientation (Fig. 1(d)), and hence we ignore this orientation in the following discussion.

#### **Results and discussion**

Now, let us consider a 30 nm diameter spherical nanoparticle. This can be distorted into several ellipsoidal shapes while keeping the volume constant. The equivalent diameters that would be inferred for the different-shaped nanoparticles using SEM/TEM (Dxy) and AFM (Dz) are then

evaluated as shown in Figures 2 (a) to (d) (the numerical values are presented in Table 2). Comparing the inferred Dxy with Dz, we can observe that the SEM measured diameters are always larger than the AFM measured diameters (except for a perfectly spherical particle). The same trend can be observed in the published literature (Table 1), indicating that the model we assume so far points in the right direction.



**Figure 2 Equivalent diameters calculation.** (a) A representative nanoparticle. This could be of any irregularly-shaped, three-dimensional nanoparticle. (b) Top-down area of projection as measured by tools such as SEM and TEM. For an IS-NP this would be an irregular area of projection. (c) Peak-height (Dz) as measured by tools such as AFM, and would be the reported diameter. (d) Calculated projection area equivalent diameter (Dxy) using the measured projection area from (b). Dxy would be the reported diameter using SEM or TEM. (e) Constructed ellipsoid by combining SEM/TEM and AFM reported diameters, where Dx = Dy = Dxy, and Dz = peak-height. (f) Calculated spherical equivalent volume diameter (Dxyz) from the constructed ellipsoid, where Dx = Dy = Dzy.

Here we propose a third method to calculate the volume of some types of NPs, especially that of the type of IS-NPs shown in Fig. 1(a). In this method, top-down projection area as measured by SEM/TEM, and peak-height as measured by AFM are combined to obtain a 3D volume from which spherical equivalent volume diameter (Dxyz) is calculated as depicted in Figs. 2(e) and 2(f). Comparing the calculated Dz, Dxy and Dxyz in Table 2, we can see that the diameters determined by SEM/TEM are always larger while the diameters determined by AFM are always smaller than the correct diameter of 30 nm (except for the perfect spherical particle), which exactly equals Dxyz. We can also infer another minor point: the difference in the SEM/TEM and AFM measurements is proportional to the deviation of a NP from the perfect spherical shape.

Dx	Dy	Dz	SEM	AFM	Dxyz	SEM	AFM	Shape ratios			
			Dia.	Dia.		(Dxy)	(Dz)				
			(Dxy)	(Dz)		Dev.	Dev.	R1	R2	R3	R4
nm	nm	nm	nm	nm	nm	%	%	Dx/Dy	Dxy/Dz	R2/R1	
30.00	30.00	30.00	30.00	30.00	30.00	0.00	0.00	1.000	1.000	1.000	1.000
31.00	29.50	29.52	30.24	29.52	30.00	0.80	-1.59	1.051	1.024	0.975	-1.976
36.00	27.00	27.78	31.18	27.78	30.00	3.92	-7.41	1.333	1.122	0.842	-1.888
50.00	20.00	27.00	31.62	27.00	30.00	5.41	-10.00	2.500	1.171	0.468	-1.849
31.00	31.00	28.10	31.00	28.10	30.00	3.33	-6.35	1.000	1.103	1.103	-1.904
35.00	35.00	22.04	35.00	22.04	30.00	16.67	-26.53	1.000	1.588	1.588	-1.592
37.00	32.00	22.80	34.41	22.80	30.00	14.70	-23.99	1.156	1.509	1.305	-1.632

**Table 2. Inferred equivalent diameter evaluation.** A 30 nm diameter nanoparticle was distorted by keeping the volume constant to produce different Dx, Dy, and Dz diameters. Inferred equivalent diameters of the distorted nanoparticles were then calculated (Figure 2) as if they were measured using SEM/TEM (using the projection area), AFM (using the peak-height), and the combined method (using projection area and peak-height) proposed here. Dz = diameter in the Z-direction as measured by AFM, Dxy = projection area equivalent diameter as measured by SEM/TEM, Dxyz = spherical equivalent volume diameter determined by combining Dxy and Dz. If the shape ratio R3 is less than one, then the nanoparticles are dominated more by elongation than by flattening. If the shape ratio R3 is more than one, then the nanoparticles are dominated more by flattening than by elongation. R4 = (DzDev./Dxydev) or (AFMDev./SEMDev.)

In this section, we employ a similar analysis concept, but use experimentally measured realistic shape variations for Au nanoparticles, as an example. The shape variations thus obtained are then applied to 10 nm, 30 nm, 60 nm and 100 nm diameter nanoparticles for which the diameters as measured by SEM/TEM and AFM were calculated. Detailed results are presented in Tables S1(a) to S1(d). A summary of the results is presented in Table 3, from which the following important inferences can be made.

- I. It appears possible to accurately obtain the spherical equivalent volume diameter of IS-NPs by combining SEM/TEM (projection-area) and AFM (peak-height) measurements.
- II. The SEM/TEM-measured mean diameters are always larger than the AFM-measured mean diameters.
- III. The SEM/TEM-measured diameters are larger, while the AFM-measured diameters are smaller than the correct equivalent diameters.

IV. The SEM/TEM measured diameters are more accurate than the AFM-measured diameters.

Nevertheless, the data presented in Tables 2 and 3 using partly modeled and partly experimental analysis cannot fully prove the accuracy of the proposed method; validation using experimental data is needed. Figure 15 in Ref. <sup>3</sup> provides this experimental data where for the first time, direct SEM and AFM measurements were provided for the same set of nanoparticles; suitable for the current type of analysis. The extracted SEM- and AFM-measured diameters and the calculated spherical equivalent volume diameters for each pair are tabulated in Supplementary Table S2. The average values obtained are presented in Table 4. AFM measured diameter is 5.48 nm smaller, while SEM measured diameter is 2.96 nm larger than the Dxyz.

Dxyz	SEM	AFM	SEM	AFM	Ratio
	(Dxy)	(Dz)	Dev.	Dev.	
nm	nm	nm	%	%	R4
10.00	10.37	9.29	3.75	-7.06	-1.884
30.00	31.13	27.89	3.75	-7.03	-1.874
60.00	61.62	56.93	2.69	-5.12	-1.901
100.00	102.38	95.44	2.38	-4.56	-1.917

Table 3. Calculated spherical equivalent volume diameters (Dxyz) using realistic shape variations of IS-NPs made of Au. A summary of statistically derived Dxyz by combining the calculated SEM/TEM (Dxy) and AFM (Dz) diameters (from Supplementary Table S1) after applying measured shape variations of Au NPs to the four diameter sizes. SEM Dev. and AFM Dev. are the percentage deviations from the inferred correct diameters (Dxyz).

Method	SEM (Dxy)	AFM (Dz)	Dxyz	SEM Dev.	AFM Dev.	Ratio
	nm	nm	nm	%	%	R4
Direct	58.17	49.73	55.21	5.36	-9.92	-1.85
Adjusted	55.00	49.73	53.18	3.41	-6.49	-1.90

**Table 4. A summary of the SEM and AFM data presented in Supplementary Table S2.** Dxyz in this table was calculated by combining the reported SEM and AFM diameters. Deviations of the reported SEM and AFM diameters from Dxyz were evaluated as a percentage. In the 'Adjusted' row, the SEM value was floated by maintaining a fixed AFM value such that the 'Ratio' magnitude equals 1.90.

The above shown experimental data are strikingly similar to those in Table 3. Identical inferences can be made using both the tables. The deviation ratios are also similar. This experimental data supports the modeled approach followed earlier in the paper and hence appears to confirm all the important inferences made above. However, we have not yet provided independent evidence for the accuracy of the proposed Dxyz. We provide this proof below.

Here we apply the combined method to evaluate the spherical equivalent volume diameter to irregularly-shaped macroparticles, such as pebbles. A key benefit of this approach is that the spherical equivalent volume diameter obtained using the projection area and the height can be corroborated by measuring the actual volume. This provides a definitive validation of the combined method presented in this paper. To perform this, we selected irregularly shaped aquarium glass pebbles weighing approximately 15 g to 21 g (Fig. 3). After following the procedure described above, we obtained a mean Dxyz of 23.8 mm using the projection area and the height of the 54 pebbles (Supplementary Table S3). The independently calculated Dxyz using the measured volume (obtained from mass and density) is 23.96 mm. As we can see, both the diameter values are statistically the same (from the statistical analysis given in Table S3) and hence confirm the accuracy of the combined method presented earlier. To obtain the true mean size, statistically-sufficient large number of particles must be measured. Of course, the particle size could range from nanoparticles to macroparticles. However, the 'Ratio' seems to depend on the nature and extent of the particle irregularities. In the case of gold NPs, the 'Ratio' varied from 1.88 to 1.92, depending on the size. However, in the case of pebbles, the measured 'Ratio' is approximately 1.63.

Since we have confirmed the modeled approach of the combined measurement presented earlier (Tables 2 and 3) with the experimental values (Supplementary Tables S2 and S3), we can now say with confidence that the most accurate Dxyz from the Au nanoparticle data presented in Ref. <sup>3</sup> is 55.21 nm (Supplementary Table S2). The SEM and AFM measured diameters are approximately 5.36 % and 9.92 % larger and smaller, respectively, than the Dxyz.



**Figure 3.** Macroscale pebbles used to compare the spherical equivalent volume diameters calculated based on Dxy and Dz; and from mass and density.

There is some degree of uncertainty associated with the SEM data presented in Ref. <sup>3</sup>. This stems from the fact that no accurate edge detection method has been developed for SEM images of Au nanoparticles. Because of this, the reported diameters of Au NPs were extracted at the full-width-half-max of the SEM intensity profiles <sup>3</sup>. This selection is a convenience based on the general nature of SEM imaging, but is not based on a rigorous model for imaging of Au NPs. Hence, we can expect an unknown measurement bias in the SEM measured diameters. However, no such wiggle room exists for the AFM measured diameters (peak-heights) and hence these can be considered fixed. Based on this, we can adjust the value of the SEM measured diameter such that the deviation ratio equals approximately 1.90 (from Table 3, based on the similarity in the size of particles to a 60 nm diameter and assuming similar shape irregularities). When the experimental data were adjusted to this ratio by keeping the AFM value constant and floating the SEM value, we obtain a new adjusted Dxyz of 53.18 nm (Table 4). The new adjusted SEM value is approximately 3.17 nm smaller than the reported value.

We can find some validation for this type of analysis from the published literature; however, we suggest caution. TEM measurements provide more accurate edge detection than SEM measurements, as it is possible to identify the end of atoms at the edges in TEM images more precisely than any non-model-based assignment of the physical edge location within the intensity profile of the SEM image. Because of this, we can expect a more accurate diameter measurement with TEM than with SEM when the measurements are performed on the same batch of NPs. The data presented in Table 1 show the measured diameters using SEM and TEM (in addition to AFM) for the same set of NPs. In these data, we can clearly see that the diameters measured using TEM are smaller than those measured using SEM, with a magnitude similar to the SEM adjusted value used in Table 4.

If the deviation in the diameter for a batch of monodispersed (in spherical equivalent volume diameters —or equivalently, volume) IS-NPs was known with certainty, then the measured diameters could be 'adjusted' by the known deviation to obtain the true diameter of Dxyz. For example, in Table 3 the AFM measured diameter (for 30 nm size) is 6.91 % smaller than the Dxyz. By adjusting the AFM measured diameter with this known deviation one could, in principle, obtain the Dxyz by using only one tool. This paper provides a basis for this type of adjustment.

For regularly-shaped, symmetric NPs such as spheres and cubes, the proposed combined method may not be of much use, as either SEM or AFM can provide a one dimensional measurement which will be same in other directions also. For perfect cylindrical-shaped NPs (with the axis along the substrate), SEM alone can provide all the necessary dimensions to calculate the volume; however, AFM cannot provide both necessary dimensions. For some other types of regularly-shaped NPs such as prisms, cuboids, and cylinders (resting on the circular end), the height information cannot be automatically obtained using SEM measurements alone, and lateral dimensions cannot be obtained using AFM alone. In such cases, again, combining SEM and AFM measurements will provide a more accurate volume determination compared to using either SEM or AFM alone; however, for these regular shapes, it is unnecessary and in fact preferable

not to assume a spheroidal shape. Because of their regular shapes, basic geometric formulas can be used to determine their volumes from the two or three dimensions measured using SEM and/or AFM.

#### Conclusions

In this paper, we have described a method to determine the volume of certain types of irregularly shaped, and regularly-shaped NPs by combining measurements obtained using top-down projection areas (as measured by tools such as SEM/TEM), and peak-heights (as measured by tools such as AFM) with more accuracy than can be measured by using either SEM or AFM alone. Both modeled and experimental evidence are provided to validate this concept. Based on the more accurate volumetric measurement of certain types of IS-NPs that this method provides, we demonstrate that for these IS-NPs, SEM/TEM measurements are more accuracy. We expect this method to facilitate the certification of accurate reference materials for IS-NPs, which will increase the reliability and enhance trust in the results of studies of size-dependent properties of NPs, such as in toxicology.

#### Acknowledgements

The author would like to thank John Kramar for the useful discussions, and Andras Vladar for providing high-quality SEM images of Au nanoparticles used in Ref. <sup>28</sup>.

#### **Competing financial interests**

The author declares no competing financial interests.

#### References

- 1 Hole, P. *et al.* Interlaboratory comparison of size measurements on nanoparticles using nanoparticle tracking analysis (NTA). *J Nanopart Res* **15**, doi:10.1007/s11051-013-2101-8 (2013).
- 2 Roco, M. C. The long view of nanotechnology development: the National Nanotechnology Initiative at 10 years. *J Nanopart Res* **13**, 427-445, doi:10.1007/s11051-010-0192-z (2011).
- 3 Delvallee, A., Feltin, N., Ducourtieux, S., Trabelsi, M. & Hochepied, J. F. Direct comparison of AFM and SEM measurements on the same set of nanoparticles. *Meas Sci Technol* **26**, doi:10.1088/0957-0233/26/8/085601 (2015).
- 4 Oberdorster, G., Oberdorster, E. & Oberdorster, J. Nanotoxicology: An emerging discipline evolving from studies of ultrafine particles. *Environ Health Persp* **113**, 823-839, doi:10.1289/ehp.7339 (2005).

- 5 Ko, S. H. *et al.* High-speed, high-purity separation of gold nanoparticle-DNA origami constructs using centrifugation. *Soft Matter* **10**, 7370-7378, doi:10.1039/c4sm01071j (2014).
- 6 Crist, R. M. *et al.* Common pitfalls in nanotechnology: lessons learned from NCI's Nanotechnology Characterization Laboratory. *Integr Biol-Uk* **5**, 66-73, doi:10.1039/c2ib20117h (2013).
- Fadeel, B., Fornara, A., Toprak, M. S. & Bhattacharya, K. Keeping it real: The importance of material characterization in nanotoxicology. *Biochem Biophys Res Commun*, doi:10.1016/j.bbrc.2015.06.178 (2015).
- 8 Krug, H. F. Nanosafety Research-Are We on the Right Track? *Angew Chem Int Edit* **53**, 12304-12319, doi:10.1002/anie.201403367 (2014).
- Hussain, S. M. *et al.* Toxicity Evaluation for Safe Use of Nanomaterials: Recent Achievements and Technical Challenges. *Adv Mater* **21**, 1549-1559, doi:10.1002/adma.200801395 (2009).
- 10 Schrurs, F. & Lison, D. Focusing the research efforts. *Nat Nanotechnol* **7**, 546-548, doi:10.1038/nnano.2012.148 (2012).
- 11 Napierska, D. *et al.* Size-Dependent Cytotoxicity of Monodisperse Silica Nanoparticles in Human Endothelial Cells. *Small* **5**, 846-853, doi:10.1002/smll.200800461 (2009).
- 12 Teeguarden, J. G., Hinderliter, P. M., Orr, G., Thrall, B. D. & Pounds, J. G. Particokinetics in vitro: Dosimetry considerations for in vitro nanoparticle toxicity assessments. *Toxicol Sci* **95**, 300-312, doi:10.1093/toxsci/kfl165 (2007).
- 13 Cabral, H. *et al.* Accumulation of sub-100 nm polymeric micelles in poorly permeable tumours depends on size. *Nature Nanotechnology* **6**, 815-823, doi:10.1038/Nnano.2011.166 (2011).
- 14 Blanco, E., Shen, H. & Ferrari, M. Principles of nanoparticle design for overcoming biological barriers to drug delivery. *Nat Biotechnol* **33**, 941-951, doi:10.1038/nbt.3330 (2015).
- 15 Linsinger, T. P. J., Roebben, G., Solans, C. & Ramsch, R. Reference materials for measuring the size of nanoparticles. *Trac-Trend Anal Chem* **30**, 18-27, doi:10.1016/j.trac.2010.09.005 (2011).
- 16 Reich, E. S. Nano rules fall foul of data gap. *Nature* **480**, 160-161 (2011).
- 17 Join the dialogue. *Nature Nanotechnology* **7**, 545-545 (2012).
- 18 Meli, F. *et al.* Traceable size determination of nanoparticles, a comparison among European metrology institutes. *Meas Sci Technol* **23**, doi:10.1088/0957-0233/23/12/125005 (2012).
- 19 Wang, C. Y., Fu, W. E., Lin, H. L. & Peng, G. S. Preliminary study on nanoparticle sizes under the APEC technology cooperative framework. *Meas Sci Technol* **18**, 487-495, doi:10.1088/0957-0233/18/2/S23 (2007).
- Motzkus, C. *et al.* Size characterization of airborne SiO2 nanoparticles with on-line and off-line measurement techniques: an interlaboratory comparison study. *J Nanopart Res* 15, doi:10.1007/s11051-013-1919-4 (2013).

- 21 MacCuspie, R. I. *et al.* Challenges for physical characterization of silver nanoparticles under pristine and environmentally relevant conditions. *Journal of Environmental Monitoring* **13**, 1212-1226, doi:10.1039/c1em10024f (2011).
- 22 Coleman, V. A., Jamting, A. K., Catchpoole, H. J., Roy, M. & Herrmann, J. Nanoparticles and metrology: a comparison of methods for the determination of particle size distributions. *Instrumentation, Metrology, and Standards for Nanomanufacturing, Optics, and Semiconductors V* **8105**, doi:10.1117/12.894297 (2011).
- 23 Cho, E. J. *et al.* Nanoparticle Characterization: State of the Art, Challenges, and Emerging Technologies. *Mol Pharmaceut* **10**, 2093-2110, doi:10.1021/mp300697h (2013).
- 24 Thomas, J. M., Midgley, P. A., Ducati, C. & Leary, R. K. Nanoscale electron tomography and atomic scale high-resolution electron microscopy of nanoparticles and nanoclusters: A short surveyNanoscale electron tomography and atomic scale high-resolution electron microscopy of nanoparticles and nanoclusters: A short surveyretain--&gt. *Progress in Natural Science: Materials International* **23**, 222-234, doi:10.1016/j.pnsc.2013.04.003 (2013).
- 25 Midgley, P. A., Ward, E. P. W., Hungria, A. B. & Thomas, J. M. Nanotomography in the chemical, biological and materials sciences. *Chemical Society Reviews* **36**, 1477-1494, doi:10.1039/B701569K (2007).
- 26 Hayashida, M., Kumagai, K. & Malac, M. Three dimensional accurate morphology measurements of polystyrene standard particles on silicon substrate by electron tomography. *Micron* **79**, 53-58, doi:10.1016/j.micron.2015.08.003 (2015).
- 27 Migunov, V. *et al.* Rapid low dose electron tomography using a direct electron detection camera. *Sci Rep-Uk* **5**, 14516, doi:10.1038/srep14516 (2015).
- 28 Vladar, A. *NIST NCL Joint Assay Protocol, PCC-15* (National Institute of Standards and Technology, DOC, 2011).
- 29 Grobelny, J. *et al.* Size measurement of nanoparticles using atomic force microscopy. *Methods in molecular biology (Clifton, N.J.)* **697**, 71-82, doi:10.1007/978-1-60327-198-1\_7 (2011).
- 30 Dagata, J. SRM 1964 Nominal 60 nm Diameter Polystyrene Spheres (2014) https://www-s.nist.gov/srmors/view\_detail.cfm?srm=1964 (Date of access:25/02/2016).
- 31 Malysheva, A., Lombi, E. & Voelcker, N. H. Bridging the divide between human and environmental nanotoxicology. *Nature Nanotechnology* **10**, 835-844, doi:10.1038/Nnano.2015.224 (2015).
- 32 Montano, M. D., Lowry, G. V., von der Kammer, F., Blue, J. & Ranville, J. F. Current status and future direction for examining engineered nanoparticles in natural systems. *Environ Chem* **11**, 351-366, doi:10.1071/En14037 (2014).
- 33 Attota, R., Kavuri, P. P., Kang, H., Kasica, R. & Chen, L. Nanoparticle size determination using optical microscopes. *Appl Phys Lett* **105**, doi:10.1063/1.4900484 (2014).
- 34 Kang, H., Attota, R., Tondare, V., Vladár, A. E. & Kavuri, P. A method to determine the number of nanoparticles in a cluster using conventional optical microscopes. *Appl Phys Lett* **107**, 103106, doi:10.1063/1.4930994 (2015).

### Supplementary information for A Better Method for Volume Determination of Regularly and Irregularly Shaped Nanoparticles with Enhanced Accuracy

Ravi Kiran Attota NSMG, EPD, PML, NIST, Gaithersburg, MD 20899, USA Ravikiran.attota@nist.gov

# Table S1. Spherical equivalent volume diameters (Dxyz) evaluated using realistic shape variations from measurements of Au nanoparticles.

In Table 2, Dx and Dy values were randomly selected to demonstrate that combining SEM/TEM and AFM measurements will result in a correct Dxyz. Here we obtain realistic variations in Dx and Dy from the measured projection areas using actual data sets. For this we analyzed high-quality SEM images of nominally 10 nm, 30 nm, 60 nm and 100 nm diameter, irregularly shaped Au nanoparticles (RM certified by NIST [16,17]) using ImageJ to obtain areas of projection (Aproj), Dx (Feret maximum diameter), and Dy (Feret minimum diameter) for isolated nanoparticles. The Dx and Dy thus obtained are assumed to be the major and the minor axes, respectively, of the ellipsoid-shaped projected areas, even though they are not always exactly perpendicular to each other. With these assumptions we obtain typical, realistic shape variations from the projected areas. The measured shape variations thus obtained are then applied to (a) 10 nm, (b) 30 nm, (c) 60 nm, and (d) 100 nm diameter nanoparticles, and the diameters as measured by SEM/TEM and AFM are calculated.

Aproj = measured area of projection; Dx = measured maximum Feret diameter; Dy = measured minimum Feret diameter; Dxy = calculated diameter of a circle that has the same Aproj;  $\Delta$ Dx = percentage deviation in Dx from Dxy;  $\Delta$ Dy = percentage deviation in Dy from Dxy; Db = base diameter; Since this is a modeled calculation, Dxyz must be equal to Db.

By first evaluating  $\Delta Dx$  and  $\Delta Dy$  and then applying them to modeled sizes (Db), size variations (typically Gaussian distribution) present in the measured data (in addition to shape/area variations) within a size group was normalized.

Table	S1 a
-------	------

	Measured data for nominally 10 nm diameter						Measure	d ∆Dx & /	Dy applie	ed for Db =	100 nm	S	nape ratio	os
S.No.	Aproj	Dx	Dy	Dxy	ΔDx	∆Dy	Dx	Dy	Dz(AFM)	Dxy(SEM)	Dxyz	R1	R2	R3
	nm <sup>2</sup>	nm	nm	nm	%	%	nm	nm	nm	nm	nm	Dx/Dy	Dxy/Dz	R2/R1
1	138.52	14.89	12.77	13.28	12.13	-3.86	11.21	9.61	9.28	10.38	10.00	1.166	1.119	0.960
2	126.71	14.34	12.13	12.70	12.95	-4.50	11.30	9.55	9.27	10.39	10.00	1.183	1.120	0.947
3	109.19	12.93	11.42	11.79	9.64	-3.12	10.96	9.69	9.41	10.31	10.00	1.132	1.095	0.967
4	105.52	13.02	10.85	11.59	12.34	-6.37	11.23	9.36	9.51	10.26	10.00	1.200	1.079	0.899
5	114.89	13.33	11.86	12.09	10.22	-1.89	11.02	9.81	9.25	10.40	10.00	1.123	1.124	1.001
6	123.04	14.70	11.49	12.51	17.43	-8.19	11.74	9.18	9.28	10.38	10.00	1.279	1.119	0.875
7	123.04	13.65	12.64	12.51	9.04	0.99	10.90	10.10	9.08	10.49	10.00	1.080	1.156	1.070
8	105.12	12.93	10.85	11.57	11.75	-6.19	11.17	9.38	9.54	10.24	10.00	1.191	1.073	0.901
9	96.15	12.93	10.21	11.06	16.84	-7.68	11.68	9.23	9.27	10.39	10.00	1.266	1.120	0.885
10	96.97	12.33	10.85	11.11	10.97	-2.32	11.10	9.77	9.23	10.41	10.00	1.136	1.129	0.993
11	112.04	12.93	11.49	11.94	8.24	-3.79	10.82	9.62	9.60	10.20	10.00	1.125	1.063	0.945
12	110.01	13.66	11.69	11.83	15.45	-1.23	11.54	9.88	8.77	10.68	10.00	1.169	1.218	1.042
13	122.64	13.57	12.13	12.49	8.62	-2.92	10.86	9.71	9.48	10.27	10.00	1.119	1.083	0.968
14	112.04	13.42	11.49	11.94	12.37	-3.79	11.24	9.62	9.25	10.40	10.00	1.168	1.124	0.962
15	131.60	14.01	12.77	12.94	8.29	-1.36	10.83	9.86	9.36	10.34	10.00	1.098	1.104	1.006
16	101.86	12.57	11.28	11.39	10.43	-0.89	11.04	9.91	9.14	10.46	10.00	1.114	1.145	1.028
17	112.45	14.00	11.27	11.96	17.02	-5.79	11.70	9.42	9.07	10.50	10.00	1.242	1.157	0.932
18	116.52	13.57	11.49	12.18	11.43	-5.66	11.14	9.43	9.51	10.25	10.00	1.181	1.078	0.913
							A	verage ->	9.29	10.37	10.00	1.165	1.117	0.961
							:	St.Dev>	0.20	0.12	0.00	0.056	0.038	0.055

#### Table S1 b.

	Measu	red data	for nom	ninally 30	) nm dia	meter	Measured ΔDx &		$\Delta Dy$ applied for Db = 30 nm		Shape ratios		S	
S.No.	Aproj	Dx	Dy	Dxy	∆Dx	∆Dy	Dx	Dy	Dz(AFM)	Dxy(SEM)	Dxyz	R1	R2	R3
	nm <sup>2</sup>	nm	nm	nm	%	%	nm	nm	nm	nm	nm	Dx/Dy	Dxy/Dz	R2/R1
1	478.23	27.00	24.54	24.67	9.45	-0.51	32.83	29.85	27.55	31.30	30.00	1.100	1.136	1.033
2	531.73	27.85	25.86	26.01	7.07	-0.59	32.12	29.82	28.18	30.95	30.00	1.077	1.098	1.020
3	541.77	29.85	25.25	26.26	13.69	-3.85	34.11	28.85	27.44	31.37	30.00	1.182	1.143	0.967
4	603.63	29.49	27.99	27.72	6.38	0.98	31.92	30.29	27.93	31.09	30.00	1.054	1.113	1.057
5	494.95	27.28	25.30	25.10	8.68	0.80	32.61	30.24	27.38	31.40	30.00	1.078	1.147	1.063
6	688.91	34.43	27.76	29.61	16.28	-6.24	34.88	28.13	27.52	31.33	30.00	1.240	1.138	0.918
7	551.80	29.32	26.52	26.50	10.63	0.06	33.19	30.02	27.10	31.56	30.00	1.106	1.165	1.053
8	595.27	32.33	24.57	27.52	17.45	-10.74	35.24	26.78	28.62	30.72	30.00	1.316	1.073	0.816
9	707.31	37.05	25.41	30.00	23.49	-15.31	37.05	25.41	28.69	30.68	30.00	1.458	1.070	0.733
10	499.96	28.36	21.98	25.23	12.43	-12.85	33.73	26.14	30.62	29.69	30.00	1.290	0.970	0.752
11	551.80	30.08	25.86	26.50	13.50	-2.41	34.05	29.28	27.09	31.57	30.00	1.163	1.166	1.002
12	/15.6/	32.41	28.45	30.18	/.3/	-5.74	32.21	28.28	29.64	30.18	30.00	1.139	1.018	0.894
13	5/1.86	32.33	23.72	26.98	19.83	-12.08	35.95	26.38	28.48	30.79	30.00	1.363	1.081	0.793
14	553.47	30.74	24.57	26.54	15.81	-7.43	34.74	27.77	27.98	31.06	30.00	1.251	1.110	0.887
15	580.22 702 E9	31.09	25.80	27.17	14.40	-4.83	34.32	28.55	27.55	31.30	30.00	1.202	1.130	0.945
17	192.30	33.31	29.74	31.70	0.02	-0.50	22.05	20.09	20.00	21.42	20.00	1.194	1.071	1.046
10	504 08	27.00	24.57	24.30	9.05	-0.00	22.95	29.90	27.55	21.45	30.00	1.099	1.130	1.040
10	/73 21	20.39	23.45	23.33	12.02	-10.42	33.55	26.87	20.08	30.05	30.00	1.115	1.192	0.804
20	496 62	27.45	23.28	25.14	11.02	-7 42	33.01	20.07	29.50	30.03	30.00	1 1 9 9	1.005	0.869
20	469.87	27.51	23.20	23.14	13 17	-2.05	33.95	29.38	27.07	31 58	30.00	1 1 1 5 5	1.042	1 010
21	511 67	28.39	23.55	25.52	11 25	-3 72	33.33	29.30	27.07	31.05	30.00	1 1 1 5 5	1.107	0.959
23	700.62	34 14	28.45	29.86	14 33	-4.73	34 30	28.58	27.54	31.31	30.00	1,200	1 137	0.935
24	570.19	29.77	25.86	26.94	10.51	-4.00	33.15	28.80	28.28	30.90	30.00	1.151	1.093	0.949
25	570.19	29.54	25.86	26.94	9.67	-4.00	32.90	28.80	28.49	30.78	30.00	1.142	1.080	0.946
26	449.80	27.19	23.28	23.93	13.62	-2.72	34.09	29.18	27.14	31.54	30.00	1.168	1.162	0.995
27	489.93	28.48	24.54	24.97	14.04	-1.71	34.21	29.49	26.76	31.76	30.00	1.160	1.187	1.023
28	598.62	33.94	22.50	27.60	22.97	-18.50	36.89	24.45	29.93	30.03	30.00	1.509	1.003	0.665
29	603.63	32.41	25.86	27.72	16.91	-6.69	35.07	27.99	27.50	31.33	30.00	1.253	1.139	0.909
30	580.22	31.49	25.77	27.17	15.88	-5.18	34.76	28.44	27.31	31.45	30.00	1.222	1.152	0.942
31	667.17	32.12	28.35	29.14	10.23	-2.71	33.07	29.19	27.97	31.07	30.00	1.133	1.111	0.980
32	596.95	31.57	26.55	27.56	14.53	-3.69	34.36	28.89	27.20	31.51	30.00	1.189	1.159	0.974
33	478.23	27.67	23.77	24.67	12.17	-3.64	33.65	28.91	27.75	31.19	30.00	1.164	1.124	0.965
34	548.45	29.77	25.86	26.42	12.67	-2.11	33.80	29.37	27.20	31.51	30.00	1.151	1.158	1.006
35	630.39	31.94	27.16	28.33	12.75	-4.13	33.83	28.76	27.75	31.19	30.00	1.176	1.124	0.956
36	573.54	29.52	25.86	27.02	9.25	-4.28	32.77	28.72	28.69	30.68	30.00	1.141	1.069	0.937
37	601.96	30.19	27.16	27.68	9.06	-1.89	32.72	29.43	28.04	31.03	30.00	1.112	1.107	0.996
38	513.34	27.91	25.60	25.56	9.21	0.16	32.76	30.05	27.43	31.38	30.00	1.090	1.144	1.049
39	553.47	29.49	25.86	26.54	11.10	-2.56	33.33	29.23	27.71	31.21	30.00	1.140	1.126	0.988
40	499.96	28.36	24.57	25.23	12.43	-2.60	33.73	29.22	27.40	31.39	30.00	1.154	1.146	0.993
41	652.13	31.46	28.45	28.81	9.21	-1.25	32.76	29.62	27.82	31.15	30.00	1.106	1.120	1.013
42	551.80	29.06	26.50	26.50	9.65	0.00	32.90	30.00	27.36	31.42	30.00	1.096	1.148	1.047
43	635.40	32.74	28.45	28.44	15.13	0.04	34.54	30.01	26.05	32.19	30.00	1.151	1.230	1.074
44	015.34 E21 72	30.30	27.10	27.98	0.27	-2.97	32.48	29.11	28.30	30.75	30.00	1.110	1.077	0.905
45	551.75 620 7E	20.92	24.57	20.01	12.01	-5.50	22.60	20.35	20.30	30.74	20.00	1.1//	1.070	1.004
40	678.89	2/ 92	27.92	20.31	12.01	-2.07	25.50	29.50	27.35	31.42	30.00	1.144	1.149	0.805
48	493 27	28 24	23.28	25.06	12 72	-7.10	33.81	27.87	28.65	30.70	30.00	1 213	1.071	0.883
40	588 59	30.27	27.18	27.37	10.60	-0.68	33.18	29.80	27.31	31.44	30.00	1.213	1.071	1 034
50	596.95	31.94	24.69	27.56	15.87	-10.43	34.76	26.87	28.91	30.56	30.00	1.294	1.057	0.817
51	615.34	35.18	22.54	27.98	25.70	-19.47	37.71	24.16	29.64	30.18	30.00	1.561	1.018	0.652
52	643.77	31.09	28.45	28.62	8.61	-0.62	32.58	29.82	27.79	31.17	30.00	1.093	1.121	1.026
53	847.76	36.78	32.00	32.85	11.97	-2.57	33.59	29.23	27.50	31.33	30.00	1.149	1.139	0.991
54	837.73	37.05	32.33	32.65	13.47	-0.99	34.04	29.70	26.70	31.80	30.00	1.146	1.191	1.039
55	608.65	30.66	27.16	27.83	10.14	-2.43	33.04	29.27	27.92	31.10	30.00	1.129	1.114	0.987
56	633.73	31.30	28.45	28.40	10.22	0.17	33.07	30.05	27.17	31.52	30.00	1.100	1.160	1.054
57	503.31	30.19	24.47	25.31	19.28	-3.34	35.78	29.00	26.02	32.21	30.00	1.234	1.238	1.003
58	810.98	35.32	31.98	32.13	9.94	-0.45	32.98	29.87	27.41	31.38	30.00	1.104	1.145	1.037
59	688.91	32.64	28.45	29.61	10.22	-3.93	33.07	28.82	28.33	30.87	30.00	1.147	1.090	0.950
60	506.65	27.76	24.57	25.39	9.34	-3.25	32.80	29.03	28.36	30.86	30.00	1.130	1.088	0.963
							A	verage ->	27.89	31.13	30.00	1.182	1.118	0.953
								St.Dev>	0.91	0.50	0.00	0.099	0.053	0.096

Та	ble	2 S1	с.

	Meas	sured dat	a for nom	inally 60	nm diam	eter	Measured ∆Dx &		$\Delta \Delta Dy$ applied for Db = 60 nm		Shape ratios		S	
S.No.	Aproj	Dx	Dy	Dxy	∆Dx	∆Dy	Dx	Dy	Dz(AFM)	Dxy(SEM)	Dxyz	R1	R2	R3
	nm <sup>2</sup>	nm	nm	nm	%	%	nm	nm	nm	nm	nm	Dx/Dy	Dxy/Dz	R2/R1
1	2150.47	57.45	51.65	52.32	9.82	-1.27	65.89	59.24	55.34	62.47	60.00	1.112	1.129	1.015
2	2266.16	57.92	52.17	53.70	7.85	-2.85	64.71	58.29	57.26	61.42	60.00	1.110	1.073	0.966
3	2027.98	58.62	48.72	50.80	15.39	-4.10	69.23	57.54	54.22	63.12	60.00	1.203	1.164	0.967
4	2395.46	60.23	52.17	55.22	9.07	-5.51	65.44	56.69	58.22	60.91	60.00	1.154	1.046	0.906
5	2837.81	74.75	49.81	60.10	24.38	-17.13	74.63	49.72	58.21	60.92	60.00	1.501	1.046	0.697
6	2368.24	61.57	53.49	54.90	12.14	-2.56	67.29	58.46	54.91	62.72	60.00	1.151	1.142	0.992
7	2204.92	58.33	52.00	52.97	10.11	-1.84	66.07	58.89	55.51	62.38	60.00	1.122	1.124	1.002
8	2504.35	65.69	49.57	56.46	16.35	-12.21	69.81	52.68	58.74	60.64	60.00	1.325	1.032	0.779
9	2361.44	62.23	53.24	54.82	13.51	-2.88	68.10	58.27	54.43	63.00	60.00	1.169	1.157	0.990
10	2259.36	61.07	49.33	53.62	13.88	-8.02	68.33	55.19	57.28	61.41	60.00	1.238	1.072	0.866
11	1932.70	56.44	48.38	49.60	13.79	-2.45	68.27	58.53	54.05	63.21	60.00	1.166	1.169	1.003
12	2170.89	56.19	52.17	52.56	6.90	-0.74	64.14	59.56	56.54	61.81	60.00	1.077	1.093	1.015
13	2096.03	60.23	46.69	51.65	16.61	-9.61	69.96	54.23	56.93	61.60	60.00	1.290	1.082	0.839
14	2041.59	55.77	49.57	50.97	9.40	-2.76	65.64	58.34	56.40	61.88	60.00	1.125	1.097	0.975
15	2640.45	63.10	58.35	57.97	8.84	0.65	65.30	60.39	54.77	62.80	60.00	1.081	1.147	1.060
16	1851.04	52.17	46.96	48.54	7.49	-3.26	64.50	58.05	57.70	61.19	60.00	1.111	1.060	0.954
17	2300.19	58.86	51.98	54.11	8.78	-3.93	65.27	57.64	57.42	61.34	60.00	1.132	1.068	0.943
18	2204.92	59.49	52.17	52.97	12.29	-1.51	67.38	59.09	54.25	63.10	60.00	1.140	1.163	1.020
19	1959.92	54.28	49.57	49.94	8.69	-0.76	65.21	59.54	55.63	62.31	60.00	1.095	1.120	1.023
20	2300.19	65.22	46.36	54.11	20.53	-14.32	72.32	51.41	58.10	60.97	60.00	1.407	1.049	0.746
21	2048.39	55.03	49.57	51.06	7.78	-2.93	64.67	58.24	57.35	61.37	60.00	1.110	1.070	0.964
22	2238.94	59.49	52.17	53.38	11.44	-2.26	66.86	58.64	55.09	62.62	60.00	1.140	1.137	0.997
23	1653.69	49.22	44.35	45.88	7.29	-3.33	64.37	58.00	57.85	61.10	60.00	1.110	1.056	0.952
24	2116.45	57.27	51.98	51.90	10.35	0.15	66.21	60.09	54.29	63.08	60.00	1.102	1.162	1.055
25	1776.18	52.82	46.96	47.55	11.10	-1.24	66.66	59.26	54.68	62.85	60.00	1.125	1.149	1.022
26	1980.34	55.03	49.57	50.20	9.61	-1.27	65.77	59.24	55.44	62.42	60.00	1.110	1.126	1.014
2/	2347.83	60.06	52.17	54.66	9.87	-4.55	65.92	57.27	57.22	61.44	60.00	1.151	1.074	0.933
28	2136.86	58.62	51.26	52.15	12.41	-1.72	67.45	58.97	54.31	63.07	60.00	1.144	1.161	1.015
29	1943.03	53.45	47.16	49.73	7.48	-5.16	64.49	56.90	58.86	60.58	60.00	1.133	1.029	0.908
30	2429.61	58.74	53.69	55.61	5.63	-3.45	63.38	57.93	58.83	60.59	60.00	1.094	1.030	0.942
31	2391.17	50.19 FC 44	53.03	55.17	9.10	-3.80	05.40	57.08	57.21	61.45	60.00	1.135	1.074	0.940
22	1706.05	50.44	47.16	17 02	0.40	-5.10	62.02	50.10	57.12	61.02	60.00	1.120	1.070	0.901
24	2251 62	50.15	47.10 E1.01	47.0Z	4.07	-1.50	64.24	59.17	50.01	50.00	60.00	1.005	1.052	0.969
25	2331.03	50.00	14 42	52 14	14.25	-0.70	68 56	55.95	61.64	59.99	60.00	1.130	0.960	0.809
35	2130.34	68.06	61 12	61.82	14.20	-14.02	66.02	50.22	54.40	62.01	60.00	1.341	1 159	1.027
37	1854.06	53.04	46 69	48 58	9.18	-3.89	65 51	57.67	57 18	61.46	60.00	1.120	1.138	0.946
38	2275 84	56.41	53 72	53.82	4 81	-0.19	62.89	59.89	57.10	61 37	60.00	1.150	1.075	1 019
39	2399 95	58.28	53.44	55.02	5 45	-3 30	63.27	58.02	58.84	60.59	60.00	1.090	1.070	0.944
40	3121 59	65 59	62 72	63.03	4.06	-0.49	62.44	59.70	57.94	61.06	60.00	1.036	1.050	1 008
41	2441 69	64 78	51 61	55 75	16 21	-7.42	69.73	55 55	55.77	62.24	60.00	1 255	1 116	0.889
42	2156 11	57 56	51.01	52 38	9.87	-2 47	65.92	58 52	55.99	62.11	60.00	1 126	1 109	0.985
43	2147.33	55.90	51.36	52.28	6.93	-1.76	64.16	58.94	57.12	61.50	60.00	1.088	1.077	0.989
44	1972.69	53.87	48.78	50.11	7.51	-2.65	64.51	58.41	57.33	61.38	60.00	1.104	1.071	0.969
45	2574.60	62.55	56.04	57.24	9.27	-2.10	65.56	58.74	56.09	62.06	60.00	1.116	1.106	0.991
46	2160.51	56.63	50.31	52.44	8.00	-4.07	64.80	57.56	57.91	61.07	60.00	1.126	1.055	0.937
47	2238.49	60.76	49.03	53.38	13.83	-8.14	68.30	55.12	57.38	61.35	60.00	1.239	1.069	0.863
48	2453.78	61.77	55.02	55.88	10.54	-1.54	66.32	59.07	55.13	62.59	60.00	1.123	1.135	1.011
49	2172.59	56.50	50.39	52.58	7.44	-4.17	64.46	57.50	58.27	60.88	60.00	1.121	1.045	0.932
50	2069.34	55.56	49.38	51.32	8.25	-3.77	64.95	57.74	57.60	61.24	60.00	1.125	1.063	0.945
51	2416.43	62.08	51.55	55.46	11.95	-7.05	67.17	55.77	57.66	61.20	60.00	1.204	1.061	0.881
52	2198.95	56.53	51.35	52.90	6.85	-2.93	64.11	58.24	57.85	61.11	60.00	1.101	1.056	0.960
53	2074.84	54.71	50.68	51.39	6.46	-1.39	63.88	59.17	57.15	61.48	60.00	1.080	1.076	0.996
54	2438.40	60.80	51.39	55.71	9.13	-7.76	65.48	55.34	59.61	60.20	60.00	1.183	1.010	0.854
55	2236.49	57.09	51.13	53.35	7.00	-4.16	64.20	57.50	58.51	60.76	60.00	1.117	1.038	0.930
56	3116.28	68.03	61.64	62.98	8.02	-2.13	64.81	58.72	56.75	61.69	60.00	1.104	1.087	0.985
57	2665.50	63.09	56.35	58.24	8.33	-3.26	65.00	58.05	57.25	61.42	60.00	1.120	1.073	0.958
58	2172.25	55.74	51.98	52.58	6.00	-1.14	63.60	59.32	57.26	61.42	60.00	1.072	1.073	1.000
59	2207.09	57.72	51.13	53.00	8.91	-3.53	65.35	57.88	57.11	61.50	60.00	1.129	1.077	0.954
60	2259.36	57.77	51.19	53.62	7.73	-4.54	64.64	57.28	58.34	60.85	60.00	1.129	1.043	0.924
							A	Average ->	56.93	61.62	60.00	1.148	1.084	0.949
1								St.Dev>	1.61	0.87	0.00	0.083	0.046	0.076

#### Table S1 d.

	Measure	d data f	or nomi	nally 100	nm dia	meter	Measured $\Delta Dx \& \Delta Dy$ applied for D			ed for Db =	= 100 nm	S	hape rati	os
S.No.	Aproj	Dx	Dy	Dxy	∆Dx	∆Dy	Dx	Dy	Dz(AFM)	Dxy(SEM)	Dxyz	R1	R2	R3
	nm <sup>2</sup>	nm	nm	nm	%	%	nm	nm	nm	nm	nm	Dx/Dy	Dxy/Dz	R2/R1
1	7957.22	108.36	97.66	100.63	7.68	-2.96	107.68	97.04	95.70	102.22	100.00	1.110	1.068	0.963
2	7480.38	104.50	96.10	97.57	7.10	-1.51	107.10	98.49	94.81	102.70	100.00	1.087	1.083	0.996
3	8577.11	123.55	88.04	104.48	18.25	-15.73	118.25	84.27	100.35	99.82	100.00	1.403	0.995	0.709
4	9584.43	116.47	109.86	110.45	5.46	-0.53	105.46	99.47	95.33	102.42	100.00	1.060	1.074	1.013
5	8761.88	112.89	104.74	105.60	6.90	-0.81	106.90	99.19	94.31	102.97	100.00	1.078	1.092	1.013
6	8589.03	112.86	102.54	104.55	7.95	-1.93	107.95	98.07	94.46	102.89	100.00	1.101	1.089	0.990
7	7629.40	103.58	97.66	98.54	5.11	-0.90	105.11	99.10	96.00	102.06	100.00	1.061	1.063	1.002
8	9530.78	121.90	105.30	110.14	10.68	-4.39	110.68	95.61	94.50	102.87	100.00	1.158	1.089	0.940
9	9155.27	113.49	107.42	107.95	5.14	-0.48	105.14	99.52	95.58	102.29	100.00	1.057	1.070	1.013
10	8052.59	109.32	99.37	101.24	7.98	-1.84	107.98	98.16	94.34	102.95	100.00	1.100	1.091	0.992
11	8845.33	121.63	97.66	106.10	14.63	-7.96	114.63	92.04	94.78	102.72	100.00	1.245	1.084	0.870
12	7885.70	108.36	97.66	100.18	8.10	-2.52	108.16	97.48	94.84	102.68	100.00	1.110	1.083	0.976
13	0854.53	117.24	83.01	93.40	25.52	-11.13	125.52	88.87	89.64	105.62	100.00	1.412	1.1/8	0.834
14	7212.10	112.41	06.69	95.81	27.33	-12.43	102.00	00 04	97.33	101.30	100.00	1.340	1.041	1.000
15	9936.09	119 73	112 73	112 /15	6.47	0.25	106.00	100.25	93.00	103.32	100.00	1.095	1.103	1.009
17	8779.76	113.75	102 54	105 71	7.74	-3.00	107.74	97.00	95.69	102.23	100.00	1 111	1.105	0.962
18	9638.07	116 70	102.54	110 75	5 37	-1.80	105.37	98.20	96.64	102.23	100.00	1.111	1.000	0.902
19	8529.43	111 83	102.54	104 19	7 33	-1 58	107.33	98.42	94 67	102.72	100.00	1.075	1.035	0.995
20	10526.18	129.42	109.55	115.75	11.81	-5.35	111.81	94.65	94.49	102.87	100.00	1.181	1.089	0.922
21	8755.92	111.83	104.98	105.56	5.93	-0.55	105.93	99.45	94.93	102.64	100.00	1.065	1.081	1.015
22	12809.04	158.24	108.12	127.68	23.93	-15.32	123.93	84.68	95.29	102.44	100.00	1.464	1.075	0.735
23	8678.44	111.19	104.31	105.10	5.79	-0.75	105.79	99.25	95.23	102.47	100.00	1.066	1.076	1.009
24	8410.22	111.19	101.67	103.46	7.47	-1.73	107.47	98.27	94.69	102.76	100.00	1.094	1.085	0.992
25	9280.44	116.83	107.03	108.68	7.50	-1.52	107.50	98.48	94.46	102.89	100.00	1.092	1.089	0.998
26	9506.94	114.77	107.42	110.00	4.34	-2.34	104.34	97.66	98.14	100.94	100.00	1.068	1.029	0.963
27	9167.19	115.06	108.08	108.02	6.52	0.06	106.52	100.06	93.82	103.24	100.00	1.065	1.100	1.034
28	10216.24	123.77	109.86	114.03	8.54	-3.65	108.54	96.35	95.62	102.26	100.00	1.127	1.069	0.949
29	9167.19	115.76	105.31	108.02	7.17	-2.51	107.17	97.49	95.71	102.21	100.00	1.099	1.068	0.972
30	10496.38	123.16	109.86	115.58	6.56	-4.95	106.56	95.05	98.73	100.64	100.00	1.121	1.019	0.909
31	11271.24	140.88	100.31	119.77	17.63	-16.25	117.63	83.75	101.51	99.25	100.00	1.405	0.978	0.696
32	9614.23	120.62	104.83	110.62	9.04	-5.23	109.04	94.77	96.77	101.66	100.00	1.151	1.051	0.913
33	9781.12	120.60	110.07	111.57	8.09	-1.35	108.09	98.65	93.78	103.26	100.00	1.096	1.101	1.005
34	9787.08	118.90	108.96	111.61	6.54	-2.38	106.54	97.62	96.15	101.98	100.00	1.091	1.061	0.972
35	9530.78	116.91	109.86	110.14	6.15	-0.25	106.15	99.75	94.44	102.90	100.00	1.064	1.090	1.024
36	9852.65	123.16	102.54	111.98	9.99	-8.43	109.99	91.57	99.29	100.36	100.00	1.201	1.011	0.841
37	0707.00	114.93	95.22	103.31	11.24	-7.84	111.24	92.10	97.54	101.25	100.00	1.207	1.038	0.800
20	7771 77	124.32	05.22	06.20	7.07	-7.52	107.07	92.00	90.80	101.01	100.00	1.202	1.049	1.009
39	9921 /0	112.26	106 55	105.06	6.80	-1.05	106.80	100 56	94.57	102.94	100.00	1.062	1.091	1.008
40	9083 75	114 36	103.63	107.52	6 36	-3.62	106.36	96 38	97.55	101.00	100.00	1 103	1.114	0.941
42	7551 91	104 30	95 22	98.04	6 38	-2.88	106.38	97 12	96 79	101.25	100.00	1.105	1.050	0.959
43	7551.91	105.32	97.09	98.04	7.43	-0.97	107.43	99.03	93.99	103.15	100.00	1.085	1.097	1.012
44	7230.04	101.28	94.95	95.93	5.58	-1.02	105.58	98.98	95.69	102.23	100.00	1.067	1.068	1.002
45	9167.19	115.47	104.84	108.02	6.90	-2.94	106.90	97.06	96.38	101.86	100.00	1.101	1.057	0.960
46	8451.94	110.70	101.61	103.72	6.73	-2.03	106.73	97.97	95.63	102.26	100.00	1.089	1.069	0.981
47	10454.66	126.98	110.49	115.35	10.08	-4.22	110.08	95.78	94.85	102.68	100.00	1.149	1.083	0.942
48	8726.12	111.56	103.42	105.38	5.86	-1.86	105.86	98.14	96.26	101.92	100.00	1.079	1.059	0.982
49	7319.45	101.55	95.22	96.52	5.21	-1.35	105.21	98.65	96.35	101.88	100.00	1.066	1.057	0.991
50	7247.93	119.13	87.40	96.04	24.04	-9.00	124.04	91.00	88.59	106.24	100.00	1.363	1.199	0.880
51	9912.25	123.16	108.71	112.32	9.66	-3.21	109.66	96.79	94.22	103.02	100.00	1.133	1.093	0.965
52	10216.24	125.85	108.37	114.03	10.36	-4.97	110.36	95.03	95.35	102.41	100.00	1.161	1.074	0.925
53	7992.98	113.60	92.08	100.86	12.63	-8.70	112.63	91.30	97.25	101.40	100.00	1.234	1.043	0.845
54	8380.41	111.83	97.20	103.28	8.28	-5.88	108.28	94.12	98.13	100.95	100.00	1.150	1.029	0.894
55	8940.70	118.53	102.66	106.67	11.11	-3.76	111.11	96.24	93.52	103.41	100.00	1.155	1.106	0.958
56	9816.89	119.33	109.86	111.78	6.76	-1.71	106.76	98.29	95.30	102.43	100.00	1.086	1.075	0.990
57	8249.28	108.88	102.37	102.47	6.26	-0.09	106.26	99.91	94.19	103.04	100.00	1.064	1.094	1.028
58	8356.57	110.70	101.55	103.13	7.34	-1.53	107.34	98.47	94.61	102.81	100.00	1.090	1.087	0.997
59	7104.87	101.16	95.22	95.09	6.38	0.13	106.38	100.13	93.88	103.21	100.00	1.062	1.099	1.035
60	10031.46	121.51	108.87	112.99	7.54	-3.65	107.54	96.35	96.51	101.79	100.00	1.116	1.055	0.945
									05.44	102.20	100.00	1 4 2 7	1.072	0.051
							A	verage ->	95.44	102.38	100.00	1.13/	1.073	0.951
1								JL.DEV>	2.00	1.11	0.00	0.099	0.035	0.080

**Table S2.** Spherical equivalent volume diameters (Dxyz) calculated by combining SEM (Dxy) and AFM (Dz) reported diameters of the Au nanoparticles extracted from Fig. 15 in Ref. [3].

S.No.	AFM(Dz)	SEM(Dxy)	Dxyz		
	nm	nm	nm		
1	43.50	42.50	42.83		
2	42.20	45.50	44.37		
3	43.20	49.00	46.98		
4	46.50	49.00	48.15		
5	49.60	51.00	50.53		
6	48.80	51.50	50.58		
7	47.90	52.00	50.60		
8	50.70	52.10	51.63		
9	46.90	52.40	50.50		
10	50.40	52.50	51.79		
11	53.80	52.50	52.93		
12	51.80	53.00	52.60		
13	46.50	53.30	50.93		
14	43.80	53.50	50.05		
15	46.80	54.00	51.48		
16	46.20	54.20	51.39		
17	50.70	54.30	53.07		
18	49.00	54.80	52.79		
19	54.10	54.80	54.57		
20	50.00	55.10	53.34		
21	54.30	55.10	54.83		
22	46.50	55.50	52.32		
23	49.20	55.50	53.32		
24	49.40	55.50	53.39		
25	50.90	55.70	54.05		
26	49.90	55.80	53.76		
27	45.20	56.00	52.14		
28	55.20	56.80	56.26		
29	47.10	57.50	53.80		
30	49.50	58.00	55.02		
31	50.10	58.20	55.36		
32	50.80	58.20	55.62		
33	53.20	58.40	56.61		
34	55.80	58.40	57.52		
35	52.80	58.50	56.53		
36	48.60	58.60	55.06		
37	49.70	58.60	55.47		
38	56.00	58.80	57.85		
39	50.80	59.00	56.13		
40	52.30	59.00	56.68		

#### Table S2 continued...

41	43.80	59.50	53.72
42	49.70	59.50	56.04
43	46.90	59.80	55.15
44	45.80	60.00	54.83
45	51.10	60.50	57.19
46	46.00	60.70	55.34
47	51.10	61.20	57.63
48	50.70	61.30	57.54
49	58.30	61.40	60.35
50	52.10	61.70	58.32
51	49.00	62.30	57.51
52	53.40	62.40	59.24
53	50.20	62.90	58.34
54	52.30	63.30	59.40
55	54.10	63.40	60.13
56	54.60	63.80	60.57
57	52.70	64.00	59.99
58	45.40	64.50	57.38
59	39.90	64.50	54.96
60	50.10	64.60	59.35
61	53.50	64.70	60.73
62	54.40	64.90	61.19
63	36.50	65.00	53.63
64	48.80	65.30	59.26
65	51.10	66.00	60.60
66	52.60	67.80	62.30
67	59.60	68.00	65.08
68	47.90	68.40	60.74
Average->	49.73	58.17	55.14
St.Dev>	4.12	5.39	4.13

## Table S3. Spherical equivalent volume diameters of the pebbles measured using the two independent methods.

(a) Spherical equivalent volume diameters of the irregularly shaped pebbles calculated by combining the projection area and the peak-height. Images of irregularly shaped, aquarium glass pebbles spread onto a smooth surface were acquired and then analyzed using ImageJ to determine Aproj, from which Dxy was calculated. The peak-height of each pebble (Dz) was measured using a digital caliper. The volume was calculated using Dxy and Dz assuming an ellipsoidal shape. From this, an experimental spherical equivalent volume diameter (Dxyz1) was then calculated.

(b) Spherical equivalent volume diameters of the irregularly shaped pebbles calculated using mass and density. The mass of each pebble was measured, from which volume was calculated using the density of the pebbles (we measured the density of the pebbles to be  $2.55 \times 10^{-3}$  g/mm<sup>3</sup>). The spherical equivalent volume diameter (Dxyz2) was then calculated based on the volume.

(c) Shape ratios. The pebbles used have a mean R3 ratio of 1.1 indicating that they are dominated relatively more by flattening than by elongation.

Averages, standard deviations, and standard deviations of the mean (SDOM =  $\frac{Std.Dev.}{n^{1/2}}$ , where n is the number of measurements) of Dxyz1 and Dxyz2 are provided at the bottom of the table. Statistical standard uncertainty ( $\sigma = \sqrt{\sum SDOM^2}$ ) of Dxyz1 and Dxyz2 as estimated using error propagation (quadrature addition) is 0.172. The fact that the absolute difference (0.16) between Dxyz1 and Dxyz2 is less than  $\sigma$  means that the difference is consistent with zero. This implies that Dxyz1 and Dxyz2 are statistically in agreement.

	a							b				С		
	Dxyz	Dxyz calculated using projection area and height							Dxyz calculated using mass and density			Shape ratios		
S.No.	Aproj	Dx	Dy	Dxy	Dz	Volume	Dxyz1	Mass	Volume	Dxyz2	R1	R2	R3	
	mm <sup>2</sup>	mm	mm	mm	mm	mm <sup>3</sup>	mm	g	mm <sup>3</sup>	mm	Dx/Dy	Dxy/Dz	R2/R1	
1	596.78	30.49	25.78	27.56	18.56	7384.17	24.16	20.43	8011.76	24.82	1.18	1.48	1.26	
2	617.50	30.43	26.28	28.03	18.03	7422.39	24.20	20.68	8109.80	24.92	1.16	1.55	1.34	
3	557.03	29.05	24.43	26.63	16.70	6201.59	22.79	16.23	6364.71	22.99	1.19	1.59	1.34	
4	511.92	28.53	24.18	25.53	18.77	6405.81	23.04	16.23	6364.71	22.99	1.18	1.36	1.15	
5	443.41	27.25	21.21	23.76	22.45	6636.37	23.31	16.39	6427.45	23.07	1.28	1.06	0.82	
6	564.64	31.51	23.28	26.81	17.86	6722.99	23.41	17.48	6854.90	23.57	1.35	1.50	1.11	
/	584.94	30.83	25.15	27.28	20.56	8017.58	24.83	21.00	8235.29	25.05	1.23	1.33	1.08	
8	567.04	28.20	26.55	20.80	19.68	7439.53	24.22	16.50	7294.12	24.06	1.06	1.37	1.28	
10	522.56	27.60	24.90	27.54	16.62	5789.99	22.14	16.75	6/11 76	23.25	1.20	1.92	1.00	
11	468 22	26.74	23.52	23.75	22.26	6948 38	23.67	17.60	6901.96	23.05	1.15	1.55	0.94	
12	538.13	30.56	22.18	26.17	22.31	8003.80	24.81	21.48	8423.53	25.24	1.38	1.17	0.85	
13	528.34	28.88	23.39	25.93	15.46	5445.44	21.82	14.75	5784.31	22.27	1.23	1.68	1.36	
14	438.27	28.17	20.15	23.62	22.66	6620.72	23.29	17.45	6843.14	23.55	1.40	1.04	0.75	
15	546.39	29.37	23.54	26.37	21.65	7886.17	24.69	20.53	8050.98	24.86	1.25	1.22	0.98	
16	509.45	27.87	23.99	25.46	16.26	5522.47	21.93	14.42	5654.90	22.10	1.16	1.57	1.35	
17	521.72	27.92	23.28	25.77	17.00	5912.79	22.43	15.88	6227.45	22.82	1.20	1.52	1.26	
18	575.98	29.74	25.25	27.08	19.71	7568.38	24.36	18.75	7352.94	24.12	1.18	1.37	1.17	
19	606.93	30.32	25.65	27.79	17.70	7161.79	23.91	19.63	7698.04	24.49	1.18	1.57	1.33	
20	596.71	30.67	24.93	27.56	19.53	7769.18	24.57	20.52	8047.06	24.86	1.23	1.41	1.15	
21	574.37	28.07	26.66	27.04	23.67	9063.53	25.86	21.10	8274.51	25.09	1.05	1.14	1.08	
22	466.39	30.19	19.39	24.36	24.30	7555.47	24.34	20.14	7898.04	24.70	1.56	1.00	0.64	
23	565.91	29.00	25.62	26.84	19.41	7322.88	24.09	20.96	8219.61	25.04	1.13	1.38	1.22	
24	464.06	27.39	21.30	24.30	25.73	7960.19	24.77	19.96	7827.45	24.63	1.29	0.94	0.73	
25	553.01	30.56	22.80	26.53	21.56	7948.61	24.76	19.71	7729.41	24.53	1.34	1.23	0.92	
20	511.92	29.69	21.50	25.53	23.30	7951.81	24.76	20.38	7992.10	24.80	1.38	1.10	0.79	
27	186 60	26.74	22.57	20.49	23.55	7180.26	25.47	19.88	7105.08	24.00	1.30	1.12	0.01	
20	499.80	20.74	23.59	24.09	22.13	7776.83	23.55	18.12	7447.06	23.85	1.13	1.12	0.33	
30	533.14	29.06	23.78	26.05	22.84	8117.87	24.93	20.13	7894.12	24.70	1.20	1.14	0.93	
31	661.70	31.28	28.16	29.02	15.04	6634.61	23.31	17.56	6886.27	23.60	1.11	1.93	1.74	
32	647.32	33.43	24.19	28.70	16.51	7124.81	23.87	20.74	8133.33	24.95	1.38	1.74	1.26	
33	486.55	29.04	20.78	24.88	22.41	7269.00	24.03	17.38	6815.69	23.52	1.40	1.11	0.79	
34	469.91	27.30	22.56	24.46	22.73	7120.73	23.87	18.15	7117.65	23.86	1.21	1.08	0.89	
35	570.98	28.74	25.77	26.96	17.42	6631.03	23.31	18.31	7180.39	23.93	1.12	1.55	1.39	
36	583.87	29.68	24.37	27.26	17.83	6940.32	23.66	18.99	7447.06	24.23	1.22	1.53	1.26	
37	501.91	27.81	23.06	25.27	19.86	6645.30	23.32	17.87	7007.84	23.74	1.21	1.27	1.06	
38	570.77	29.27	24.47	26.95	16.28	6194.79	22.78	16.85	6607.84	23.28	1.20	1.66	1.38	
39	527.64	29.70	22.20	25.91	22.11	7777.37	24.58	19.91	7807.84	24.61	1.34	1.17	0.88	
40	619.41	33.05	25.30	28.08	16.66	6879.54	23.59	17.78	6972.55	23.70	1.31	1.69	1.29	
41	595.30	29.84	25.72	27.53	17.89	7099.96	23.84	19.36	/592.16	24.38	1.16	1.54	1.33	
42	440 50	20.32	22.4/	23.44	15.20	3133.79	21.43	14.08	5/ 50.00	22.23	1.20	1.07	1.33	
43	552 26	20.04	21.59	23.08	24.00	6832 20	23.98	17.72	6949.02	23.07	1.23	1 /2	U./8 1 2/	
44	687 92	30.89	24.41	29.54	16.52	7727 59	23.34	20.07	7870 59	23.07	1.10	1.45	1.59	
46	563.37	29.58	24.78	26.78	15.06	5656 25	22.10	14.98	5874.51	22.38	1.19	1.78	1.49	
47	600.31	31.02	24.40	27.64	20.48	8196.16	25.01	20.81	8160.78	24.98	1.27	1.35	1.06	
48	463.00	29.26	18.66	24.27	21.65	6682.58	23.37	17.47	6850.98	23.56	1.57	1.12	0.72	
49	504.38	29.65	22.21	25.34	23.27	7824.58	24.63	19.12	7498.04	24.28	1.34	1.09	0.82	
50	508.47	29.85	21.66	25.44	20.62	6989.71	23.72	17.40	6823.53	23.53	1.38	1.23	0.90	
51	582.54	30.93	23.81	27.23	21.20	8233.29	25.05	19.18	7521.57	24.31	1.30	1.28	0.99	
52	615.46	31.88	24.30	27.99	16.66	6835.70	23.54	19.33	7580.39	24.37	1.31	1.68	1.28	
53	495.14	27.83	23.44	25.10	22.21	7331.43	24.10	17.50	6862.75	23.57	1.19	1.13	0.95	
54	403.09	26.08	20.71	22.65	23.50	6315.12	22.93	16.32	6400.00	23.03	1.26	0.96	0.77	
		Average ->								23.96	1.25	1.36	1.10	
				Stan	dard de	viation ->	0.97			0.82	0.11	0.26	0.26	
						SDOM ->	0.13		o /-	0.11				
								σ =	0.17					
							Uxyz	1 - DXYZZ =	0.16					