On-the-fly, Sample-tailored Optimisation of NMR Experiments

Jonathan R. J. Yong and Mohammadali Foroozandeh*

Chemistry Research Laboratory, Department of Chemistry, University of Oxford, 12 Mansfield Road, Oxford, OX1 3TA, U.K.

* mohammadali.foroozandeh@chem.ox.ac.uk



NMR parameters must often be adjusted for maximum performance on different samples or instruments. We show that, using our implementation of POISE (Parameter Optimisation by Iterative Spectral Evaluation), a wide range of experiments may be optimised on the fly in a fully automated fashion, using only a few spectrum acquisitions.

Abstract

NMR experiments, indispensable to chemists in many areas of research, are often run with generic, unoptimised experimental parameters. This approach makes robust and automated acquisition on different samples and instruments extremely challenging. Here, we introduce NMR-POISE (Parameter Optimisation by Iterative Spectral Evaluation), the first demonstration of on-the-fly, sample-tailored, and fully automated optimisation of a wide range of NMR experiments. We illustrate how POISE maximises spectral sensitivity and quality with a diverse set of 1D and 2D examples, ranging from HSQC and NOESY experiments to ultrafast and pure shift techniques. Our Python implementation of POISE has an interface integrated into Bruker's TopSpin software, one of the most widely used platforms for NMR acquisition and automation, allowing NMR optimisations to be run without direct user supervision. We predict that POISE will find widespread usage in academia and industry, where sample-specific and automated experiment optimisation is mandatory.

Keywords: NMR spectroscopy \cdot Analytical methods \cdot NMR-POISE \cdot Optimisation \cdot Automation

Nuclear magnetic resonance (NMR) spectroscopy is firmly established as one of the most important analytical methods in chemistry.^{1,2} Although many NMR experiments are designed to be as general as possible, the sheer diversity of molecular structure means that experimental parameters must often be varied in order to achieve optimal performance for a particular sample. Using poorly set parameters can lead to inferior spectra, as measured by (for example) the attainable signal-to-noise ratio (SNR) per unit time, or the presence of artefacts arising from unwanted coherences. A typical way around this is to use a "compromise" value, chosen to provide acceptable (albeit suboptimal) performance across a range of samples.

In contrast, the direct optimisation of experimental parameters on a per-sample basis is rarely carried out. When it is, it typically relies on trial-and-error or otherwise inefficient methods, leading to highly tedious and time-consuming processes. For example, the current best routines (e.g. **popt** in Bruker's TopSpin software) rely on systematic evaluation of parameter values using an exhaustive grid search. A more desirable approach would be to utilise an optimisation algorithm³ which iteratively monitors the quality of a series of spectra acquired with different parameter values, eventually converging to optimal values for the sample under study. Although there exist isolated examples of such iterative optimisations in laser,⁴ electron spin resonance,⁵ and nuclear quadrupole resonance⁶ spectroscopies, such a procedure has only been used in NMR by Emsley and co-workers in the design of dipolar decoupling pulses.^{7–9} Note that the proposed method is distinct from other model-based "optimisations", e.g. for sparse sampling in multidimensional NMR^{10,11} or relaxation measurements,¹² where data from multiple acquisitions are *aggregated* until a certain threshold, such as a sufficient confidence in peak locations or extracted parameters, is reached.

Despite its apparent simplicity, iterative optimisation of NMR experiments clearly enjoys little popularity. This can perhaps be attributed to a perception that such a process would require too much time, as well as the challenge of implementing it in a manner that can be suitably generalised for multiple disparate contexts.

In this work, we dispel these notions by introducing NMR-POISE (*Parameter Optimisation by Iterative Spectral Evaluation*), a method which enables the efficient optimisation of a wide range of NMR experimental parameters and yields demonstrably better spectra in only a matter of minutes. NMR-POISE is provided as a free Python package, with a powerful user interface allowing optimisations to be run directly from the command line of Bruker's TopSpin software. This means that optimisations can proceed in a fully automated fashion without any user input, a benefit which is especially relevant in high-throughput settings or when non-expert users are present. The optimisation process (Figure 1) requires a *cost function*, a user-defined function which measures how good a spectrum is; POISE provides a number of built-in cost functions for typical use cases, e.g. for maximising or minimising signal intensity. Currently, three optimisation algorithms are provided: Py-BOBYQA, a trust-region interpolation method¹³ (enabled by default); the Nelder–Mead (NM) simplex method;¹⁴ and the multidimensional search (MDS) method.¹⁵



Figure 1: Flowchart depicting the main steps in a POISE optimisation.

Experimental section

All NMR spectra in this work were acquired on Bruker AV III or NEO spectrometers, with field strengths ranging from 400 to 700 MHz. POISE was executed on the spectrometer terminals using Python 3.8.1, 3.8.5, and 3.9.2 on Windows 10. In general, POISE can be run on any Bruker spectrometer; the minimum version of Python needed is 3.6. Experimental details of individual optimisations are given in the *Supporting Information*.

Results and discussion

The first example of a POISE control optimisation is the calibration of a 90° 1 H pulse, a quantity which is fundamental to practically every NMR experiment.¹⁶ In a simple 1 H pulse-acquire experiment, a perfect 360° pulse would generate no detectable magnetisation. The intensity of the resulting signal can thus be used as a cost function which we seek to minimise.



Figure 2: Optimisation trajectory for pulse width calibration on a sample of 45 mM ferulic-acid in DMSO- d_6 using: (a) TopSpin's popt, which uses a grid search; and (b) POISE, using the BOBYQA algorithm. The grey line indicates the "true" objective function. Each coloured arrow indicates a single function evaluation, unfilled dots represent the initial point, and filled dots indicate optima at which the cost function is minimised.

The benefits of direct optimisation, as compared to the grid search of popt, can be easily appreciated upon consideration of Figure 2. By manually acquiring an array of spectra with different pulse widths, we first ascertained that the spectral intensity was minimised at 48.4 µs (Figure S5); knowing this allows us to evaluate the accuracy of the optimisation routines. **popt** faithfully reproduces this optimum, but must evaluate a total of 41 spectra to search the entire parameter space, which takes 5 minutes (Figure 2a). Although the search may in principle be terminated immediately after locating the minimum, this requires constant monitoring by the user, since the optimum point is not known a priori. On the other hand, POISE quickly homes in on the minimum without any user intervention: for example, using the BOBYQA algorithm, it requires only 70 seconds and 9 spectrum acquisitions to identify an accurate minimum of 48.3 us, even when given a suboptimal starting point of 43 µs (Figure 2b). Convergence is even faster when a good initial point is provided (Tables S2 to S4). As a "rapid" alternative for pulse width calibration, the Bruker pulsecal routine (which uses a nutation experiment¹⁷) completes in 37 seconds; however, it underestimates the 360° null by almost 2 µs (Table S2). For routine 1D ¹H NMR spectra, the exact 90° pulse width is often not of interest. However, by instead searching for a *maximum* in the spectral intensity, POISE can likewise be used to identify the pulse width corresponding to the Ernst angle¹⁸ (Section S5), which vields the maximal sensitivity for a specific repetition time.

As long as a suitable cost function can be defined, POISE is capable of providing rapid results with no compromise in accuracy. Our POISE interface is thus designed to grant users great flexibility in defining their own cost functions, which allows the optimisation of a wide range of parameters which simply cannot be carried out using currently available tools. One such example is the mixing time in nuclear Overhauser effect spectroscopy (NOESY) experiments, which typically varies from hundreds of milliseconds to several seconds. By finding the optimal mixing time for a given sample, the sensitivity of the resulting NOESY spectrum can be maximised. Here, we perform this optimisation using a selective 1D NOESY experiment (Figure 3a), using the intensity of the spectrum excluding the selectively irradiated peak as the cost function. Starting from an initial value of 0.5 s, a typical value for many organic molecules, POISE is capable of finding the optimum mixing time for the sample under study (3.5 s). This optimisation converges in around 5 minutes (Table S8), leading to an approximate twofold increase in NOESY crosspeak intensity for this particular case (Figure 4).



Figure 3: Pulse sequences used for POISE optimisations; (a) 1D selective NOESY, (b) ASAP-HSQC,¹⁹ (c) echo-planar spectroscopic imaging (EPSI) acquisition scheme used in ultrafast NMR,^{20,21} and (d) spin echo sequence used for evaluation of PSYCHE²² refocusing element, (e) 1D NOESY sequence with presaturation, used for solvent suppression.²³ Details of the optimisation processes, along with complete lists of experimental parameters used, can be found in the *Supporting Information*.

Where possible, POISE optimisations should be performed using 1D pulse sequences, in order to minimise the optimisation duration. The NOE mixing time optimisation shown here, for example, takes only a fraction of the time needed for a 2D NOESY. With more dilute samples, a larger number of transients per function evaluation may be required in order for a reliable cost function to be measured; consequently, the optimisation as a whole will take a longer time to run. However, since the same can be said of the final 2D spectrum that is acquired, the benefits that can be reaped are correspondingly larger, making the "return on investment" the same. A similar strategy can be used to search for the "null" in an inversion-recovery (IR) experiment: by running the optimisation on a 1D IR (i.e. a single increment of the full experiment), T_1 relaxation times can be measured in a far shorter time than the full 2D IR (Section S7).

An example of a 2D experiment which can be optimised in a reasonable time is the ASAP-HSQC experiment (Figure 3b),¹⁹ where the recovery delay between transients is dramatically shortened. Decreasing the duration of the INEPT delay leads to greater sensitivity through Ernst angle-style excitation of ¹³C-bound proton magnetisation. Using the projection of the full 2D spectrum onto the f_2 dimension as a cost function, POISE is able to optimise this delay to approximately 60% of its original length. This takes approximately 3 minutes for the same sample as above (Table S12)



Figure 4: 2D NOESY spectra of a sample of 120 mM 3-fluorophenylboronic acid in DMSO- d_6 , obtained (a) before and (b) after optimising the mixing time. Optimisation leads to an approximate $2 \times$ improvement in crosspeak sensitivities. In order to accomplish this sensitivity improvement, one would need to acquire the unoptimized experiment with four times as many transients.

and provides an 18–27% increase in sensitivity (Figure 5). This procedure can be extended to the optimisation of other delays in 2D spectra, particularly with "fast" variants such as the ASAP method used here.²⁴



Figure 5: (a) 2D ASAP-HSQC spectrum of a sample of 120 mM 3-fluorophenylboronic acid in DMSO- d_6 . (b) f_2 projections of ASAP-HSQC spectra before (grey dotted line) and after POISE optimisation of the INEPT delay (green solid line). The sensitivity increases attained are indicated above each peak.

The logical extension of the preceding example is ultrafast 2D NMR experiments, in which the combination of spatially encoded t_1 evolution and the echo-planar spectroscopic imaging (EPSI) detection method enables single-scan acquisition of 2D spectra.^{20,21} Despite the substantial efforts of Giraudeau and coworkers in the development of an accessible protocol,²¹ obtaining high-quality ultrafast spectra remains a nontrivial task.

We illustrate here how POISE can be used to aid in the setup of ultrafast experiments, in this case by compensating for systematic imbalances in the positive and negative gradients used during EPSI acquisition. Such an imbalance leads to a progressive shifting of echoes in the spatial frequency (k)domain (Figure 6a), ultimately causing degradation of lineshapes in the resulting spectra.²⁰ The pulse programme used for optimisation consists of 90° excitation, immediately followed by EPSI acquisition (Figure 3c); the amplitudes of the negative gradients are multiplied by a factor α in order to compensate for any imbalance. Here, the cost function is the slope generated by the echoes in a (k, t_2) data matrix (Figures 6a and 6b), which directly reflects the extent of gradient imbalance. On our system, POISE finds the optimal value of α to be 1.00038 in under a minute (Table S13). Such a small imbalance may appear inconsequential, but in fact leads to marked improvements in the resulting lineshapes of ultrafast spectra, such as the TOCSY spectra shown in Figures 6c and 6d.



Figure 6: (a) and (b): Magnitude-mode (k, t_2) data matrix obtained from the excitation-EPSI pulse sequence (Figure 3c), before and after optimisation of the gradient imbalance factor α . Data points acquired during the negative gradients have been discarded. (c) and (d): The corresponding 2D ultrafast TOCSY spectra of a sample of 45 mM ferulic acid in DMSO- d_6 (10 ms MLEV-16 mixing), before and after optimisation.

Finally, another major advantage of POISE is that its algorithms are capable of simultaneously optimising multiple parameters in an efficient manner. Although this is possible with a grid search, the time required scales exponentially with the number of parameters, quickly rendering such a procedure unfeasible. We show how this capability of POISE can be used to search for optimal performance in PSYCHE broadband homonuclear decoupling.²² Here, homodecoupling is accomplished by the application of a J-refocusing element (JRE), comprised of two saltire chirped pulses²⁵ and a weak gradient. These can be expressed in terms of several parameters, namely the flip angle β , gradient amplitude g, bandwidth ΔF , and duration τ_p .

We used a 1D experiment for optimisation of these parameters, which consists of a 90° excitation pulse followed by a spin echo with the JRE in the middle (Figure 3d). Any non-idealities in the JRE can be captured in a "spectral difference" cost function (Section S10), which compares this against a standard pulse-acquire experiment. The simultaneous optimisation of all four parameters can be completed in around 15 minutes and leads to noticeably better decoupling performance for the sample under study, especially for the peak at 1.36 ppm (Figure 7a). Furthermore, the optimised spectrum exhibits improved suppression of the artefacts at 2.5 ppm (Figure 7b), which arise from a challenging case of strong coupling. Although the optimised spectrum has a decreased SNR, the cost function used here also penalises sensitivity losses, thereby ensuring that any SNR loss is not excessive.

In a similar fashion, we applied POISE to the 1D NOESY / presaturation experiment commonly used for water suppression (Figure 3e),²³ where up to four parameters (transmitter offset, presaturation power, mixing time, and presaturation duration) are optimised. This procedure takes between 4 and 25 minutes, depending on the number of parameters under optimization (Section S11). By applying appropriate upper bounds on the presaturation power and duration, POISE greatly reduces the intensity of the residual water peak, while also ensuring that nearby peaks of interest are not lost (Figure 8). Finally, a further example of multiple-parameter optimisation is presented in the context of diffusion NMR, enabling a fully automated diffusion ordered spectroscopy (DOSY) experiment (Section S12).



Figure 7: Insets from spectra of a sample of 45 mM andrographolide in DMSO- d_6 ; artefacts in PSYCHE spectra are indicated with asterisks. Left column: Ordinary coupled ¹H spectrum; middle column: unoptimised PSYCHE spectrum; right column: PSYCHE spectrum after simultaneous optimisation of four parameters. (a) Peak at 1.36 ppm, where recoupling artefacts are eliminated through optimisation. (b) Peak at 2.5 ppm, where optimisation leads to a substantial reduction in strong coupling artefacts.

(a)	М			M			Μ	
~~~	www.	mm	www	mon	•∿~	mm	mar h	mm
	5.2	5.1	5.2	5.1		5.3	2	5.1
(b)		/9			/9			/9
4.8	4.7	4.6	4.8	4.7	4.6	4.8	4.7	4.6
(c)	M			V				
~~~		····			<del></del>	- www	$\sqrt{4}$	
	4.6	4.5	4.6	4.5	;	4	.6	4.5
	δ ¹ Η / μ	opm	δ ¹ Η	l / ppm		δ	¹ H/p	pm

Figure 8: Spectra of a sample of rodent urine. *Left column*: Unoptimised spectrum; *middle column*: spectrum following optimisation of two parameters (transmitter offset and presaturation power); *right column*: spectrum after optimisation of all four parameters. (a) Peak at 5.15 ppm. (b) Residual water peak. (c) Peak at 4.56 ppm. Note that the flanking peaks in (a) and (c) are largely unaffected by the improvement in water suppression.

Through a series of examples, we have shown how POISE can be used to great effect in the optimisation of NMR parameters. Each optimisation provides parameter values that are tailored to the instrument and sample under study, typically in only a few minutes and in a fully automated fashion through the TopSpin command line or acquisition scripts. Importantly, POISE allows users to define arbitrary cost functions, making it a highly versatile tool that can be adapted to the myriad uses of NMR in modern chemistry. POISE is fully compatible with, and will be immediately usable in, other fields such as solid-state and biomolecular NMR. More generally, we also note that the iterative optimisation process demonstrated here can be extended to other forms of spectroscopy.

The source code of NMR-POISE is publicly hosted on GitHub: https://github.com/foroozand ehgroup/nmrpoise. Detailed instructions on installation, cost function definition, and optimisation setup may be found in the user documentation, which is included in the *Supporting Information* and can also be found online at https://foroozandehgroup.github.io/nmrpoise.

Acknowledgements

We thank Coralia Cartis, Tim Claridge, Jean-Nicolas Dumez, David Goodwin, Iain Swan, and Chris Waudby for helpful discussions; as well as Fay Probert and Abi Yates for providing the rodent urine sample. J.R.J.Y. thanks the Clarendon Fund (University of Oxford) and the EPSRC Centre for Doctoral Training in Synthesis for Biology and Medicine (EP/L015838/1) for a studentship, generously supported by AstraZeneca, Diamond Light Source, Defence Science and Technology Laboratory, Evotec, GlaxoSmithKline, Janssen, Novartis, Pfizer, Syngenta, Takeda, UCB, and Vertex. M.F. thanks the Royal Society for a University Research Fellowship and a University Research Fellow Enhancement Award (Grant Nos. URF\R1\180233 and RGF\EA\181018).

References

- 1. Cavanagh, J., *Protein NMR Spectroscopy: Principles and Practice*, 2nd ed.; Academic Press: Burlington, Mass., 2007.
- 2. Claridge, T. D. W., *High-Resolution NMR Techniques in Organic Chemistry*, 3rd ed.; Elsevier: Amsterdam, 2016.
- 3. Nocedal, J.; Wright, S. J., Numerical Optimization, 2nd ed.; Springer: New York, 2006.
- Bardeen, C. J.; Yakovlev, V. V.; Wilson, K. R.; Carpenter, S. D.; Weber, P. M.; Warren, W. S. Chem. Phys. Lett. 1997, 280, 151–158.
- 5. Goodwin, D. L.; Myers, W. K.; Timmel, C. R.; Kuprov, I. J. Magn. Reson. 2018, 297, 9–16.
- 6. Monea, C. J. Magn. Reson. 2020, 321, 106858.
- 7. De Paëpe, G.; Hodgkinson, P.; Emsley, L. Chem. Phys. Lett. 2003, 376, 259–267.
- 8. Elena, B.; de Paëpe, G.; Emsley, L. Chem. Phys. Lett. 2004, 398, 532–538.
- Salager, E.; Dumez, J.-N.; Stein, R. S.; Steuernagel, S.; Lesage, A.; Elena-Herrmann, B.; Emsley, L. Chem. Phys. Lett. 2010, 498, 214–220.
- Eghbalnia, H. R.; Bahrami, A.; Tonelli, M.; Hallenga, K.; Markley, J. L. J. Am. Chem. Soc. 2005, 127, 12528–12536.
- 11. Hansen, A. L.; Brüschweiler, R. Angew. Chem. Int. Ed. 2016, 55, 14169–14172.
- 12. Song, Y.-Q.; Tang, Y.; Hürlimann, M. D.; Cory, D. G. J. Magn. Reson. 2018, 289, 72-78.
- 13. Cartis, C.; Fiala, J.; Marteau, B.; Roberts, L. ACM Trans. Math. Softw. 2019, 45, 1–41.
- 14. Nelder, J. A.; Mead, R. The Computer Journal 1965, 7, 308–313.
- 15. Dennis Jr., J. E.; Torczon, V. SIAM J. Optim. 1991, 1, 448-474.
- 16. Keifer, P. A. Concepts Magn. Reson. 1999, 11, 165–180.
- 17. Wu, P. S. C.; Otting, G. J. Magn. Reson. 2005, 176, 115–119.
- 18. Ernst, R. R.; Anderson, W. A. Rev. Sci. Instrum. 1966, 37, 93–102.
- 19. Schulze-Sünninghausen, D.; Becker, J.; Luy, B. J. Am. Chem. Soc. 2014, 136, 1242–1245.
- 20. Frydman, L.; Lupulescu, A.; Scherf, T. J. Am. Chem. Soc. 2003, 125, 9204–9217.
- 21. Gouilleux, B.; Rouger, L.; Giraudeau, P. Annu. Rep. NMR Spectrosc. 2018, 75–144.
- Foroozandeh, M.; Adams, R. W.; Meharry, N. J.; Jeannerat, D.; Nilsson, M.; Morris, G. A. Angew. Chem. Int. Ed. 2014, 53, 6990–6992.
- 23. Mckay, R. T. Concepts Magn. Reson. 2011, 38A, 197–220.
- 24. Schanda, P. Prog. Nucl. Magn. Reson. Spectrosc. 2009, 55, 238–265.
- 25. Foroozandeh, M. J. Magn. Reson. 2020, 318, 106768.

Supporting Information

On-the-fly, Sample-tailored Optimization of NMR Experiments

Jonathan R. J. Yong and Mohammadali Foroozandeh*

Chemistry Research Laboratory, Department of Chemistry, University of Oxford, 12 Mansfield Road, Oxford, OX1 3TA, U.K.

* mohammadali.foroozandeh@chem.ox.ac.uk

Contents

S 1	Raw data and reproducibility S4
S1.1	Datasets
S1.2	Figures
$\mathbf{S2}$	Routines S6
S2.1	poise_1d AU programme
S2.2	poise_1d_noapk AU programme S7
S2.3	poise_2d AU programme
$\mathbf{S3}$	1D ¹ H spectra for all samples S8
$\mathbf{S4}$	90° pulse width calibration S10
S4.1	Cost function
S4.2	Reference grid search
S4.3	Optimisation
S4.4	Using POISE under automation
S4.4	4.1 AU programme
S4.4	4.2 Python script
S 5	Franst angle S16
SJ S5 1	Cost function S16
55.1 S5.2	$\begin{array}{c} \text{Cost function} & \dots & $
55.2 S5.2	$\begin{array}{c} \text{Inconclusion} \\ \text{Optimization} \\ \text{S16} \end{array}$
55.5 S5.4	Variable fin angle pulse acquire pulse programme
55.4	Variable inplangie puise-acquire puise programme
$\mathbf{S6}$	NOE mixing time S19
S6.1	Cost function
S6.2	Reference grid search
S6.3	Optimisation
S6.4	1D NOE pulse programme
S 7	Inversion-recovery S24
S7.1	Cost function
S7.2	2D inversion-recovery
S7.3	Optimisation
S7.4	1D inversion-recovery pulse programme
S 8	ASAP-HSQC S27
S8.1	Cost function
S8.2	Reference grid search
S8.3	Optimisation
S8.4	ASAP-HSQC pulse programme
S 9	EPSI S33
S9.1	Cost function

S9.2	Reference grid search
S9.3	Optimisation
S9.4	Excitation–EPSI pulse programme
S10	PSYCHE S38
S10.1	Cost function
S10.2	Flip angle reference grid search
S10.3	Optimisations
S10.4	AU programme used for optimisations
S10.5	Python script for generating double saltire
S10.6	1D spin echo pulse programme
S 11	1D NOESY-based solvent suppression S48
S11.1	Cost function
S11.2	Transmitter offset reference grid search
S11.3	Optimisations
S11.4	1D NOESY pulse programme
$\mathbf{S12}$	Diffusion NMR S55
S12.1	Sequential optimisation strategy
S12.2	Python script for automation
S12.3	Sequential optimisation details
S12.4	Simultaneous optimisation strategy
S12.5	Oneshot DOSY pulse programmes
S12	5.1 2D Oneshot $\ldots \ldots \ldots$
S12	5.2 1D Oneshot
$\mathbf{S13}$	Documentation S69

S13**Documentation**

Note: The Supporting Information sections mainly contain summaries of the results presented in this paper, as well as technical details for reproducing these results. It does not explain how to set up POISE, or to how to create your own optimisation routines. For step-by-step instructions on these, please refer to the documentation. A copy of it is attached at the end of the Supporting Information, but the most updated version will be accessible at https://foroozandehgroup.github.io/nmr poise (a PDF copy is located at https://foroozandehgroup.github.io/nmrpoise/poise.pdf).

S1 Raw data and reproducibility

All raw data used for this paper, as well as the Python scripts to generate the plots in the main paper and SI, can be downloaded from Zenodo at the DOI 10.5281/zenodo.4698423.

S1.1 Datasets

The datasets that were used for each section are tabulated in Table S1.

Section	Description	Folder name	Expnos
90° pulse width	Reference grid search	P1_Scan	1 - 161
	Optimisation (initial value 48 µs)	P1_Opt	4
	Optimisation (initial value $43 \mu s$)	P1_Opt	5
	Optimisation (initial value $53 \mu s$)	P1_Opt	6
Ernst angle	2D inversion-recovery	$\texttt{Invrec_Opt}$	2
	Optimisation $(6 \text{ to } 8 \text{ ppm})$	$\texttt{Ernst_Opt}$	1001
	Optimisation (6.79 ppm only)	$\texttt{Ernst_Opt}$	1002
NOE mixing time	Reference grid search	NOESY_OptScan	101 - 159
	Optimisation	NOESY_OptScan	3
Inversion-recovery	2D inversion-recovery	Invrec_Opt	2
	Optimisation $(6 \text{ to } 8 \text{ ppm})$	$\texttt{Invrec_Opt}$	3
	Optimisation (7.08 ppm only)	$\texttt{Invrec_Opt}$	4
ASAP-HSQC	Reference grid search	$\tt HSQC_OptScan$	101 - 133
	Optimisation	$\tt HSQC_OptScan$	2
EPSI	Reference grid search	EPSI_Scan	95 - 116
	Optimisation	\texttt{EPSI}_{Opt}	1
	Ultrafast 2D TOCSY	$\texttt{EPSI_Opt}$	10 - 11
PSYCHE	Flip angle reference grid search	PSYCHE_FAScan	301 - 331
	Optimisation (1 parameter)	$PSYCHE_Opt$	101
	Optimisation (2 parameters)	$PSYCHE_Opt$	102
	Optimisation (3 parameters)	$PSYCHE_Opt$	103
	Optimisation (4 parameters)	$PSYCHE_Opt$	104
	Pseudo-2D homodecoupled spectra	$PSYCHE_Opt$	1400 - 1404
Solvent suppression	Offset reference grid search	Solvsupp_O1Scan	101 - 141
	Optimisation (1 parameter)	Solvsupp_Opt	11
	Optimisation (2 parameters)	Solvsupp_Opt	14
	Optimisation (3 parameters)	Solvsupp_Opt	15
	Optimisation (4 parameters)	$\texttt{Solvsupp_Opt}$	16
	Representative optimised spectra	$\texttt{Solvsupp}_\texttt{Opt}$	31, 35 37
DOSY	Optimisation (sequential)	$DOSY_Opt$	101 - 102
	Optimised Oneshot DOSY	$DOSY_Opt$	$1,\ 10011021$
	Optimisation (simultaneous)	DOSY_Opt	201 - 202

Table S1: Datasets used for each section of the Supporting Information.

All of the POISE optimisation data was logged in poise.log files, which can be found in the

corresponding raw data folders. These log files detail the optimisation trajectory and the cost function at each point, and can be parsed using the parse_log() function (please see the documentation for more information). However, some of these date from an older development version of POISE in which less information was recorded, so some details (e.g. of the AU programme used for acquisition) may be blank.

Full details of acquisition and processing parameters for each optimisation can be found in these folders. For this reason, except for a few key values, they are not tabulated in the main text or here.

S1.2 Figures

All figures are done in Python 3 and use the **penguins** package (written by J.R.J.Y.) to parse and plot NMR data. As of the time of writing, this package is still in development; consequently, the interface and functionality may change in a backward-incompatible manner. In order to ensure full reproducibility, please install version 0.4.1 of **penguins**, which was released specifically for this paper. This can be done using the command

pip install penguins=0.4.1

After downloading the raw data, the NMR datasets themselves should be found in a directory called datasets, and the scripts in a directory called figures. The figures themselves are provided in this folder, but can be regenerated at any point in time by re-running the scripts (the scripts also provide a "trace" of which data is being used to plot the data).

The scripts will be executable as long as this directory structure is preserved, as they use a relative path to locate the datasets. If the directories are moved for any reason, the environment variable POISE_DATA should be set such that it points to the datasets directory. For example, in bash and similar shells, this can be done using (replace the path accordingly)

export POISE_DATA=/path/to/nmrpoise_data/datasets

S2 Routines

In POISE, routines are stored in the human-readable JSON format, and are thus provided in the sections which follow. These JSON files should be placed in the TopSpin directory .../exp/stan/nmr/ py/user/poise_backend/routines. Each routine consists of several entries:

- name: The name of the routine, which can be arbitrarily chosen (although it should match the filename: thus a routine with name my_routine should be stored in my_routine.json).
- pars: A list of string parameters which are to be optimised. These correspond exactly to parameter names in TopSpin.
- 1b: A list of numbers acting as the lower bounds for each parameter. Units are given in TopSpin's native units, i.e. pulses are in microseconds, delays in seconds, etc.
- ub: A list of numbers acting as the upper bounds for each parameter.
- init: A list of numbers acting as the initial values for each parameter.
- tol: A list of numbers acting as the tolerances for each parameter.
- cf: The name of the cost function. A function with this name must be defined in either of the files poise_backend/costfunctions.py or poise_backend/costfunctions_user.py. For a list of pre-installed cost functions and descriptions, please refer to the POISE user documentation.
- au: The name of the AU programme used for acquisition and processing. This can be left blank (see below).

If no AU programme is specified by the user upon setup, the **poise_1d** or **poise_2d** AU programmes will automatically be used depending on the dimensionality of the dataset. These default AU programmes are used for many of the optimisations described in the *Supporting Information*, and are detailed below.

All the routines in the *Supporting Information* can be found in the .../exp/stan/nmr/py/user/ poise_backend/example_routines directory, and can be directly copied to the routines directory if desired.

S2.1 poise_1d AU programme

```
ZG // acquire
EFP // exponential apodisation, Fourier transform
APBK // automatic phase and baseline correction
QUIT
```

S2.2 poise_1d_noapk AU programme

This AU programme is not used by default, but comes installed together with POISE, used in some of the optimisations described below.

```
ZG // acquire
EFP // exponential apodisation, Fourier transform
XCMD("abs n") // automatic baseline correction
QUIT
```

S2.3 poise_2d AU programme

```
ZG // acquire
XFB // 2D Fourier transform
XCMD("apk2d") // automatic phase correction
ABS2 // f2 baseline correction
QUIT
```

S3 1D ¹H spectra for all samples

As a reference, we provide here the ¹H spectra of all samples used in this work. These can also be found in the raw data (Section S1).



Figure S1: ¹H spectrum of 45 mM ferulic acid in DMSO- d_6 (500 MHz resonance frequency).



Figure S2: ¹H spectrum of 120 mM 3-fluorophenylboronic acid in DMSO- d_6 (700 MHz resonance frequency).



Figure S3: ¹H spectrum of 45 mM and rographolide in DMSO- d_6 (600 MHz resonance frequency).



Figure S4: ¹H spectrum of rodent urine in D_2O (400 MHz resonance frequency). The water peak has been suppressed using the 1D NOESY / presaturation sequence.

S4 90° pulse width calibration

All spectra in this section were acquired on a 600 MHz Bruker AVIII spectrometer, equipped with a Prodigy N₂ broadband cryoprobe with a nominal maximum z-gradient strength of $66 \,\mathrm{G \, cm^{-1}}$. The sample used was $45 \,\mathrm{mM}$ ferulic acid in DMSO- d_6 .

S4.1 Cost function

The cost function used in this section, minabsint, comes pre-installed with POISE. It can be mathematically expressed as

$$f = \sum_{i} \sqrt{\mathbf{S}_{\text{real},i}^2 + \mathbf{S}_{\text{imag},i}^2} \tag{1}$$

where \mathbf{S}_{real} and \mathbf{S}_{imag} are vectors of length SI corresponding to the real and imaginary parts of the spectrum respectively. The elements \mathbf{S}_i are the intensities of the spectrum at each point along the spectral width.

S4.2 Reference grid search

The grid search, as used in the builtin **popt**, is a highly time-consuming method of locating the optimum. It is used here only as a reference, i.e. to prove that the faster optimisation algorithms are indeed finding the true optimum.



Figure S5: Dependence of the cost function on the pulse width (p1). The minimum value of the cost function is 3.34×10^7 and occurs at p1 = 48.3 µs.

S4.3 Optimisation

Table S2: 360° pulse width optimisations run with an initial point of 48 µs. popt optimisations were run between values of 40 µs and 56 µs, with a linear increment of 0.4 µs. The POISE routine is: {"name": "p1cal", "pars": ["p1"], "lb": [40.0], "ub": [56.0], "init": [48.0], "tol": [0.2], "cf": "minabsint", "au": "poise_1d"}. The pulse programme used was the Bruker standard zg. Key acquisition parameters: NS=1, DS=0, SW=12.0, TD=16384, D1=1.5. n_{fev} refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm.

Entry	Method	Description	Optimum found (μs)	$n_{\rm fev}$	Time taken (s)
1	popt	MAGMIN cost function	48.40	41	299
2	pulsecal	—	46.64	—	37
3a	POISE	Nelder-Mead	48.38	10	76 - 79
3b		MDS	48.38	10	77 - 80
3c		BOBYQA	48.29-48.41	6 - 7	46 - 54

The BOBYQA algorithm takes the *exact value* of the cost function into consideration when deciding what points to sample. On the other hand, the simplex-based methods only consider how the various points are ordered according to the value of the cost function at each point. This is partly why BOBYQA tends to converge faster than the other two algorithms, which in turn explains why we have set it as the default algorithm for POISE; it also accounts for the slight variability in the optimum found (note that this range is still well within the specified tolerance of $0.2 \,\mu$ s, so it does not reflect a failing). However, in edge cases when the cost function is extremely unbalanced or noisy, we find that BOBYQA can be unduly biased by the resulting cost function values (this is not the case here). We find that the simplex-based NM and MDS methods, although slower on average, can be slightly more reliable in such situations.

The next two tables describe optimisations starting from an inaccurate initial point.

Table S3: Pulse width optimisations with an initial point of 43 µs. The POISE routine is the same as in Table S2, except with "init": [43.0]. All other acquisition parameters are identical. The results are aggregated from five separate runs for each optimisation algorithm.

Entry	Method	Description	Optimum found (μs)	$n_{\rm fev}$	Time taken (s)
1a	POISE	Nelder-Mead	48.38	14	109-114
1b		MDS	48.38	14	108 - 112
1c		BOBYQA	48.27 - 48.33	9	70

 $\frac{1}{1a} \begin{array}{c} \text{POISE} & \text{Nelder-Mead} \end{array} \begin{array}{c} \text{Algorithm.} \\ \hline \text{Method} \\ \hline 1a \end{array} \begin{array}{c} \text{POISE} & \text{Nelder-Mead} \end{array} \begin{array}{c} \text{Algorithm.} \\ \hline 48.38 \end{array} \begin{array}{c} n_{\text{fev}} \\ \hline 14 \\ \hline 110-114 \end{array} \end{array}$

48.25 - 48.38

48.26 - 48.33

16

9

123 - 126

69 - 70

Table S4: Pulse width optimisations with an initial point of 53 µs. The POISE routine is the same as in Table S2, except with "init": [53.0]. All other acquisittion parameters are identical. The results are aggregated from five separate runs for each optimisation algorithm.

As can be seen, POISE is capable of "recovering" from a bad start. Although this requires	more
time, even in the worst-case scenario, the time taken is still less than half of that needed f	or the
full grid search (299 seconds, cf. Entry 1, Table S2).	

S4.4 Using POISE under automation

1b

1c

POISE provides a command-line interface which allows it to be called from within AU and Python scripts once the appropriate routines and cost functions have been set up. Thus, for example, a typical optimisation 'procedure' would involve the following steps:

1. Create a new dataset specially for the optimisation.

MDS

BOBYQA

- 2. Read in a parameter set, where the appropriate pulse programme for optimisation, as well as all other unoptimised parameters, have been chosen.
- 3. Run the POISE optimisation with poise <routine_name> <options>. (For the available options please see the user documentation.)
- 4. Transfer the optimised parameter values back to the original dataset, and run the optimised experiment.

All of this can be easily automated using either AU or Python scripts. Here we provide examples of both AU and Python scripts, called **poisecal**, which perform a similar role to the existing **pulsecal** AU programme in TopSpin. Running either of these will carry out a POISE optimisation in expno 99999, then set the optimised value of **p1** to the currently active dataset. (These examples are also included in the documentation with extra explanation.) Please note that both of these can be downloaded from GitHub; the links are near the top of the source code.

In general, the AU programme should be preferred as it contains some extra functionality (namely the option to dynamically specify the region of the spectrum to optimise on: please see the header of the AU programme for more information).

S4.4.1 AU programme

```
* pulsecal: optimises the value of p1 in expno 99999, then sets the optimised
 * value to the current dataset and shows a message.
 * To optimise on only one particular part of the spectrum, either:
 *
   1) Use the command `dpl` to specify a particular region, then run
       `poisecal`; or
 *
   2) use the command
          poisecal LEFT RIGHT
 *
 *
      where LEFT and RIGHT are the edges of the desired spectral window in
 *
      ppm. The order in which these are specified does not matter, so
 *
       `poisecal 1 4` and `poisecal 4 1` behave exactly the same way.
 * To optimise only one end of a spectrum (e.g. above 4 ppm only), use option 2
 * and pass a large value as the other argument, e.g. `poisecal 4 1000`.
 * SPDX-License-Identifier: GPL-3.0-or-later
 GETCURDATA
int old_expno = expno;
double f1p, f2p;
// Argument parsing
if (i_argc == 4) { // argv[0] is always the programme name, argv[1] is 'exec'
   double d1, d2;
   char *cp;
   d1 = strtod(i_argv[2], &cp);
   if (*cp) STOPMSG("Invalid arguments. Usage: poisecal [left right]")
   d2 = strtod(i_argv[3], &cp);
   if (*cp) STOPMSG("Invalid arguments. Usage: poisecal [left right]")
   f1p = d1 > d2 ? d1 : d2;
   f2p = d1 > d2 ? d2 : d1;
}
else if (i_argc == 2) {
   FETCHPAR("F1P", &f1p)
   FETCHPAR("F2P", &f2p)
}
else STOPMSG("Invalid arguments. Usage: poisecal [left right]")
// Use EXPNO 99999 in the current folder for optimisation.
DATASET(name, 99999, procno, disk, user)
// Set some key parameters. Notice that these lines can be substantially cut
// if an appropriate parameter set is set up beforehand.
RPAR("PROTON", "all")
GETPROSOL
STOREPAR("PULPROG", "zg")
STOREPAR("NS", 1)
STOREPAR("DS", 0)
STOREPAR("D 1", 1.0)
STOREPAR("RG", 1)
```

```
STOREPAR("F1P", f1p)
STOREPAR("F2P", f2p)
STOREPARS("F1P", f1p)
STOREPARS("F2P", f2p)
// Run optimisation (uses BOBYQA by default).
XCMD("sendgui xpy poise p1cal -q")
// POISE stores the optimised value in p1 after it's done. We can retrieve it
// here. Don't try to get the *status* parameter, since that is not the
// optimised value (it is the value used for the last function evaluation!)
float plopt;
FETCHPAR("P 1", &plopt)
plopt = plopt/4;
// Move back to old dataset and set p1 to optimised value.
DATASET(name, old_expno, procno, disk, user)
VIEWDATA_SAMEWIN // not strictly necessary, just re-focuses the original spectrum
STOREPAR("P 1", plopt)
Proc_err(INFO_OPT, "Optimised value of p1: %.3f", p1opt);
// (Optional) Run acquisition.
// ZG
ΩUIT
// vim: ft=c
```

S4.4.2 Python script

""" poisecalpy.py

Download from: https://git.io/JUHOU

Python script which automates p1 optimisation using POISE. Behaves similarly to pulsecal: optimises the value of p1 in expno 99999, then sets the optimised value to the current dataset and shows a message.

This does not have the command-line argument functionality of the poisecal AU programme, and tends to suffer from some edge-case bugs which are not present in poisecal. This is generally a result of TopSpin's Python programmes being more buggy than AU programmes. Nonetheless, as a result of this, we recommend using the poisecal AU programme instead of this (which is retained mainly to serve as a illustrative example in the documentation).

```
SPDX-License-Identifier: GPL-3.0-or-later
```

```
# Read in F1P and F2P from current dataset.
f1p = GETPAR("F1P")
f2p = GETPAR("F2P")
# Use EXPNO 99999 in the current folder for optimisation.
old_dataset = CURDATA()
opt_dataset = CURDATA()
opt_dataset[1] = "99999"
```

```
# Create new dataset and move to it.
NEWDATASET(opt_dataset, None, "PROTON")
RE(opt_dataset)
# Set some key parameters. Notice that these lines can be cut if an appropriate
# parameter set is set up beforehand (and loaded using NEWDATASET()).
XCMD("getprosol")
PUTPAR("PULPROG", "zg")
PUTPAR("NS", "1")
PUTPAR("DS", "0")
PUTPAR("D 1", "1")
PUTPAR("RG", "1")
XCMD("s f1p {}".format(f1p)) # PUTPAR("status F1P") doesn't work.
XCMD("s f2p {}".format(f2p)) # Even though the documentation says it should.
# Run optimisation (uses BOBYQA by default).
XCMD("poise p1cal -q")
# POISE stores the optimised value in p1 after it's done. We can retrieve it
# here. Don't try to get the *status* parameter, since that is not the
# optimised value (it is the value used for the last function evaluation!)
plopt = float(GETPAR("P 1"))/4
# Move back to old dataset and set p1 to optimised value.
RE(old_dataset)
PUTPAR("P 1", str(plopt))
ERRMSG("Optimised value of p1: {:.3f}".format(p1opt))
# A TopSpin quirk: using MSG() will block subsequent commands until user hits
# "OK". You can use ERRMSG(), as is done here. There are other ways around
# this, see MSG_nonmodal() in the POISE frontend script.
# (Optional) Run acquisition.
# ZG()
```

S5 Ernst angle

All spectra in this section were acquired on a 500 MHz Bruker AVIII spectrometer, equipped with a triple-resonance TBO probe with a nominal maximum z-gradient strength of $50 \,\mathrm{G \, cm^{-1}}$. The sample used was $45 \,\mathrm{mM}$ ferulic acid in DMSO- d_6 .

S5.1 Cost function

The cost function used in this section is **maxrealint**, which corresponds simply to a summation of all peaks in the real part of the spectrum:

$$f = -\sum_{i} \mathbf{S}_{\mathrm{real},i} \tag{2}$$

The minus sign is necessary because optimisation algorithms always seek to *minimise* the cost function. Thus, in this case, a larger peak intensity corresponds to a more negative cost function, i.e. a "better" spectrum.

S5.2 Theoretical optima via T_1

In lieu of a reference grid search, we sought to measure the longitudinal relaxation times (T_1) of each peak in the sample. The Ernst angle θ which maximises the sensitivity for a given repetition time τ_r can then be calculated from this via

$$\cos\theta = \exp\left(-\frac{\tau_{\rm r}}{T_1}\right) \tag{3}$$

thus giving us the theoretical optimal value, without having to go through a grid search. The T_1 values extracted from a gradient-enhanced inversion-recovery experiment, and the corresponding Ernst angles (for a repetition time of $\tau_r = 1.20 \text{ s}$, as used in the optimisation), are summarised in Table S5.

Table S5: T_1 (at 500 MHz) and Ernst angles for each peak in ferulic acid, calculated for a repetition time of 1.20 s.

Peak	¹ H chemical shift (ppm)	T_1 (s)	θ (°)
1	7.49	1.750	59.8
2	7.27	0.977	73.0
3	7.08	1.279	67.0
4	6.79	1.615	61.6
5	6.36	1.415	64.6
6	3.49	0.949	73.6

S5.3 Optimisation

It is possible to show, both theoretically and experimentally, that in order to maximise the sensitivity per unit time, the relaxation delay should be set to zero and the Ernst angle applied with τ_r simply equal to the acquisition time AQ. See, e.g., Concepts Magn. Reson. **1992**, 4 (2), 153–160 (DOI: 10.1002/cmr.1820040204) and J. Mol. Spectrosc. **1970**, 35 (2), 298–305 (DOI: 10.1016/0022-2852(70)90205-5).

Furthermore, in order to make an Ernst angle optimisation worthwhile, the repetition time must be on the order of T_1 . With longer repetition times on the order of $3T_1$ or so, the Ernst angle is very close to 90°, which transforms the issue into the 90° calibration already described in the previous section. In the event, we elected to use a repetition time of 1.20 s, which is on the order of the T_1 times shown in Table S5.

In the first optimisation, we set the region under optimisation to be just the aromatic and olefinic peaks between 6 and 8 ppm (i.e. peaks 1–5). Instead of changing the actual spectral window, this is most easily done using the dpl command in TopSpin; for more information, please refer to the user documentation.

Table S6: Ernst angle optimisations on all aromatic and olefinic peaks in ferulic acid (between 6 and 8 ppm). The POISE routine used here is: {"name": "ernst", "pars": ["cnst20"], "lb": [10.0], "ub": [90.0], "init": [30.0], "tol": [3.0], "cf": "maxrealint", "au": "poise_1d"}. The pulse programme used was zgvfa_jy (Section S5.4). Key acquisition parameters: NS=1, DS=4, D1=0, SW=14, TD=16384, AQ=1.17. n_{fev} refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm.

Entry	Algorithm	Optimum found (°)	$n_{\rm fev}$	Time taken (s)
1	Nelder-Mead	67.5 - 73.1	9–13	91 - 132
2	MDS	67.5 - 73.1	9	90 - 92
3	BOBYQA	70.1 - 70.7	7	70 - 71

Of course, the user may select the spectral bounds themselves in order to maximise the overall sensitivity for the entire spectrum, a particular region, or even a specific peak of interest. In the second case, we ran an optimisation on only the peak at 6.79 ppm.

Table S7: Ernst angle optimisations on the peak at 6.79 ppm in ferulic acid. The POISE routine and all acquisition parameters are the same as in Table S6, but the spectral region under optimisation was set to be 6.71–6.87 ppm via the dpl TopSpin command. n_{fev} refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm.

Entry	Algorithm	Optimum found (s)	$n_{\rm fev}$	Time taken (s)
1	Nelder-Mead	60.0-67.5	9–11	91–111
2	MDS	65.6 - 67.5	11	110 - 111
3	BOBYQA	60.0 - 65.2	6 - 7	59 - 71

S5.4 Variable flip angle pulse-acquire pulse programme

Note that the following pulse programme, which requires p1 to be the calibrated 90° pulse width and sets cnst20 to the flip angle, should only be used if the *actual flip angle* is of interest.

If the sole aim is to maximise sensitivity, then the standard Bruker **zg** pulse programme suffices. The routine should be modified to optimise **p1** directly instead of the flip angle itself. In essence, this means that we do not care what the flip angle is, we merely want the pulse width which maximises the sensitivity. This approach has an advantage in that **p1** does not need to be calibrated beforehand.

```
;zgvfa_jy
;download from: https://git.io/JnokS
;pulse-acquire with variable flip angle (set by cnst20)
;set p1 to be calibrated 90 degree pulse width
; Jonathan Yong - University of Oxford - 20 June 2021
;$CLASS=HighRes
;$DIM=1D
;$TYPE=
;$SUBTYPE=
;$COMMENT=
#include <Avance.incl>
"p20=p1*cnst20/90"
"acqt0=-p20*2/3.1416"
1 ze
2 30m
  d1
  p20 ph1
  go=2 ph31
  30m mc #0 to 2 F0(zd)
exit
ph1=0 2 2 0 1 3 3 1
ph31=0 2 2 0 1 3 3 1
;pl1 : f1 channel - power level for pulse (default)
;p1 : f1 channel - high power pulse
;d1 : relaxation delay; 1-5 * T1
;ns: 1 * n, total number of scans: NS * TDO
;cnst20: excitation flip angle in degrees
```

S6 NOE mixing time

All spectra in this section were acquired on a 700 MHz Bruker AVIII spectrometer, equipped with a TCI H/C/N cryoprobe with a nominal maximum z-gradient strength of $53 \,\mathrm{G \, cm^{-1}}$. The sample used was $120 \,\mathrm{mM}$ 3-fluorophenylboronic acid in DMSO- d_6 .

S6.1 Cost function

The cost function used for the optimisation here, named noe_1d , is pre-installed with POISE. It is only applicable to 1D selective NOE spectra. In short, it takes the processed spectrum, removes a region of 50 Hz around the selectively excited peak, and calculates the intensity of the remaining parts of the spectrum. Since we aim to *maximise* this intensity, the cost function must return the *negative* of this intensity. Further details are best explained through the comments in the source code, which is reproduced here:

```
def noe_1d():
   ......
   Measures the intensity of peaks in the spectrum, *except* for anything
   within 25 Hz of the selectively excited peak.
    ппп
    # SPOFFS2 plus 01 is the frequency of the selective pulse, in Hz.
    # We can get the values of both using getpar(), which automatically
    # converts numeric parameters to floats.
   f = getpar("SPOFFS2") + getpar("01")
    # Define the (rough) bandwidth of the selective pulse. We will ignore the
    # region of the spectrum that is bw Hz wide and centred on the selective
    # pulse.
   bw = 50.0
    # To get the spectrum *without* the ignored region, we split up the
    # spectrum into two "portions". We need to construct the 'bounds' parameter
    # for both portions.
    sfo1 = getpar("SF01")
    lowerhalf = f''(f + bw/2)/sfo1:...f' # frequencies (f + bw)/2 and above
    upperhalf = f"...{(f - bw/2)/sfo1:..3f}" # frequencies up to (f - bw)/2
    # We use get1d_real() to get both portions of the spectrum, passing the
    # bounds as appropriate. The intensities of the spectrum are obtained by
    # summing all points (np.sum).
    # Note that manually passing the bounds parameter will override the F1P/F2P
    # parameters if set by the user (this is probably a good thing).
   upper_integral = np.sum(get1d_real(bounds=upperhalf))
    lower_integral = np.sum(get1d_real(bounds=lowerhalf))
    # Lastly, we take the negative absolute value of the resulting sum, in
    # order to make sure that the resulting cost function is negative
    # regardless of how the NOE peaks are phased (positive or negative).
   return -abs(upper_integral + lower_integral)
```

S6.2 Reference grid search

The grid search was carried out using actual 2D NOESY spectra; this has the advantage of proving that the mixing time optimised using 1D spectra does in fact accurately represent the best mixing time for 2D spectra. The intensity of each crosspeak here refers to the intensity of the most negative point within a 0.12 ppm by 0.12 ppm box centred on the crosspeak. A broad minimum around a mixing time of 3–4 seconds can be seen.



Figure S6: Dependence of crosspeak intensity on the NOE mixing time (d8).

S6.3 Optimisation

Table S8: NOE mixing time optimisations. The POISE routine used here is: {"name": "1dnoe", "pars": ["d8"], "lb": [0.2], "ub": [6.0], "init": [0.5], "tol": [0.1], "cf": "noe_1d", "au": "poise_1d"}. The pulse programme used here is spfgsenoezs2.1_jy (Section S6.4). Key acquisition parameters: NS=2, DS=1, D1=1. n_{fev} refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm.

Entry	Algorithm	Optimum found (s)	$n_{\rm fev}$	Time taken (s)
1	Nelder-Mead	3.25 - 3.88	16 - 18	268 - 312
2	MDS	3.63 - 3.75	16 - 18	269 - 305
3	BOBYQA	3.38 - 3.80	6 - 10	88 - 162

S6.4 1D NOE pulse programme

The 1D NOE spectra were acquired with a modified pulse programme with two inversion pulses during the NOE mixing time to minimise artefacts arising from relaxation. (In principle, these artefacts do not contribute to the cost function as they ideally sum to zero due to their antiphase nature. However, it is almost always a good idea to try to minimise noise in the cost function where possible.) The pulse programme is attached below, and can also be found in the main POISE repository. For most compounds, the existing family of Bruker pulse sequences selnogpzs should suffice.

```
;spfgsenoezs2.1_jy
;download from: https://git.io/JUEOP
;Tim Claridge & Jonathan Yong - University of Oxford
;1D NOESY using selective refocussing with two shaped pulses
;dipolar coupling may be due to noe or chemical exchange.
;H. Kessler, H. Oschkinat, C. Griesinger & W. Bermel,
    J. Magn. Reson. 70, 106 (1986)
;J. Stonehouse, P. Adell, J. Keeler & A.J. Shaka, J. Am. Chem. Soc 116,
   6037 (1994)
;K. Stott, J. Stonehouse, J. Keeler, T.L. Hwang & A.J. Shaka,
   J. Am. Chem. Soc 117, 4199-4200 (1995)
;Tim C 2/2008 Using two nulling 180s and ZQC filter at end of mixing time
;Includes adiabatic pulse/gradient combination for ZQC suppression
;Thrippleton & Keeler, Angew. Chem. Int. Ed., 2003, 42, 3938-3941.
;Tim C 10/2007 updated to match Bruker format
;Tim C 11/2008 matches Bruker ZQS nomenclature
; Jon Y 21/07/2020 added CALC_SPOFFS acquisition flag
;$CLASS=HighRes
;$DIM=1D
;$TYPE=
;$SUBTYPE=
;$COMMENT=
#include <Avance.incl>
#include <Grad.incl>
"p2=p1*2"
"d21=d8*0.3-p16-d16-3u"
"d22=d8*0.5-p16*2-d16*2-6u"
"d23=d8*0.2-p16-d16*2-p32-53u"
#ifdef CALC_SPOFFS
"spoff2=bf1*(cnst21/1000000)-o1"
#else
#endif /*CALC_SPOFFS*/
1 ze
2 30m
  20u BLKGRAD
  d1
  50u UNBLKGRAD
```

(p1 ph1):f1 3u p16:gp1 d16 pl0:f1 3u p12:sp2:f1 ph2:r 3u p16:gp1 d16 pl1:f1 3u (p1 ph4):f1 d21 p16:gp3 d16 3u (p2 ph5):f1 Зu p16:gp3*-1 d16 d22 p16:gp4 d16 Зu (p2 ph5):f1 Зu p16:gp4*-1 d16 d23 20u pl0:f1 10u gron0 ;ZQC suppression (p32:sp29 ph7):f1 20u groff d16 pl1:f1 (p1 ph6):f1 go=2 ph31 30m mc #0 to 2 F0(zd) 20u BLKGRAD exit ph1=0 2 ph2=0 0 1 1 2 2 3 3 ph4=0 ph5=0 ph6=0 ph7=0 ph31=0 2 2 0 0 2 2 0 ;pl0 : 120dB ;pl1 : f1 channel - power level for pulse (default)

```
;sp2: f1 channel - shaped pulse
;sp29: f1 channel - shaped pulse (adiabatic sweep)
;p1 : f1 channel - 90 degree high power pulse
;p2 : f1 channel - 180 degree high power pulse
;p12: f1 channel - 180 degree shaped pulse
                                                          [1 msec]
;p16: homospoil/gradient pulse
;p32: f1 channel - 180 shaped adiabatic pulse
                 smoothed CHIRP
;d1 : relaxation delay; 1-5 * T1
;d8 : mixing time
;d16: delay for homospoil/gradient recovery
;d21: d8*0.3-p16-d16-3u
;d22: d8*0.5-p16*2-d16*2-6u
;d23: d8*0.2-p16-d16*2-p32-53u
;NS: 2 * n
;DS: 4
;cnst21: chemical shift for selective pulse (offset, in ppm)
; choose p12 according to desired selectivity
;the flip-angle is determined by the amplitude
;set O1 on resonance on the multiplet to be excited or use spoffs
;use gradient ratio:
                     gp 0 : gp 1 : gp 3 : gp 4
                       ~10 : 17 : 40 : 25
;
; for z-only gradients:
;gpz0: 10% or so
;gpz1: 17%
;gpz3: 40%
;gpz4: 25%
;use gradient files:
;gpnam1: SINE.100
;gpnam2: SINE.100
;gpnam3: SINE.100
;gpnam4: SINE.100
```

S7 Inversion-recovery

All spectra in this section were acquired on a 500 MHz Bruker AVIII spectrometer, equipped with a triple-resonance TBO probe with a nominal maximum z-gradient strength of $50 \,\mathrm{G \, cm^{-1}}$. The sample used was $45 \,\mathrm{mM}$ ferulic acid in DMSO- d_6 .

S7.1 Cost function

In this optimisation, we sought to make the real part of the spectrum as close to zero as possible: this indicates the "null" in the inversion-recovery experiment where the delay τ is exactly $T_1 \ln 2$. This is accomplished by the **zerorealint** cost function, which comes installed with POISE:

$$f = \sum_{i} |\mathbf{S}_{\text{real},i}| \tag{4}$$

S7.2 2D inversion-recovery

 T_1 for each peak was measured via a gradient-selected inversion-recovery experiment. The data presented here is the same as in Table S5, but additionally includes the theoretical optimum of $T_1 \ln 2$ for each peak.

Table S9: T_1 (at 500 MHz) and the corresponding theoretical null in the inversion-recovery profile, given by $T_1 \ln 2$, for each peak in ferulic acid.

Peak	$^{1}\mathrm{H}$ chemical shift (ppm)	T_1 (s)	$T_1 \ln 2$ (s)
1	7.49	1.750	1.213
2	7.27	0.977	0.677
3	7.08	1.279	0.887
4	6.79	1.615	1.119
5	6.36	1.415	0.981
6	3.49	0.949	0.658

S7.3 Optimisation

As previously mentioned, it is possible to use the dpl TopSpin command to restrict the spectral region that the cost function acts on. In the first optimisation, we set this region to be between 6 and 8 ppm. This provides an "averaged" value of $T_1 \ln 2$, which can be useful in e.g. choosing relaxation delays for 2D experiments, which invariably need to be a compromise between multiple peaks.

Table S10: Inversion-recovery optimisations on all aromatic and olefinic peaks in ferulic acid (between 6 and 8 ppm. The POISE routine used here is: {"name": "invrec", "pars": ["d27"], "lb": [0.35], "ub": [1.75], "init": [0.6], "tol": [0.01], "cf": "zerorealint", "au": "poise_1d"}. The pulse programme used was t1irgp1d_jy (Section S7.4). Key acquisition parameters: NS=1, DS=0, D1=5, AQ=3.6 (note that there is also an additional ca. 5 seconds between experiments due to pulse programme compilation, etc.). n_{fev} refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm.

Entry	Algorithm	Optimum found (s)	$n_{\rm fev}$	Time taken (s)
1	Nelder-Mead	0.938 – 0.969	14 - 16	204 - 235
2	MDS	0.956 – 0.975	16	233 - 235
3	BOBYQA	0.953 – 0.971	9 - 11	130 - 160

In the second optimisation we focused on only the peak at 7.08 ppm. This provides a relatively "fast" method of determining T_1 for one specific peak in a spectrum.

Table S11: Inversion-recovery optimisations on only the peak at 7.08 ppm. The POISE routine and all acquisition parameters are the same as in Table S10, but the spectral region under optimisation was set to be 7.02–7.15 ppm via the dpl TopSpin command. n_{fev} refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm.

Entry	Algorithm	Optimum found (s)	$n_{\rm fev}$	Time taken (s)
1	Nelder-Mead	0.863 - 0.875	14	202 - 205
2	MDS	0.863 - 0.869	14	203 - 204
3	BOBYQA	0.862 - 0.873	9–10	128 - 145

S7.4 1D inversion-recovery pulse programme

```
;tlirgp1d_jy
;
;download from: https://git.io/Jnok6
;
;T1 measurement using inversion recovery
;using broadband 1H inversion pulse and field gradients
;1D version with only one increment
;set inversion-recovery delay to d27
;
;Tim Claridge & Jonathan Yong - University of Oxford - 20 June 2021
;
;$CLASS=HighRes
;$DIM=1D
;$TYPE=
;$SUBTYPE=
;$COMMENT=
#include <Avance.incl>
#include <Delay.incl>
```

```
#include <Grad.incl>
"acqt0=-p1*2/3.1416"
1 ze
2 50m BLKGRAD
  d1
  "TAU=(d27-p19-d16-4u)"
  if "TAU > 0.96"
  {
  "DELTA=TAU-0.45"
  }
  50u UNBLKGRAD
  4u pl0:f1
  (p44:sp30 ph1):f1
  p19:gp1
  d16
  if "TAU > 0.96"
  {
  0.4s
  50m BLKGRAD
  DELTA pl1:f1
  }
  else
  {
  TAU pl1:f1
  }
 p1 ph2
  go=2 ph31
  50m mc #0 to 2 F0(zd)
50m BLKGRAD
exit
ph1=0 2
ph2=0 0 2 2 1 1 3 3
ph31=0 0 2 2 1 1 3 3
;pl1 : f1 channel - power level for pulse (default)
;p1 : f1 channel - 90 degree high power pulse
;d1 : relaxation delay; 1-5 * T1
;ns: 1
;ds: 0
;spnam30: 1H inversion [Bip720,50,20.1]
;p44: broadband inversion pulse [240 us]
;d27: inversion-recovery delay
```

S8 ASAP-HSQC

All spectra in this section were acquired on a 700 MHz Bruker AVIII spectrometer, equipped with a TCI H/C/N cryoprobe with a nominal maximum z-gradient strength of $53 \,\mathrm{G \, cm^{-1}}$. The sample used was 120 mM 3-fluorophenylboronic acid in DMSO- d_6 .

S8.1 Cost function

The cost function used here, asaphsqc, simply takes the f_2 projection of the acquired 2D spectrum and calculates the total intensity by summing up all points. Unlike the previous cost functions, it is not builtin: to enable this cost function, please uncomment the asaphsqc() function in $poise_backend/costfunctions.py$.

S8.2 Reference grid search

Crosspeak intensities here refer to the maximum peak height, i.e. the intensity of the highest point found within a $3 \text{ ppm}(f_1)$ by $0.12 \text{ ppm}(f_2)$ box centred on the crosspeak.



Figure S7: Dependence of crosspeak intensity on the INEPT delay $\Delta = (1/4 \cdot \texttt{cnst3})$.

The exact point of maximum intensity is not obvious (and not really useful in this case, as essentially *any* point near the maximum will give essentially the same result). However, broadly speaking, it can be seen that the intensity reaches a maximum at around 230 Hz, after which it plateaus off.

S8.3 Optimisation

Table S12: ASAP-HSQC INEPT delay optimisations. The POISE routine used was: {"name": "asaphsqc", "pars": ["cnst3"], "lb": [120.0], "ub": [280.0], "init": [150.0], "tol": [10.0], "cf": "asaphsqc", "au": "poise_2d"}. The pulse programme used here is asaphsqc_jy (Section S8.4). Key acquisition parameters: NS=1, DS=2, TD1=36, TD2=512, D1=0.1. n_{fev} refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm.

Entry	Algorithm	Optimum found	$n_{\rm fev}$	Time taken (s)
1	Nelder-Mead	250.0 - 256.3	8 - 9	169 - 195
2	MDS	237.5 - 243.8	8	171 - 179
3	BOBYQA	229.8 - 245.6	4 - 7	114 - 157

S8.4 ASAP-HSQC pulse programme

The pulse programme used here is almost identical to the published version but includes several minor changes for convenience.

```
;asaphsqc_jy
;download from: https://git.io/JUHOo
;modified from version by Luy et al. to include internal power calculations for
;DIPSI-2 mixing, as well as URPs
;order of gradients also switched round to correct frequencies
;added extra delay d30 before experiment starts
;Jonathan Yong - University of Oxford - 26 August 2020
;ASAP-HSQC
;2D H-1/X correlation via double inept transfer
;phase sensitive using Echo/Antiecho-TPPI gradient selection
;with decoupling during acquisition
;with DIPSI-2 mixing after acquisition for ASAP
;using shaped pulses on f1 and f2.
;David Schulze-Sunninghausen, Johanna Becker, Burkhard Luy
;"Rapid Heteronuclear Single Quantum Correlation NMR Spectra at Natural Abundance"
; J. Am. Chem. Soc. 2014, 136 (4), 1242-1245
#include <Avance.incl>
#include <Grad.incl>
#include <Delay.incl>
"cnst12=1/(1+sfo2/sfo1)"
"cnst13=(1-sfo2/sfo1)"
define list<gradient> EA1 = { 1 cnst12 }
define list<gradient> EA2 = { cnst13 1 }
"p2=p1*2"
"d0=3u"
"d11=30m"
"d16=200u"
```
```
"d20=1m"
"FACTOR1=(d9/(p6*115.112))/2"
"l1=FACTOR1*2" ; TOCSY loop counter
"d0=3u"
"in0=inf1/2"
"d4=1s/(4*cnst4)" ; 1/(4J)XH
"d3=1s/(4*cnst3)" ; reduced d4
"DELTA=0.5*(2*d3+p21-p1)"
"DELTA1=0.5*(2*d3-p21+p1)"
"DELTA4=d4-p16-d16-3u"
"DELTA3=p16+3u+d16+p21-p25+p1"
"DELTA2=0.5*(DELTA3+p16+d16+3u-2*d0)"
"cnst21=20000"
"spw21=plw1*(cnst21/(1000000/(4*p1)))*(cnst21/(1000000/(4*p1)))"
"spw22=spw21"
"cnst23=10000"
"spw23=plw2*(cnst23/(1000000/(4*p3)))*(cnst23/(1000000/(4*p3)))"
"spw24=spw23"
"spw25=spw23"
"spw26=spw23"
; calculate mixing length
"cnst9=p6*115.112*11/1000"
"cnst9=cnst9" ; so that it's displayed in ased
; calculate 90deg pulse duration & power level for mixing
"p6=(1000000/(4*cnst10))"
"plw10=plw1*(p1/p6)*(p1/p6)"
1 ze
  d30
 d11 pl12:f2
2 d1 do:f2
3 d20 UNBLKGRAD
 p16:gp2
  d16 pl10:f1
                                                 ;begin DIPSI2
4 p6*3.556 ph23
 p6*4.556 ph25
 p6*3.222 ph23
 p6*3.167 ph25
 p6*0.333 ph23
 p6*2.722 ph25
 p6*4.167 ph23
 p6*2.944 ph25
 p6*4.111 ph23
 p6*3.556 ph25
  p6*4.556 ph23
 p6*3.222 ph25
  p6*3.167 ph23
```

```
S29
```

p6*0.333 ph25 p6*2.722 ph23 p6*4.167 ph25 p6*2.944 ph23 p6*4.111 ph25 p6*3.556 ph25 p6*4.556 ph23 p6*3.222 ph25 p6*3.167 ph23 p6*0.333 ph25 p6*2.722 ph23 p6*4.167 ph25 p6*2.944 ph23 p6*4.111 ph25 p6*3.556 ph23 p6*4.556 ph25 p6*3.222 ph23 p6*3.167 ph25 p6*0.333 ph23 p6*2.722 ph25 p6*4.167 ph23 p6*2.944 ph25 p6*4.111 ph23 lo to 4 times 11 ;end DIPSI2 p16:gp3 d16 (p21:sp21 ph1):f1 DELTA (center (p22:sp22 ph1):f1 (p22:sp23 ph6):f2) DELTA1 pl1:f1 (ralign (p1 ph1):f1 (p21:sp24 ph9):f2) DELTA2 d0 (p22:sp22 ph5):f1 d0 DELTA2 pl1:f1 (lalign (p1 ph7):f1 (p25:sp25 ph10):f2) DELTA3 (p22:sp22 ph1):f1 Зu p16:gp1*EA1 d16 (p21:sp26 ph10):f2 d4 (center (p22:sp22 ph1):f1 (p22:sp23 ph6):f2) Зu p16:gp1*EA2 d16 pl12:f2

```
DELTA4 BLKGRAD
  go=2 ph31 cpd2:f2
  d1 do:f2 mc #0 to 2
  F1EA(igrad EA1 & igrad EA2, id0 & ip9*2 & ip6*2 & ip31*2)
exit
ph1=0
ph2=1
ph5=0 0 2 2
ph6=0
ph7=3
ph8=3
ph9=2 0
ph10=0 0 0 0 2 2 2 2
ph23=3
ph25=1
ph31=0 2 0 2 2 0 2 0
;pl0 : OW
;pl1 : f1 channel - power level for pulse (default)
;pl2 : f2 channel - power level for pulse (default)
;pl10 : f1 channel - power level for mixing sequence (DIPSI-2)
;pl12: f2 channel - power level for CPD/BB decoupling
;p1 : f1 channel - 90 degree high power pulse
;p2 : f1 channel - 180 degree high power pulse
;p6 : f1 channel - 90 degree pulse for mixing sequence (DIPSI-2)
;p21 : 550u excitation pulses
;p22 : 600u UR and Inversion pulses
;p25 : 1100u 13C refocusing pulse
;p16: homospoil/gradient pulse
;spnam21 : jc01_BEBOP_zx_550u_BW10_RF20_pm20_Hc0.99997119.pul
;spnam22 : jc02_BURBOP_x_600u_BW10_RF20_pm20_matched.pul
;spnam23 : jc03_BIBOP_600u_BW37.5_RF10_pm5_matched.pul
;spnam24 : jc05_BEBOP_zy_550u_BW37.5_RF10_pm5_matched.pul
;spnam25 : jc07_BURBOP_y_1100u_BW37.5_RF10_pm5_Hc0.999876221.pul
;spnam26 : jc09_BEBOP_-yz_550u_BW37.5_RF10_pm5_matched.pul
;sp21 : 20 kHz Rf Amplitude
;sp22 : 20 kHz Rf Amplitude
;sp23 : 10 kHz Rf Amplitude
;sp24 : 10 kHz Rf Amplitude
;sp25 : 10 kHz Rf Amplitude
;sp26 : 10 kHz Rf Amplitude
;d0 : incremented delay (2D)
                                             [3 usec]
;d1 : relaxation delay; 1-5 * T1
;d4 : 1/(4J)XH
;d3 : reduced d4
;d9 : TOCSY mixing time
;d11: delay for disk I/O
                                             [30 msec]
;d16: delay for homospoil/gradient recovery [200ms]
;d30: delay before experiment starts (use if running extended multizgs)
; cnst4: = J(XH)
```

;cnst3: reduced cnst4 ('ernst angle') ; inf1: 1/SW(X) = 2 * DW(X);in0: 1/(2 * SW(X)) = DW(X);nd0: 2 ;NS: 1 * n ;DS: >= 16 ;td1: number of experiments ;FnMODE: echo-antiecho ;cpd2: decoupling according to sequence defined by cpdprg2 ;pcpd2: f2 channel - 90 degree pulse for decoupling sequence ;use gradient ratio: gp 1 : gp 2 : gp 3 80 : 33 : 43 ; ;for z-only gradients: ;gpz1: 80% ;gpz2: 33% ; gradient before mixing ;gpz3: 43% ; gradient after mixing ;use gradient files: ;gpnam1: SMSQ10.100 ;gpnam2: SMSQ10.100 ;gpnam3: SMSQ10.100

S9 EPSI

All spectra in this section were acquired on a 600 MHz Bruker AVIII spectrometer, equipped with a Prodigy N₂ broadband cryoprobe with a nominal maximum z-gradient strength of $66 \,\mathrm{G \, cm^{-1}}$. The sample used was $45 \,\mathrm{mM}$ ferulic acid in DMSO- d_6 .

S9.1 Cost function

The cost function used for EPSI optimisations (called epsi_gradient_drift) comes pre-installed with POISE. It is more involved than usual, as the entirety of the EPSI processing must be carried out within the cost function (POISE itself only providing a function to get the FID). Briefly, this processing consists of the following steps:

- Removal of the "group delay" (by circularly shifting the FID).
- Reshaping of the 1D FID into a 2D (k, t_2) matrix.
- Removing parts of the FID acquired under negative gradients, as well as during the delay between gradients.
- Absolute-value processing.
- Apodisation.
- Removing any rows (i.e. any values of t_2) for which the maximum intensity is below a certain threshold. This ensures that only parts of the (k, t_2) matrix with strong signal are used for the calculation of the slope.

The slope of k versus t_2 can then be calculated by finding, for each row, the value of k where the intensity is maximum, and performing a least-squares fitting. Technically k/k_{max} , not k, is used instead of k, but k_{max} is merely a constant factor which has no influence on the results.

For more details, please refer to the source code which is attached below:

```
def epsi_gradient_drift():
    .....
    Calculates the amount of 'gradient drift' seen in a 1D EPSI acquisition, as
    reflected by the position of the 'echo' moving over time. This can be
    caused by imbalanced positive and negative gradients.
    .....
    from numpy.polynomial.polynomial import Polynomial
    # --- Calculate key parameters ---
    fid = get1d_fid(remove_grpdly=True)
    # Number of complex points in k-space, i.e. number of points per EPSI
    # gradient. TD2 is total number of real & imag points in FID. L3 is number
    # of EPSI loops (one loop includes both pos + neg gradient). So TD2 / (2 *
    # L3) is the number of points in one EPSI gradient (i.e. only positive
    # gradient). The extra factor of 2 is because we're interested only in
    # complex points.
    td_k = int(getpar("TD") / (2 * 2 * getpar("L3")))
```

```
# Time between consecutive EPSI positive gradients, in seconds
td_t2 = int(getpar("L3"))
dw_eff = getpar("AQ") / td_t2
# --- Perform EPSI processing -----
# Reshape into 2D matrix
ser = fid.reshape((-1, td_k))
# Discard the part acquired with negative gradients
ser = ser[0::2, :]
# Discard any part of the spectrum that was not acquired during an EPSI
# gradient, i.e. if the delay D6 was nonzero.
if getpar("D6") > 0:
   ineligible_epsi_points = int(1e6 * getpar("D6") / (getpar("DW") * 2))
   td_k = td_k - ineligible_epsi_points
   ser = ser[:, :td_k]
# --- Apodisation -----
                                 _____
# Along k-dimension (Hamming window)
alpha_0 = 0.54
ks = np.linspace(0, 1, td_k)
k_winfunc = (alpha_0)
            - (1 - alpha_0) * np.cos(2 * np.pi * np.linspace(0, 1, td_k)))
ser = ser * k_winfunc[np.newaxis, :]
# Along direct dimension
t2_winfunc = np.sin(np.pi * np.linspace(0, 1, td_t2))
ser = ser * t2_winfunc[:, np.newaxis]
abs_ser = np.abs(ser)
# Calculate k- and t2-axes
t2_values = np.arange(td_t2) * dw_eff
k_values = np.linspace(-0.5, 0.5, td_k)
# Drop all rows (i.e. all values of t2) for which the maximum is less than
# 20% of the overall maximum.
threshold_amp = 0.2 * np.max(abs_ser)
maxima_along_rows = np.max(abs_ser, axis=1)
indices_to_use = np.nonzero(maxima_along_rows > threshold_amp)
# Locate the maxima
maximal_indices_along_k = np.argmax(abs_ser, axis=1)
# Get the value of k at each maximum
k_max_values = k_values[maximal_indices_along_k]
poly = Polynomial.fit(x=t2_values[indices_to_use],
                     y=k_max_values[indices_to_use], deg=1)
return np.abs(poly.coef[1])
```

S9.2 Reference grid search

The parameter optimised here is cnst16, which represents the ratio of nominal negative gradient amplitude to nominal positive gradient amplitude during the EPSI acquisition. In the main text, this is referred to as α . On an ideal instrument, the optimised value of cnst16 would simply be 1.



Figure S8: Dependence of the cost function on the value of cnst16. The minimum occurs at around cnst16 = 1.0004.

S9.3 Optimisation

Table S13: Results of EPSI gradient imbalance optimisations. The POISE routine used here is: {"name": "epsi", "pars": ["cnst16"], "lb": [0.99], "ub": [1.01], "init": [1.0], "tol": [0.0001], "cf": "epsi_gradient_drift", "au": "poise_1d"}. The pulse programme used was zg_epsi.2_jy (Section S9.4). Key acquisition parameters: NS=1, DS=0, L4=0, D1=0.5, DW=0.7 (note units of DW are µs), AQ=0.092, GPZ15=35.644 (the maximum gradient amplitude on this machine is 67.5 G cm⁻¹). To avoid hardware damage AQ should always be kept short, ideally under 100 ms. Note that such a short value of D1 is only possible because there is always an extra delay between successive experiments required for pulse programme compilation, etc. A short recovery delay should not be used when DS > 0 or NS > 1. n_{fev} refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm.

Entry	Algorithm	Optimum found	$n_{\rm fev}$	Time taken (s)
1	Nelder–Mead	1.000375	10	49 - 52
2	MDS	1.000375	10	49 - 50
3	BOBYQA	1.000383 - 1.000388	7	34 - 35

S9.4 Excitation–EPSI pulse programme

```
;zg_epsi.2_jy
;
```

```
;download from: https://git.io/JYGEa
;
;90 pulse followed by EPSI acquisition, with prephasing gradient
;modified from uftocsy, Giraudeau et al.
;https://madoc.univ-nantes.fr/course/view.php?id=24710&section=1
;IMPORTANT: set p15 + d6 = DW x TD(F2)/(2xL3)
;$CLASS=HighRes
;$DIM=1D
;$TYPE=
;$SUBTYPE=
;$COMMENT=
#include <Avance.incl>
#include <Grad.incl>
#include <De.incl>
#include <Delay.incl>
"d0=3u"
"d11=30m"
"p15=(td*dw/(2*13))-d6"
"cnst17=cnst16" ; so that cnst16 shows up in ased
1 ze
  "cnst17=cnst17" ; so that cnst16 shows up in ased
  d11
  30m zd
  ; dummy scans
4 30m
  10u
  4u
  d1
  4u
  100u UNBLKGRAD
  20u pl1:f1
  (p1 ph1):f1
  10u
  p15:gp25
  10u
5 p15:gp15
  d6
  p15:gp16
  d6
  lo to 5 times 13
  p27:gp27
  100u BLKGRAD
  lo to 4 times 14
  ; end of dummy scans
2 30m
  10u
```

4u d1 4u 100u UNBLKGRAD 20u pl1:f1 ; excitation (p1 ph1):f1 ; prephasing gradient 10u p15:gp25 10u ; EPSI acquisition ACQ_START(ph30,ph31) 1u DWELL_GEN:f1 3 p15:gp15 d6 p15:gp16*cnst16 d6 lo to 3 times 13 100u BLKGRAD rcyc=2 100u UNBLKGRAD p27:gp27 100u BLKGRAD 30m mc #0 to 2 F0(zd) d17 exit ph1=0 ph2=0 1 ph3=0 ph29=0 ph30=0 ph31=0 2 ;pl1: f1 channel - power level for pulse (default) ;p1: f1 channel - 90 degree high power pulse ;d1: relaxation delay; 1-5 * T1 ;d6: delay between EPSI gradients ;d17: delay after experiment (not necessary) ;p15: EPSI gradient duration ;gpz15: EPSI positive gradient ;gpz16: EPSI negative gradient ;gpz25: prephasing gradient ;gpz27: purge gradient ;ns: 1 ;13: number of loops for EPSI acquisition

S10 PSYCHE

All spectra in this section were acquired on a 600 MHz Bruker AVIII spectrometer, equipped with a Prodigy N₂ broadband cryoprobe with a nominal maximum z-gradient strength of $66 \,\mathrm{G \, cm^{-1}}$. The sample used was $45 \,\mathrm{mM}$ and rographolide in DMSO- d_6 .

S10.1 Cost function

The cost function, named specdiff, is given by:

$$f = \left\| \frac{\mathbf{S}}{||\mathbf{S}||} - \frac{\mathbf{T}}{||\mathbf{T}||} \right\|$$
(5)

where **S** and **T** are vectors corresponding to the spectrum under optimisation and a target spectrum respectively. The target spectrum was acquired using a proton pulse-acquire sequence (zg). The optimisation spectrum was acquired using a 90°-delay-refocusing element-delay-acquire sequence, as described in the main text. (This pulse programme is attached further below.) The spectral region used in the optimisation was restricted to 1.15 ppm and above, in order to exclude methyl singlets. This is not strictly necessary (roughly the same results will be obtained if the entire spectral window is used), but helps to reduce noise in the cost function.

specdiff is not pre-installed with POISE: to enable it, uncomment the specdiff() function in poise_backend/costfunctions.py.

S10.2 Flip angle reference grid search

Due to the large amount of time required for multiple-parameter grid searches, we chose to only run this over a single parameter, namely the flip angle β .



Figure S9: Dependence of cost function on PSYCHE flip angle (cnst20).

The cost function is minimised at around 15–18°. We can use this opportunity to have a closer look at exactly how this cost function works. It penalises the low signal-to-noise seen at smaller flip angles: in this case $||\mathbf{S}||$ is small and consequently the noise in $\mathbf{S}/||\mathbf{S}||$ has a larger intensity. This contributes to a larger value of $\mathbf{S}/||\mathbf{S}|| - \mathbf{T}/||\mathbf{T}||$ and hence a larger value of the cost function. On the other hand, at larger flip angles this is not a problem: instead one starts to observe phase distortions in the peaks. This causes $\mathbf{S}/||\mathbf{S}|| - \mathbf{T}/||\mathbf{T}||$ to be larger in the vicinity of the peaks, which also pushes f up. The cost function is therefore minimised in a "compromise" region which has acceptable signal-to-noise but also not too much phase distortion. This is depicted in Figure S10.



Figure S10: A closer look at the cost function at different flip angles. Left column: The target spectrum \mathbf{T} (this is the same for all three rows). Middle column: The 1D test spectrum \mathbf{S} , acquired using different flip angles (the optimum is around 15°). Right column: The vector $\mathbf{S}/||\mathbf{S}|| - \mathbf{T}/||\mathbf{T}||$. One can see how at low flip angles, the main contribution to the cost function is noise. At higher flip angles, the noise contributes less and the increases in cost function are due to phase distortion.

S10.3 Optimisations

In this section we present the overall results of all the PSYCHE optimisations ran (with 1, 2, 3, and 4 parameters being optimised at once). From the aggregated results (Table S14), it is clear that increasing the number of parameters leads to larger decreases in the cost function, i.e. better performance as a refocusing element, although there is a drawback in that more time is required for convergence. We note here that an increase in the number of parameters leads to greater complexity in the cost function. This increases the likelihood of converging to a local minimum which may be suboptimal. The tradeoff between various factors such as potentially improved performance, the time required for optimisation, and the chances of being caught in a shallow local minimum must be carefully weighed when deciding how complicated an optimisation to run. In general, our opinion is that optimising two parameters (the flip angle and gradient amplitude) provides a good balance between all of these competing factors, and is most suitable for routine use.

First, we present an overview of all the optimisations run.

Table S14: Summary of optimisations on PSYCHE refocusing element. Details of the routines used are described in Tables S15 to S18. Parameters in parentheses were not optimised (they are inherited from the initial point). n_{fev} refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm. The optimised parameters in each row are taken from the best run, as judged by the lowest value of the cost function found.

	Aggregated results from all runs		Optimised parameters from best run			best run	
Description	$f_{ m best}$	$n_{\rm fev}$	Time (s)	β (°)	g~(%)	ΔF (kHz)	$\tau_{\rm p}~({\rm ms})$
Initial point	0.353	—	_	(25.0)	(2.00)	(10.0)	(30.0)
1 parameter	0.340 – 0.343	5 - 10	84 - 168	17.5	(2.00)	(10.0)	(30.0)
2 parameters	0.325 – 0.338	12 - 33	205 - 565	16.5	1.00	(10.0)	(30.0)
3 parameters	0.320 – 0.328	18 - 77	315 - 1344	17.4	1.73	15.9	(30.0)
4 parameters	0.316 – 0.333	39 - 85	705 - 1504	13.9	1.45	15.9	36.0

We follow this with details of the individual classes of optimisation. For the 1-parameter optimisation, the results are summarised in a similar format to previous tables. However, for the remainder of the optimisations, we present only summaries of n_{fev} and the time taken. For details of the optima found, we point the interested reader in the direction of the **poise.log** files, as well as the **parse_log()** function to help with parsing these (please refer to the documentation for more details).

Table S15: PSYCHE 1-parameter (flip angle) optimisations. The POISE routine used was: {"name": "psyche1", "pars": ["cnst20"], "lb": [10.0], "ub": [35.0], "init": [25.0], "tol": [2.0], "cf": "specdiff", "au": "poise_1d"}. Key acquisition parameters: NS=2, DS=1, D1=1.5.

Time taken (s)
135 - 168
133 - 136
84 - 85

Table S16: PSYCHE 2-parameter (flip angle and gradient amplitude) optimisations. The POISE routine used was: {"name": "psyche2", "pars": ["cnst20", "gpz10"], "lb": [10.0, 0.2], "ub": [35.0, 5.0], "init": [25.0, 2.0], "tol": [2.0, 0.2], "cf": "specdiff", "au": "poise_1d"}. Key acquisition parameters: NS=2, DS=1, D1=1.5.

Entry	Algorithm $n_{\rm fev}$		Time taken (s)		
1	Nelder-Mead	18 - 24	307 - 410		
2	MDS	25 - 33	424 - 565		
3	BOBYQA	12 - 16	205 - 271		

Table S17: PSYCHE 3-parameter (flip angle, gradient amplitude, and bandwidth) optimisations. The POISE routine used was: {"name": "psyche3", "pars": ["cnst20", "gpz10", "cnst21"], "lb": [10.0, 0.2, 1000.0], "ub": [35.0, 5.0, 20000.0], "init": [25.0, 2.0, 10000.0], "tol": [2.0, 0.2, 500.0], "cf": "specdiff", "au": "poise_psyche"}. Key acquisition parameters: NS=2, DS=1, D1=1.5.

Entry	Algorithm $n_{\rm fev}$		Time taken (s)
1	Nelder-Mead	33-43	576 - 770
2	MDS	46 - 77	804 - 1344
3	BOBYQA	18 - 31	315 - 540

Table S18: PSYCHE 4-parameter (flip angle, gradient amplitude, bandwidth, and pulse duration) optimisations. The POISE routine used was: {"name": "psyche4", "pars": ["cnst20", "gpz10", "cnst21", "p40"], "lb": [10.0, 0.2, 1000.0, 5000.0], "ub": [35.0, 5.0, 20000.0, 75000.0], "init": [25.0, 2.0, 10000.0, 30000.0], "tol": [2.0, 0.2, 500.0, 2000.0], "cf": "specdiff", "au": "poise_psyche"}. Key acquisition parameters: NS=2, DS=1, D1=1.5.

Entry	Algorithm	$n_{\rm fev}$	Time taken (s)		
1	Nelder–Mead	40 - 47	733 - 845		
2	MDS	57 - 85	1006 - 1504		
3	BOBYQA	39 - 62	705 - 1130		

The best of each of these optimisations are shown in the next figure.



Figure S11: Insets of PSYCHE spectra of a sample of 45 mM and rographolide in DMSO- d_6 , obtained before and after optimisation of up to 4 parameters. The original coupled ¹H spectrum is shown as a reference. The time required for each optimisation is indicated on each spectrum.

S10.4 AU programme used for optimisations

The poise_psyche AU programme was used for all the optimisations above. It is the same as poise_1d, except that it calls on a separate Python script to create a new double saltire pulse with the specified parameters prior to acquisition (similar to TopSpin's WaveMaker software):

```
XCMD("sendgui xpy make_double_saltire")
ZG
EFP
APBK
QUIT
```

The Python script make_double_saltire is described next.

S10.5 Python script for generating double saltire

Please note that this is a Python script to be run *inside* TopSpin. It can be downloaded from the link in the docstring.

```
.....
make_double_saltire.py
Download from: https://git.io/JUHOZ
This script reads in a pulse duration (stored in us as P40) and bandwidth
(stored in Hz as CNST21), then creates a double saltire shaped file called
"auto_double_saltire" which can be directly used as the shaped file in
PSYCHE experiments.
This can be included in an AU programme using the line:
   XCMD("sendgui xpy make_double_saltire")
SPDX-License-Identifier: GPL-3.0-or-later
11 11 11
from __future__ import division
from math import sin, cos, pi, sqrt, degrees, atan2
def make_double_saltire(duration, bandwidth):
   11.11.11
   Parameters
    _____
   duration : float
        Time in units of seconds.
   bandwidth : float
       Bandwidth in units of Hz.
   Returns
```

```
A : list of float
   Amplitudes, scaled between 0 and 100.
phi : list of float
   Phases, in degrees between -180 and 180.
These lists can be directly passed to the Bruker SAVE_SHAPE routine, as
demonstrated in this script. The amplitude of the pulse should be
automatically calculated in the pulse programme.
.....
# Our strategy is to generate the first saltire first.
# Number of points in half of the pulse. We use 1 point per us.
npoints = int(duration * 1000000 / 2)
# 20% smoothing.
s = 0.2
# This is equivalent to ls = np.linspace(0, 1, npoints), but you can't use
# numpy in TopSpin...
ts = [i/(npoints - 1) for i in range(0, npoints)]
# Construct the amplitude piecewise function for a single chirp.
A_{chirp} = [sin((pi * t) / (2 * s)) if t < s]
           else 1 if s <= t and t <= 1 - s
           else sin((pi * (1 - t)) / (2 * s))
           for t in ts]
# These are phases for a single chirp, not a saltire.
phi_chirp = [pi * (duration/2) * bandwidth * ((t - 0.5) ** 2) for t in ts]
# Calculate x-coefficients (Cx) for the saltire.
Cx = [100 * a * cos(phi) for a, phi in zip(A_chirp, phi_chirp)]
# The Cy's are 100 * a * sin(phi) ..., but that's only for one chirp. For a
# saltire, the two component chirps have equal and opposite Cy's (because
# sin(phi) = -sin(-phi)), so the total Cy is simply zero for every point.
# In principle we could make a list full of zeroes, like Cy = [0 for _ in
# range(npoints)], but that's a waste of memory and time.
# Incidentally, the Cx's for both chirps add up since cos(phi) = cos(-phi),
# but because the amplitudes are scaled to [0, 100], there's no point
# actually adding them up.
# Convert back to polar coordinates. Ordinarily A = sqrt(xi * xi + yi *
# yi), but since all the yi's are zero, we just leave it out of the
# equation.
A = [sqrt(xi * xi) for xi in Cx]
# Likewise here we use atan2(0, xi) instead of atan2(yi, xi).
phi = [degrees(atan2(0, xi)) for xi in Cx]
# Finally, we double the lists (because it's a double saltire).
return A * 2, phi * 2
```

```
if __name__ == "__main__":
    # Let the user know what's going on
    SHOW_STATUS("make_double_saltire: generating saltire...")
```

```
# Pulses are by default in us, so must convert first
duration = float(GETPAR("P 40"))/1000000
# Bandwidth is specified in Hz.
bandwidth = float(GETPAR("CNST 21"))
# Generate the amplitudes and phases.
A, phi = make_double_saltire(duration, bandwidth)
# Save it to disk.
shape_name = "auto_double_saltire"
SAVE_SHAPE(shape_name, "Excitation", A, phi)
# Use it as SPNAM40. It's too easy to forget this!
PUTPAR("SPNAM 40", shape_name)
```

S10.6 1D spin echo pulse programme

```
;psyche_1dopt_jy
;
;download from: https://git.io/JUHOS
;
;90 - 15m delay - double saltire - 15m delay - acquire
;Jonathan Yong - University of Oxford - 26 Aug 2020
;Optimisation parameters:
;cnst20 : flip angle (degrees)
;gpz10 : gradient amplitude (%)
;cnst21 : bandwidth (Hz)
      : saltire + gradient duration (us)
;p40
#include <Avance.incl>
#include <Grad.incl>
"d2=15m"
"cnst10=(cnst20/360)*sqrt((2*cnst21)/(p40/2000000))" ; rf amplitude of saltire
"cnst11=1000000/(p1*4)"
                                                      ; rf amplitude hard pulse
"spw40=plw1*(cnst10/cnst11)*(cnst10/cnst11)"
                                                      ; power level of saltire
"acqt0=0"
baseopt_echo
1 ze
2 30m
  50u BLKGRAD
  d1 pl1:f1
  50u UNBLKGRAD
  p1 ph1
  d2
  d16
  p16:gp1
  d16
  (center (p40:sp40 ph3):f1 (p40:gp10) )
  d16
  p16:gp1
  d16
  d2
  go=2 ph31
  30m mc #0 to 2 F0(zd)
50u BLKGRAD
exit
```

```
ph1 = 0
ph3 = 0 1 2 3
ph31= 0 2 0 2
;pl1 : f1 channel - power level for pulse (default)
;p1 : f1 channel - high power pulse
;p40: shaped pulse in PS element
;p21: gradient during PS element
;d1 : relaxation delay; 1-5 * T1
;d2 : delay on either side of double saltire
;d16: gradient recovery delay
;p16: CTP gradient duration
;ns: 1 * n, total number of scans: NS * TDO
;cnst20: flip angle (degrees)
;cnst21: bandwidth of PSYCHE element (Hz) (ca. 10000)
;gpnam10: RECT.1
;gpz10 : 2%
;gpnam1 : SINE.100
;gpz1 : 50%
```

S11 1D NOESY-based solvent suppression

All spectra in this section were acquired on a 400 MHz Bruker Avance NEO spectrometer, equipped with a broadband Smart probe with a nominal maximum z-gradient strength of $50 \,\mathrm{G\,cm^{-1}}$. A sample of rodent urine in D₂O was used, kindly provided by Fay Probert and Abi Yates (Department of Chemistry, University of Oxford).

S11.1 Cost function

The cost function used for all optimisations here was zerorealint_squared, which attempts to minimise the squared intensity of the real spectrum:

$$f = \sum_{i} \mathbf{S}_{\mathrm{real},i}^2 \tag{6}$$

The aim of this is to strongly penalise large deviations from the baseline, which are more likely to potentially obscure nearby peaks. This cost function comes pre-installed with POISE. The region of the spectrum under consideration was restricted to between 4.65 and 4.75 ppm using the dpl function in TopSpin.

S11.2 Transmitter offset reference grid search

It is not feasible to run a grid search over four parameters at the same time. As with the PSYCHE section, we therefore only ran one grid search for the first parameter, i.e. the transmitter offset 01.



Figure S12: Insets of the residual water peak, and the value of the cost function thus calculated, as a function of transmitter offset (01).

It is evident here that there are multiple possible minima which POISE could optimise to. POISE is not a *global* minimiser; it will only converge to a *local* minimum. In fact, a similar situation was observed when we ran grid searches over the other three parameters: multiple potential minima were observed, due to an oscillatory dependence on the parameters under optimisation.

S11.3 Optimisations

It is likely that there are multiple different sets of optimised parameters which give similar minima in the cost function. We assume that the user does not mind exactly *which* minimum is obtained, as long as there is a sufficient decrease in the cost function. In other words, we assume that the user is not concerned with exactly which combination of parameters provides the best solvent suppression (as long as they are within reasonable bounds, e.g. presaturation power is not too high; these should be set appropriately in the POISE routine).

As with the PSYCHE optimisations (Section S10), we note that optimising more parameters simultaneously takes a longer time and risks being caught in local minima. However, it is also clear that this has the potential of leading to a greater decrease in the cost function. In our experience, optimising two parameters (O1 and CNST20) is usually sufficient to give a good decrease in the cost function (i.e. good water suppression). For routine usage this two-parameter optimisation should be considered.

The detailed optimisation results are presented below, in a similar format to the PSYCHE section (Section S10.3). On the 400 MHz spectrometer we used, the default **01** value of 1880.61 Hz corresponds to exactly 4.70 ppm.

Table S19: Summary of optimisations on 1D NOESY / presaturation pulse sequence for solvent suppression.. Details of the routines used are described in Tables S20 to S23. Parameters in parentheses were not optimised (they are inherited from the initial point). $n_{\rm fev}$ refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm. The optimised parameters in each row are taken from the best run, as judged by the lowest value of the cost function found.

	Aggregated results from all runs			Optimised parameters from best run			
Description	$f_{\rm best}/10^{18}$	$n_{\rm fev}$	Time (s)	01 (Hz)	CNST20 (Hz)	D8 (s)	D1 (s)
Initial point	14.7	_	—	(1880.61)	(50.0)	(0.100)	(2.00)
1 parameter	1.85 - 2.49	6 - 12	259 - 520	1880.41	(50.0)	(0.100)	(2.00)
2 parameters	1.21 - 9.68	9 - 21	390 - 911	1880.20	51.94	(0.100)	(2.00)
3 parameters	1.24 - 4.20	19 - 26	825 - 1128	1879.90	47.79	0.118	(2.00)
4 parameters	0.165 - 1.65	25 - 53	1143 - 2314	1881.10	53.28	0.150	(3.00)

Table S20: Solvent suppression 1-parameter (transmitter offset) optimisations. The POISE routine used was: {"name": "solvsupp1", "pars": ["o1"], "lb": [1870.61], "ub": [1890.61], "init": [1880.61], "tol": [0.2], "cf": "zerorealint_squared", "au": "poise_1d_noapk"}. Key acquisition parameters: NS=2, DS=4, AQ=4.19.

Entry	Algorithm	Optimum found (Hz)	$n_{\rm fev}$	Time taken (s)
1	Nelder–Mead	1880.24 - 1880.49	10 - 12	436 - 519
2	MDS	1880.24 - 1880.36	12	518 - 520
3	BOBYQA	1880.34 - 1880.47	6 - 7	259 - 303

Table S21: Solvent suppression 2-parameter (transmitter offset and presaturation power) optimisations. The POISE routine used was: {"name": "solvsupp2", "pars": ["o1", "cnst20"], "lb": [1870.61, 10.0], "ub": [1890.61, 55.0], "init": [1880.61, 50.0], "tol": [0.2, 2.5], "cf": "zerorealint_squared", "au": "poise_1d_noapk"}. Key acquisition parameters: NS=2, DS=4, AQ=4.19.

Entry	Algorithm $n_{\rm fev}$		Time taken (s)		
1	Nelder-Mead	17 - 21	737 - 911		
2	MDS	15	649 - 652		
3	BOBYQA	9 - 14	390 - 608		

Table S22: Solvent suppression 3-parameter (transmitter offset, presaturation power, and mixing time) optimisations. The POISE routine used was: {"name": "solvsupp3", "pars": ["o1", "cnst20", "d8"], "lb": [1870.61, 10.0, 0.010], "ub": [1890.61, 55.0, 0.150], "init": [1880.61, 50.0, 0.100], "tol": [0.2, 2.5, 0.010], "cf": "zerorealint_squared", "au": "poise_1d_noapk"}. Key acquisition parameters: NS=2, DS=4, AQ=4.19.

Entry	Algorithm	Algorithm $n_{\rm fev}$	
1	Nelder–Mead	23 - 26	1002 - 1128
2	MDS	24 - 25	1034 - 1080
3	BOBYQA	19 - 26	825 - 1133

Table S23: Solvent suppression 4-parameter (transmitter offset, presaturation power, mixing time, and presaturation duration) optimisations. The POISE routine used was: {"name": "solvsupp4", "pars": ["o1", "cnst20", "d8", "d1"], "lb": [1870.61, 10.0, 0.010, 1.0], "ub": [1890.61, 55.0, 0.150, 3.0], "init": [1880.61, 50.0, 0.100, 2.0], "tol": [0.2, 2.5, 0.010, 0.1], "cf": "zerorealint_squared", "au": "poise_1d_noapk"}. Key acquisition parameters: NS=2, DS=4, AQ=4.19.

Entry	Algorithm	$n_{\rm fev}$	Time taken (s)
1	Nelder–Mead	29-40	1250 - 1726
2	MDS	34 - 53	1487 - 2314
3	BOBYQA	25 - 40	1143 - 1843

The best of each of these optimisations are shown in the next figure.



Figure S13: Insets of 1D NOESY / presaturation spectra of a sample of rodent urine, obtained before and after optimisation of up to 4 parameters. The time required for each optimisation is indicated on each spectrum.

S11.4 1D NOESY pulse programme

This pulse programme is almost the same as the Bruker default noesygppr1d, except that the presaturation power is controlled via the CNST27 parameter (which represents the RF amplitude in Hz).

```
;noesygppr1d_jy
;
;download from: https://git.io/Jnokc
;modified from noesygppr1d to use CNST20 as presaturation RF amplitude (in Hz)
; Jonathan Yong - University of Oxford - 20 June 2021
;$CLASS=HighRes
;$DIM=1D
;$TYPE=
;$SUBTYPE=
;$COMMENT=
;$RECOMMEND=
#include <Avance.incl>
#include <Grad.incl>
"d12=20u"
"plw9=plw1*(cnst20/(250000/p1))*(cnst20/(250000/p1))"
"acqt0=-p0*2/3.1416"
1 ze
2 30m
    ifdef FLAG_BLK
#
  4u BLKGRAD
    else
#
  4u
#
    endif /*FLAG_BLK*/
 d12 pl9:f1
 d1 cw:f1 ph29
  4u do:f1
 50u UNBLKGRAD
 p16:gp1
 d16 pl1:f1
 p1 ph1
 4u
 p1 ph2
 d12 pl9:f1
 d8 cw:f1
  4u do:f1
 p16:gp2
 d16 pl1:f1
```

ifdef FLAG_BLK # 4u # else 4u BLKGRAD # endif /*FLAG_BLK*/ p0 ph3 go=2 ph31 30m mc #0 to 2 F0(zd) ifdef FLAG_BLK # 4u BLKGRAD else # 4u # endif /*FLAG_BLK*/ exit ph1=0 2 ph3=0 0 2 2 1 1 3 3 ph29=0 ph31=0 2 2 0 1 3 3 1 2 0 0 2 3 1 1 3 ;pl1 : f1 channel - power level for pulse (default) ;pl9 : f1 channel - power level for presaturation ;p0 : for any flip angle ;p1 : f1 channel - 90 degree high power pulse ;p16: homospoil/gradient pulse ;d1 : relaxation delay; 1-5 * T1 ;d8 : mixing time ;d12: delay for power switching [20 usec] ;d16: delay for homospoil/gradient recovery ;ns: 8 * n, total number of scans: NS * TDO ;ds: 4 ;cnst20: presaturation power (Hz) ;use gradient ratio: gp 1 : gp 2 50 : -10 ; ;for z-only gradients: ;gpz1: 50% ;gpz2: -10% ;use gradient files: ;gpnam1: SMSQ10.100 ;gpnam2: SMSQ10.100

;preprocessor-flags-start ;FLAG_BLK: for BLKGRAD before d1 rather than go ; option -DFLAG_BLK: (eda: ZGOPTNS) ;preprocessor-flags-end

S12 Diffusion NMR

All spectra in this section were acquired on a 600 MHz Bruker AVIII spectrometer, equipped with a Prodigy N₂ broadband cryoprobe with a nominal maximum z-gradient strength of $66 \,\mathrm{G \, cm^{-1}}$. The sample used was $45 \,\mathrm{mM}$ and rographolide in DMSO- d_6 .

S12.1 Sequential optimisation strategy

Traditionally, setting up a diffusion experiment is considered to be a tedious process. At the very least, one must search for a range of gradient amplitudes (G_{\min} through G_{\max}) which provide an appropriate range of signal attenuation: thus, for example, the increment acquired with G_{\max} must exhibit a certain level of attenuation with respect to G_{\min} . If G_{\max} is too large or too small, then the resulting fitting will be less accurate. In more complicated cases, the diffusion delay Δ and/or the gradient pulse duration δ must also be tweaked.

Here, we use the oneshot DOSY experiment (*Magn. Reson. Chem.* **2002**, 40 (13), S147–S152, DOI: 10.1002/mrc.1107); the process can easily be adapted to all other DOSY variants. Throughout the optimisation, the script maintains a "reference" 1D spectrum acquired with $G_{\min} = 10\%$. The entire optimisation can be run using a single wrapper Python script, dosy_opt.py.

We seek to first optimise Δ by performing stepwise increases from a small initial value, until a separate spectrum, obtained with $G_{\max} = 80\%$, displays "sufficient" attenuation with respect to the reference spectrum. Here, we define "sufficient" as 75% attenuation. To determine this, the script leverages the dosy_aux routine and cost function. Since this is not a traditional "optimisation" in the sense that the value can only go upwards, we use the wrapper script to check the value of the cost function: thus we set the maxfev option to 1, such that on every function evaluation, control is returned to the wrapper script. The routine is (note that 1b, ub, and tol have no use in this routine and can be set to any valid value):

{"name": "dosy_aux", "pars": ["gpz1"], "lb": [79.0], "ub": [81.0], "init": → [80.0], "tol": [0.01], "cf": "dosy_aux", "au": "poise_1d"}

The dosy_aux cost function is defined by

$$f_{\text{att}} = \frac{\sum_{i} \mathbf{S}_{i}}{\sum_{i} \mathbf{R}_{i}} - 0.25 \tag{7}$$

where **R** is the reference spectrum (acquired with a minimum gradient amplitude of 10%) and **S** is the maximally attenuated spectrum (acquired with a gradient amplitude of 80%).

Once a suitable diffusion delay has been found, we then carry out a simple POISE optimisation to find the actual value of G_{max} that provides 75% attenuation. This utilises the dosy routine and cost function. The routine is as follows:

```
{"name": "dosy", "pars": ["gpz1"], "lb": [20.0], "ub": [80.0], "init":

→ [50.0], "tol": [2.0], "cf": "dosy", "au": "poise_1d"}
```

and the dosy cost function is simply the absolute value of f_{att} .

The total time taken for the entire process depends on how many times Δ must be increased from its initial value (roughly 1 additional minute per increase). In order to minimise this, it is recommended to start with a sensible value of Δ based on known precedent. The second step takes approximately 3 minutes in our hands regardless of optimisation algorithm, for an overall total time of ca. 5 minutes (using the parameters DS=2, NS=4).

In Figure S14, we plot the typical decay profiles obtained for CH and OH protons in andrographolide, obtained from a oneshot DOSY experiment using parameters optimised according to the procedure described above. The two profiles are slightly different because of chemical exchange between the OH protons and water present in the sample. This data is clearly amenable towards further analysis using processing tools such as Bruker's Dynamic Centre software or the DOSY Toolbox.



Figure S14: Example diffusion decay profiles obtained using optimised values of $G_{\text{max}} = 75.6\%$ (equivalent to 49.7 G/cm on the spectrometer we used) and $\Delta = 100 \text{ ms}$ in the oneshot DOSY experiment. G_{min} is 10%, or 6.57 G/cm. Dashed lines are the Gaussian profiles which have been fitted to the observed intensities. The sample used was 45 mM andrographolide in DMSO- d_6 .

S12.2 Python script for automation

The steps to use this script are:

- 1. Make sure that you have enabled the DOSY routines and cost functions. The routines can be found in the .../py/user/poise_backend/example_routines directory and should be copied to the .../py/user/poise_backend/routines directory; the cost functions are in .../py/user/poise_backend/costfunctions.py and must be uncommented before they can be used.
- 2. Download this Python script. There is a GitHub link in the source code below.
- 3. Set up a 2D diffusion experiment (most easily done by reading in a parameter set). The initial gradient amplitudes do not matter, but the diffusion delay should be set to something

that is not too large (since the script only *increases* the diffusion delay if it is found to be too short).

- 4. If the particular DOSY experiment has not been registered in this script yet, then do so by modifying the dictionary near the top of the script. This only needs to be done once. (If the DOSY experiment has never been set up before, then you need to find (or create) a 1D version of the diffusion experiment. It will basically be the same as the 2D version, except that the encoding gradient is not incremented.)
- 5. Run this Python script (just type dosy_opt in the TopSpin command line).

dosy_opt.py

11.11.11

Download from: https://git.io/JUHOY

Script that optimises parameters for a diffusion experiment. Currently works with oneshot DOSY (Pelta et al., Magn. Reson. Chem. 2002, 40 (13), S147-S152, DOI: 10.1002/mrc.1107), as well as the stegp1s and ledgp2s DOSY variants which are standard TopSpin pulse programmes.

This is meant to be called from the TopSpin command line on a 2D diffusion experiment.

Note that this script does not tweak the gradient lengths ("little delta"), so these should be set to a sensible value (e.g. 1 ms) before running this script.

The algorithm is as follows:

- Create 1D versions of the diffusion experiment, one with the minimum gradient (10%, the "reference spectrum", expno 99998) and one with the maximum gradient (80%, "optimisation spectrum", expno 99999).
- 2. Acquires both spectra.
- Check whether the opt. spectrum has at least 75% signal attenuation relative to the ref. spectrum. If not, increases the diffusion delay ("big Delta") and repeats steps 1 and 2.
- 4. If the attenuation is sufficient, it then runs a POISE optimisation on the opt. spectrum to find the gradient strength at which 75% signal attenuation is achieved.
- 5. Once this gradient strength has been found, returns to the 2D experiment and runs the `dosy` Bruker AU programme using the new parameter.

SPDX-License-Identifier: GPL-3.0-or-later

Dictionary of DOSY pulse programmes, to be entered in the form

```
"2d_pulprog": ["1d_pulprog", "GRAD_AMPLITUDE_PARAM", "DIFF_DELAY_PARAM"]
#
# In more detail:
  - 2d_pulprog is the name of the 2D DOSY pulse sequence.
#
# - 1d_pulprog is the name of the corresponding 1D DOSY pulse sequence, which
      should be just a single increment of the 2D sequence, but with Difframp
#
      removed. See the doneshot_nd_jy pulse programmes for an example.
#
#
  - GRAD_AMPLITUDE_PARAM is the parameter name in the 2D DOSY pulse sequence
      that is multiplied by Difframp. In the 1D DOSY pulse sequence it should
#
      not be multiplied by Difframp, so simply corresponds to the gradient
#
      strength. Note that you need a space between GPZ and the number.
#
# - DIFF_DELAY_PARAM is the parameter name in both 1D and 2D sequences that
#
      corresponds to the diffusion delay. Again, a space between D and the
#
      number is required.
# Other DOSY variants can be entered here to "register" them and make them
# available for use with this script.
pulprog_dict = {
    # doneshot_nd_jy can be downloaded from https://git.io/JUHOX and
    # https://qit.io/JUHOP.
    "doneshot_2d_jy": ["doneshot_1d_jy", "GPZ 1", "D 20"],
    # the following two entries are Bruker standard sequences and already come
    # with TopSpin.
    "stegp1s": ["stegp1s1d", "GPZ 6", "D 20"],
    "ledgp2s": ["ledgp2s1d", "GPZ 6", "D 20"],
}
def main():
    # Look up the 1D pulse programme and associated parameters based on the
    # current 2D pulse programme.
   pp2d = GETPAR("PULPROG")
    if pp2d in pulprog_dict:
       pp1d, gpzparam, dparam = pulprog_dict[pp2d]
   else:
       MSG("dosy_opt: unsupported pulse programme. New pulse programmes can"
            " be added in the dosy_opt.py script.")
       EXIT()
    # Create new routines. We need to do this within this script since the
    # routine depends on the GPZ parameter being optimised; we can't hardcode
    # that in a predefined routine.
    gpzparam_short = gpzparam.replace(" ", "").lower()
   XCMD("poise --create dosy {} 20 80 50 2"
         " dosy poise_1d".format(gpzparam_short))
   XCMD("poise --create dosy_aux {} 79 81 80 0.1"
         " dosy_aux poise_1d".format(gpzparam_short))
```

```
S58
```

```
# Get the expno of the 2D dataset, and store the parameter set. This is
# important because we want the optimisation to inherit most parameters,
# e.g. D1, SW, etc. from the parent experiment (which has presumably been
# set up by the user).
original_2d_expno = int(CURDATA()[1])
original_2d_dataset = CURDATA()
XCMD("wpar dosytemp all")
# Set up the reference experiment.
reference_expno = 99998
reference_dataset = list(original_2d_dataset)
reference_dataset[1] = str(reference_expno)
NEWDATASET(reference_dataset, None, "dosytemp")
RE(reference_dataset)
PUTPAR("PULPROG", pp1d)
PUTPAR("PARMODE", "1D")
PUTPAR("PPARMOD", "1D")
PUTPAR(gpzparam, "10")
XCMD("sendgui browse_update_tree")
# Acquire reference dataset.
msg_nonmodal("dosy_opt: optimising diffusion delay...")
XCMD("poise_1d", wait=WAIT_TILL_DONE)
# Set up optimisation expno.
optimisation_expno = 99999
optimisation_dataset = list(original_2d_dataset) # make a copy
optimisation_dataset[1] = str(optimisation_expno)
NEWDATASET(optimisation_dataset, None, "dosytemp")
RE(optimisation_dataset)
PUTPAR("PULPROG", pp1d)
PUTPAR("PARMODE", "1D")
PUTPAR("PPARMOD", "1D")
PUTPAR(gpzparam, "80")
XCMD("sendgui browse_update_tree")
# Optimise the diffusion delay first.
while True:
    # Here we use the `--maxfev 1` trick to evaluate a cost function at 80%
    # gradient. The cost function is stored as the `TI` parameter after the
    # experiment has been recorded.
    RE(optimisation_dataset)
    PUTPAR("TI", "") # has to be one space, not an empty string (try it)
   XCMD("poise dosy_aux -a nm --maxfev 1 -q", wait=WAIT_TILL_DONE)
    # If TI isn't a valid float, something went wrong with POISE.
```

```
try:
        cf = float(GETPAR("TI"))
    except ValueError:
       EXIT()
    # The dosy_aux cost function is (opt. intensity/ref. intensity) minus
    # 0.25. If the opt. spectrum has been attenuated too little, then this
    # value will be positive.
    if cf > 0:
        # Increase delay, i.e. big Delta, by 100 ms
        current_delay = float(GETPAR(dparam))
        new_delay = current_delay + 0.1
       PUTPAR(dparam, str(new_delay))
        # Set it in the reference spectrum too
       RE(reference_dataset)
       PUTPAR(dparam, str(new_delay))
        # Re-acquire the reference spectrum.
        XCMD("xau poise_1d", wait=WAIT_TILL_DONE)
    # If the cost function is negative, then that means it's sufficiently
    # attenuated, and we can go on to the actual optimisation.
    else:
       break
# At this point, we know that 80% gradient strength provides at least 75%
# attenuation. So we can run an optimisation to find the actual value that
# does provide 75% attenuation (which is somewhere between 10% and 80%).
RE(optimisation_dataset)
best_delay = GETPAR(dparam)
msg_nonmodal("dosy_opt: diffusion delay set to {} seconds.\n"
                 now optimising gradient amplitude...".format(best_delay))
PUTPAR("TI", " ")
XCMD("poise dosy -a bobyqa -q", wait=WAIT_TILL_DONE)
# Again, check to make sure that TI is a valid float.
try:
   float(GETPAR("TI"))
except ValueError:
   EXIT()
best_gpz = GETPAR(gpzparam)
# Jump back to original 2D expt. Set dparam, but not gpzparam as that
# should always be 100% (the actual gradient amplitude is controlled by
# Difframp, which will be generated by the `dosy` AU programme).
RE(original_2d_dataset)
PUTPAR(dparam, best_delay)
PUTPAR(gpzparam, "100")
```

```
# Generate Difframp, etc. and then run the AU programme.
```

```
msg_nonmodal("dosy_opt: maximum gradient set to {}%.\n"
                 " starting 2D experiment with optimised"
                 " parameters...".format(best_gpz))
   XCMD("xau dosy 10 {} 16 l y".format(best_gpz))
def msg_nonmodal(s):
    ппп
   Show a message box that does not force the user to click OK before the
   programme continues. See also the same function in poise.py (the frontend
   script).
    ппп
   dialog = mfw.BInfo()
   dialog.setMessage(s)
   dialog.setTitle("dosy_opt")
   dialog.setBlocking(0) # in ordinary MSG() this is 1
   dialog.show()
if __name__ == "__main__":
```

```
main()
```

S12.3 Sequential optimisation details

The second part of the optimisation (G_{max} optimisation) is the only place where using different algorithms can have an effect. In order to evaluate the effect of using different algorithms, we ran the second optimisation using each of the algorithms, using the diffusion delay $\Delta = 100 \text{ ms}$ obtained from the first step. The results are summarised in Table S24.

Table S24: Optimisations of maximum gradient amplitude (G_{max}) for diffusion NMR spectroscopy. The POISE routine used here is: {"name": "dosy", "pars": ["gpz1"], "lb": [20.0], "ub": [80.0], "init": [50.0], "tol": [2.0], "cf": "dosy", "au": "poise_1d"}. The pulse programme used was doneshot_1d_jy (Section S12.5). Key acquisition parameters: NS=4, DS=2, D1=1.5, AQ=1.14. n_{fev} refers to the number of function evaluations, i.e. the number of spectra acquired in the course of the optimisation. The results are aggregated from five separate runs for each optimisation algorithm.

Entry	Algorithm	Optimum found $(\%)$	$n_{\rm fev}$	Time taken (s)
1	Nelder-Mead	75.00	9	195 - 197
2	MDS	75.00	9	196 - 197
3	BOBYQA	75.00 - 75.60	8–9	175 - 197

S12.4 Simultaneous optimisation strategy

We also tried directly performing a simultaneous two-parameter optimisation of both the gradient amplitude and the diffusion delay, instead of the entire procedure above. The POISE routine and its associated cost function are called dosy_2p. The routine is:

{"name": "dosy_2p", "pars": ["gpz1", "d20"], "lb": [20.0, 0.05], "ub": → [80.0, 0.4], "init": [50.0, 0.05], "tol": [2.0, 0.1], "cf": "dosy_2p", → "au": "poise_2pdosy"}

In order to bias the optimisation towards smaller values of Δ (which is desirable as it minimises T_1 losses), we added a term to the cost function that is proportional to Δ . This cost function can be found in .../py/user/poise_backend/costfunctions.py as well, although it is commented out by default.

$$f = |f_{\text{att}}| + (\Delta/s) \tag{8}$$

Unfortunately, although this optimisation did work, it was far slower: for example, using the Nelder-Mead algorithm, this took 36 function evaluations and 1572 seconds (over 26 minutes). The optimum converged to was $G_{\text{max}} = 79.3\%$ and $\Delta = 91.6 \text{ ms}$, which is very similar to that found above. Part of the reason why this is so slow is that every time the diffusion delay Δ is changed, the reference spectrum must be reacquired: this means that every function evaluation in this two-parameter optimisation takes twice as long as every function evaluation in the individual optimisation of G_{max} .

S12.5 Oneshot DOSY pulse programmes

S12.5.1 2D Oneshot

Some modifications from the published version were made so that it could be run on our spectrometers.

```
;doneshot_2d_jy
;download from: https://git.io/JUHOX
; Jonathan Yong - University of Oxford - 26 Aug 2020
;minor modifications to original oneshot DOSY sequence to make it play better with
\hookrightarrow TopSpin.
;modified from:
;Ralph Adams, Juan Aguilar, Robert Evans, Mathias Nilsson and Gareth Morris
;University of Manchester
;Release 1.0c (27Mar2012)
;Source citation:
;M.D. Pelta, G.A. Morris, M.J. Stchedroff, S.J. Hammond, Magn. Reson. Chem. 40 (2002)
\hookrightarrow S147-S152.
;Other relevant papers that could be of use include:
;A. Botana, J.A. Aguilar, M. Nilsson, G.A. Morris, J. Magn. Reson. 208 (2011) 270-278.
;2D Doneshot DOSY pulse sequence
;$CLASS=HighRes
;$DIM=2D
;$TYPE=
```

```
;$SUBTYPE=
;$COMMENT=
#include <Avance.incl>
#include <Grad.incl>
#include <Delay.incl>
define list<gradient> diff=<Difframp>
;JY - alpha is hardcoded as 0.2
"cnst14=0.2"
"cnst17=(2*p1+d16+p30)*0.000001"; Dtau
;Assuming half-sine gradient pulses [most common on Bruker systems]
"cnst18=0.000001*p30*2*0.000001*p30*2*(d20 - (2*0.000001*p30*(5-3*cnst14*cnst14)/16) -
_{\leftrightarrow} (cnst17*(1-cnst14*cnst14)/2) )" ; Dosytimecubed
"cnst15=1+cnst14" ; 1 + alpha
"cnst16=1-cnst14" ; 1 - alpha
"p2=p1*2"
"DELTA1=d20-4.0*p1-4.0*p30-5.0*d16-p19"
"acqt0 = -p1*2/3.1416"
1 ze
2 d1
3 50u UNBLKGRAD
                           ;Spoiler gradient balancing pulse
  p19:gp2*-1
 d16
                           :1st 90
  p1 ph1
  p30:gp1*diff*cnst16
                           ;1 - alpha
  d16
  p2 ph2
                           ;First 180
  p30:gp1*diff*cnst15*-1
                           ;1 + alpha
  d16
  p1 ph3
                           ;2nd 90
                           ;Lock refocusing pulse pulse
  p30:gp1*diff*cnst14*2
  d16
  p19:gp2
                           ;Spoiler gradient balancing pulse
  d16
  DELTA1
  p30:gp1*diff*cnst14*2
                           ;Lock refocusing pulse pulse
  d16
                           ;3rd 90
  p1 ph4
  p30:gp1*diff*cnst16
                           ;1 - alpha
  d16
  p2 ph5
  p30:gp1*diff*cnst15*-1 ;1 + alpha
  d16
  4u BLKGRAD
  go=2 ph31
  d1 mc #0 to 2 F1QF(calgrad(diff))
```

;d20 : diffusion time (big DELTA)

;gpnam1: SMSQ10.100 ;gpnam2: SMSQ10.100

;gpz1 : diffusion encoding gradient (100%)

;gpz2 : Oneshot CTP gradients

exit
```
;cnst14 : alpha, typically 0.2
;cnst15 : 1+alpha
;cnst16 : 1-alpha
;cnst17 : Dtau
;cnst18 : dosytimecubed
;NS : 1 * n
;DS : 1 * m
;td1: number of experiments
;FnMODE: QF
;use AU-program dosy to calculate gradient ramp-file Difframp
         use xf2 and DOSY processing
         use "setdiffparmUoM" if 'setdiffparm' does not work
;
     or use "setdiffparm STEbp" but this gives slighly distorted diffusion coefficients
;
S12.5.2 1D Oneshot
; doneshot_1d_jy
:
; download from: https://git.io/JUHOP
; 1D version of Doneshot
; Jonathan Yong - University of Oxford - 26 Aug 2020
; modified from
; Doneshot
; DOSY pulse sequence
; Ralph Adams, Juan Aguilar, Robert Evans, Mathias Nilsson and Gareth Morris
; University of Manchester
; Release 1.0c (27Mar2012)
; Source citation:
; M.D. Pelta, G.A. Morris, M.J. Stchedroff, S.J. Hammond, Magn. Reson. Chem. 40 (2002)
\hookrightarrow S147-S152.
; Other relevant papers that could be of use include:
; A. Botana, J.A. Aguilar, M. Nilsson, G.A. Morris, J. Magn. Reson. 208 (2011) 270-278.
;$CLASS=HighRes
;$DIM=1D
;$TYPE=
;$SUBTYPE=
;$COMMENT=
#include <Avance.incl>
#include <Grad.incl>
#include <Delay.incl>
"cnst17=(2*p1+d16+p30)*0.000001" ; Dtau
```

```
;JY - alpha is hardcoded as 0.2
"cnst14=0.2"
;Assuming half-sine gradient pulses [most common on Bruker systems]
"cnst18=0.000001*p30*2*0.000001*p30*2*(d20 - (2*0.000001*p30*(5-3*cnst14*cnst14)/16) -
\rightarrow (cnst17*(1-cnst14*cnst14)/2))"; Dosytimecubed
"cnst15=1+cnst14" ; 1 + alpha
"cnst16=1-cnst14" ; 1 - alpha
"p2=p1*2"
"DELTA1=d20-4.0*p1-4.0*p30-5.0*d16-p19"
;The below line should be commented out for Topspin < 2.0
"acqt0 = -p1*2/3.1416"
1 ze
2 d1
3 50u UNBLKGRAD
 p19:gp2*-1
                       ;Spoiler gradient balancing pulse
 d16
 p1 ph1
                         ;1st 90
 p30:gp1*cnst16
                           ;1 - alpha
 d16
                         ;First 180
 p2 ph2
 p30:gp1*cnst15*-1
                     ;1 + alpha
 d16
 p1 ph3
                         ; 2nd 90
 p30:gp1*cnst14*2
                     ;Lock refocusing pulse pulse
 d16
                          ;Spoiler gradient balancing pulse
 p19:gp2
 d16
 DELTA1
 p30:gp1*cnst14*2
                      ;Lock refocusing pulse pulse
 d16
                         ; 3rd 90
 p1 ph4
 p30:gp1*cnst16
                           ; 1 - alpha
 d16
 p2 ph5
 p30:gp1*cnst15*-1
                           ;1 + alpha
 d16
 4u BLKGRAD
 go=2 ph31
 d1 mc #0 to 2 F0(zd)
evit
```

;cnst16: 1-alpha ;cnst17: Dtau ;cnst18: dosytimecubed ;ns: 1*n ;ds: 1*m ;td1: number of experiments ;gpnam1: SMSQ10.100

;gpnam2: SMSQ10.100

;gpz1: DOSY gradient strength. DO NOT EXCEED 80 ;gpz2: Oneshot CTP gradients

S13 Documentation

Finally, the current version of the POISE documentation is attached here. This is meant to act as a user guide, with specific instructions on how to set up and run optimisations.

Please note that the latest version of the documentation can always be found online at

https://foroozandehgroup.github.io/nmrpoise

and a PDF version of this can be downloaded at

https://foroozandehgroup.github.io/nmrpoise/poise.pdf

Contents

1	Installation1.1Updating POISE	S72 S72 S73 S73 S73 S73
2	Setting up a Routine	S 75
3	Running an optimisation3.1 Parsing the log3.2 Errors	S78 . S79 . S80
4	Frontend options4.1General options4.2Options for running optimisations	S81 . S81 . S82
5	POISE under automation5.1 Basics5.2 A simple(ish) example5.3 Terminating scripts which call POISE	S83 S83 S83 S85
6	Builtin cost functions6.1minabsint6.2maxabsint6.3minrealint6.4maxrealint6.5zerorealint6.6noe_1d6.7epsi_gradient_drift	S87 S87 S87 S87 S88 S88 S88 S88 S88 S88
7	Builtin AU programmes7.1 poise_1d7.2 poise_2d	S89 . S89 . S89
8	Custom cost functions8.1The rules for cost functions8.2Accessing spectra and parameters8.3Logging8.4Premature termination8.5Examples	S90 . S90 . S91 . S94 . S95 . S96
9	Developer notes 9.1 Testing	S97 . S98

POISE (*Parameter Optimisation by Iterative Spectral Evaluation*) is a Python package for on-the-fly optimisation of NMR parameters in Bruker's TopSpin software.

The code is licensed under the GNU General Public License v3 and can be found on GitHub.

Note: The documentation you are currently reading is for version 1.2.0 of POISE. To check your current version of POISE, type poise --version in TopSpin.

The most recent HTML version of the documentation can always be found at

https://foroozandehgroup.github.io/nmrpoise.

A PDF version can be found at

https://foroozandehgroup.github.io/nmrpoise/poise.pdf

In here you will find guides on setting POISE up and using it in routine NMR applications. This guide can largely be read in sequence. However, depending on your level of interaction with the software, you may not need to read all of it. For example, if somebody else has already set up POISE for you, you can probably skip to *Running an optimisation*.

Installation

POISE comprises a *frontend* script, which is accessible from within TopSpin itself, as well as a *backend* script, which has to be run using Python 3 (i.e. not TopSpin's native Python interpreter).

This means that Python 3 must be installed separately first. In particular, POISE requires a minimum version of **Python 3.6.** For Windows, your best bet is to download an installer from the Python website. For Unix machines, we suggest using a package manager to do so (such as Homebrew for macOS or apt/yum & their equivalents on Linux), although the installers are also fine.

Once that's done, you can install POISE using pip (or pip3):

```
pip install nmrpoise
```

The package requirements are numpy, scipy, pandas, and Py-BOBYQA; these will be automatically downloaded if necessary.

pip tries to take care of installing the scripts to your TopSpin directory. To do so, it checks for TopSpin installations in standard directories (/opt on Unix and C:\Bruker on Windows). If pip exits without errors, this should have succeeded; you can test it by typing poise -h into TopSpin's command-line, which should spawn a popup. If that is the case, congratulations — you can move on to the next chapter, *Setting up a Routine*.

1.1 Updating POISE

Simply use:

pip install --upgrade nmrpoise

All other steps (including troubleshooting, if necessary) are the same.

1.2 Troubleshooting

The installation can occasionally fail if TopSpin is installed to a non-standard location. To solve this issue, you can specify the TopSpin installation directory as an environment variable TSDIR before installing POISE. The way to do this depends on what operating system (and shell) you use.

On Windows PowerShell, run the following command:

\$env:TSDIR = "C:\Bruker\TopSpinX.Y.Z\exp\stan\nmr"

replacing the part in quotes with your actual TopSpin installation directory (it must point to the exp/stan/ nmr folder).

On old-school Windows cmd, use:

set TSDIR="C:\Bruker\TopSpinX.Y.Z\exp\stan\nmr"

On Unix systems, use:

export TSDIR="/opt/topspinX.Y.Z/exp/stan/nmr"

(Unless you're using csh or the like, in which case you use setenv, although you probably didn't need to be told that!)

After running the appropriate command for your operating system, pip install nmrpoise should be able to detect the TSDIR environment variable and install the scripts accordingly.

1.3 From source

If you obtained the source code (e.g. from git clone or a GitHub release) and want to install from there, simply cd into the top-level nmrpoise directory and run:

pip install .
or equivalently:
python setup.py install

The installation to the TopSpin directory is subject to the same considerations as above.

1.4 Without an internet connection

If your spectrometer does not have an Internet connection, then the installation becomes a bit more protracted. On a computer that *does* have an Internet connection:

- 1. Make a new folder.
- 2. Download the CPython installer for the spectrometer operating system. Place it in that folder.
- 3. cd to the folder and run these commands:

```
pip download Py-BOBYQA nmrpoise --no-deps --no-binary=:all:
pip download numpy pandas scipy --only-binary=:all: --python-version <3.X> --platform
$\to <PLATFORM>$
```

where <3.X> is the version of Python to be installed on the spectrometer, and <PLATFORM> is one of win32, win_amd64, macosx_10_9_x86_64, or manylinux1_x86_64 depending on the spectrometer operating system. (n.b. These options are not well-documented; I have figured them out by trawling Stack Overflow. Any additions are welcome.)

4. Copy the whole folder over to your spectrometer. It should contain a bunch of .whl files and two .tar.gz files.

Now, on the spectrometer:

- 1. Install CPython with the installer.
- 2. cd to the folder and run

pip install ./nmrpoise-<VERSION>.tar.gz --no-index --find-links .

(replace <VERSION> with whichever version you downloaded).

It should then install properly, unless your TopSpin installation location is non-standard: in that case, set the \$TSDIR environment variable (described above) before retrying Step 2.

Setting up a Routine

If you're coming here from the *Installation*, you should make sure that POISE has been installed correctly. A simple check is to type in poise -1 into the TopSpin command line: if it shows a text box, then you should be good to go.

Each optimisation in POISE is controlled by a **routine**, which contains all information necessary for an optimisation. The ingredients of a routine are:

- A name
- The parameters to be optimised
- Lower bounds, upper bounds, initial values, and tolerances for each parameter
- A cost function which determines the 'badness' of a spectrum
- (Optional) The name of an AU programme for acquisition and processing

We will now walk through how to set a routine up, elaborating on each of these ingredients as we come to them. To get this process started, type poise into the TopSpin command line.

Note: POISE has a number of command-line options. If you're interested in finding out more about these, poise -h will give you a short summary of each of them, and *Frontend options* has additional info.

The routine we will set up now is one for the calibration of the 360° pulse width. The first ingredient we need to provide is a **name**. I've used p1cal for this, but of course you can choose anything you prefer:

icuse ei	oose a name to sa	re this routine	
Name:	p1cal		
			OK Canc

After clicking OK, you will be prompted to provide the **parameters** that are being optimised. Here we are just optimising one parameter, p1.

Please select e	xperimental parameters to be optimised (separated by commas or spaces
Parameters:	p1
	OK Cancel

Note: Only parameters that take on float values can be optimised (pulses p, delays d, constants cnst...) Integer values, like td or loop counters 1, will not work.

At this stage, you will be prompted to enter the **bounds**, **initial value**, and **tolerances**. The lower and upper bounds simply reflect a range within which the optimum can reasonably be assumed to lie within, and the initial value should be your best guess at where the optimum is.

On the spectrometer we're currently using, the Prosol value for a 90° pulse is 12 μ s, so we'll go ahead and set the initial guess for the 360° pulse to be 48 μ s. (*If your Prosol value differs, you should adjust these values accordingly.*) The lower and initial bounds can be 40 and 56 μ s respectively, corresponding to a 90° pulse of between 10 and 14 μ s. The tolerance, on the other hand, roughly reflects the degree of accuracy that you want in the answer. Here we've used a value of 0.2 μ s.

lease choose t	he minimum value, maximum value, initial value, and tolerance for p
Minimum:	40
Maximum:	56
Initial value:	48
Tolerance:	0.2
	OK

Choosing tolerances can be tricky sometimes. Too large a tolerance can lead to inaccurate answers (as the optimisation converges before it's really found the minimum); and too small a tolerance is meaningless, as often the resulting spectra are barely different. Generally, it's a good idea to choose the smallest value where going in either direction will give you an appreciable difference in the spectrum. However, it doesn't have to be *too* precise: as long as you aren't off by an order of magnitude POISE will still work reasonably well.

Note: The default TopSpin units for pulse lengths are microseconds, so the unitless 48 is equal to 48 μ s. However, for delays the default units are seconds.

POISE also allows you to specify units using the suffixes 'u', 'm', and 's' for microseconds, milliseconds, and seconds respectively. This is designed to mimic TopSpin's parameter settings, where 30m means 30 ms (for example). So you can enter 48u in this screen as well, or indeed 0.048m.

Finally, we have to choose a **cost function**, as well as (optionally) an **AU programme**.

The cost function is a Python function which reads the spectrum and returns a 'cost', i.e. how bad the spectrum is. All optimisations in POISE seek to minimise the cost function. In our case, the best value of p1 is one for which the intensity of a pulse-acquire spectrum (zg) is minimised, i.e. magnetisation is returned to the positive z-axis. So, we can conveniently use the intensity of the magnitude-mode spectrum as the cost function. This cost function is also bundled with POISE, and is called minabsint. (For those who are familiar with TopSpin's built-in popt, this is equivalent to the MAGMIN criterion.)

Note: For more information about the built-in cost functions, check out *Builtin cost functions*.

The AU programme controls spectrum acquisition and processing, and can be left blank in this case. All we need to do for this routine is to acquire the spectrum, Fourier transform, then perform phase and baseline correction. For 1D and 2D datasets, if the AU programme option is left blank, POISE will automatically do exactly these steps. Therefore, there is no need to specify an AU programme unless you want to customise this process.

Cost function:	minabsint	
AU programme (optional):		
	ОК	Cancel

That's it — congratulations, you've set up a POISE routine! If you type poise -1 now, you should now see the p1cal routine (or whatever you named it) appear in the text box:

e	Edit Search	
	Available cost functions (and their descriptions)	
	<pre>* maxabsint Cost function which maximises the absolute (magnitude-mode) intensity of the spectrum. * minobcit</pre>	
	* miniousin Cost function which minimises the absolute (magnitude-mode) intensity of the spectrum. This is probably the easiest cost function. :-)	
	* maxrealint None	
	★ noe_ld Measures the intensity of peaks in the spectrum, ★except★ for anything within 25 Hz of the selectively excited peak.	
	* zerorealint None	
	* minrealint None	
	Available routines	
	slcal au : cf : minabsint	
	init : [48.0] lb : [40.0] name : plcal	
	pars : [u'p1'] tol : [0.2] ub : [56.0]	

Running an optimisation

Assuming you or someone else has already created a routine (see *Setting up a Routine* if not), this page will show you how to run the optimisation. We'll use the same p1 calibration routine that we described on that page, but the principles apply equally to all routines.

The first thing to do is to set up the NMR experiment. Use edc or new to create a new proton pulse-acquire experiment. You should use the pulse programme zg (not zg30 or zg60!). Set the other experimental parameters, such as the spectral width SW, relaxation delay d1, etc. as desired for your compound of interest.

All these steps can in principle be done most easily by loading a parameter set (rpar). On Bruker systems, there should already be a builtin PROTON parameter set. After loading this parameter set you will have to run getprosol, then change the pulse programme to zg. Alternatively, there might be a different parameter set that has been set up by a member of the NMR staff for simple proton spectra. As long as you make sure the pulse programme is zg the optimisation will work fine.

3 poise_demo 1 1 C:	NMR\data\mfgroup\nm							
Spectrum ProcPa	ars AcquPars Title	PulseProg P	eaks Integrals	Sample	Structure	Plot	Fid	Acqu
ю Л S 📙 🖾	1,2, ▼ C 🆓		Probe:	Prodig	y BBO			
Experiment Width Receiver Nucleus Durations Power Probe Lists Wobble Lock Automation Miscellaneous	Experiment PULPROG AQ_mod TD DS NS TD0 Width SW (ppm) SW (ppm)	zg DQD 16384 0 1 1 1	· · · · · · · · · · · · · · · · · · ·		E	Curr Acq Size Num Loo	rent p uisitio of fio nber o p cou	bulse program on mode d of dummy scans of scans int for 'td0' width
User Routing	SWH [HZ] AQ [sec] FIDRES [HZ] FW [HZ]	7211.539 1.1359574 0.880315 4032000.0	00			Spe Acq Fid Filte	ctral uisitic resol er wid	wiath on time ution th

Note: Apart from the pulse programme, basically every other parameter can be set to whatever you like. However, to reduce the overall time taken, it's generally a good idea to try to make each experiment as short as possible. In this case, even with very dilute samples, 1 scan (NS=1) will suffice. We find that dummy scans are not needed to obtain accurate results, so it is permissible to set DS=0 (although see next paragraph for a caveat). You could cut this even further by lowering TD to 8192 (for a given SW, this translates to a shorter acquisition time AQ).

For routine usage we recommend using at least 1 dummy scan. This (p1cal) is the only optimisation example in which we have used 0 dummy scans. Skipping dummy scans altogether can lead to inaccuracies

in the cost functions (as the system has not reached a steady state).

Lock, shim, and tune as usual (if you haven't already). Once that's done, simply enter into the TopSpin command line:

```
poise p1cal
```

Sit back and watch it run! You should get a result in 1–2 minutes.



The best value(s) will automatically be stored in the corresponding parameter so that any subsequent acquisition will use the optimised parameters. Don't forget that for this particular optimisation, you will have to divide the optimised value by 4 to get back the 90° pulse.

3.1 Parsing the log

If you are interested in analysing data from (possibly multiple) optimisation runs, all information is logged in a poise.log file. This log file is stored inside the expno folder (the post-optimisation popup also tells you where it can be found — see above for an example). Some of the key information inside the poise.log file can be extracted in Python:

```
>>> from nmrpoise import parse_log
>>> parse_log("C:/NMR/data/mfgroup/nmr/poise_demo/1")
 routine initial param
                           1b
                                 ub tol algorithm
                                                       costfn
                                                                 auprog
                                                                         optimum
→fbest nfev time
0
   p1cal
             48.0
                   [p1]
                         40.0 56.0 0.2
                                                nm minabsint poise_1d
                                                                          48.125 6.
                     77
→849146e+06
               10
```

parse_log returns a pandas.DataFrame object which contains most of the information in the log file. However, this object does not include details of individual cost function evaluations (even though these are fully logged). If you want to analyse that data, you will have to write your own function!

The full documentation for parse_log (which really doesn't say much more than the previous example) is as follows:

```
nmrpoise.parse_log(fname='.')
Parse a poise.log file.
```

Parameters

fname [Path or str or int, optional] Path to poise.log file, or the folder containing it (this would be the TopSpin EXPNO folder). If an int is passed, it is interpreted as the string "./<fname>" (i.e. expno X in the current working directory). If not specified, defaults to the current working directory.

Returns

log_df [DataFrame] DataFrame with rows corresponding to optimisations which successfully terminated. The time taken is given in seconds.

3.2 Errors

POISE tries its best to exit gracefully from errors, and often you won't need to care about any of them. However, if something *does* go wrong during an optimisation, errors will be logged to the two files poise_err_frontend.log and poise_err_backend.log, depending on which script runs into an error. These files reside in the same folder as poise.log, i.e. the "expno folder". We welcome bug reports — please submit an issue on GitHub or drop us an email (addresses can be found on the paper).

Frontend options

The frontend script (i.e. the script which you run in TopSpin) has a variety of flags which control its behaviour. These are fully described here. You can also access short descriptions of these by typing poise -h in TopSpin.

4.1 General options

```
--create NAME PARAM MIN MAX INIT TOL CF AU
```

Create a new named routine. This can only handle one-parameter routines; if you want a multiple-parameter routine, please use the GUI. The AU programme must be specified, it can't be left blank like in the GUI. Also, short forms such as 4u for 4 μ s are not allowed here.

--delete DELETE

Delete a named routine. For example, run poise --delete p1cal to delete a routine named p1cal.

If you want to edit a routine, you can just re-create a new routine with the same name: the old one will be overwritten.

```
-h, --help
```

Show a help message then exit.

--kill

Kill POISE backends that may still be running.

Running poise --kill should be the first course of action if you find unusual behaviour after terminating a POISE optimisation (e.g. being unable to delete a log file as it is still in use). If this does not work (very rare), then you may need to manually kill the Python processes: for example, on Windows PowerShell, run:

Stop-Process -name python

or on Unix systems:

killall -9 python

(replace python with python3 as appropriate for your system). If you find that you need to do this, and can reproduce the error, please do submit a bug report or drop us an email.

-l, --list

List all available cost functions, as well as all available routines and their parameters.

4.2 Options for running optimisations

-a ALG, --algorithm ALG

Use the algorithm ALG for the optimisation. ALG can be one of nm (for Nelder–Mead), mds (for multidirectional search), or bobyqa (for Py-BOBYQA). The default is bobyqa.

--maxfev MAXFEV

Maximum function evaluations to allow (i.e. maximum number of spectra to acquire during the optimisation run). If the optimisation reaches the limit, it will terminate, reporting the best value so far as the 'optimum'.

This is useful for enforcing an upper limit on the time taken to perform an optimisation. Since by far the majority of the time is spent on acquiring the NMR spectra, MAXFEV evaluations will simply take roughly MAXFEV * t time to run (where t is the time taken for one spectrum — you can find this out using TopSpin's expt command).

If you don't want to have a limit on function evaluations, just don't use this flag, or pass the value of 0. Technically, there is always a hard limit on the number of function evaluations (which is 500 times the number of parameters being optimised). However, it is probably almost impossible to run into that hard limit.

-q, --quiet

Don't display the final popup at the end of the optimisation informing the user that the optimisation is done. This is mostly a matter of taste, as the final popup does not block any subsequent commands from being executed.

-s, --separate

Use a separate expno for each function evaluation. Note that if POISE runs into an expno which already exists, it will terminate with an error!

POISE under automation

The command-line interface that POISE offers (see *Frontend options*) allows us to wrap POISE within a larger script. Here we present some examples of how this can be accomplished: after this, the extension to automation is largely straightforward. For example, if POISE is used inside an AU programme, then set the the AUNM parameter in TopSpin appropriately.

5.1 Basics

To incorporate POISE in an AU programme, you can use the syntax:

XCMD("sendgui xpy poise <routine_name> [options]")

inside the AU script.

(Optional: To suppress POISE's final popup (telling the user that the optimum has been found), you can add the -q flag in the options. The popup won't stop TopSpin from running whatever it was going to run, though, so it's completely safe to show the message.)

Alternatively, you can wrap POISE within a Python script, which is arguably easier to write. The corresponding syntax for running an optimisation would be:

XCMD("xpy poise <routine_name> [options]")

As you can see, it is the same except that sendgui isn't needed.

5.2 A simple(ish) example

Here's an example of how the p1 optimisation (shown in *Setting up a Routine*) can be incorporated into an AU script (download from here). This AU script performs a very similar task to the existing pulsecal script: it finds the best value of p1 and plugs it back into the current experiment. However, as we wrote in the paper, it tends to provide a much more accurate result. In practice, we've already used it many times to calibrate p1 before running other experiments.

Note: This AU programme comes installed with POISE. However, you must create the p1cal routine before you can use this. Please see *Setting up a Routine* for a full walkthrough.

```
GETCURDATA
int old_expno = expno;
// Use EXPNO 99999 in the current folder for optimisation.
DATASET(name, 99999, procno, disk, user)
// Set some key parameters. Notice that these lines can be substantially cut
// if an appropriate parameter set is set up beforehand.
RPAR("PROTON", "all")
GETPROSOL
STOREPAR("PULPROG", "zg")
STOREPAR("NS", 1)
STOREPAR("DS", 0)
STOREPAR("D 1", 1.0)
STOREPAR("RG", 1)
STOREPAR("F1P", f1p)
STOREPAR("F2P", f2p)
STOREPARS("F1P", f1p)
STOREPARS("F2P", f2p)
// Run optimisation (uses BOBYQA by default).
XCMD("sendgui xpy poise p1cal -q")
// POISE stores the optimised value in p1 after it's done. We can retrieve it
// here. Don't try to get the *status* parameter, since that is not the
// optimised value (it is the value used for the last function evaluation!)
float plopt;
FETCHPAR("P 1", &plopt)
p1opt = p1opt/4;
// Move back to old dataset and set p1 to optimised value.
DATASET(name, old_expno, procno, disk, user)
VIEWDATA_SAMEWIN // not strictly necessary, just re-focuses the original spectrum
STOREPAR("P 1", plopt)
Proc_err(INFO_OPT, "Optimised value of p1: %.3f", p1opt);
// (Optional) Run acquisition.
// ZG
QUIT
```

Note that the six lines underneath "set some key params" can be collapsed to one line if an appropriate parameter set is set up beforehand.

Here's the Python equivalent of the AU programme above (download from here):

```
# Read in F1P and F2P from current dataset.
f1p = GETPAR("F1P")
f2p = GETPAR("F2P")
# Use EXPNO 99999 in the current folder for optimisation.
old_dataset = CURDATA()
opt_dataset = CURDATA()
opt_dataset[1] = "99999"
# Create new dataset and move to it.
NEWDATASET(opt_dataset, None, "PROTON")
RE(opt_dataset)
# Set some key parameters. Notice that these lines can be cut if an appropriate
# parameter set is set up beforehand (and loaded using NEWDATASET()).
XCMD("getprosol")
PUTPAR("PULPROG", "zg")
```

(continues on next page)

(continued from previous page)

```
PUTPAR("NS", "1")
PUTPAR("DS", "0")
PUTPAR("D 1", "1")
PUTPAR("RG", "1")
XCMD("s f1p {}".format(f1p)) # PUTPAR("status F1P") doesn't work.
XCMD("s f2p {}".format(f2p)) # Even though the documentation says it should.
# Run optimisation (uses BOBYQA by default).
XCMD("poise p1cal -q")
# POISE stores the optimised value in p1 after it's done. We can retrieve it
# here. Don't try to get the *status* parameter. since that is not the
# optimised value (it is the value used for the last function evaluation!)
plopt = float(GETPAR("P 1"))/4
# Move back to old dataset and set p1 to optimised value.
RE(old_dataset)
PUTPAR("P 1", str(plopt))
ERRMSG("Optimised value of p1: {:.3f}".format(p1opt))
# A TopSpin quirk: using MSG() will block subsequent commands until user hits
# "OK". You can use ERRMSG(), as is done here. There are other ways around
# this, see MSG_nonmodal() in the POISE frontend script.
# (Optional) Run acquisition.
# ZG()
```

5.3 Terminating scripts which call POISE

So far we have seen how POISE can be included inside an AU programme or a Python script in TopSpin (a "parent script"). One problem that we haven't dealt with yet is how to kill the parent script when POISE errors out. In general, if POISE is called via XCMD(poise ...), then even if POISE fails, the parent script will continue running.

One way to deal with this is to use a trick involving the TI TopSpin parameter, which can be set to any arbitrary string. POISE, upon successful termination, will store the value of the cost function in the TI parameter. If it doesn't successfully run (for example if a requested routine or cost function is not found, or some other error), then TI will be left untouched.

In order to detect when POISE fails from the top-level script, we therefore:

1. Set TI to be equal to some sentinel value, i.e. any string whose exact value is just used as a marker. Note that this should not be a numeric value. You can set it to be blank if you like, but please read the note below.

Note: In an AU or Python programme, to set a parameter to a blank value, you have to set it to be a non-empty string that contains only whitespace. For example, in a Python script:

```
PUTPAR("TI", " ") # this will work
PUTPAR("TI", "") # this will NOT work!
```

TopSpin mangles empty strings: instead of putting an empty string in, it puts the string "0" in. On the other hand, if the string is not empty but contains whitespace, TopSpin automatically trims it to an empty string after it's been put in. I don't know why. The same applies to the STOREPAR macro in AU programmes.

2. Run POISE.

3. Check if TI is equal to that sentinel value. If it is, then quit unceremoniously with an error message of your choice.

Briefly, here is an example of this strategy in action. We show a Python script here, but the AU script is essentially the same, just with different function names.

```
PUTPAR("TI", "poise") # here "poise" serves as the sentinel value.
XCMD("poise p1cal -q")
# Here POISE should be done. If it succeeded then TI will no longer be "poise".
if GETPAR("TI") == "poise":
    raise RuntimeError("POISE failed!")
# After this you can continue with whatever you wanted to do.
```

Builtin cost functions

POISE comes with a few, very basic, builtin cost functions. These largely mirror those that are in TopSpin's native popt screen.

Note: Just like in popt, it is possible to use the dpl command in TopSpin to select a *portion* of the spectrum to be optimised. This stores the left and right region of the currently active view to the parameters F1P and F2P respectively. This works for all the builtin cost functions except for noe_1d.

More generally, any cost function that uses any of the get1d or get2d functions will respect the bounds placed in F1P and F2P. See *Custom cost functions* for a more in-depth explanation.

6.1 minabsint

Seeks to minimise the intensity of the magnitude-mode spectrum. The intensity is measured by integration of the entire spectral region, i.e. summation of every point.

Note that this is different from the MAGMIN criterion in popt, which (from what the Bruker documentation suggests) seeks to minimise the highest point in the magnitude-mode spectrum.

6.2 maxabsint

Seeks to maximise the intensity of the magnitude-mode spectrum.

6.3 minrealint

Seeks to minimise the intensity of the real spectrum (this is probably equivalent to INTMIN in popt).

Note that this does *not* behave in the same way as minabsint. Because the real spectrum can have negative peaks, this essentially tries to maximise the intensity of negative peaks.

6.4 maxrealint

Seeks to maximise the intensity of the real spectrum (equivalent to INTMAX).

6.5 zerorealint

Seeks to make the intensity of the real spectrum as close as possible to zero (equivalent to ZERO).

6.6 noe_1d

Seeks to minimise the intensity of the spectrum, *except* for a region of 50 Hz centred on the parameter SP0FFS2 (which corresponds to the frequency of the selective pulse).

Since NOE crosspeaks are typically negative (and apk typically phases them to be so), this essentially seeks to maximise the intensity of the crosspeaks.

Please see the POISE paper for example usage.

6.7 epsi_gradient_drift

Performs EPSI processing on a 1D FID, and seeks to minimise the 'drift' seen in the echo positions. This is only valid when the echo locations are *supposed* to be constant, e.g. when there is no indirect-dimension evolution period (otherwise, echoes will be observed at indirect-dimension frequencies). This is to be used for optimising the positive/negative gradient balance in EPSI acquisitions.

Please see the POISE paper for example usage.

Builtin AU programmes

POISE is bundled with two basic AU programmes for 1D and 2D spectra respectively. These are mostly selfexplanatory, so the text is just given below. Of course, you can write your own AU programmes to be used with POISE routines.

These are the "default" AU programmes that POISE uses (for 1D and 2D spectra respectively) if an AU programme is not specified with the routine.

7.1 poise_1d

ZG		
EFP		
APBK		
QUIT		

The APBK command on older versions of TopSpin falls back to APK then ABS, so can be safely used.

7.2 poise_2d

ZG	
XFB	
XCMD("apk2d")	
ABS2	
QUIT	

Custom cost functions

All user-defined cost functions are stored inside the file:

\$TS/exp/stan/nmr/py/user/poise_backend/costfunctions_user.py

where \$TS is your TopSpin installation path. In order to modify or add cost functions, you will need to edit this file (with your favourite text editor or IDE).

The corresponding file containing builtin cost functions is costfunctions.py. You *can* edit this file directly: if you add a cost function there, it will work. However, there are two risks with this. Firstly, if you ever reinstall POISE, this file will be reset to the default (whereas costfunctions_user.py will not). Secondly, any cost functions defined in costfunctions_user.py will shadow (i.e. take priority over) the cost functions defined in costfunctions.py if they have the same name.

8.1 The rules for cost functions

Cost functions are defined as a standard Python 3 function which takes no parameters and returns a float (the value of the cost function).

- 1. Do write a useful docstring if possible: this docstring will be shown to the user when they type poise -1 into TopSpin (which lists all available cost functions and routines).
- 2. The spectrum under optimisation, as well as acquisition parameters, can be accessed via helper functions. These are described more fully below.
- 3. Never print anything inside a cost function directly to stdout. This will cause the optimisation to stop. If you want to perform debugging, use the log function described below.
- 4. To terminate the optimisation prematurely and return the best point found so far, raise Cost-FunctionError(). See below for more information.

8.2 Accessing spectra and parameters

The most primitive way of accessing "outside" information is through the class _g, which is imported from shared.py and contains a series of global variables reflecting the current optimisation. For example, _g. p_spectrum is the path to the procno folder: you can read and parse the 1r file inside this to get the real spectrum as a numpy.ndarray (for example).

class nmrpoise.poise_backend.shared.**_g** Class to store the "global" variables.

Attributes

optimiser [str from {'nm', 'mds', 'bobyqa'}] The optimiser being used.

routine_id [str] The name of the routine being used.

- p_spectrum [Path] The path to the procno folder of the spectrum just acquired. (e.g. /path/to/data/1/pdata/1)
- **p_optlog** [Path] The path to the currently active poise.log file.
- **p_errlog** [Path] The path to the currently active poise_err_backend.log file.
- **maxfev** [int] The maximum number of function evaluations specified by the user. Can be zero, indicating no limit (beyond the hard limit of 500 times the number of parameters).
- **p_poise** [Path] The path to the \$TS/exp/stan/nmr/py/user/poise_backend folder.
- spec_f1p [float or tuple of float] The F1P parameter. For a 1D spectrum this is a float. For a 2D spectrum this is a tuple of floats (indirect, direct) corresponding to the values of F1P in both spectral dimensions.
- **spec_f2p** [float or tuple of float] The F2P parameter.
- **xvals** [list of ndarray] The points sampled during the optimisation, in chronological order. These are ndarrays which contain the values of the parameters being optimised at each spectrum acquisition. The parameters are ordered in the same way as specified in the routine.
- **fvals** [ndarray] The values of the cost functions calculated at each stage of the optimisation.

However, this is quite tedious and error-prone, so there are a number of helper methods which use these primitives. All the existing cost functions (inside costfunctions.py) only use these helper methods. All of these methods are stored inside cfhelpers.py and are already imported by default.

The ones you are likely to use are the following:

nmrpoise.poise_backend.cfhelpers.make_p_spec(path=None, expno=None, procno=None)

Constructs a Path object corresponding to the procno folder <path>/<expno>/pdata/<procno>. If parameters are not passed, they are inherited from the currently active spectrum (_g.p_spectrum).

Thus, for example, make_p_spec(expno=1, procno=1) returns a path to the spectrum with EXPNO 1 and PROCNO 1, but with the same name as the currently active spectrum.

Parameters

path [str or Path, optional] Path to the main folder of the spectrum (one level above the expno folders).

expno [int, optional]

procno [int, optional]

Returns

p_spec [Path] Path pointing to the requested spectrum.

nmrpoise.poise_backend.cfhelpers.get1d_fid(remove_grpdly=True, p_spec=None) Returns the FID as a ndarray.

Parameters

- **remove_grpdly** [bool, optional] Whether to remove the group delay (to be precise, it is shifted to the end of the FID). Defaults to True.
- p_spec [Path, optional] Path to the procno folder of interest. (The FID is taken from the expno folder two levels up.) Defaults to the currently active spectrum (i.e. _g. p_spectrum).

Returns

ndarray Complex-valued array containing the FID.

```
nmrpoise.poise_backend.cfhelpers.get1d_real(bounds=", p_spec=None)
```

Return the real spectrum as a ndarray. This function accounts for TopSpin's NC_PROC variable, scaling the spectrum intensity accordingly.

Note that this function only works for 1D spectra. It does *not* work for 1D projections of 2D spectra. If you want to work with projections, you can use get2d_rr to get the full 2D spectrum, then manipulate it using numpy functions as appropriate. A documented example can be found in the asaphsqc() function in costfunctions.py (commented out by default).

The *bounds* parameter may be specified in the following formats:

- between 5 and 8 ppm: bounds="5..8" OR bounds=(5, 8)
- greater than 9.3 ppm: bounds="9.3.." OR bounds=(9.3, None)
- less than -2 ppm: bounds="..-2" OR bounds=(None, -2)

Parameters

- **bounds** [str or tuple, optional] String or tuple describing the region of interest. See above for examples. If no bounds are provided, uses the F1P and F2P processing parameters, which can be specified via dpl. If these are not specified, defaults to the whole spectrum.
- **p_spec** [Path, optional] Path to the procno folder of interest. Defaults to the currently active spectrum (i.e. _g.p_spectrum).

Returns

ndarray Array containing the spectrum or the desired section of it (if bounds were specified).

```
nmrpoise.poise_backend.cfhelpers.get1d_imag(bounds=", p_spec=None)
Same as get1d_real, except that it reads the imaginary spectrum.
```

nmrpoise.poise_backend.cfhelpers.get2d_rr(f1_bounds=", f2_bounds=", p_spec=None)

Return the real part of the 2D spectrum (the "RR" quadrant) as a 2D ndarray. This function takes into account the NC_PROC value in TopSpin's processing parameters.

The *f1_bounds* and *f2_bounds* parameters may be specified in the following formats:

- between 5 and 8 ppm: f1_bounds="5..8" OR f1_bounds=(5, 8)
- greater than 9.3 ppm: f1_bounds="9.3.." OR f1_bounds=(9.3, None)
- less than -2 ppm: f1_bounds="..-2" OR f1_bounds=(None, -2)

Parameters

- **f1_bounds** [str or tuple, optional] String or tuple describing the indirect-dimension region of interest. See above for examples. If no bounds are provided, uses the 1 F1P and 1 F2P processing parameters, which can be specified via dp1. If these are not specified, defaults to the whole spectrum.
- **f2_bounds** [str or tuple, optional] String or tuple describing the direct-dimension region of interest. See above for examples. If no bounds are provided, uses the 2 F1P and 2 F2P processing parameters, which can be specified via dpl. If these are not specified, defaults to the whole spectrum.

Returns

ndarray 2D array containing the spectrum or the desired section of it (if *f1_bounds* or *f2_bounds* were specified).

nmrpoise.poise_backend.cfhelpers.get2d_ri(f1_bounds=", f2_bounds=", p_spec=None)
Same as get2d_rr, except that it reads the '2ri' file.

nmrpoise.poise_backend.cfhelpers.get2d_ir(f1_bounds=", f2_bounds=", p_spec=None)
Same as get2d_rr, except that it reads the '2ir' file.

nmrpoise.poise_backend.cfhelpers.get2d_ii(f1_bounds=", f2_bounds=", p_spec=None)
Same as get2d_rr, except that it reads the '2ii' file.

nmrpoise.poise_backend.cfhelpers.getpar(par, p_spec=None)

Obtains the value of a numeric (acquisition or processing) parameter. Non-numeric parameters (i.e. strings) are not currently accessible! Works for both 1D and 2D spectra (see return type below), but nothing higher.

Parameters

par [str] Name of the parameter.

Returns

float or ndarray Value(s) of the requested parameter. None if the given parameter was not found.

For parameters that exist for both dimensions of 2D spectra, getpar() returns an ndarray consisting of (f1_value, f2_value). Otherwise (for 1D spectra, or for 2D parameters which only apply to the direct dimension), getpar() returns a float.

nmrpoise.poise_backend.cfhelpers.getndim(p_spec=None)

Obtains the dimensionality of the spectrum, i.e. the status value of PARMODE, plus one (becaus PARMODE is 0 for 1D spectra, etc.)

Parameters

Returns

int Dimensionality of the spectrum.

nmrpoise.poise_backend.cfhelpers.getnfev()

Returns the number of NMR spectra evaluated so far. This will be equal to 1 (not 0) the first time the cost function is called, since the cost function is only called after the NMR spectrum has been acquired.

8.3 Logging

As noted above, printing anything to stdout will cause the optimisation to crash. The reason for this is because stdout is reserved for communication between the POISE backend (the Python 3 component) and frontend (which runs in TopSpin). Please use log() instead, which will print to the poise.log file in the expno folder. It works in exactly the same way as the familiar print(), and accepts the same kind of arguments.

nmrpoise.poise_backend.cfhelpers.log(*args)

Prints something to the poise.log file.

If this is called from inside a cost function, the text is printed *before* the cost function is evaluated, so will appear above the corresponding function evaluation.

Parameters

args The arguments that log() takes are exactly the same as those of print().

Returns

None

8.4 Premature termination

In order to terminate an optimisation prematurely, you can raise any exception you like, for example with raise ValueError. POISE will stop the optimisation, and the error will be displayed in TopSpin. The drawback of this naive approach is that *no information from the incomplete optimisation will be retained*. That means that even if you have found a point that is substantially better than the initial point, it will not be saved.

If you want to terminate *and* return the best point found so far, please raise a CostFunctionError instead of any other exception (such as ValueError). CostFunctionError takes up to 2 arguments: the number of arguments to supply depends on whether you want to include the final point which caused the error to be raised.

If you want to discard the last point, i.e. just stop the optimisation right away, then raise a CostFunctionError with just 1 argument. The argument should be the message that you want to display to the user:

It is possible, and probably advisable, to use this string to show the user helpful information (further steps to take, or the value of the cost function, for example).

Alternatively, you may want the current point (and the corresponding cost function value) to be saved as part of the optimisation. For example, it may be the case that a certain threshold is "good enough" for the cost function and any value below that is acceptable. In that situation, you would want to raise CostFunctionError once the cost function goes below that threshold, but *also* save that point as the best value. To do so, pass the value of the cost function as the *second* parameter when raising CostFunctionError:

The value passed as the second argument can in general be any number. If you want to make sure that the final point (which raised the CostFunctionError) is *always* chosen to be the optimal point, then you can do this:

raise CostFunctionError("A message here", -np.inf)

Because -np.inf is smaller than any other possible number, it will always be picked as the "best" point.

8.5 Examples

There are a number of cost functions which ship with POISE. These can all be found inside the costfunctions.py file referred to above. This file also contains a number of more specialised cost functions, which were used for the examples in the POISE paper. (Instead of opening the file, you can also find the source code on GitHub.)

A number of these are thoroughly commented with detailed explanations; do consider checking these out if you want more guidance on how to write your own cost function.

Developer notes

POISE does not contain all that much code. However, if you are thinking of modifying it, it helps to have an understanding of the distinction between the *frontend* and the *backend*. This graphic (taken from an old presentation) is a slightly more technical description of the flowchart presented in the POISE paper, and shows which part is responsible for which step. It's a good starting point for understanding how POISE works internally.



Beyond this, we recommend reading the source code of POISE: start at the frontend script (poise.py, which is run from inside TopSpin), then go to the backend script (poise_backend/backend.py) at the appropriate time (when the frontend launches it as a subprocess). The source code is quite thoroughly commented.

There are two main things to point out. Probably the most important thing worth mentioning is the location of the relevant files. The backend script is always ran from *TS/py/user/poise_backend* (putting it here allows the frontend to access it much more easily). The entire *TS/py/user/poise_backend* folder is treated as if it is a Python package, by virtue of some code near the top of poise_backend.py:

```
if __name__ == "__main__" and __package__ is None:
    __package__ = "poise_backend"
    sys.path.insert(1, str(Path(__file__).parents[1].resolve()))
    __import_(__package__)
```

This allows relative imports of the other files in the same directory, such as costfunctions.py, in which the user-defined cost functions reside.

Note: These backend files *also* reside inside the Python 3 site-packages directory, where all packages are installed to. However, these files will *never* be used by POISE. So, there is nothing to be gained by modifying these at all.

The only file inside the site-packages directory which has any effect is nmrpoise/__init__.py, where the parse_log() function is defined. This allows you to (for example) run:

>>> from nmrpoise import parse_log

Python will look inside the site-packages directory, not \$TS/py/user, to find this function.

The second thing is that you should never, ever, do anything with the backend's stdin and stdout, because these are exclusively reserved for communication with the frontend. So you should never print anything from the backend, since the frontend will just interpret that as an error. This applies to all files inside the \$TS/py/user/poise_backend directory, including cost functions, which is why custom cost functions should always use POISE's log function instead of plain old print.

Note: If you really just want to do some quick-and-dirty debugging, you *can* actually use this behaviour to your advantage. The frontend will echo any "invalid" message it receives from the backend, so if you print some unexpected text from the backend (on purpose), you should see it pop up as a TopSpin message when you run an optimisation. This is slightly less hassle than printing to a file and opening the file.

9.1 Testing

Tests are carried out using the excellent pytest and tox tools. To run all tests, simply run:

```
pip install tox
tox
```

from anywhere inside the nmrpoise directory. This runs tests on Python 3.6, 3.7, and 3.8. If you only have one of these versions, use:

tox -e py38 # or py36 or py37

To build the Sphinx documentation, use:

tox -e docs

The HTML documentation will be built in docs/dirhtml, and the PDF documentation in docs/latex (this assumes you have a working installation of pdflatex on your system).