Advanced technique for measuring relative length changes under control of temperature and helium-gas pressure

Y. Agarmani,¹ S. Hartmann,¹ J. Zimmermann,¹ E. Gati,^{1,2} C. Delleske,¹ U. Tutsch,¹ B. Wolf,¹ and M. Lang¹ ¹⁾Institute of Physics, Goethe University Frankfurt, Max-von-Laue-Straße 1, 60428 Fourtheat on Main Conversion

²⁾Max-Planck-Institute for Chemical Physics of Solids, Nöthnitzer Straße 40,

01187 Dresden, Germany

(*The author to whom correspondence may be addressed: agarmani@physik.uni-frankfurt.de)

(Dated: December 12, 2022)

We report the realization of an advanced technique for measuring relative length changes $\Delta L/L$ of mm-sized samples under control of temperature (*T*) and helium-gas pressure (*P*). The system, which is an extension of the apparatus described in Manna *et al.*, Rev. Sci. Instrum. **83**, 085111 (2012), consists of two ⁴He-bath cryostats each of which houses a pressure cell and a capacitive dilatometer. The interconnection of the pressure cells, the temperature of which can be controlled individually, opens up various modes of operation to perform measurements of $\Delta L/L$ under variation of temperature and pressure. Special features of this apparatus include the possibilities (1) to increase the pressure to values much in excess of the external pressure reservoir, (2) to substantially improve the pressure stability during temperature sweeps, (3) to enable continuous pressure sweeps both with decreasing and increasing pressure, and (4) to simultaneously measure the dielectric constant of the pressure-transmitting medium helium, $\varepsilon_r^{\text{He}}(T, P)$, along the same T-P trajectory as used for taking the $\Delta L(T, P)/L$ data. The performance of the setup is demonstrated by measurements of relative length changes ($\Delta L/L$)_T at T = 180 K of single crystalline NaCl upon continuously varying the pressure in the range 6 MPa $\leq P \leq 40$ MPa.

I. INTRODUCTION

The coefficient of thermal expansion, describing how the size of a material changes upon changing the temperature, is a fundamental thermodynamic quantity. In fact, measurements of relative length changes, $\Delta L_i(T)/L_i = [L_i(T) - L_i)/L_i$ $L_i(T_0)]/L_i(T_0)$, with L_i the sample length along the axis i and T_0 a reference temperature, have proven a powerful tool for exploring the electronic, magnetic and lattice properties and their directional dependencies. These investigations have been of great importance for studying cooperative phenomena such as phase transformations or crossovers and their variations with external parameters such as temperature and magnetic fields. Among the various methods for measuring length changes of a sample, capacitive dilatometry stands out due to the extraordinarily high resolution of $\Delta L/L > 10^{-10}$ which exceeds other techniques, including X-ray diffraction¹, optical interferometers² and strain gauges³, by at least one order of magnitude. A first implementation for dilatometry under high pressure using a capacitive technique was presented in Ref. 4. In a recent work, we reported on the development of an apparatus, which enables to expand the range of application of capacitive dilatometry to finite hydrostatic pressure⁵, providing a more sensitive alternative to commonly used methods for dilatometric studies under pressure such as X-ray diffraction and strain gauges 6,7 . The latter two techniques along with ultrasonic studies⁸ are also used to measure the compressibility of solids (see Ref. 9 for a review), albeit with reduced sensitivity. The setup described in Ref. 5 is based on the use of a high-resolution capacitive dilatometer in combination with a variable helium-gas pressure environment for measuring relative length changes under control of temperature and hydrostatic pressure. This technique was successfully applied in measurements of $\Delta L(T, P)/L$ on an organic chargetransfer salt, which is located close to the Mott metal-insulator transition, for investigating the coupling of the critical electronic system to the lattice degrees of freedom¹⁰. The results obtained by Gati et al., disclosing a pronounced non-linear strain-stress relation on approaching the critical endpoint of the Mott transition¹⁰, has demonstrated the high potential of dilatometric studies under variable helium-gas pressure. This technique is particularly useful for those investigations where a fine and *in situ* pressure tuning is demanded. In general, two modes of operation are of interest including (A) temperature sweeps over a sufficiently wide range of temperatures at $P \approx \text{const.}$ conditions, as well as (B) pressure sweeps at $T \approx \text{const.}$ conditions. Depending on the problem under investigation, these experiments may pose challenges on the process control. This includes the stabilization of pressure while ramping the temperature up and down (A) as well as a sufficiently good temperature control in studies where the pressure is varied (B). The latter also comprises the need to change the pressure in a continuous way at a desired and sufficiently low rate to ensure thermal equilibrium. Moreover, in specific cases, measurements at both decreasing and increasing pressure at $T \approx \text{const.}$ could be of interest, in particular for studying pressure-induced phase changes. The setup described here, consisting of two interconnected pressure cells each of which is equipped with its own temperature-control unit, considerably improves the range of applications and the quality of the pressure- and temperature control as compared to a system with a single measuring cell and a gas reservoir controlled at room temperature (cf. Ref. [5]). Here, the following challenges have been addressed successfully: 1) The maximum accessible pressure, P_{max} , has been increased significantly by drastically reducing the gas volume at room temperature. In particular, P_{max} is no longer limited by the maximum pressure of the external pressure reservoir - a value

⁶⁰⁴³⁸ Frankfurt am Main, Germany



Figure 1. Sketch of the setup with all relevant components. A gas bottle (a), loaded with helium gas at $P \le 30$ MPa, serves as a pressure reservoir. The gas is injected into a two-stage compressor unit (b). A digital manometer (c) is used to determine the pressure. The compressor is connected via CuBe capillaries to the sample- (d) and reference-setup (e) each of which consists of a ⁴He-bath cryostat and an insert including a pressure cell and a dilatometer cell. The pressure loading of the pressure cells in (d) and (e) can be controlled individually through the valves (IV) and (V). The pressure can be released into the recovery line (f) via valve (VI).

which, due to safety regulations, strictly limits gas-containing parts at room temperature by placing a limit on the product of volume × pressure. 2) The unwanted change of pressure during a temperature sweep – an inevitable consequence of working with a closed system – has been drastically reduced. 3) Continuous temperature-driven pressure sweeps can be performed both with decreasing and increasing pressure for a wide range of sweep rates with a quality unmatched by any mechanical pressure control. 4) As the capacitive dilatometer is surrounded by the pressure-transmitting medium helium, a precise knowledge of its dielectric constant $\varepsilon_r^{He}(T,P)$ is required in the data analysis, since ε_r^{He} varies significantly with temperature and pressure at low temperatures. Our setup allows for a simultaneous measurement of $\varepsilon_r^{He}(T,P)$ along the same *T-P* trajectory as used for taking the $\Delta L(T,P)/L$ data.

II. INTERCONNECTION OF TWO PRESSURE CELLS

Figure 1 shows a sketch of the setup with all relevant components. The system consists of two He-bath cryostats (d) and (e) both of which are equipped with a pressure cell and a capacitive dilatometer. The cryostat (d) and its insert, referred to as sample setup from here on, are described in detail in Ref. 5. This setup has been complemented by adding a second, almost identical system, referred to as a reference setup. Both pressure cells are connected via capillaries to a compressor unit (b) and a helium reservoir (a). The compressor unit is equipped with high-pressure valves (I) - (VI) for controlling the gas flow.

The helium gas, used as pressure-transmitting medium, is provided by utilizing a standard gas bottle (a) as a pressure reservoir. The gas bottle can be loaded with a maximum allowed pressure of 30 MPa. The purity of the helium used is typically 99.999 %. The bottle is connected via valve (I) to the two-stage helium-gas compressor U11 300 MPa (b), which is commercially available (for more details see Ref. 5) and was developed in cooperation with the Institute of High Pressure Physics, Polish Academy of Sciences, Unipress Equipment Division, in the following abbreviated as Unipress. In order to reach pressures P > 30 MPa, valve (I) needs to be closed and the two stages of the compressor pressurize the helium gas to the desired value. Whereas the first (low-pressure) stage enables pressures P < 70 MPa, the second (high-pressure) stage reaches pressures of 70 MPa $\leq P \leq 250$ MPa. The highpressure stage is connected via valve (III) to two CuBe capillaries of similar length (inner diameter 0.3 mm, outer diameter 3 mm), one of which is connected to the sample setup, the other one is connected to the reference setup. The pressure cells and capillaries are manufactured by Unipress. Valves (IV) and (V) allow for an individual pressure loading of the pressure cells. Valve (VI) enables to relief the pressure into the helium recovery system.

The two identical pressure cells, used in the sample- and reference setups, are made of CuBe (see Refs. 5 and 10 for details) with an inner/outer diameter of 36 mm/56 mm and a maximum permissible pressure of 250 MPa. Each pres-

sure cell houses a capacitive dilatometer cell¹¹ with the same design, manufactured by Kuechler Innovative Measurement Technology¹². The pressure cell of the sample setup (d) is opened regularly for mounting the sample to be investigated. On the other hand, the reference-pressure cell (e) is kept closed and is opened only for maintenance. A highpurity (99.99%) aluminum sample of 4.5 mm length is installed in the reference-dilatometer cell and serves as a reference material. The criteria used for selecting this material are (i) a known and well-reproducible thermal expansion and (ii) a small compressibility ensuring a negligibly small effect of pressure in the pressure range used here. A digital manometer LEO5 (KELLER AG) (c) is used to determine the pressure with a resolution of $|\Delta P| \leq 0.1$ MPa for pressures up to 100 MPa and $|\Delta P| \leq 0.03$ MPa for pressures up to 30 MPa. The home-made inserts, used for the ⁴He cryostats (d) and (e), were designed with special focus placed on ensuring a good temperature control of the pressure cells. The essential construction elements, described in Ref. 5, include the pressure cell, the capillary, the capillary heaters, the thermometers (Cernox®, LakeShore Cryotronics), and the inner vacuum can (IVC) made of copper (diameter 66 mm). This can is wrapped with heating foil (Kapton[®], thermofoil, Minco company) to ensure a homogeneous heat input into the pressure cells. The IVC is surrounded by an outer vacuum can made of stainless steel (diameter 75 mm). Both cans are connected to pump lines that can be vented independently with ⁴He exchange gas. The IVCs of both inserts are filled with ⁴He gas at low pressure of typically $0 \le P \le 0.01$ MPa (at 300 K) for ensuring a good thermal contact between the heater and the pressure cell, see Ref. 5 for details of the inserts including thermometry and temperature control unit.

The accessible temperature range for the inserts is $1.4 \text{ K} \le T \le 300 \text{ K}$. The maximum accessible pressure P_{max} is constrained by the operating principle of the capacitive dilatometer, requiring a freely movable upper capacitor plate. Thus, P_{max} is limited by the solidification line $T_{\text{sol}}^{\text{He}}(T,P)$ of the pressure-transmitting medium helium, see Ref. 13 and references cited therein, which is 60 MPa at 10 K and 1.45 GPa at 77 K, for example.

III. ACCESSIBLE MODES OF OPERATION

Depending on the physical problem under investigation and the desired temperature- and pressure range of the experiment, different modes of operation are accessible. For experiments under finite pressure, where the variations in temperature and pressure imply significant changes of the dielectric constant of helium, simultaneous measurements of $\varepsilon_r^{\text{He}}(T, P)$ can be performed. This procedure will be described in section III A. If temperature sweeps at $P \approx \text{const.}$ conditions are demanded, special measures can be taken to keep temperatureinduced pressure changes small, see section III B. In section III C we describe a mode of operation, which allows for pressure sweeps with both decreasing and increasing pressure. In the description of the various modes of operation, we use a simplified sketch of the relevant parts of the inserts of the sample- and reference setups as displayed on the right side of figure 1. While the sample setup (Fig. 1 (d)) contains the dilatometer cell with the sample of interest, the reference setup (Fig. 1 (e)) contains the dilatometer cell with the aluminum sample installed. The valves (III) to (VI) allow for an individual control of the gas flow. Both setups are supplied with their own temperature-control units enabling to run independent temperature-time (*T*-*t*) profiles.

A. Simultaneous determination of relative length change and the dielectric constant of the pressure-transmitting medium

The pressure-transmitting medium helium in the pressure cell acts as dielectric between the capacitor plates of the dilatometer cell. It therefore affects the capacitance *C* of the dilatometer and thus has to be taken into account in converting changes of the capacitance into relative length changes, see Eq. 2 below and Ref. 5. As there is only limited data on $\varepsilon_r^{\text{He}}(T,P)$ available in the literature (see Refs. 5 and 14), a simultaneous determination of both the relative length change of the sample $\Delta L(T,P)/L$ and $\varepsilon_r^{\text{He}}(T,P)$ along the same *T*-*P* trajectory brings major benefits. As will be shown below, $\varepsilon_r^{\text{He}}(T,P)$ can be derived from measuring capacitance changes in the reference setup. In vacuum, the capacitance *C* in the reference setup is described by:

$$C(T, P=0) = \varepsilon_0 \frac{A}{d(T, P=0)} , \qquad (1)$$

with the vacuum permittivity ε_0 , the area of the capacitor plates *A*, and the distance between the plates *d*. If the system is under helium pressure *P*, this expression changes into:

$$C(T,P) = \varepsilon_0 \varepsilon_{\rm r}^{\rm He}(T,P) \frac{A}{d(T,P)} , \qquad (2)$$

with the dielectric constant of helium $\varepsilon_{\rm r}^{\rm He}(T, P)$. It is reasonable to assume that in the pressure range discussed here, i.e., $P \leq 250$ MPa, the area *A* does not change significantly with pressure, see below for explanation. By combining Eq. 1 and 2 we obtain:

$$\varepsilon_{\rm r}^{\rm He}(T,P) = \frac{C(T,P)}{C(T,P=0)} \cdot \frac{d(T,P)}{d(T,P=0)} .$$
 (3)

The compressibility of the dilatometer cell (made of a copper beryllium alloy with a small beryllium concentration of $1.84 \%^{12}$) and the reference sample (aluminium sample installed in the reference setup) is rather small, with $\kappa_{CuBe}(300 \text{ K}) = 7.9 \cdot 10^{-6} \text{ MPa}^{-1.15}$ and $\kappa_{Al}(300 \text{ K}) = 13.85 \cdot 10^{-6} \text{ MPa}^{-1.16}$, respectively. Thus, the corresponding pressure-induced changes in both *d* (caused by length changes of aluminum) and the dimensions of the dilatometer cell give rise to changes in the capacitance of $\mathcal{O}(10^{-3} \text{ pF})$ for capacitance of typically $C \approx 13 \text{ pF}$. This has to be compared with changes in $\varepsilon_{r}^{He}(T, P)$ throughout the ranges of temperature $(1.4 \text{ K} \leq T \leq 300 \text{ K})$ and pressure ($P \leq 250 \text{ MPa}$) covered in the experiment, which cause changes in the capacitance of



Figure 2. Sketch of the mode of operation for a simultaneous determination of length changes of a sample $\Delta L(T, P)/L$ and the dielectric constant of the pressure-transmitting medium helium $\varepsilon_r^{\text{He}}(T, P)$. The temperature control units for both systems are operated in a way so that the temperature is the same for both setups, illustrated by the same color. The valves (IV) and (V) are open so that the pressure is the same in both pressure cells. For *T*-sweeps (a), the temperature in both setups is changed with the same rate $q_1 = q_2$, while *P* is approximately held constant. For *P*-sweeps (b), the pressure of both setups is changed with the same rate, while *T* is held constant, ensuring the same conditions for the sample- and reference-measurement.

 $\mathcal{O}(10^{-1} \text{ pF})$. It is therefore reasonable to set d(T, P = 0) = d(T, P) in Eq. 3, yielding:

$$\varepsilon_{\rm r}^{\rm He}(T,P) = \frac{C(T,P)}{C(T,P=0)} . \tag{4}$$

Thus, results obtained in measurements of C(T,P) in the reference setup, combined with C(T,P=0) data collected in an independent run at P = 0, allow to determine $\varepsilon_{\rm r}^{\rm He}(T,P)$ simultaneously with the measurements of $\Delta L(T,P)/L$.

A sketch of the corresponding idealized modes of operation is shown in figure 2. The valves (IV) and (V) are open so that the pressure equilibrates in both pressure cells. For *T*-sweeps (a), both pressure cells are stabilized at the same starting temperature T_{start} . Then the temperature is varied in the same way for both setups, normally by applying a linear *T*-*t* profile as sketched in figure 2(a). A typical ramp rate for such thermal expansion measurements is 1.5 K/h. Ideally, the pressure should stay constant during a temperature sweep. In reality, however, the pressure varies by an amount ΔP depending on the *T*-*P* range of operation, see section III B. For *P*-sweeps, shown in Fig. 2 (b), both pressure-cells are hold at the same constant temperature and the pressure is released [increased] by opening valve (VI) [(III)]. The change in pressure is the same for both pressure cells. No matter which option is chosen, the C(T,P) data for the reference setup are collected simultaneously with the $\Delta L(T,P)/L$ data for the sample under investigation.

B. Temperature sweeps with improved stabilization of pressure

The mode of operation described in this section reduces unwanted changes of pressure during temperature sweeps. For a closed system, where the volume and the number of Heparticles are kept constant, the ideal gas law implies that a change in the temperature ΔT gives rise to a change in pressure $\Delta P \propto \Delta T$. The volumes involved include the gas bottle (a) (cf. Fig. 1) $V_{\text{bottle}} = 5 \cdot 10^4 \text{ cm}^3$, the first- and second stage of the compressor (b) with maximum volumes of $V_{1^{\text{st}}, \text{max}} =$ 720 cm^3 and $V_{2^{\text{nd}}, \text{max}} = 103 \text{ cm}^3$, respectively, and the pressure cells $V_{\text{cell}} = 81.4 \text{ cm}^3$, each of which houses a dilatometer cell $V_{\text{dil}} \sim 20 \text{ cm}^3$. Both pressure cells are connected to the compressor by a capillary of length $l \sim 10 \text{ m}$ with an in-



Figure 3. Sketch of the mode of operation for improving the pressure stability during temperature sweeps. While the temperature is increased from T_1 to T_2 in the sample-pressure cell (illustrated by the red color in the IVC) at a rate q_1 , the temperature in the reference-pressure cell is decreased from T_2 to T_1 with $q_2 = -q_1$ (illustrated by the blue color in the IVC).

ner diameter of 0.3 mm. For measurements at $P \le 30$ MPa, the gas bottle remains connected to the system and serves as a large-volume room-temperature gas reservoir. The entire multi-component system can be considered to be in equilibrium ensuring the same pressure in all components of the system. As the volume of the gas bottle is much bigger than the volume of the pressure cell, $V_{\text{bottle}} \approx 600 \cdot V_{\text{cell}}$, the pressure changes within the cell, induced by changes of its temperature, are negligibly small.

For measurements at P > 30 MPa, the maximum allowed pressure for the gas bottle, valve (I) is closed and the two stages of the compressor are being activated. The higher the target pressure is, the smaller becomes the remaining volume of the pressurized system, i.e., $V_{1^{st}} \leq V_{1^{st}, \max}$ and $V_{2^{nd}} \leq$ $V_{2^{nd},max}$. Thus, with increasing pressure, the temperatureinduced pressure changes become more and more significant. This effect is further reinforced with decreasing the temperature. When high-precision measurements of $\Delta L(T, P)/L$ over a limited range of temperature $T_1 \leq T \leq T_2$ are demanded, where only small variations in pressure can be tolerated, the reference setup can be used for stabilizing the pressure. The corresponding mode of operation is sketched in figure 3. Prior to the measurements, the sample setup is stabilized at a temperature $T_{\text{start}} = T_1$ whereas the reference setup is stabilized at a temperature $T_{end} = T_2$. The pressure cells are kept connected via the open valves (IV) and (V) to ensure an exchange of gas particles. On starting the measurement sequence, the temperature of the sample setup is increased in a linear T-tfashion from T_1 to T_2 with the desired ramp rate $q_1 = \Delta T / \Delta t$. This process is accompanied by simultaneously decreasing the temperature of the reference system from T_2 down to T_1 at a rate $q_2 = -q_1$. With this mode of operation, the temperatureinduced changes in the sample pressure cell can be significantly reduced. Simulations of this effect based on the ideal gas law are shown in figure 4. In the calculations, the volume of the pressure cells and the maximum available volume of the two stages of the compressor are taken into account, whereas the volume of the capillaries has been neglected.

Figure 4(a) shows the calculated change in pressure re-

sulting from sweeping the temperature from 25 K to 35 K at two different starting pressures of 10 MPa (lower panel) and 30 MPa (upper panel). The red curves correspond to the standard operation mode where only one pressure cell (single) containing the sample is connected to the compressor giving rise to a significant increase of $\Delta P \sim 2$ MPa (at an initial pressure of 10 MPa) and $\Delta P \sim 6$ MPa (at 30 MPa). This can be compared to simulations of ΔP (blue curve) for the operation mode described above, where both cells (coupled) are connected to the compressor unit and their temperatures are swept at opposite rates $q_2 = -q_1$ as indicated in figure 3. Here, the corresponding P-T curves (blue curves) are characterized by



Figure 4. Simulations of the variation of pressure during a temperature sweep. For details of the simulations see main text. The red lines represent the change in pressure when only one pressure cell (single) is connected to the compressor unit, corresponding to the standard setup. The blue lines correspond to pressure changes when both pressure cells (coupled) are coupled and connected to the compressor unit while running a process shown in figure 3. (a) Variation in pressure for two different target pressures: The upper panel presents the calculation for P = 30 MPa while the lower panel shows the pressure changes for P = 10 MPa. (b) Variation in pressure for a target pressure of 10 MPa while sweeping from $T_1 = 25$ K to $T_2 = 45$ K for different operation modes, i.e., traversing the temperature interval in a single step of width $\Delta T = 20$ K, or in two steps of width $\Delta T = 10$ K (in green), as indicated in the figure.



Figure 5. (a) Sketch of the mode of operation for performing pressure sweeps (with increasing pressure) at T = const. in the sample-pressure cell. To this end the reference-pressure cell is warmed up at a rate q_2 resulting in an overall increase of pressure in both pressure-cells. (b) Reversed process: The sample-pressure cell is held at constant T while the reference-pressure cell is cooled down at a rate $-q_2$ resulting in an overall decrease of pressure in both pressure cells.

a dome shape with $P(T_1) = P(T_2)$ and a maximum pressure change of only $\Delta P \sim 0.2$ MPa (at 10 MPa) and $\Delta P \sim 0.6$ MPa (at 30 MPa). This mode of operation thus indicates a massive improvement in the stabilization of pressure during temperature sweeps. Figure 4(b) illustrates how the performance can be even further improved. The figure compares the variations in pressure, intended to be stabilized at 10 MPa, for three different modes of operation while sweeping the temperature from $T_1 = 25$ K to $T_2 = 45$ K. The red line indicates the pressure increase by using a single pressure cell, yielding an increase of $\Delta P \sim 3$ MPa. This effect can be reduced to $\Delta P \sim 0.6$ MPa by the simultaneous use of two pressure cells as described above. The ΔP can be even further reduced by dividing the temperature interval (T_1, T_2) in half, and performing consecutive runs in both subsections (dotted green line). The corresponding maximum pressure change amounts to only $\Delta P \sim 0.2$ MPa.

C. Pressure sweeps at an adjustable sweep rate

In the previous version of the setup described in Ref. 5, consisting of a single pressure cell connected to the compressor unit and gas bottle, continuous pressure sweeps were lim-

ited to a maximum pressure of $P_{\text{max}} = 30 \text{ MPa}$, the maximum allowed pressure for the gas bottle. In addition, measurements could be performed only on decreasing the pressure (by slightly opening valve (VI)), which implied a nonlinear P-t profile. As sketched in figure 5, the use of two interconnected setups opens up new possibilities for performing pressure sweeps at an adjustable sweep rate including an increased P_{max} . In the initial stage of the process, the pressure cells are connected via the open valves (IV) and (V) and are loaded with an initial pressure P_1 from the compressor unit. The connection to the compressor and backline, valves (III) and (VI), is then closed. While the sample-pressure cell is held at a constant temperature, the reference-pressure cell is warmed up with a certain rate q_2 . This process is accompanied by a smooth increase of pressure in the coupled system. Likewise, decreasing the temperature in the referencepressure cell at a rate of $-q_2$ induces a continuous decrease in pressure. The choice of the process parameters strongly depends on the targeted temperature and pressure range. In general, the lower the temperature and the higher the pressure of the sample are, the wider is the temperature change the reference-pressure cell has to pass. The possibility of performing pressure sweeps with increasing and decreasing pressure enables for example detailed investigations of phase transitions, in particular studying potential hysteresis effects. In addition, the interconnected two-pressure-cell system, which can be run independently from the compressor unit and gas bottle, enables controlled pressure sweeps up to pressures significantly higher than 30 MPa.

IV. TEST MEASUREMENTS OF THE VARIOUS MODES OF OPERATION

In the previous paragraph III we have described various modes of operations offering unique possibilities for performing experiments for which special and well-defined conditions are demanded. As the scenarios described there consider idealized cases, the performance under real conditions may show some deviations the size of which depend on the chosen parameter range. In this paragraph we describe the results of test measurements for the different modes of operation. For selected cases we compare the performance with model calculations for the idealized system.

A. Temperature sweeps at $P \approx \text{const.}$ conditions

This process in its idealized form is shown in Fig. 2(a). In figure 6 we show the variation of pressure, P(T), during temperature sweeps with the initial pressure set to P_0 as indicated in the figure. The conditions apply to typical thermal expansion experiments aiming at measuring relative length changes as a function of temperature over an extended range of temperatures, while attempting to maintain $P \approx \text{const.}$ conditions. The warming rate chosen of $q_1 = q_2 = 1.5$ K/h for the sample- and reference system ensures thermal equilibrium of all components inside the IVCs. Both pressure cells were interconnected by opening valves (IV) and (V), and connected to the compressor unit via an open valve (III). The figure demonstrates that for the experiment at lower pressure around 11 MPa (green) the increase in temperature from $T_1 \sim 140 \,\mathrm{K}$ to $T_2 \sim 175$ K results in a moderate change of pressure with $|\Delta P|/P \sim 5\%$. However, on increasing the target pressure to around 46 MPa (red), the pressure change is significantly higher reaching values of $|\Delta P|/P \sim 10\%$ upon traversing the same temperature window. As mentioned in section III B, this is due to the fact that at pressures P > 30 MPa the compressor stages have to be activated, which leads to a reduction in the remaining volume of the pressurized system. These two examples indicate that the performance in terms of pressure stability depends strongly on the chosen process parameters. In this mode of operation, corresponding to section III A, $P \approx$ const. conditions can only be realized at a sufficiently low pressure (due to the use of a larger gas volume) and/or for traversing a sufficiently narrow temperature interval.

B. Pressure sweeps at $T \approx \text{const.}$ conditions

In order to perform pressure sweeps at $T \approx \text{const.}$ conditions, corresponding to a realization of the process shown in



Figure 6. Experimental data showing the variation of pressure during *T*-sweeps for starting values $P_0 = 11.2$ MPa (green) and $P_0 = 46.4$ MPa (red). In these experiments both setups were interconnected and connected to the compressor unit as described in chapter III A and Fig. 2. The *T*-sweeps were performed in both setups simultaneously, using heating rates $q_1 = q_2 = 1.5$ K/h, while measuring the resulting variations in P(T).

Fig. 2 (b), special measures would be required for ensuring a smooth, approximately t-linear change of pressure. In principle, this could be realized to some extent by using a controllable metering valve (VI) which would allow for a precise control of the flow rate by which the gas is released from the pressure cells. In the present system, however, which lacks such a metering valve, the flow rate can be adjusted by slightly opening valve (VI) manually. An example of the resulting pressure change as a function of time is depicted in Fig. 7. It shows a relatively rapid initial drop in P(t), which turns into an approximately exponential decrease on a longer time scale. Apart from the initial phase of this process (t < 3 h), where the change in pressure is rather large $|\Delta P(t)|/\Delta t \ge 2$ MPa/h and accompanied by some temperature instabilities, a smooth, approximately exponential variation in P(t) can be realized in this way while keeping the variations in temperature small $|\Delta T|/T < 0.03$ %. Thus, this mode of operation is useful for experiments where decreasing pressure is desired/sufficient and the pressure range to be covered is not too large, typically $|\Delta P| < 25$ MPa. As mentioned above, however, this mode of operation cannot be used for pressure sweeps with increasing pressure.

C. Temperature sweeps with improved stabilization of pressure

Figure 8 shows the results of a temperature sweep experiment at $P \approx \text{const.}$, as sketched in Fig. 3, aiming at a high degree of pressure stabilization. In the experiment the sample

cell was warmed from $T_1 = 25$ K to $T_2 = 35$ K with a heating rate $q_1 = 1.5$ K/h while simultaneously cooling the reference cell with $q_2 = -q_1$. The target pressure was ~ 10 MPa. The observed variation of pressure, which qualitatively conforms to the simulations shown in Fig. 4(a), reveal a high degree of stability of $|\Delta P|/P < 0.7 \%$. For P(t) we observe a domelike shape, however, with some asymmetry. We assign this fact to small differences in the volumes of the pressure components for the sample (V_s) - and reference (V_r) setups. Indeed, accurate measurements of the involved volumes reveal $V_{\rm s}/V_{\rm r} \approx 0.93$, which can account for this asymmetry.

D. Pressure sweeps at an adjustable sweep rate

In Fig. 9 we show results for a pressure sweep experiment performed as sketched in Fig. 5 (a). In this experiment, it was intended that the temperature of the sample setup, T_s , remains constant while the temperature of the reference system, T_r , was varied by a rate of $q_2 = 5$ K/h to induce a smooth increase of pressure in the whole system. Prior to the pressure sweep the system was loaded with a starting pressure of $P_1 \sim 10$ MPa. By raising T_r from 15 K to 100 K the pressure increases smoothly from $P_1 \sim 10$ MPa to $P_2 \sim 33$ MPa. While in the initial phase of the sweep the pressure increase occurred at a rate $\Delta P / \Delta t \approx 2.4$ MPa/h, the rate decreased with time and reached a value of $\Delta P / \Delta t \approx 1$ MPa/h towards the end of the sweep. Apart from a small instability of T_s in the initial phase of the process (t < 1 h), T_s remained constant within $|\Delta T| < 0.025$ K, corresponding to $|\Delta T|/T < 0.02$ %.

The induced pressure changes depend sensitively on the pa-





Figure 9. Experimental data for the sample temperature T_s (black, left scale) and pressure P (red, right scale) as a function of time for a pressure sweep with intended $T \approx \text{const.}$ conditions. In this experiment both setups were interconnected and connected to the compressor unit as described in Fig. 6 (a). The pressure sweep was performed through heating the reference setup at a rate $q_2 = 5$ K/h (inset), while measuring the accompanied variations in P(t) in the sample pressure cell. Prior to the sweep the system was loaded with a starting pressure of ~ 10 MPa.

30 Ρ T (K) P (MPa) t (h) 0 4 8 12 16 151.9 100 T_r 75 T(K) 50 151.8 10 25 0 4 12 16 8 t (h)

Figure 7. Experimental data taken in the sample setup for temperature T_s (black, left scale) and pressure P (red, right scale) as a function of time during a pressure sweep with intended $T \approx \text{const. con-}$ ditions. In this experiment both setups were interconnected and connected to the compressor unit as described in Fig. 2 (b). The pressure sweeps were performed simultaneously in both setups by slightly opening valve (VI) and releasing the pressure into the recovery line while measuring the associated temperature changes T(t).



time during a temperature sweep with intended $P \approx \text{const.}$ conditions. In this experiment, both setups were interconnected and connected to the compressor unit as described in Fig. 3. The temperature sweeps were performed simultaneously in both setups using rates $q_1 = -q_2 = 1.5$ K/h, while the associated pressure changes P(t)

Figure 8. Experimental data for the evolution of the sample tem-

rameters chosen in the reference system. In order to generate a high increase in pressure, a low starting temperature should be selected. Since a *t*-linear heating in the reference cell does not result in a *t*-linear increase of pressure in the sample system,



Figure 10. (a) Relative length change $(\Delta L(P)/L)_T$ as a function of pressure *P* at T = 180 K for single crystalline NaCl with L(300 K) = 6.8 mm. (b) Results of the dielectric constant of helium $\varepsilon_r^{\text{He}}(180 \text{ K}, P)$ derived from simultaneous measurements along the same *T*-*P* trajectory in the reference setup. (c) Uniaxial isothermal compressibility, κ_T^{uni} , of NaCl, derived as described in the text. Representative error bars resulting from the statistical treatment of the raw data and the estimated systematic uncertainties.

a suitable, non-linear T_r -t profile can be chosen to generate an approximately t-linear increase in P(t) if required.

V. MEASUREMENTS OF RELATIVE LENGTH CHANGES ($\Delta L/L$) $_{\mathcal{T}}$ AS A FUNCTION OF PRESSURE ON SODIUM CHLORIDE

In order to test the performance of the setup, measurements of relative length changes $(\Delta L/L)_{T\sim\text{const}}$ were conducted as a function of pressure on a single crystal of NaCl. This material was selected for these test studies due to its relatively high compressibility and its simple (cubic) structure. These facts along with the lack of any isomorphic transition up to 30 GPa^{17,18} make this material a frequently used pressurecalibration standard^{19–21}. Due to the large amount of experimental data on this material, there have been various approaches to derive the equation of states, see Refs. 19–23.

While the majority of studies were focusing on high pressures $P \ge 100$ MPa and temperatures $T \ge 300 \text{ K}^{24-30}$, much less is known for the *P*-*T* range of interest here, i.e., for pressures P < 100 MPa and temperatures T < 300 K. Nevertheless, the available data allow for an extrapolation into the low-*P* range with sufficient accuracy, yielding an isothermal volume compressibility at T = 300 K of $\kappa_T = (4.20 \pm 0.01) \cdot 10^{-5} \text{ MPa}^{-1}$, see Refs. 9 and 28 and references cited therein. In addition, an extrapolation of the isothermal bulk modulus for lower temperatures can be obtained from Refs. 27, 30, and 31, which is useful for our purposes.

In Fig. 10(a) we show results of relative length changes

 $(\Delta L/L)_T$ as a function of pressure for 6 MPa $\leq P \leq$ 40 MPa at a constant temperature $T = 180 \,\mathrm{K}$ taken on a single crystalline sample of NaCl. The sample length at 300 K was L = 6.8 mm. The experiment was performed in a mode of operation as described above in paragraph IV B, i.e., both setups were interconnected and connected to the compressor unit as described in Fig. 2 (b). *P*-sweeps were performed simultaneously in both setups by slightly opening valve (VI) and releasing the pressure into the recovery line. This resulted in an initial sweep rate of $|\Delta P|/\Delta t \sim 2$ MPa/hour for $P \leq 40$ MPa, cf. Fig. 8. During the *P*-sweep the variation in temperature was $|\Delta T| \leq 1$ K. The experiment allowed for measurements of $(\Delta L/L)_{T\approx 180 \text{ K}}$ in the sample setup while simultaneously measuring the dielectric constant of helium $\varepsilon_r^{\text{He}}(T, P)$ for the same T-P trajectory. The obtained results for the dielectric constant of helium $\varepsilon_r^{\text{He}}(T, P)$ are shown in Fig. 10 (b) for the corresponding pressure range. These $\varepsilon_r^{\text{He}}(T,P)$ data were used to transform the measured capacitance changes $(\Delta C/C)_{T\approx 180 \text{ K}}$ into relative length changes $(\Delta L/L)_{T\approx 180 \text{ K}}$ by using Eq. 2, see also Ref. 5. The data in Fig. 10 (a) yield a smooth decrease in $(\Delta L(P)/L)_T$ the slope of which gradually decreases with increasing pressure, indicating a hardening of the material. The high relative resolution of the measurement becomes even more clear by looking at the uniaxial compressibility $\kappa_T^{\text{uni}} = -1/L \cdot (\partial L/\partial P)_T$, approximated by $\kappa_T^{\text{uni}} \approx -1/L \cdot (\Delta L/\Delta P)_T$, shown in Fig. 10 (c). For determining the differential quotient, the $\Delta L/L$ -data were divided into equal pressure intervals of width 2.5 MPa in each of which a mean slope is determined by linear regression. The data, in the limit $P \rightarrow 0$, correspond to a volume compressibility of $\kappa_T^{\text{vol}} = 3 \cdot \kappa_T^{\text{uni}} = (3.60 \pm 0.24) \cdot 10^{-5} \text{ MPa}^{-1}$. This value can be compared with the extrapolated value (for T = 180 K) of the isothermal bulk modulus $B_T = 1/\kappa_T^{\text{vol}}$ of 25.6 GPa reported in Refs. 27 and 30, corresponding to $\kappa_T^{\text{vol}} = 3.9 \cdot 10^{-5} \text{ MPa}^{-1}$. Note that the calculation of the equation of state, which forms the basis for the latter value, is determined from experimental quantities associated with relative errors of typically a few percent, cf. Ref. 27. Even though the present apparatus has been designed aiming at a high relative resolution, the above comparison of compressibility data for NaCl demonstrates a satisfactory absolute accuracy.

VI. CONCLUSION

We have realized an advanced technique for measuring relative length changes of mm-sized samples under control of temperature and helium-gas pressure. The key elements of the apparatus include a room-temperature helium-gas reservoir, which is connected to two He-bath cryostats each of which is equipped with a pressure cell and a capacitive dilatometer. The system enables to run different, well-controlled T(t) and P(t) profiles which can be identical for both setups or individually customized for generating a desired T-P trajectory in the sample setup. As compared with an earlier version of this apparatus⁵, consisting of a single cryostat housing the sample setup, this extended version offers new possibilities for performing high-resolution measurements of relative length changes under well-controlled conditions. The system is distinguished by the following features: 1) It enables to significantly increase the maximum accessible pressure, P_{max} , being no longer limited by the external pressure reservoir. 2) The unwanted change of pressure during a temperature sweep has been drastically reduced. 3) Pressure sweeps can be performed both with decreasing and increasing pressure for a wide range of sweeping rates. 4) The setup allows for a simultaneous measurement of $\varepsilon_{\rm r}^{\rm He}(T,P)$ along the same T-P trajectory as used for taking the $\Delta L(T, P)/L$ data. This setup is well suited for detailed investigations of lattice effects and their coupling to other degrees of freedom in certain ranges of the T-P phase diagram.

ACKNOWLEDGMENTS

We acknowledge support by the Deutsche Forschungsgemeinschaft (DFG, German Research Foundation) through TRR 288 - 422213477 (projects A01 and A06) and K.D. Luther for providing us with the NaCl crystal.

AUTHOR DECLARATIONS

Conflict of Interest

The authors have no conflicts to disclose.

DATA AVAILABILITY STATEMENT

The data that support the findings of this study are available from the corresponding author upon reasonable request.

REFERENCES

- ¹T.H.K. Barron and J.G. Collins and G.K. White, "Thermal expansion of solids at low temperatures," Advances in Physics **29**, 609–730 (1980).
- ²I. Hamann, J. Sanjuan, R. Spannagel, M. Gohlke, G. Wanner, S. Schuster, F. Guzman, and C. Braxmaier, "Laser-dilatometer calibration using a single-crystal silicon sample," International Journal of Optomechatronics 13, 18–29 (2019).
- ³N. Kabeya, K. Imura, K. Deguchi, and N. K. Sato, "Thermal expansion measurements using the strain gauge technique with kelvin double bridge," Journal of the Physical Society of Japan 80, SA098 (2011).
- ⁴W. H. Fietz, K. Grube, and H. Leibrock, "Dilatometry under high pressure," High Pressure Research **19**, 373–378 (2000).
- ⁵R. S. Manna, B. Wolf, M. de Souza, and M. Lang, "High-resolution thermal expansion measurements under helium-gas pressure," Review of Scientific Instruments 83, 085111 (2012).
- ⁶T. Sakai, T. Kagayama, and G. Oomi, "Application of the strain gauge method to thermal expansion measurement under high pressure and high magnetic field," Journal of Materials Processing Technology 85, 224–228 (1999).
- ⁷E. Gati, J. M. Wilde, R. Khasanov, L. Xiang, S. Dissanayake, R. Gupta, M. Matsuda, F. Ye, B. Haberl, U. Kaluarachchi, R. J. McQueeney, A. Kreyssig, S. L. Bud'ko, and P. C. Canfield, "Formation of short-range magnetic order and avoided ferromagnetic quantum criticality in pressurized LaCrGe₃," Physical Review B **103**, 075111 (2021).
- ⁸O. L. Anderson, "The use of ultrasonic measurements under modest pressure to estimate compression at high pressure," Journal of Physics and Chemistry of Solids **27**, 547–565 (1966).
- ⁹D. L. Decker, W. A. Bassett, L. Merrill, H. T. Hall, and J. D. Barnett, "Highpressure calibration: A critical review," Journal of Physical and Chemical Reference Data 1, 773–836 (1972).
- ¹⁰E. Gati, M. Garst, R. S. Manna, U. Tutsch, B. Wolf, L. Bartosch, H. Schubert, T. Sasaki, J. A. Schlueter, and M. Lang, "Breakdown of Hooke's law of elasticity at the Mott critical endpoint in an organic conductor," Science Advances 2, e1601646 (2016).
- ¹¹R. Pott and R. Schefzyk, "Apparatus for measuring the thermal expansion of solids between 1.5 and 380k," Journal of Physics E: Scientific Instruments 16, 444–449 (1983).
- ¹²R. Küchler, T. Bauer, M. Brando, and F. Steglich, "A compact and miniaturized high resolution capacitance dilatometer for measuring thermal expansion and magnetostriction," Review of Scientific Instruments 83, 095102 (2012).
- ¹³R. L. Mills, D. H. Liebenberg, and J. C. Bronson, "Equation of state and melting properties of ⁴He from measurements to 20 kbar," Physical Review B **21**, 5137–5148 (1980).
- ¹⁴J. W. Schmidt and M. R. Moldover, "Dielectric permittivity of eight gases measured with cross capacitors," International Journal of Thermophysics 24, 375–403 (2003).
- ¹⁵J. Kamarád, Z. Machátová, and Z. Arnold, "High pressure cells for magnetic measurements—Destruction and functional tests," Review of Scientific Instruments **75**, 5022–5025 (2004).
- ¹⁶C. Kittel, *Introduction to Solid State Physics* (John Wiley and Sons Ltd, 2004).
- ¹⁷Q. Johnson, "High-pressure polymorphism in sodium chloride: A reinvestigation," Science 153, 419–420 (1966).
- ¹⁸W. A. Bassett, T. Takahashi, H.-K. Mao, and J. S. Weaver, "Pressureinduced phase transformation in NaCl," Journal of Applied Physics 39, 319–325 (1968).
- ¹⁹D. L. Decker, "Equation of state of NaCl and its use as a pressure gauge in high-pressure research," Journal of Applied Physics **36**, 157–161 (1965).
- ²⁰D. L. Decker, "High-pressure equation of state for NaCl, KCl, and CsCl," Journal of Applied Physics **42**, 3239–3244 (1971).

- ²¹J. M. Brown, "The NaCl pressure standard," Journal of Applied Physics 86, 5801–5808 (1999).
- ²²P. Dorogokupets, "Critical analysis of equations of state for NaCl," Geochemistry INTERNATIONAL C/C OF GEOKHIMIIA **40**, S132–S144 (2002).
- ²³T. Sumita and A. Yoneda, "Anharmonic effect on the equation of state (EoS) for NaCl," Physics and Chemistry of Minerals **41**, 91–103 (2013).
- ²⁴P. Bridgman, "The compression of twenty-one halogen compounds and eleven other simple substances to 100,000 kg/cm²," in *Proc. Am. Acad. Arts Sci*, Vol. 76 (1945).
- ²⁵G. R. Barsch and Z. P. Chang, "Adiabatic, isothermal, and intermediate pressure derivatives of the elastic constants for cubic symmetry. II. numerical results for 25 materials," physica status solidi (b) **19**, 139–151 (1967).
- ²⁶H. Spetzler, C. Sammis, and R. O'Connell, "Equation of state of NaCI: Ultrasonic measurements to 8 kbar and 800°c and static lattice theory," Journal

of Physics and Chemistry of Solids 33, 1727–1750 (1972).

- ²⁷S. Yamamoto, I. Ohno, and O. L. Anderson, "High temperature elasticity of sodium chloride," Journal of Physics and Chemistry of Solids 48, 143–151 (1987).
- ²⁸L. C. Chhabildas and A. L. Ruoff, "Isothermal equation of state for sodium chloride by the length-change-measurement technique," Journal of Applied Physics 47, 4182–4187 (1976).
- ²⁹R. Boehler and G. C. Kennedy, "Equation of state of sodium chloride up to 32 kbar and 500°c[†]," Journal of Physics and Chemistry of Solids 41, 517–523 (1980).
- ³⁰K. Wang and R. Reeber, "Thermal expansion of alkali halides at high pressure: NaCl as an example," Physics and Chemistry of Minerals 23, 354–360 (1996).
- ³¹P. Vinet, J. R. Smith, J. Ferrante, and J. H. Rose, "Temperature effects on the universal equation of state of solids," Physical Review B 35, 1945–1953 (1987).