## Practical Crystallography with a Transmission Electron Microscope

Benjamin L Weare,<sup>1,2</sup> Kayleigh L Y Fung,<sup>3</sup> Ian Cardillo-Zallo,<sup>1,4</sup> William J Cull,<sup>2</sup> Michael W Fay,<sup>1,5</sup> Stephen P Argent,<sup>2</sup> Paul D Brown<sup>5</sup>

1 - Nanoscale and Microscale Research Centre, University of Nottingham, Nottingham NG7 2RD, United Kingdom.

2 - School of Chemistry, University of Nottingham, Nottingham NG7 2RD, United Kingdom.

3 - Department of Computer Science, Nottingham Trent University, Nottingham NG11 8NS, United Kingdom.

4 – School of Pharmacy, University of Nottingham, Nottingham NG7 2RD, United Kingdom.

5 - Department of Mechanical, Materials, & Manufacturing Engineering, Faculty of Engineering, University of Nottingham, Nottingham NG7 2RD, United Kingdom.

## Abstract

Three-dimensional electron diffraction (3DED) is a powerful technique providing for crystal structure solutions of sub-micron sized crystals too small for structure determination via Xray techniques. The entry requirement, however, of a transmission electron microscope (TEM) adapted with bespoke software for coordinated sample stage rotation and continuous electron diffraction data acquisition has generally inhibited the wider uptake of 3DED. To address this limitation, we present novel software *GiveMeED* appropriate for controlled 3DED data acquisition. The collection of useable reflections beyond 0.8 Å makes 3DED crystallographic processing effectively routine, using standard software and workflows derived from single-crystal X-ray diffraction (SCXRD) techniques. A full experimental workflow for 3DED on a conventional TEM is described in practical terms, in combination with direct imaging, and energy dispersive X-ray spectroscopy (EDS) and electron energy loss spectroscopy (EELS), for the return of comprehensive correlative descriptions of crystal morphologies and sample compositions, with due regard for the quantification of electron flux at each stage of the characterisation process. The accuracy and effectiveness of *GiveMeED* is demonstrated through structure solutions for case study paracetamol, copper(II) phthalocyanine, and percholorocoronene samples, characterised in their nearnative states under controlled low dose conditions at either room or cryogenic temperatures, with determined unit cell parameters and atomic connectivity matching accepted literature X-ray structures for these compounds. To promote the wider adoption of 3DED, we make GiveMeED freely available for use and modification, in support of greater uptake and utilisation of structure solution procedures via electron diffraction.

#### 1. Introduction



**3D Electron Diffraction** 

Figure 1: Illustration of how three-dimensional electron diffraction (3DED) on a transmission electron microscope (TEM) may be used to solve crystal structures of small molecules. Diffraction patterns recorded whilst a specimen is rotated are used to solve the structure of the compound. Complementary microscopy data (images and spectra), used to inform the solution/refinement process, also provide a wealth of additional information about the sample. The recording of high quality 3DED data is rendered facile, with standard operation protocols developed from practice and fostered by open knowledge outlined in Supplementary Information (SI).

Three-dimensional electron diffraction (3DED) provides for the analysis of nano- and microcrystals that are too small for single-crystal X-ray diffraction (SCXRD) or powder X-ray diffraction (PXRD) techniques (Figure 1). 3DED has been used successfully for structure determination of small molecules,<sup>1-4</sup> zeolites,<sup>5-7</sup> proteins,<sup>8-15</sup> coordination polymers,<sup>16-18</sup> coordination complexes,<sup>19,20</sup> and other compounds.<sup>21-23</sup> In addition to routine structure determination, 3DED has also been shown to be capable of advanced crystallographic analysis such as refining hydrogen positions.<sup>24</sup>

3DED methodology has undergone several iterations.<sup>25–32</sup> Rotation electron diffraction (RED)<sup>33</sup> utilised stepped rotation of both the microscope stage and electron beam to collect diffraction tomography series. Precession-assisted electron diffraction tomography (PEDT)<sup>34</sup> coupled stepped stage rotation with beam precession to reduce multiple electron scattering events and increase crystallographic resolution, whilst the later innovation of continuous rotation electron diffraction (CRED) adopted continuous stage rotation to minimise the time needed for data collection whilst maximising the volume of reciprocal space interogated.<sup>35</sup> Data reduction in 3DED is currently performed by a small number of programs derived from SCXRD software: *CrysAlisPro,<sup>2</sup> PETS2,<sup>36</sup> Scipion-ED,<sup>37</sup> MicroED suite,<sup>38</sup> DIALS,<sup>39</sup>* and *XDS.*<sup>40</sup>

3DED has been performed using conventional TEM<sup>3-5,17,21,33,41-44</sup> operating at a range

of accelerating voltages, using both field-emission gun and thermionic sources, and typically without aberration correction. TEMs enabled with charged-couple device (CCD),<sup>5,45</sup> complementary metal-oxide-semiconductor (CMOS),<sup>5</sup> or direct-electron detecting<sup>3</sup> cameras have all been used for data collection, under conditions of low electron flux to help maintain samples in their near-native state. Accordingly, 3DED may be used to investigate delicate crystal structures before they become compromised by electron beam damage.<sup>46,47</sup> This approach in TEM offers the additional advantage, compared to commercial electron diffractometers,<sup>1,2</sup> of detailed morphological investigation via high-resolution imaging combined with elemental analysis from the same volume of material, delivered via energydispersive X-ray spectroscopy (EDS) or electron energy loss spectroscopy (EELS). It is recognised that a lack of generalised 3DED setups means implementation requires competence in both crystallography and electron microscopy,<sup>48</sup> rendering it "...a highly specialised technique for microscopy experts".<sup>1</sup> This reservation, however, may be overcome by the creation of more accessible methods for performing 3DED, and by establishing strong bridges between research communities, allowing crystal growers and crystallographers to make their data requirements clear and microscopists to design and carefully perform 3DED experiments.

The vital missing piece to turn a TEM into an electron diffractometer is software to synchronise stage tilt and ED data acquisition. Some 3DED software is commercially available<sup>49</sup> or open-source, *e.g. Instamatic*,<sup>50</sup> *InsteaDMatic*,<sup>5</sup> *Fast-ADT*, <sup>43</sup> and others,<sup>51</sup> but may not be applicable to local hardware/software setups. Hence, it's important that ED data collection software is easy to modify, so that it can be adapted and updated to future proof local instruments.

Here, we present a novel script *GiveMeED* for recording continuous rotation 3DED datasets, with general applicability illustrated via three case study small molecule crystal structure solutions. For microscopists, we describe how to adapt TEMs to collect high quality diffraction data using low electron fluence; for crystallographers, three crystal structures solved using data collected with *GiveMeED* are described, demonstrating the effectiveness of 3DED compared to X-ray diffraction approaches; and for the community of researchers growing crystals, we emphasise that 3DED on a conventional TEM can be used effectively to provide crystal structure solutions for sub-micrometre sized samples not suitable for X-ray investigation.

#### 2. Materials and Methods

#### 2.1. 3DED and TEM Data Acquisition Conditions



Flow 1: The workflow of a 3DED experiment. The steps executed by *GiveMeED* during 3DED data collection are shown in coloured boxes. Once the user presses "start", the processes in blue boxes are automated – the microscope stage is rotated whilst the camera records diffraction data, which is then saved alongside necessary metadata. Automation gave a streamlined, repeatable data collection procedure, helping to reduce operator error.

3DED was performed using *GiveMeED* on a JEOL 2100Plus at 200 kV with a Gatan OneView camera at 175 fps, with image size 1024x1024. The camera was operated in both image ("I") and video ("In-Situ") modes for 3DED. Diffraction patterns were captured at 600 mm camera length (pixel size =  $0.029 \text{ nm}^{-1}$ ) using selected area aperture 1 (100 µm diameter). The tiltaxis of the stage was calibrated prior to recording 3DED,<sup>44</sup> and could be recalibrated perdataset. The stage was rotated at eucentric height using the lowest speed setting, around 9.1° s<sup>-1</sup>. Specimens were located by low flux bright-field TEM imaging with defocussed objective lens (~100 nm).

The pre-sample optics and emission current were used to control electron flux at the sample and the size of the illuminated area: *i.e.* brightness (C3 lens strength), alpha (convergence angle via condenser mini lens strength), spot size (minimum beam radius via condenser system), and emission current ("beam target"). Searching for samples was performed with 50  $\mu$ m diameter condenser aperture, spot size 5, largest convergence angle, fully spread beam, and 104  $\mu$ A emission current. Once a suitable sample was identified, the beam blank was used to minimise total electron fluence. 3DED was then performed with 100  $\mu$ m diameter selected area aperture, spot size 1 or 2, and 108  $\mu$ A emission current. The corresponding electron flux measured at the sample during 3DED is indicated per-specimen in Section 3.1. To minimise beam damage the electron flux while searching for samples was maintained around 0.01 e<sup>-</sup> Å<sup>-2</sup> s<sup>-1</sup>. The electron beam stability of each sample was qualitatively judged by exposing a sacrificial region to an intense beam and observing how quickly the diffraction pattern faded.

During 3DED *GiveMeED* was used to control stage position, camera acquisition and metadata collection (Flow 1). Data was collected over the largest possible angle range at eucentric height, typically taking 12-14 seconds for 120° stage rotation. To judge whether samples diffracted beyond 0.8 Å the script *AutoResRings* was used to draw resolution rings on the camera live view. Full details of all scripts are included in SI. To achieve high-quality

crystal structures, data was recorded for samples that showed at least diffraction at 0.8 Å crystallographic resolution.

Complementary TEM images and associated spectra from the same sample were acquired following 3DED. Bright-field TEM imaging was performed with the OneView camera in both "I" and "Normal" modes, with 4K imaging performed at typically 1 s exposure and either 12 kX or 15 kX magnification (pixel size 0.024 nm or 0.77 nm). EEELS was performed using a Gatan Quefina EL spectrometer, with typically 5 s acquisition time. EDS was performed with an Oxford Instruments X-Mantle detector, typically with 1 minute live time. Samples were loaded at room temperature (RT) into the TEM column. JEOL single-tilt or HTR holders were used for 3DED investigations performed at room RT. A Gatan Elsa cryotomography holder was used for cryoTEM investigations, with samples cooled to 120 K *in situ*.

# 2.2. TEM Data Processing

*Digital Micrograph* (ver. 3.60.4435.0) was used for processing 3DED datasets into a format that could be imported into *PETS2* (ver. 2.2.20231114.1332), *CrysAlisPro* (ver. 171.43) and *DIALS* (ver. 3.23). The EDS spectra and images presented were prepared with *Digital Micrograph* (ver. 3.60). EELS spectra were background subtracted using a first order power law and plotted with *MATLAB* R2021a via the literature method.<sup>52</sup> Crystal structure models were prepared using *Vesta* (ver. 3.5.8).<sup>53</sup>

## 2.3. Data Reduction and Structure Solution/Refinement

Crystallographic processing of 3DED datasets was achieved with routine procedures and standard SCXRD software. Full details are given in SI. Data reduction was performed with *CrysAlisPro* and *PETS2*, then the output from *CrysAlisPro* was used for solution and refinement with *SHELXT* and *SHELXL* in *Olex2* (ver. 1.5)<sup>54–56</sup> Complete refinement information was retained in corresponding crystallographic information files (CIFs) (Table S1).

# 2.4. Electron Flux Measurements

The camera gain was calibrated by measuring beam diameter for a range of beam currents, then plotting intensity versus flux and taking the gradient (Figure S1). Beam intensity was measured with the Gatan OneView camera; beam currents were measured at the maximum current value with a Keithley 6485 Picoammeter connected to a Faraday cup integrated into a Gatan Model 646 Double Tilt Analytical Holder. Electron flux measurements were performed according to the method in Fung *et al.*<sup>57,58</sup> An image was taken with a 1 s exposure, the mean intensity per pixel was measured for a region of vacuum, then converted to electron flux using the calibrated camera gain (Table S2). From this the applied fluence for each stage of the 3DED analytical procedure described, and the cumulative fluence over the entire analysis, was calculated for each sample (Table S3). Beam exposure time was measured using how long a detector was running, so fluence applied during scouting for samples and aligning the microscope for different analytical techniques was calculated using an estimated beam exposure time. For EDS and EELS, an estimated electron flux based on typical illumination conditions was used.

# 2.6. Sample Preparation

Methods of sample preparation have been discussed in literature.<sup>48</sup> Perchlorocoronene and copper(II) phthalocyanine were supplied pre-prepared on lacy carbon on copper 200 mesh TEM grids. Paracetamol was recrystallised from hot water prior to use, then gently ground and drop-cast from n-hexane onto a lacy carbon on copper 200 mesh grid.

# 3. Results

# 3.1. Solving crystal structures using GiveMeED

The crystal structures of paracetamol, perchlorocoronene (PCC) and copper(II) phthalocyanine (CuPC) were determined from data collected with *GiveMeED*. Complementary bright-field TEM images and EDS/EELS spectra were also recorded for all samples (Figures 2 - 4).



Figure 2: Graphical overview of how the structure of (a) paracetamol was solved via 3DED.
(b) Low dose, cryogenic, bright-field TEM imaging located a suitable crystal (yellow circle indicates diffraction aperture position) then (c) 3DED data was recorded using *GiveMeED*.
The dashed red circle indicates reflections beyond 0.8 Å resolution could be captured during 3DED. (d) EDS and (e) EELS confirmed the specimen contained nitrogen, oxygen and carbon, before (f) the structure was solved and refined from electron diffraction data.

The structure of paracetamol was solved from a single crystal grain at 120 K (Figure 2). The specimen was irregular, measuring 1.1  $\mu$ m by 0.89  $\mu$ m in projection. The target

crystal was rotated 121° over 13.3 s with a mean electron flux of 1.90 e- Å<sup>-2</sup> s<sup>-1</sup> for a total electron fluence of 24.9 e- Å<sup>-2</sup>. 1152 unique reflections were collected for 75.8 % completeness to 0.8 Å. The lattice parameters were determined to be: a = 7.06 Å, b = 9.19 Å, c = 11.47 Å,  $\alpha = 90$ °,  $\beta = 96.2$ °,  $\gamma = 90$ °, with space group  $P2_1/c$  (#14). For the anisotropic structure, R<sub>1</sub> was 23.8 % and R<sub>int</sub> was 13.3 %.



Figure 3: Graphical overview of how the structure of (a) perchlorocoronene (PCC) was solved via 3DED. (b) Target crystal found in bright-field TEM mode at cryogenic temperatures (yellow circle indicates diffraction aperture), then 3DED data was acquired using *GiveMeED*. Dashed red circle indicates 0.8 Å resolution. (d) EDS and (e) EELS confirmed the presence of chlorine and sp<sup>2</sup> carbon in the sample prior to (f) structure solution and refinement from electron diffraction data.

The structure of PCC was solved from a single crystal grain at 120 K rotated 121° over 13.2 s (Figure 3). The specimen comprised the tip of a crystal needle 0.98 µm at its widest point and 1.0 µm long. The mean electron flux was 1.69 e- Å<sup>-2</sup> s<sup>-1</sup> and total electron fluence was 22.3 e- Å<sup>-2</sup>. 1213 unique reflections were collected for 99.0 % completeness to 0.8 Å. The lattice parameters determined were: *a* = 22.37 Å, *b* = 8.16 Å, *c* = 12.81 Å,  $\alpha = \beta = \gamma = 90$ °, with space group *Cmce* (#64). R<sub>1</sub> for the refinement 33.3 % and R<sub>int</sub> was 13.1 % for the anisotropic structure.



Figure 4: Structure solution for (a) copper(II) phthalocyanine (CuPC) from 3DED data. (b) A crystal was located at room temperature via bright-field TEM (yellow circle indicates diffraction aperture); and 3DED data was collected using *GiveMeED*. Dashed red circle indicates 0.8 Å resolution. (d) EDS and (e) EELS confirmed the presence of copper, nitrogen, and sp<sup>2</sup> carbon prior to structure solution. (f) The structure was solved and refined from electron diffraction data.

CuPC was solved from three crystal grains at room temperature (Figure 4). The crystals were approximately 1.5  $\mu$ m in width and several micrometres long, from which a small region was defined by the SAED aperture. 2137 unique reflections were collected for 94.5 % completeness to 0.8 Å resolution. Over all three grains, the mean flux was 0.62 e- Å<sup>-2</sup> s<sup>-1</sup> and the total electron fluence was 7.68 e- Å<sup>-2</sup>. The lattice parameters determined were: *a* = 14.29 Å, *b* = 4.75 Å, *c* = 17.0 Å,  $\alpha$  = 90 °,  $\beta$  = 104.8 °,  $\gamma$  = 90 °, with space group *P2*<sub>1</sub>/*c* (#14). For the anisotropic structure, R<sub>1</sub> was 33.4 % and R<sub>int</sub> was 12.0 %.

#### 4. Discussion

#### 4.1 Overview

A novel script *GiveMeED* has been developed to collect 3DED datasets from sub-micron sized crystalline specimens within a conventional TEM, under low dose conditions at room or cryogenic temperatures, as appropriate, for the purpose of structure solution and refinement of samples maintained near their native state. In particular, 3DED combined with high spatial resolution TEM imaging and chemical analysis provided for the comprehensive

investigation of individual specimens that would otherwise be intractable to identify and characterise using light microscopy and X-ray diffraction techniques. Three case study crystal structures, in the form of paracetamol, perchlorocoronene (PCC) and copper(II) phthalocyanine (CuPC) have been determined here, by way of illustration of the general applicability of 3DED.

In each case, initial low dose, bright-field TEM imaging during the scoping stage provided for quick assessment of sample crystallinity, via the presence of diffraction contrast observed with the objective lens defocused. Thereafter, the sequential workflow of 3DED data acquisition, mediated by *GiveMeED*, to provide datasets for the determination of sample structure, followed by correlative imaging and spectroscopy to confirm each specimen contained the elements expected, ensured time examining spurious contaminants was avoided, in advance of the more labour-intensive stage of crystallographic data analysis.

Data reduction, *i.e.*, the conversion of experimental diffraction patterns into a matrix of indexed reflections with measured intensity, required the development of protocols to convert 3DED datasets into a format that could be imported into standard X-ray software. Standard data reduction and structure solution/refinement procedures could then be used to solve the crystal structures from 3DED data; *i.e.*, the creation of atomic models describing experimental electron potential in the unit cell. Comparison with literature X-ray data<sup>59–61</sup> confirmed the structures returned by 3DED to be reliable, demonstrating the suitability of this workflow to solve molecular crystal structures, in particular crystals that cannot be grown with sufficient size suitable for characterisation via X-ray techniques.

It is recognised that both microscopy and crystallography are rich in technical terminology, hence plain language clarifications have added to the text and a list of abbreviations has been included in Supporting Information (SI Table S5).

#### 4.2. GiveMeED

*GiveMeED* is a 3DED data collection script for continuous-rotation electron diffraction measurements. For routine operation, a user interface has been provided. It was designed to be flexible via user modification, and adaptable to future developments in electron crystallography methodology. Hence, the code for *GiveMeED* is freely available and full operational details are included in SI.

*GiveMeED* has been designed to be understood easily and modified. Each process in the script has been compartmentalised into a function, letting users more easily adjust individual aspects of the data collection software, and to allow quick transplantation of functionality to new software by copying entire functions. These design choices are intended to make development of new software built on *GiveMeED* relatively straightforward.

*GiveMeED* records 3DED metadata at the point of data capture. This contains information including electron wavelength and the rotational relationship between each diffraction pattern; *i.e.*, all of the information required to perform crystallographic analysis. Pathways have been developed to export data from *GiveMeED* into common data reduction software (*DIALS*, *PETS2* and *CrysAlisPro*) facilitated by scripts written in DM-script and Python. Flexibility is key, hence data recorded with *GiveMeED* can converted into arbitrary image file formats, suitable for the individual needs of users and their preferred data reduction software. Such design choices will allow *GiveMeED* to be modified as standard practices for electron diffractometry develop, empowering users to collect the exact data considered necessary for their scientific endeavours.

## 4.3. Data collection workflow

Accurate crystal structures required collection of high angle reflections from crystal planes with small d-spacing, giving sub-Ångstrom precision of atomic positions in solved crystal structures. We aimed to collect reflections of at least 0.8 Å crystallographic resolution, judged using resolution rings on the camera live view, and routinely collected reflections corresponding to smaller d-spacing. This allowed confidence that the collected data would be useable for solving crystal structures.

The microscope was operated under conditions of very low electron flux to limit beam damage whilst searching for samples, which meant objective lens defocus was required to give enough contrast via Fresnel effects to see the sample on the camera live view. Once a suitable specimen was located, recording data was as simple as pressing "Start 3DED" on the *GiveMeED* user interface (Flow 1). Thereafter, qualitative TEM spectroscopy confirmed the presence of elements expected within the sample (Figure S2, Table S5), sufficient to support atom assignments performed during structure solution and refinement, without need for quantitative determination of elemental concentrations. In particular, if a crystalline contaminant was located in the sample, it could be disregarded quickly, before performing any labour-intensive crystallographic data analysis. Overall, data collection was routine and the most time-consuming step was locating suitable crystals in each sample.

## 4.4. Electron fluence and beam damage



Figure 5: Electron fluence applied at each step of analysis (bar chart) and the cumulative electron fluence at each step (line chart) for paracetamol (purple), PCC (blue), and CuPC (green). Inset shows the first four steps in greater detail. The majority of electron fluence was applied during spectroscopy, while 3DED and imaging involved a much smaller amount of fluence. Between each technique, "Setup" indicates the estimated fluence applied while the microscope was adjusted for a new analysis mode. Fields marked with asterisk (\*) were calculated estimated electron flux, and rows with dagger (†) used estimated time, to

calculate electron fluence.

Beam damage causes reflections to be lost from the diffraction pattern as the crystal lattice becomes amorphous – starting with high-angle reflections – and this is typically proportional to the applied electron fluence. Hence, prior to and during 3DED measurements we aimed to use minimum electron flux to preserve the crystal lattice at near native-state.

By measuring the electron flux and the beam exposure time the applied electron fluence during 3DED, imaging, and spectroscopy could be calculated; the electron fluence applied while changing the microscope between analytical modes was then estimated to allow calculation of the cumulative electron fluence over the entire analysis (Figure 5, Table S3). This allowed tracking of which techniques applied the most fluence to the specimen; it clearly showed the importance of performing 3DED before BF-TEM imaging or spectroscopy as each technique increased the cumulative fluence, taking the specimen further from its native state by damaging the crystal structure. Spectroscopy applied the most fluence, causing the most damage to the sample. This was borne out in the diffraction pattern; after 3DED and TEM imaging, we observed reflections above 0.8 Å that indicated the structure was still useful for crystallographic analysis, while after EDS and EELS no reflections were observed in the diffraction pattern. A high electron flux was used for EDS and EELS as it allowed rapid collection of spectra. While it would be possible to collect spectra with a lower flux to preserve the crystal structure this would require a proportional increase in time to collect spectra with acceptable signal-to-noise, so sample throughput was prioritised. Degradation of the crystal lattice during spectroscopy was not expected to affect the qualitative elemental information provided by EDS and EELS. Overall, fluence control allowed confidence in the quality of the crystallographic data acquired by limiting the extent of sample damage.

## 4.5. Electron crystal structure determination

The structures of three small molecules were solved from sub-micron sized specimens, which it would not have been possible to examine using SCXRD techniques. As the data collected contained useable reflections above 0.8 Å, crystallographic processing was routine, using standard software and workflows derived from SCXRD. Despite the known practical challenges of electron crystallography,<sup>26</sup> we determined unit cell parameters and atomic connectivity that matched the literature X-ray structures for the compounds in this study<sup>59–61</sup> providing confidence in the utility of data recorded with *GiveMeED* to generate useful crystal structures in future studies.

## 5. Conclusions

3DED is primed to become a major tool for crystallographers and microscopists, for the determination of accurate crystal structures for compounds that cannot be investigated using X-ray techniques. The major roadblock to 3DED uptake is availability of instrumentation – to this end *GiveMeED* software was developed to allow adaptation of conventional TEMs to perform electron diffractometry. A workflow for collecting 3DED data sets to 0.8 Å resolution at low electron flux has been described so that electron microscopists unfamiliar with crystallography can collect high-quality data using *GiveMeED*. For crystallographers we presented the structures of paracetamol, perchlorocoronene, and copper(II) phthalocyanine solved from electron diffraction data, demonstrating that crystal structures matching to literature X-ray structures can be produced from data recorded with *GiveMeED* on a TEM. A full workflow from sample preparation to data analysis was also described, and showing that elemental information from EDS and EELS can be acquired from the same specimen volume as 3DED. We anticipate that 3DED will continue to develop as more facilities acquire capacity to collect electron diffraction data for structure solution, with *GiveMeED* protocols serving as a foundation for 3DED data collection.

# Acknowledgements

*GiveMeED* is freely available online at https://github.com/benweare/GiveMeED. The structures of the compounds in this work are available from The Cambridge Crystallographic Data Centre (www.ccdc.cam.ac.uk/structures) with reference numbers 2464563, 24564, and 24565. Diffraction patterns used to solve crystal structures are available at https://doi.org/10.17639/nott.7567. The authors would like to thank the Nanoscale and Microscale Research Centre and staff for technical support and useful conversations; the

School of Chemistry at the University of Nottingham and staff for technical support; and Dr J. P. Tidey for useful conversations.

## **Author Contributions**

**B. L. Weare:** Conceptualization, Methodology, Software, Validation, Formal analysis, Investigation, Data Curation, Writing – Original Draft, Visualisation, Project administration. **K. L. Y. Fung:** Visualisation, Writing – Review and Editing. **I. Cardillo-Zallo:** Investigation, Data Curation, Resources, Writing – Review and Editing. **W. J. Cull:** Visualisation, Writing – Review and Editing. **M. W. Fay:** Methodology, Software, Validation, Project administration. **S. P. Argent:** Supervision, Formal analysis, Validation, Writing – Review and Editing. **P. D. Brown:** Conceptualization, Methodology, Supervision, Project administration, Funding acquisition, Writing – Review and Editing.

## **Funding Information**

Funding for this research was provided by: Engineering and Physical Sciences Research Council (grant No. EP/W006413/1 and EP/L022494/1); Leverhulme Trust (RPG-2022-300).

#### References

- 1 P. Simoncic, E. Romeijn, E. Hovestreydt, G. Steinfeld, G. Santiso-Quiñones and J. Merkelbach, *Acta Crystallogr. Sect. E Crystallogr. Commun.*, 2023, **79**, 410–422.
- 2 S. Ito, F. J. White, E. Okunishi, Y. Aoyama, A. Yamano, H. Sato, J. D. Ferrara, M. Jasnowski and M. Meyer, *CrystEngComm*, 2021, **23**, 8622–8630.
- 3 K. Takaba, S. Maki-Yonekura, I. Inoue, K. Tono, T. Hamaguchi, K. Kawakami, H. Naitow, T. Ishikawa, M. Yabashi and K. Yonekura, *Nat. Chem.*, 2023, **15**, 491–497.
- 4T. E. Gorelik, A. Ulmer, A. L. Schleper and A. J. C. Kuehne, *Z. Für Krist. Cryst. Mater.*, 2023, **238**, 253–260.
- 5 M. Roslova, S. Smeets, B. Wang, T. Thersleff, H. Xu and X. Zou, *J. Appl. Crystallogr.*, 2020, **53**, 1217–1224.
- 6 A. Laštovičková, D. N. Rainer and M. Mazur, *Microporous Mesoporous Mater.*, 2025, **387**, 113514.
- 7 J. I. Tirado, A. Sala, A. Bordes, P. Pratim Das, L. Palatinus, S. Nicolopoulos, J. L. Jordá, A. Vidal-Moya, T. Blasco, G. Sastre, S. Valencia and F. Rey, *Angew. Chem.*, 2025, **137**, e202416515.
- 8 J. A. Rodriguez, M. I. Ivanova, M. R. Sawaya, D. Cascio, F. E. Reyes, D. Shi, S. Sangwan, E. L. Guenther, L. M. Johnson, M. Zhang, L. Jiang, M. A. Arbing, B. L. Nannenga, J. Hattne, J. Whitelegge, A. S. Brewster, M. Messerschmidt, S. Boutet, N. K. Sauter, T. Gonen and D. S. Eisenberg, *Nature*, 2015, 525, 486–490.
- 9 H. M. E. Duyvesteyn, A. Kotecha, H. M. Ginn, C. W. Hecksel, E. V. Beale, F. De Haas, G. Evans, P. Zhang, W. Chiu and D. I. Stuart, *Proc. Natl. Acad. Sci.*, 2018, **115**, 9569–9573.
- 10 B. L. Nannenga, D. Shi, J. Hattne, F. E. Reyes and T. Gonen, *eLife*, 2014, **3**, e03600.
- 11 M. J. De La Cruz, J. Hattne, D. Shi, P. Seidler, J. Rodriguez, F. E. Reyes, M. R. Sawaya, D. Cascio, S. C. Weiss, S. K. Kim, C. S. Hinck, A. P. Hinck, G. Calero, D. Eisenberg and T. Gonen, *Nat. Methods*, 2017, **14**, 399–402.
- 12 D. Shi, B. L. Nannenga, M. J. De La Cruz, J. Liu, S. Sawtelle, G. Calero, F. E. Reyes, J. Hattne and T. Gonen, *Nat. Protoc.*, 2016, **11**, 895–904.
- 13 M. Gallagher-Jones, C. Glynn, D. R. Boyer, M. W. Martynowycz, E. Hernandez, J. Miao, C.-T. Zee, I. V. Novikova, L. Goldschmidt, H. T. McFarlane, G. F. Helguera, J. E. Evans, M. R. Sawaya, D. Cascio, D. S. Eisenberg, T. Gonen and J. A. Rodriguez, *Nat. Struct. Mol. Biol.*, 2018, **25**, 131–134.
- 14 M. W. Martynowycz, M. T. B. Clabbers, J. Hattne and T. Gonen, *Nat. Methods*, 2022, **19**, 724–729.
- 15 S. Liu and T. Gonen, *Commun. Biol.*, 2018, **1**, 38.
- A. K. Ramana, J. P. Tidey, G. M. De Lima and R. I. Walton, *Small Methods*, 2024, 2301703.
- 17 M. Quintelier, A. Hajizadeh, A. Zintler, B. F. Gonçalves, R. Fernández De Luis, L. Esrafili Dizaji,
- C. M. L. Vande Velde, S. Wuttke and J. Hadermann, *Chem. Mater.*, 2024, acs.chemmater.4c01153.
  H. Bunzen, B. Sertcan Gökmen, A. Kalytta-Mewes, M. Grzywa, J. Wojciechowski, J. Hutter, A.
- Hehn and D. Volkmer, *Chem. Eur. J.*, 2025, **31**, e202403574.
- 19 C. G. Jones, M. Asay, L. J. Kim, J. F. Kleinsasser, A. Saha, T. J. Fulton, K. R. Berkley, D. Cascio, A. G. Malyutin, M. P. Conley, B. M. Stoltz, V. Lavallo, J. A. Rodríguez and H. M. Nelson, *ACS Cent. Sci.*, 2019, **5**, 1507–1513.
- 20 V. Shvalagin, Y. Markushyna, D. Piankova, E. S. Grape, N. V. Tarakina and O. Savateev, *ChemPhotoChem*, 2025, **9**, e202400255.
- 21 S. Syniakina, B. Grushko and L. Meshi, J. Alloys Compd., 2023, 963, 171260.
- 22 P. B. Klar, Y. Krysiak, H. Xu, G. Steciuk, J. Cho, X. Zou and L. Palatinus, *Nat. Chem.*, 2023, **15**, 848–855.
- 23 S. Vergara, D. A. Lukes, M. W. Martynowycz, U. Santiago, G. Plascencia-Villa, S. C. Weiss, M. J. De La Cruz, D. M. Black, M. M. Alvarez, X. López-Lozano, C. O. Barnes, G. Lin, H.-C. Weissker, R. L. Whetten, T. Gonen, M. J. Yacaman and G. Calero, *J. Phys. Chem. Lett.*, 2017, **8**, 5523–5530.
- L. Palatinus, P. Brázda, P. Boullay, O. Perez, M. Klementová, S. Petit, V. Eigner, M. Zaarour and S. Mintova, *Science*, 2017, **355**, 166–169.

- 25 M. Gemmi, E. Mugnaioli, T. E. Gorelik, U. Kolb, L. Palatinus, P. Boullay, S. Hovmöller and J. P. Abrahams, *ACS Cent. Sci.*, 2019, **5**, 1315–1329.
- 26 M. Gemmi and A. E. Lanza, *Acta Crystallogr. Sect. B Struct. Sci. Cryst. Eng. Mater.*, 2019, **75**, 495–504.
- 27 B. L. Nannenga, *Struct. Dyn.*, 2020, **7**, 014304.
- 28 B. L. Nannenga and T. Gonen, *Curr. Opin. Struct. Biol.*, 2016, **40**, 128–135.
- 29 M. J. Gallenito and T. Gonen, *Biochem. Soc. Trans.*, 2022, **50**, 231–239.
- 30 D. Shi, B. L. Nannenga, M. G. Iadanza and T. Gonen, *eLife*, 2013, **2**, e01345.
- 31 B. L. Nannenga, G. Bu and D. Shi, *Front. Mol. Biosci.*, 2018, **5**, 114.
- 32 M. T. B. Clabbers, A. Shiriaeva and T. Gonen, *IUCrJ*, 2022, **9**, 169–179.
- 33 D. Zhang, P. Oleynikov, S. Hovmöller and X. Zou, *Z. Für Krist.*, 2010, **225**, 94–102.
- U. Kolb, E. Mugnaioli and T. E. Gorelik, Cryst. Res. Technol., 2011, 46, 542–554.
- 35 M. Gemmi, M. G. I. La Placa, A. S. Galanis, E. F. Rauch and S. Nicolopoulos, *J. Appl. Crystallogr.*, 2015, **48**, 718–727.
- 36 L. Palatinus, P. Brázda, M. Jelínek, J. Hrdá, G. Steciuk and M. Klementová, *Acta Crystallogr. Sect. B Struct. Sci. Cryst. Eng. Mater.*, 2019, **75**, 512–522.
- V. E. G. Bengtsson, L. Pacoste, J. M. De La Rosa-Trevin, G. Hofer, X. Zou and H. Xu, *J. Appl. Crystallogr.*, 2022, **55**, 638–646.
- 38 M. G. Iadanza and T. Gonen, J. Appl. Crystallogr., 2014, 47, 1140–1145.
- 39 M. T. B. Clabbers, T. Gruene, J. M. Parkhurst, J. P. Abrahams and D. G. Waterman, *Acta Crystallogr. Sect. Struct. Biol.*, 2018, **74**, 506–518.
- 40 W. Kabsch, Acta Crystallogr. D Biol. Crystallogr., 2010, 66, 125–132.
- 41 B. Wang, X. Zou and S. Smeets, *IUCrJ*, 2019, **6**, 854–867.
- 42 W. Wan, J. Sun, J. Su, S. Hovmöller and X. Zou, *J. Appl. Crystallogr.*, 2013, **46**, 1863–1873.
- 43 S. Plana-Ruiz, Y. Krysiak, J. Portillo, E. Alig, S. Estradé, F. Peiró and U. Kolb, *Ultramicroscopy*, 2020, **211**, 112951.
- J. Heidler, R. Pantelic, J. T. C. Wennmacher, C. Zaubitzer, A. Fecteau-Lefebvre, K. N. Goldie, E. Müller, J. J. Holstein, E. Van Genderen, S. De Carlo and T. Gruene, *Acta Crystallogr. Sect. Struct. Biol.*, 2019, **75**, 458–466.
- 45 A. R. Faruqi and G. McMullan, *Nucl. Instrum. Methods Phys. Res. Sect. Accel. Spectrometers Detect. Assoc. Equip.*, 2018, **878**, 180–190.
- 46 R. Henderson, *Q. Rev. Biophys.*, 1995, **28**, 171–193.
- 47 N. Vlahakis, A. Clauss and J. A. Rodriguez, *IUCrJ*, DOI:10.1107/S2052252524012132.
- 48 T. Gruene and E. Mugnaioli, *Chem. Rev.*, 2021, **121**, 11823–11834.
- 49 J. Roqué Rosell, J. Portillo Serra, T. H. Aiglsperger, S. Plana Ruiz, P. P. Das, J. Mendoza Gonzalvez, T. Trifonov and J. A. Proenza, *Microchem. J.*, 2019, **148**, 130–140.
- 50 S. Smeets, B. Wang, M. O. Cichocka, J. Ångström and W. Wan, Instamatic (version 1.0.0) Zenodo 2018.
- 51 M. J. De La Cruz, M. W. Martynowycz, J. Hattne and T. Gonen, *Ultramicroscopy*, 2019, **201**, 77–80.
- 52 K. L. Y. Fung, M. W. Fay, S. M. Collins, D. M. Kepaptsoglou, S. T. Skowron, Q. M. Ramasse and A. N. Khlobystov, *Ultramicroscopy*, 2020, **217**, 113052.
- 53 K. Momma and F. Izumi, J Appl Crystallogr 44 1272-1276 2011, 2011, 44, 1271–1276.
- 54 G. M. Sheldrick, Acta Crystallogr. Sect. Found. Adv., 2015, **71**, 3–8.
- 55 G. M. Sheldrick, Acta Crystallogr. Sect. C Struct. Chem., 2015, **71**, 3–8.
- 56 O. V. Dolomanov, L. J. Bourhis, R. J. Gildea, J. A. K. Howard and H. Puschmann, J. Appl. Crystallogr., 2009, **42**, 339–341.
- 57 K. L. Y. Fung, B. L. Weare, M. W. Fay, S. P. Argent and A. N. Khlobystov, *Micron*, 2022, 103395.
- 58 K. L. Y. Fung, S. Skowron, R. Hayter, S. E. Mason, B. L. Weare, Q. Ramasse, C. Allen and A.
- Khlobystov, Phys. Chem. Chem. Phys., 2023, 10.1039.D2CP05183D.
- 59 D. Yu. Naumov, M. A. Vasilchenko and J. A. K. Howard, 1998.

- 60 H. Jiang, P. Hu, J. Ye, R. Ganguly, Y. Li, Y. Long, D. Fichou, W. Hu and C. Kloc, 2018.
- 61 T. Baird, J. H. Gall, D. D. MacNicol, P. R. Mallinson and C. R. Michie, J. Chem. Soc. Chem. Commun., 1988, 1471.

## Supporting Information for Practical Crystallography with a Transmission Electron Microscope

Benjamin L Weare,<sup>1,2</sup> Kayleigh L Y Fung,<sup>3</sup> Ian Cardillo-Zallo,<sup>1,4</sup> William J Cull,<sup>2</sup> Michael W Fay,<sup>1,5</sup> Stephen P Argent,<sup>2</sup> Paul D Brown<sup>5</sup>

1 - Nanoscale and Microscale Research Centre, University of Nottingham, Nottingham NG7 2RD, United Kingdom.

2 - School of Chemistry, University of Nottingham, Nottingham NG7 2RD, United Kingdom.

3 - Department of Computer Science, Nottingham Trent University, Nottingham NG11 8NS, United Kingdom.

4 – School of Pharmacy, University of Nottingham, Nottingham NG7 2RD, United Kingdom.

5 - Department of Mechanical, Materials, & Manufacturing Engineering, Faculty of Engineering, University of Nottingham, Nottingham NG7 2RD, United Kingdom.

## Contents

1. Overview	2
2. Crystallographic Data Processing	2
3. Measuring Electron Flux and Calibrating Camera Gain	4
4. Calculation of Electron Fluence Errors	7
5. Exporting 3DED data from Digital Micrograph	7
6. Energy Dispersive X-ray and Electron Energy Loss Spectroscopy	8
7. Instructions for <i>GiveMeED</i> and other scripts	9
8. List of Abbreviations	11
8. References	11

## 1. Overview

This Supporting Information contains in-depth instructions for using the script *GiveMeED* to collect 3DED data, and details of the scripts used to support 3DED data collection *AutoResRings*, and *Go2Alpha*. Details are given of how to import data recorded with *GiveMeED* into 3DED data processing software. A complete description of how to measure electron flux using a CMOS camera via determining it's gain are provided, applicable to any CMOS or CCD camera. Supporting EDS and EELS spectra are presented.

	Paracetamol	PCC	CuPC
Formula	C <sub>8</sub> H <sub>9</sub> NO <sub>2</sub>	C <sub>24</sub> Cl <sub>12</sub>	$C_{32}H_{16}CuN_8$
Merged grains	1	1	3
Tilt range / °	120.9	120.9	121.0; 120.5; 96.7
Angle per frame / °	0.052	0.052	0.055; 0.055; 0.057
Frames	2333	2328	2180; 2188; 1685
Mean flux / e- Å <sup>-2</sup> s <sup>-1</sup>	1.88	1.27	0.47; 0.62; 0.91
Total time / s	13.3	13.2	12.5; 12.5; 9.68
Temp / K	120	120	Room temp.
Resolution cutoff / Å	0.8	0.8	0.8
Reflections	2784	3091	7774
Unique Reflections	1147	1213	2122
Unit cell	7.06(1), 9.19 (13),	22.37(16), 8.16(14),	14.29(11), 4.75(8), 17.0(1),
	11.47(9) <i>,</i> 96.2(2)	12.81(15)	104.8(3)
Space group	P21/n (#14)	Cmce (#64)	P21/n (#14)
Completeness / %	75.8	99.0	93.4
I/sig(I)	8.2	7.7	9.7
Rint / %	13.3	13.1	12.0
R1/%	23.8	33.2	33.4
wR2 / %	58.7	73.8	66.3
GooF	2.26	3.20	2.51

#### 2. Crystallographic Data Processing

Table S1: Data collection and crystallographic details for the case study structures of paracetamol, percholorocoronene (PCC), and copper(II) phthalocyanine (CuPC) presented in this work.

Data reduction of 3DED datasets was performed with *CrysAlisPro*. Peak hunting used the following options: "3D peak extraction" with "Single frame background subtraction" and "Weak features extraction". A resolution limit of 0.6 Å was applied for integration. Data finalisation was performed under standard (kinematic) settings, with a resolution limit of 0.8 Å. Empirical absorption correction was applied using *SCALE3 ABSPACK* (1.0.7).

For all structures, kinematic solution and refinement was performed with *SHELXT* and *SHELXL* in *Olex2* (ver. 1.5) (Table S1).<sup>1–3</sup> Applied *SHELXL* instructions are indicated in parentheses. Hydrogen atoms were observed in the electron potential map then geometrically placed using neutron bond distances and refined with a riding model. <sup>4</sup> Rigid bond (*SHELXL* RIGU) and similarity (*SHELXL* SIMU) restraints were applied to the anisotropic displacement parameters of all atoms in the structures. The values of refinement statistics

(R1, wR2, goodness-of-fit, extinction coefficient, and R<sub>int</sub>) were high due to the structure being solved and refined from electron diffraction data using a kinematic model that did not account for the dynamical diffraction of electrons. This lead to greater uncertainty in atomic positions and unit cell parameters, *i.e.* structures confirmed atomic connectivity but were not suitable for extracting detailed geometric information due to in large uncertainties on bond lengths and unit cell parameters. Experimental limitations of electron diffraction data collection method (*i.e.* limited goniometer rotation range) meant the entire diffraction sphere was not sampled. Hence, a number of reflections were not present in the data used for structure solution and refinement. Some reflections had observed intensity (I<sub>obs</sub>) much smaller than calculated intensity (I<sub>calc</sub>) due to measurement discrepancies arising from the effect of dynamical diffraction or reflections impinging on the beam stop. *Platon ADDSYM* was used to test all structures for missing symmetry.<sup>5</sup>

For paracetamol, the structure was solved in space group P2<sub>1</sub>/m and reset in P2<sub>1</sub>/n using *Platon ADDSYM*. The anisotropic displacement parameter of amide carbon atom C8 was restrained to have more isotropic character (*SHELXL* ISOR). Geometric similarity restraints were applied symmetry related bond distances with the phenyl ring moiety (*SHELXL* SAME), where 1,2-bond distances used estimated standard deviation (e.s.d) of 0.02 Å and 1,3-bond distances used e.s.d of 0.04 Å. For percholorcoronene, the anisotropic character (*SHELXL* ISOR). For copper(II) phthalocyanine, the structure was solved in space group Pm and reset in P2<sub>1</sub>/n using *Platon ADDSYM*. Data from three crystal grains was merged into a single dataset to increase data completeness using Proffit merge (*CrysAlisPro*). Same distance restraints (*SHELXL* SADI) were applied to chemically identical environments in the two benzene moieties of the asymmetric unit, where 1,2-bond distances used e.s.d of 0.04 Å.

#### 3. Measuring Electron Flux and Calibrating Camera Gain



Figure S1: Electron micrographs of vacuum illustrating how electron beam diameter was measured using (a) the whole beam and (b) an arc of the beam via the intersection chords theorem (Equation S2). Scale bars are 500 nm. (c) Linear fitting of flux versus camera intensity used to measure camera gain. The gradient of the black fit line is camera gain (*G*) used to convert from intensity to flux:  $G = 36.3 \pm 0.8$  counts electron<sup>-1</sup> (3 $\sigma$  standard deviation). Dashed grey lines indicate 99% confidence intervals. Red crosses are data points, and the blue square is an outlier (R<sup>2</sup> = 0.9987).

For accurate electron flux measurement on a CCD or CMOS cameras it was necessary to calibrate the conversion of camera counts to electrons (*i.e.* camera gain). Flux was measured using the method described by Egerton,<sup>6</sup> where flux *j* is the ratio of the beam current  $I_b$  and beam area  $A_b$ 

$$j = \frac{I_b}{A_b} = \frac{I_b}{e\pi r^2} \tag{1}$$

Where *r* was the radius of the beam in pixels and *e* was elementary charge, giving *j* in electrons pixel<sup>-1</sup> s<sup>-1</sup>. Electron beam radius was measured using the camera. When the entire beam fit on the camera sensor the beam radius was measured directly, taking the edge of the electron beam as the full width tenth maximum. To measure the radius of beams larger than the camera sensor, the intersecting chords theorem was applied

$$r = \frac{L^2}{8H} + \frac{H}{2} \tag{2}$$

Where *L* was chord length and *H* was chord height (Figure S1a,b). This allowed measurement of large radius electron beams.

For each electron flux, the electron beam intensity at the camera (Q) was measured by taking the average intensity over a region of vacuum. A plot of Q versus j yielded the camera gain (G, in counts electron<sup>-1</sup>) as the gradient (Figure S1c)

$$G = \frac{Q}{j} = \left(\frac{Q\pi r^2}{I_b}\right) \tag{2}$$

At 200 kV,  $G = 36.3 \pm 0.8$  counts electron<sup>-1</sup>. To measure the electron flux in experimental data the electron beam intensity of vacuum was measured at a region of vacuum and applying the following relation

$$j_E = \frac{Q_E}{d^2 G} \tag{4}$$

Where  $j_E$  was the experimental electron flux,  $Q_E$  was the experimental intensity, and d was pixel size in nanometres per pixel. A look up table of  $d^2G$  was then produced for a range of magnifications (Table S2). Experimental electron flux was used to calculate applied electron fluence ( $F_E$ ) (Table S3)

$$F_E = j_e t \tag{5}$$

Where t was the length of time the specimen was exposed to flux  $j_e$ . The cumulative fluence ( $F_c$ ) was then the sum of each electron fluence

$$F_C = \sum F_E \tag{6}$$

Magnification /	<i>d /</i> nm pixel <sup>-1</sup>	<i>d</i> ² <i>G</i> / e nm²
kX		count <sup>-1</sup> pixel <sup>-2</sup>
100	0.11	3.98
80	0.14	4.99
60	0.18	6.67
50	0.22	8.03
40	0.28	10.1
30	0.38	13.7
25	0.47	16.9
20	0.61	22.3
15	0.78	28.4
12	0.94	34.0
10	1.2	40.3

Table S2: Pixel length (*d*) and intensity-to-flux conversion factors ( $d^2G$ ) for a range of magnifications at 200 kV. The conversion factors are the product of the gain (G = 36.3 counts electron<sup>-1</sup>) and the square of the pixel length.

	Fluence / e- Å <sup>-2</sup>		Cumulative Fluence / e- Å <sup>-2</sup>			
	Paracetamol	PCC	CuPC	Paracetamol	PCC	CuPC
<sup>†</sup> Scouting	$0.22 \pm 1.1$	$0.22 \pm 1.1$	0.22 <u>+</u> 1.1	$0.22 \pm 1.1$	$0.22 \pm 1.1$	$0.22 \pm 1.1$
3DED	25 ± 0.076	$22\pm0.076$	7.6 <u>±</u> 0.088	25 ± 1.1	22 <u>+</u> 1.1	7.9 ± 1.1
<sup>+</sup> Image setup	9.4 ± 0.2	$6.4 \pm 0.2$	$3.4 \pm 0.2$	35 ± 1.1	23 ± 1.1	11 ± 1.1
BF-TEM	3.8 ± 0.5	$2.5 \pm 0.5$	1.3 <u>+</u> 0.017	39 ± 1.2	26 <u>+</u> 1.1	12 <u>+</u> 1.2
* <sup>+</sup> EDS setup		$(4.1 \times 10^4) \pm 0.1$		$(4.1 \times 10^4) \pm 1.2$	$(4.1 \times 10^4) \pm 1.2$	$(2.1 \times 10^4) \pm 1.3$
*EDS	$(2.5  imes 10^5) \pm 0.017$		$(2.9 \times 10^5) \pm 1.2$	$(2.9 \times 10^5) \pm 1.2$	$(2.7 \times 10^5) \pm 1.3$	
*†EELS setup	$(2.9 \times 10^6) \pm 0.1$		$(3.2 \times 10^6) \pm 1.2$	$(3.2 \times 10^6) \pm 1.2$	$(3.2 \times 10^6) \pm 1.3$	
*EELS	$(1.5  imes 10^6) \pm 0.2$		$(4.8 \times 10^6) \pm 1.3$	$(4.8 \times 10^6) \pm 1.3$	$(4.7 \times 10^6) \pm 1.3$	

Table S3: Electron fluence for each step in the 3DED workflow, and the cumulative electron fluence at each step to two significant figures. Graphical representation of the tabulated data is included in the main text Figure 5. Errors were calculated according to Section 4. Rows marked with asterisk (\*) used estimated electron flux, and rows with dagger (†) used estimated time, to calculate the electron fluence.

#### 4. Calculation of Electron Fluence Errors

Calculation of errors in electron fluence used the following method. All standard deviations were to 99% confidence intervals unless otherwise stated, and error in time was taken as  $\pm 1$  second. The error in electron flux ( $\varepsilon_j$ ) was calculated using the standard deviation of the measured pixel intensity ( $\sigma_o$ ) and the standard deviation of the camera gain ( $\sigma_G$ )

$$\varepsilon_j = \sqrt{\frac{\sigma_G^2}{G^2} + \frac{\sigma_Q^2}{Q^2}} \tag{7}$$

The error in electron fluence ( $\varepsilon_F$ ) was calculated from the error in electron flux and the error in time ( $\varepsilon_t$ )

$$\varepsilon_F = \sqrt{\frac{\varepsilon_j^2}{j^2} + \frac{\varepsilon_t^2}{t^2}} \tag{8}$$

The error in cumulative electron fluence ( $\varepsilon_{CF}$ ) was calculated from the sum of fluence errors

$$\varepsilon_{CF} = \sqrt{\Sigma \varepsilon_F^2} \tag{9}$$

#### 5. Exporting 3DED data from Digital Micrograph

3DED was recorded as .dm4 image files in the "In-Situ" dataset format, which was exported using the "In-Situ Editor" in *DigitalMicrograph* using the "Files to Folder" option to save each frame in the dataset as an individual 32-bit real image. To import data into *DIALS*, the images were saved as a single .dm4 stack. Experimental parameters were passed to *DIALS* via the command line as discussed in the literature; data reduction was then possible following the literature procedure.<sup>7</sup>

To import data into *PETS2* the "In-Situ" datasets were first exported in .tif format. The python script *write\_PETS2\_import.py* was used to write a .pts2 project file for the dataset, following the guidelines in the *PETS2* documentation; data reduction was then completed as per the literature procedure.<sup>8</sup>

To import data into *CrysAlisPro*, "In-Situ" datasets were exported in .dat format using the "Esperanto Importer" feature. The image size, byte type, rotation axis, and electron wavelength were not read from image headers so were supplied manually. A beam stop mask was added after images were imported.



6. Energy Dispersive X-ray and Electron Energy Loss Spectroscopy

Figure S2: Energy dispersive X-ray (EDS, left) and electron energy core-loss (EELS, right) spectra for (a, b) percholorcoronene, (c, d) paracetamol, and (e, f) copper(II) phthalocyanine. Relevant elements have been labelled.

Copper(II) Phthalocyanine		Percholorocoronene		Paracetamol	
Element	Atomic%	Element	Atomic%	Element	Atomic%
С	76.03	С	92.07	С	72.82
N	15.79	Ν	-	Ν	8.17
0	2.13	0	4.02	0	17.42
Cl	-	Cl	1.27	Cl	-
Cu	5.85	Cu	2.11	Cu	1.53

Table S4: Elemental abundancies from EDS of crystals in this study. Each sample exhibits copper and carbon signals due to the TEM grids used, and oxygen is present due to water ice in cryoTEM.

# 7. Instructions for *GiveMeED* and other scripts 7.1 *GiveMeED*

GiveMeED	×
Save Data Pat Path: Sample name:	h X:\TargetFolder SampleName Browse files Open work directory
Variables	0.00251
Frame rate: Notes:	175 e.g. cryo 120.0(5) K
Tilt Range	-60 60 Go to Start Go to End Tilt Neutral
Data Collection	Start 3DED Abort 3DED M may freeze during collection

Figure S3: GiveMeED user interface.

*GiveMeED* has been used successfully with a JEOL 2100Plus TEM with Gatan OneView Camera and JEOL 2100F TEM with Gatan K3 camera. It should work with any TEM using a Gatan camera that has "In-Situ" (video) mode and an installation of *DigitalMicrograph*.

When "Start 3DED" is pressed (Figure S3), *GiveMeED* rotates the TEM stage to the value entered in the "Start Angle" field. The beam blank is turned off and the camera begins recording data while the stage is rotated to the value entered in the "End Angle" field. Once rotation ends the beam is blanked and all data is saved. Pressing "Abort 3DED" will blank the beam and stop the camera recording. 3DED metadata is saved with the name entered in the "Sample name" field at the path entered in the "Path" field. The values in the "Variables" container are for later reference and do not affect data collection. The recorded metadata is a mixture of information pulled from the microscope during data collection and standard values that can be defined in the script. Modifying the metadata file to include/exclude information can be done by adding/removing lines in the "log\_message" string.

The "Abort 3DED" button will blank the beam and stop the camera recording, and the script can also be stopped using the *DigitalMicrograph* kill script shortcut (default "ctrl + numlock"). *GiveMeED* can also be operated without a user interface in which case the variables discussed above are set in the script before executing it. The maximum tilt range of the microscope stage should be checked before attempting 3DED.

Step-by-step instructions for *GiveMeED*:

- 1. Insert the camera in "In-Situ" mode and start the live view with the frame rate set to a suitable value.
- 2. In the camera control pane, check the following variables are set correctly: "save data path", and "experiment name". Ensure that the experiment has a unique name.

In the *GiveMeED* UI the "Path" field should match "save data path", and the "Sample Name" field should match the "experiment name".

- 3. Confirm that the specimen is at eucentric height and remains visible over the tilt range entered into the "Tilt Range" container. Pressing the "Go to Start" and "Go to End" buttons will rotate the stage to the value in the corresponding box. "Tilt Neutral" sets the stage angle to 0 degrees.
- 4. Put the microscope into diffraction mode and press "Start 3DED" to perform data collection. Data and metadata is automatically saved to the file path in Step 2.

# 7.2 Go2Alpha

Go2Alpha	×
- Stage Alph	na —————
Alpha2	60
Go to A1	Go to A2
Tilt Neutra	·

Figure S4: *Go2Alpha* user interface.

*Go2Alpha* is used to quickly tilt the stage between two angles (Figure S4). It has the same functionality as the "Tilt Range" container in the *GiveMeED* UI. Two fields can be set: "Alpha1" and "Alpha2" are angles in degrees, where the "Go to A1" button rotates the stage to the angle in the "Alpha1" field (and *vice versa* for "Alpha2" and "Go to A2"). The "Tilt Neutral" button sets the stage angle to 0 degrees.

# 7.3 AutoResRings

AutoResRings draws a set of resolution rings on an image in Digital Micrograph. It works on whichever image was last selected, including the camera live view. By default it will draw rings at 0.8, 1, 1.4, 2, 4 and 100 Å but different ring radii can be defined in the script. Scaling the rings is done at creation, so the live view should be in diffraction mode for correct calibration (and the rings will not change size if the camera length is changed). Rings can be deleted from an image by clicking on them and pressing the "delete" key.

#### 8. List of Abbreviations

Abbreviation	Full
3DED	Three-dimensional electron diffraction
CCD	Charge-coupled device
CIF	Crystallographic information file
CMOS	Complementary metal-oxide semiconductor
CRED	Continuous rotation electron diffraction
CuPC	Copper(II) phthalocyanine
ED	Electron diffraction
EDS	Energy dispersive X-ray spectroscopy
EELS	Electron energy loss spectroscopy
PCC	Perchlorocoronene
PEDT	Precession-assisted electron diffraction tomography
PXRD	Powder X-ray diffraction
RED	Rotation electron diffraction
SAED	Selected area electron diffraction
SCXRD	Single-crystal X-ray diffraction
TEM	Transmission electron microscope
UI	User interface

Table S5: Abbreviations used in this work.

### 8. References

- 1G. M. Sheldrick, Acta Crystallogr. Sect. Found. Adv., 2015, 71, 3–8.
- 2G. M. Sheldrick, Acta Crystallogr. Sect. C Struct. Chem., 2015, 71, 3–8.
- 3O. V. Dolomanov, L. J. Bourhis, R. J. Gildea, J. A. K. Howard and H. Puschmann, *J. Appl. Crystallogr.*, 2009, **42**, 339–341.
- 4E. Myatt, S. Lata, J. Pitcairn, D. Daisenberger, S. M. Kronawitter, S. A. Hallweger, G. Kieslich, S. P. Argent, J. P. Tidey and M. J. Cliffe, *Chem. Commun.*, 2024, **60**, 14427–14430.
- 5A. L. Spek, J. Appl. Crystallogr., 2003, 36, 7–13.
- 6R. F. Egerton, Ultramicroscopy, 2021, 229, 113363.
- 7A. Vypritskaia, X. Zou, T. Yang and D. G. Waterman, *Acta Crystallogr. Sect. C Struct. Chem.*, 2025, **81**, 1–13.
- 8L. Palatinus, P. Brázda, M. Jelínek, J. Hrdá, G. Steciuk and M. Klementová, Acta Crystallogr. Sect. B Struct. Sci. Cryst. Eng. Mater., 2019, **75**, 512–522.