# A foundation model with multi-variate parallel attention to generate neuronal activity

Francesco S. Carzaniga<sup>1, 2</sup>, Michael Hersche<sup>1</sup>, Abu Sebastian<sup>1</sup>, Kaspar Schindler<sup>2</sup>, Abbas Rahimi<sup>1</sup>

<sup>1</sup>IBM Research – Zurich, Rüschlikon, Switzerland <sup>2</sup>Department of Neurology, Inselspital, Sleep-Wake-Epilepsy-Center, Bern University Hospital, Bern University, Bern, Switzerland frc@zurich.ibm.com

### Abstract

Learning from multi-variate time-series with heterogeneous channel configurations remains a fundamental challenge for deep neural networks (DNNs). particularly in clinical domains such as intracranial electroencephalography (iEEG), where channel setups vary widely across subjects. In this work, we introduce multi-variate parallel attention (MVPA), a novel self-attention mechanism that disentangles content, temporal, and spatial attention, enabling flexible, generalizable, and efficient modeling of time-series data with varying channel counts and configurations. We use MVPA to build MVP-Former, a generative foundation model for human electrophysiology, trained to predict the evolution of iEEG signals across diverse subjects. To support this and future effort by the community, we release the SWEC iEEG dataset, the largest publicly available iEEG dataset to date, comprising nearly 10,000 hours of recordings from heterogeneous clinical sources. MVPFormer leverages MVPA to achieve strong generalization across subjects, demonstrating expert-level performance in seizure detection and outperforming state-ofthe-art Transformer baselines on our SWEC, the MAYO, and the FNUSA dataset. We further validate MVPA on standard time-series forecasting and classification tasks, where it matches or exceeds existing attentionbased models. Together, our contributions establish MVPA as a generalpurpose attention mechanism for heterogeneous time-series and MVPFormer as the first open-source, open-weights, and open-data iEEG foundation model with state-of-the-art clinical performance. The code is available at https://github.com/IBM/multi-variate-parallel-transformer. The SWEC iEEG dataset is available at https://mb-neuro.medical-blocks. ch/public\_access/databases/ieeg/swec\_ieeg.

### 1 Introduction

The increasing availability of multi-variate time-series data across domains, from financial data to sensor networks to clinical recordings, has driven demand for general-purpose neural architectures capable of learning from such data [1, 2, 3]. A fundamental challenge in this setting is channel heterogeneity: different sensors (or channels) often carry information that is both structurally and semantically non-uniform, while the number and the location of channels may vary across instances. This is especially pronounced in intracranial electroencephalography (iEEG) [4], where each subject's electrode layout is unique and tailored to clinical needs. iEEG models [5, 6, 7, 8, 9] often require subject-specific adaptation to account

for new setups, yet they still struggle to generalize. Consequently, effective learning from multi-variate time-series demands models that are flexible and channel-agnostic, without sacrificing locality or the ability to generalize.

In this work, we introduce multi-variate parallel attention (MVPA, Figure 1), a novel selfattention mechanism addressing the structural challenges of channel heterogeneity. MVPA decomposes attention into three components: content-, time-, and channel-based components. Thus, it allows the model to separately learn the semantics of the signal, its temporal dynamics, and spatial (inter-channel) structure. MVPA enables flexible and efficient processing of time-series data, without relying on fixed channel positions or global positional encodings.

To highlight MVPA's ability to handle heterogeneous and clinically relevant time-series, we apply it to the particularly challenging domain of iEEG. Indeed, as mentioned above iEEG recordings present an ideal testbed for models designed to handle variable multichannel structure. We use MVPA to build MVPFormer, a foundation model for human electrophysiology trained via generative pre-training to predict the evolution of brain signals.

MVPFormer is trained on the SWEC iEEG dataset, the largest available iEEG corpus to date with nearly 10,000 hours of multi-channel recordings (or 540,000 channel-hours), collected over a decade in clinical settings and made publicly available as part of this work. We show that MVPFormer not only models neuronal activity during both normal and ictal states, but also supports downstream tasks such as long-term seizure detection, generalizing zero-shot across patients and outperforming previous approaches on clinically relevant benchmarks. Remarkably, MVPFormer surpasses an equivalent purely discriminative version which has not undergone generative pre-training, strengthening the validity of foundation models in iEEG.

We further evaluate MVPA on classical time-series benchmarks, including ETTh and Weather [10, 11] for forecasting and EthanolConcentration, FaceDetection, and others [12] for classification. Here, MVPA matches or outperforms state-of-the-art (SOTA) models such as PatchTST [3]. These results establish MVPA as a competitive attention mechanism for general-purpose time-series beyond iEEG.

Our contributions are:

- Multi-variate parallel attention (MVPA), a novel self-attention mechanism that separately attends to content, temporal, and spatial structure, enabling strong generalization across multi-variate time-series with heterogeneous channels.
- **MVPFormer**, a foundation model for human electrophysiology, powered by MVPA and trained on the largest iEEG corpus available, showing superior generalization across subjects and clinical tasks compared to models which use vanilla attention like Brant-2 [13].
- The release of the **SWEC iEEG dataset**, the largest iEEG dataset publicly available to date, with almost 10,000 hours of highly curated and labeled iEEG recordings.

Moreover, we make all our contributions open source, realizing the first open-data, open-code, and open-weights iEEG foundation model.

## 2 Multi-variate parallel attention (MVPA)

This section introduces multi-variate parallel attention (MVPA), the first main contribution of this work, as shown in Figure 1. We start with vanilla attention for 1D sequences, then present dual-coded attention for 2D sequences, which has higher computational costs. Building on this, we derive MVPA, which *efficiently* attends to both temporal and spatial aspects of multi-variate time-series data.



Figure 1: Multi-variate parallel attention (MVPA). (a) The input signal is divided into temporal and spatial segments. Each query-key interaction is computed for all keys within a local window. (b) MVPA decomposes attention into three components: content-based, computed per segment without positional encoding; time-based, shared across channels and dependent only on temporal distance; and channel-based, shared across time steps and dependent only on spatial distance. (c) The final attention is the sum of these three components, each capturing a distinct aspect of the data.

#### 2.1 MVPA overview

Vanilla attention [14] operates on 1D sequences of embeddings  $(\boldsymbol{x}_1, \boldsymbol{x}_2, ..., \boldsymbol{x}_T)$  of dimension d  $(\boldsymbol{x}_k \in \mathbb{R}^d)$ . It computes the attention between two tokens at positions (i, j) as follows<sup>1</sup>:

$$\boldsymbol{a}_{i,j}^{\mathrm{vanilla}} = \left(\boldsymbol{x}_i + \mathcal{S}_i\right)^T W_q^T W_k \left(\boldsymbol{x}_j + \mathcal{S}_j\right),$$

where  $x_i$  is the query token and  $x_j$  the key token. S is the positional encoding, a vector with the same dimensionality (d) that helps to distinguish between different positions in the sequence.  $W_q$  and  $W_k$  are the learnable query and key matrices.

While vanilla attention has been successfully applied to 1D sequences, its extension to multi-variate time series (i.e., 2D sequences) is not obvious. Specifically, we aim to process sequences of the form  $(\boldsymbol{x}_{1,1}, \boldsymbol{x}_{1,2}, ..., \boldsymbol{x}_{C,T})$ , where *c* indicates the space and *t* the time dimension, while maintaining the embedding dimensionality (i.e.,  $\boldsymbol{x}_{c,t} \in \mathbb{R}^d$ ). One approach is to flatten the 2D data to a 1D sequence (e.g., as done by the Vision Transformer [15]); however, this will yield a loss in spatial structure. Instead, we introduce two separate learnable positional codebooks, representing space ( $\mathcal{C}$ ) and time ( $\mathcal{T}$ ). By equipping self-attention with this dual encoding we can treat the two dimensions individually, which is fundamental in recovering their interplay and would not be possible with vanilla attention:

$$\boldsymbol{a}_{c,t,c',t'}^{\text{dual}} = (\boldsymbol{x}_{c,t} + \mathcal{T}_t + \mathcal{C}_c)^T W_q^T W_k (\boldsymbol{x}_{c',t'} + \mathcal{T}_{t'} + \mathcal{C}_{c'}).$$

Dual attention allows us to exploit the relationship between time and space at the attention level, the most basic computational unit of a Transformer. We believe this allows the architecture to model the time-series at a lower level, and hence more effectively. However, the dual attention mechanism is computationally expensive, as it computes second-order correlations between time and space.

For this reason, we want to squash these cross-correlations. Specifically, we want to push as much of the spatio-temporal computation as possible to the lower levels of processing without

<sup>&</sup>lt;sup>1</sup>For better readability, we describe the attention computation for a single head without activation. In practice, we generalize it to multi-head attention and apply a consecutive softmax non-linearity.

overwhelming it. In contrast, all Transformer models equipped with vanilla attention require ancillary structures to process any relation between time and space [3, 16, 17]. Inspired by Transformer-XL [18], we encode the relative distance in the two dimensions between the segments separately and introduce new learnable bias terms (u, v, w). Finally, we remove the correlation terms and rearrange the expanded equation into three related groups:

$$\begin{split} {}^{\text{MVPA}}_{c,t,c',t'} &= \boldsymbol{x}_{c,t}^T W_q^T W_{k_e} \boldsymbol{x}_{c',t'} + u^T W_{k_e} \boldsymbol{x}_{c',t'} & \text{Content-based attention} \quad (1) \\ &+ \boldsymbol{x}_{c,t}^T W_q^T W_{k_t} \mathcal{T}_{t-t'} + v^T W_{k_t} \mathcal{T}_{t-t'} & \text{Time-based attention} \quad (2) \\ &+ \boldsymbol{x}_{c,t}^T W_q^T W_{k_c} \mathcal{C}_{c-c'} + w^T W_{k_c} \mathcal{C}_{c-c'} & \text{Channel-based attention} \quad (3) \end{split}$$

The three terms above are the attentional components of MVPA. Content-based attention only attends to the content of query and key, without any positional encoding. In this component we compute the relationship between two raw segment embeddings, so we modulate the final attention output without considering any structure of the signal. Timebased attention only attends to the query and the distance in time with the key. In this component, only the relative distance in time is considered, allowing for arbitrary signal lengths without loss of generality. Finally, channel-based attention only attends to the query and the distance in space with the key. Similarly to the time-based component, also the distance in the channel-based component is relative. This feature is particularly interesting for the channel-based component, given the heterogeneity of possible channel setups. In fact, the absolute position of a segment carries no information and is of no interest for final classification. In contrast, the relative encoding scheme allows the channel component to uncover the hidden connection map between the spatial locations, as shown in Appendix G.8.

#### 2.2 Efficient implementation of MVPA

a

As MVPA's computational cost is still quadratic in space and time, we employ several techniques to further reduce the complexity and enable the efficient processing of very large signals. We present here the main techniques, while the details can be found in Appendix A.

Computing the time- and channel-based terms efficiently requires two main techniques. First, we recognize it is unnecessary to compute the full attention matrix, which is quadratic in the context length (i.e., both time and space). By design, all elements of the time-based attention are the same for each channel (see Figure 1b, the green components are all equal), and all elements of the channel-based attention are the same for each time point (see Figure 1c, the blue components are all equal). Hence, complexity is quadratic in one dimension and constant in the other. We then simply repeat the elements along the right dimension at no additional cost. Second, we employ the shifting operation [18] to compute all relative embeddings in one pass.

Content attention remains the most expensive component. To further reduce the cost, we make use of a local attention window [19] that focuses on the most recent L (in our case 10 segments, or 50 seconds) time points. Since time-based attention is not limited, the lookup window still spans the entire context. Thus, for  $L \ll T$ , the total complexity of MVPA is  $O(T^2 \times C + T \times C^2)$ , quadratic in each dimension but subquadratic in the context length. Combining all techniques, MVPA pushes the effective total context length on an NVIDIA A100-80GB GPU to over 10,000 (e.g., 100 channels and 100 time segments).

Additionally, we use grouped query attention [20] to reduce the number of heads without loss of performance. Moreover, we develop FlashMVPA based on FlashAttention [21, 22], implemented in the OpenAI Triton language, providing us with lower-level access to CUDA primitives. FlashMVPA reaches 20 TFlops on an A100.

### 3 MVPFormer

MVPFormeris our novel Transformer-based predictive foundational model equipped with MVPA, that processes heterogeneous multi-variate iEEG data (see Figure 2). While it is customary for language-based Transformer models to employ a finite vocabulary of words, such a definition is non-trivial for iEEG. At the same time, recent works have challenged this discrete paradigm in favor of continuous latent representations [23, 24, 25] and continuous



Figure 2: **MVPFormer architecture and forward pass.** iEEG signals are segmented in time and space, encoded via a wavelet-based encoder, and arranged into a 2D embedding grid. These continuous embeddings are processed by MVPA to model temporal, spatial, and content-based dependencies. MVPFormer predicts the next-in-time embedding while reducing similarity to confounders from the same or other subjects. Notched in the bottom right is the resulting cosine similarity with the true target and the confounders after training.

chain-of-thought mechanisms [26, 27]. In a similar vein, MVPFormer predicts the development of neuronal activity in a *continuous* embedding space governed by a wavelet encoder. We build MVPFormer following the foundational paradigm, with a pre-training dedicated to predicting the future iEEG embedding using a contrastive loss function. Moreover, we show that a successive fine-tuning using LoRA [28] and a simple classification head allows MVPFormer to perform downstream classification tasks. In particular, our results indicate that a model trained in this fashion surpasses an equivalent purely discriminative model (i.e., without generative pre-training), strengthening the validity of foundation models in the iEEG domain as well.

#### 3.1 Architecture

**Wavelet encoder** The first processing step maps the raw iEEG signal to continuous embeddings. We begin by partitioning the raw data into segments of five seconds. Each segment passes independently through a db4 wavelet decomposition, which, depending on the model's overall size, is then linearly projected onto a smaller latent space. This projection, or feature vector, is the embedding. Our method is inspired by wav2vec [29], though we use learnable embeddings. We apply the encoding channel-wise, meaning each segment remains one-dimensional.

**Decoder** MVPFormer is based on the Llama2 architecture [30] with parallel attention and MLP blocks inspired by Megatron-LM [31]. This choice was informed by the selection of a

generative model powerful enough to process brain iEEG signals and computationally light enough to enable extensive testing. We provide two models to evaluate the scaling of our foundational model: MVPFormer-S with 75M parameters and MVPFormer-M (or simply MVPFormer) with 1.2B parameters. Appendix B provides more details on the architectural hyperparameters.

#### 3.2 Training

**Generative pre-training** MVPFormer is pre-trained to generate neuronal activity by predicting successive input segments in time. During pre-training, random input segments from batched windows are used as confounding targets ( $Z = \{z_1, ..., z_n\}$ ), which are plausible but different from the true target. We compute the contrastive loss as follows

$$\mathcal{L}_{c,t} = -\log \frac{\exp(\sin(\boldsymbol{o}_{c,t}, \boldsymbol{e}_{c,t+1})/\tau)}{\sum_{\boldsymbol{z}_k \in Z} \exp(\sin(\boldsymbol{o}_{c,t}, \boldsymbol{z}_k)/\tau)},\tag{4}$$

where  $o_{c,t}$  is the model's output embedding and  $e_{c,t+1}$  the ground-truth next-state embedding. Finally summing over every c, t gives us the optimization target for the generative task. The temperature  $\tau$  is 0.1 and the number of confounders n is 30. See Appendix C.1 for more details on the pre-training.

Validation of pre-training Given the architecture of MVPFormer, we need to ensure that the true target and the confounders are sufficiently well-separated in cosine similarity. We evaluate MVPFormer's ability to predict embeddings of future iEEG signals by comparing the predicted embedding at time t to the ground truth embedding at t + 1 using the cosine similarity. To do so, we introduce two references: (1) the embedding at t + 2, which is highly correlated with t + 1, and (2) a random future segment sampled within the next two minutes. In fact, given the high auto-correlation of iEEG signals, a naïve prediction model could simply predict again time t, and be moderately successful due to its similarity with t + 1. Our results (see Figure 2 and Appendix F) show that the wavelet-based encoder ensures signal features are well preserved, mapping even mildly similar signals to distinct embeddings.

**LoRA fine-tuning for seizure classification** For seizure classification, we use a small classification head (i.e., a linear layer). This layer has input size equal to the decoder's block output size, and output size equal to the dimensionality of the classification task (i.e., 2 for seizure classification). The input to this classification head is the channel-averaged output of the last signal segment in time. The output of the classification head is then passed through a softmax to compute the binary cross-entropy loss. We further adopt LoRA [28] to perform parameter-efficient fine-tuning. We only fine-tune the q and v layers of the self-attention in the base MVPFormer model, with a LoRA rank of 8 and and alpha of 16. This leads to a number of trainable parameters during fine-tuning of approximately 0.1% of the base model.

#### 4 SWEC iEEG dataset

The lack of publicly available large-scale iEEG datasets has been a significant obstacle to the development of foundation models for this modality. In fact, while EEG datasets are abundant [32, 33] and large [34], with tens of thousands of recording hours, such resources are lacking in the iEEG domain. Due to significant barriers tied to data collection and privacy, available iEEG datasets cover few hours [35] (35 subjects and 290 hours) and subjects [36] (10 subjects and 43 hours).

In an effort towards addressing this issue, together with this work we open-source the SWEC iEEG dataset, a large-scale iEEG dataset consisting of a total of 68 subjects, 9328 hours of recording, and 704 ictal events. To our knowledge, the SWEC iEEG dataset is the largest publicly available iEEG dataset, fully curated and labelled by experienced clinicians. Due to institutional data privacy concerns, the dataset does not contain information about the location of the channels in the brain. Appendix D reports more details and illustrates two example recordings.

The iEEG signals were recorded intracranially for clinical purposes with a sampling rate of either  $512 \,\text{Hz}$  or  $1024 \,\text{Hz}$ , which was then normalized to  $512 \,\text{Hz}$  before training MVPFormer.

The signals were median-referenced and digitally band-pass filtered between 0.5 and 120 Hz using a fourth-order Butterworth filter, both in a forward and backward pass to minimize phase distortions. All the recordings were inspected by an expert for identification of seizure onsets and offsets, and to remove channels corrupted by artifacts.

#### 5 Experiments

#### 5.1 Setup

**Preprocessing** We divide each recording into windows of 500 seconds each, with an overlap of 99% to increase the size of the training dataset and the number of training tokens. Then, each window is divided into 100 segments (each five seconds long), yielding a total of 39B segments of which 390M are unique.

**Pre-training** We pre-train MVPFormer on 18 subjects, leaving the remaining 50 subjects for testing. MVPFormerM is pre-trained for 1.2M steps on a single node with 8 NVIDIA A100-80GB GPUs for two weeks. The chosen optimizer is FusedAdam with 0.1 weight decay, from the Deepspeed library compiled on the specific machine. The training used bf16-mixed DeepSpeed stage 2 without activation checkpointing. The learning rate is fixed to  $10^{-4}$ . The training environment includes PyTorch 2.0, PyTorch Lightning 2.0, and Triton 2.1.0.

**Fine-tuning** After pre-training, we further fine-tune the MVPFormer on the same 18 subjects on the seizure detection task, and then test in a zero-shot manner on nearly 7,000 hours of iEEG data from 50 unseen subjects. All subjects suffer from epilepsy. This allows us to assemble a rigorous set of results which are relevant to a real-world application. Moreover, to keep computational cost moderate, we use a subset of the channels of each subject for testing. We select the channels based on a visual inspection by a non-expert to include only the least noisy ones. In a real-world clinical scenario selection and validation would comprise a minimal additional burden for the expert. The number of the chosen channels and their position varies for each subject, from 30 up to 50 channels (see Appendix G.6).

#### 5.2 Seizure detection

We evaluate MVPFormer on seizure detection. First, we consider a clinically realistic setup that compares model predictions to a board-certified neurologist annotations using Cohen's Kappa. To do so, the predictions are post-processed to yield episodic results (see Appendix E). Cohen's Kappa [37, 38, 39] is widely used to quantify inter-rater reliability in seizure classification. The Landis and Koch criteria [40] (see Appendix E) are often used in practice to evaluate human performance. Expert-level performance in the seizure classification task varies considerably, from 0.58 [41] to 0.53 [42] in EEG, to 0.57 [43] in iEEG. We consider Kappa values above 0.53 to be expert-level. We



Figure 3: **MVPFormer's inter-rater agreement on seizure detection.** The average Kappa score between MVPFormer-M and the human expert is 0.57, competitive with the values obtained between human experts.

must also consider that our evaluation setup involves many more subjects and ictal events than are reported for human experts, making this task more challenging for MVPFormer. MVPFormer achieves an average Kappa of 0.57 across 50 unseen subjects from the SWEC iEEG dataset, matching human expert performance (see Figure 3). Importantly, agreement varies by subject, reflecting the clinical reality that seizure presentation complexity strongly affects classification (see Appendix G.7). Overall, MVPFormer demonstrates expert-level seizure classification across a large, heterogeneous cohort. This performance, combined with its low false positive rate (0.17 fp/h), positions MVPFormer as a promising clinical assistant for real-world iEEG analysis.

Table 1: **Results on seizure detection.** We compare MVPFormer with Brant-2, the current SOTA Transformer model for iEEG, and MV-Llama, our vanilla attention-based baseline.

|              |           |       | SWEC     |      |      |      |      | MAYO |      | FNUSA |      |      |
|--------------|-----------|-------|----------|------|------|------|------|------|------|-------|------|------|
|              |           |       | Episodic |      |      | Raw  | Raw  |      |      | Raw   |      |      |
| Model        | Attention | Kappa | f1       | sens | fp/h | f1   | f1   | sens | spec | f1    | sens | spec |
| MVPFormer    | MVPA      | 0.57  | 0.56     | 0.73 | 0.17 | 0.50 | 0.36 | 0.38 | 0.91 | 0.46  | 0.94 | 0.10 |
| MVPFormer-S  | MVPA      | 0.54  | 0.51     | 0.69 | 0.11 | 0.49 | 0.35 | 0.41 | 0.88 | 0.46  | 0.99 | 0.03 |
| Brant-2 [13] | Vanilla   | 0.08  | 0.00     | 0.00 | 0.06 | 0.00 | 0.19 | 1.00 | 0.18 | 0.46  | 0.99 | 0.02 |
| MV-Llama     | Vanilla   | 0.05  | 0.02     | 0.01 | 0.00 | 0.01 | /    | /    | /    | /     | /    | /    |

Second, we consider a conventional evaluation (see Appendix G.1) based on F1-score, sensitivity, specificity, and number of false positives per hour (fp/h). We compare MVPFormer against two strong baselines: (1) Brant-2 [13], a SOTA Transformer for iEEG, is fine-tuned with published weights. Brant-2 requires all subjects to have the same number of channels for classification; hence, we were not able to test all subjects with it. (2) MV-Llama, an ablation of MVPFormer-S that replaces MVPA with vanilla attention, is trained identically to MVPFormer. We test the considered models zero-shot on our SWEC iEEG dataset and apply them on the MAYO and FNUSA datasets [35] (see Appendix G.9).

As shown in Table 1, both baselines fail to generalize on our SWEC iEEG dataset, achieving Kappa scores of just 0.08 and 0.05, while MVPFormer achieves a 0.57 and 0.54 in medium and small configurations, respectively. Moreover, MVPFormer outperforms the baselines on MAYO (highest f1-score of 0.36) and is on par on FNUSA. As a further baseline — to validate our choice of pre-training — we also compare MVPFormer-S with an equivalent model built without the generative base task (see Appendix G.11). With this setup, the purely discriminative model only reaches a Kappa score of 0.52, inferior to the equivalent MVPFormer-S which reached 0.54.

#### 5.3 Time-series forecasting and classification

To provide a baseline evaluation of MVPA against established alternatives, we compare MVPFormer with existing SOTA architectures on classical long-term forecasting and classification tasks. Table 2 reports the results of MVPFormer, the vanilla Transformer [14], PatchTST [3], and TimeMixer [2] on the ETTh1, ETTh2, and Weather datasets [10, 11]. MVPFormer always equals or surpasses the baselines (see Appendix G.12 for the full results).

Table 2: **Results on the time-series forecasting.** We report the mean-squared error (MSE) and mean-absolute error (MAE) averaged over all forecasting lengths.

| Model   | Model  |      | MVPFormer |      | former | Patch | $\Gamma ST [3]$ | TimeMixer [2] |      |  |
|---------|--------|------|-----------|------|--------|-------|-----------------|---------------|------|--|
| Dataset | Length | MSE  | MAE       | MSE  | MAE    | MSE   | MAE             | MSE           | MAE  |  |
| ETTh1   | Avg.   | 0.45 | 0.45      | 1.00 | 0.80   | 0.45  | 0.45            | 0.45          | 0.44 |  |
| ETTh2   | Avg.   | 0.38 | 0.41      | 3.37 | 1.48   | 0.39  | 0.41            | 0.39          | 0.41 |  |
| Weather | Avg.   | 0.25 | 0.28      | 0.59 | 0.53   | 0.26  | 0.28            | 0.25          | 0.28 |  |

Moreover, we evaluate MVPA on common classification tasks, against the vanilla Transformer and PatchTST on the EthanolConcentration (EtCo), FaceDetection (FaDe), HandWriting (HaWr), Heartbeat (HaBe), JapaneseVowels (JaVo), PEMS-SF (PEMS), SCP1, SCP2, SpokenArabic (SpAr), and Uwave datasets [12].

TimeMixer is a forecasting-only architecture, so we could not test it. Table 3 shows that MVPFormer achieves SOTA results on general classification tasks as well. At the same time, these results highlight the generalization capability of MVPFormer that, in contrast to other models such as TimeMixer, is effective in both forecasting and classification.

### 6 Related Works

Single-channel data has been treated as 1D sequences for tasks like speech recognition, where the signal is divided into patches that serve as tokens [29, 44]. Extending vanilla attention to multi-dimensional data, such as images, is more complex. The Vision Transformer [15] processes images by flattening 2D patches into a 1D sequence, losing spatial structure in the process. However, this approach is inflexible and unsuitable for generalizing to images with different heights and widths.

For multi-variate time-series, such as EEG, Transformers face challenges due to the need to preserve both time and channel information [17]. Channel-independent approaches [3] reuse vanilla self-attention and discard all information content in the time dimension, while channel-mixing promises to preserve it by either fusing the channels [45] or processing them sequentially [16]. Specifically for iEEG and EEG, there exist few Transformer-based solutions [46, 13]. Since electrode placements vary widely across subjects, these models struggle with the heterogeneous nature of the data. The complex interplay between time and space, where distant brain regions may be more strongly connected than nearby ones, makes it difficult for conventional attention mechanisms to effectively process iEEG signals.

### 7 Conclusion

We introduce MVPA, a novel attention mechanism designed to effectively process multivariate time-series data, exemplified by its application to iEEG signal analysis. MVPA enables MVPFormer, a foundation model trained on our novel SWEC iEEG dataset, to capture complex interactions between time and spatial dimensions in multi-variate time-series. We also contribute the SWEC iEEG dataset itself, as the largest iEEG dataset currently publicly available. MVPFormer is trained following the foundational paradigm to predict the next brain states, and then further fine-tuned on a seizure detection task. MVPA ensures robust performance, reaching high inter-rater agreement (0.57 Kappa score) on our large scale and challenging SWEC iEEG dataset, notably surpassing the SOTA Brant-2 [13] (0.08). Moreover, MVPA equals or surpasses the SOTA also in classical time-series forecasting and classification tasks. Overall, our results show that MVPA affords MVPFormer superior generalization capabilities while maintaining computational efficiency and scalability, marking a significant advancement in the analysis of time-series data and iEEG in particular.

Limitations While our SWEC iEEG dataset represents a notable step towards increasing the size of iEEG datasets, large language models have clearly shown [47] that ever larger amounts of data are necessary to further develop a foundation model. Indeed, even our MVPFormer-S is already oversized when adapting scaling laws from natural language processing [48, 49] (see Appendix G.2), and would benefit from a larger pre-training dataset. Moreover, recent works have highlighted the relevance of high-quality data to the performance of the final model [50], also in the (i)EEG domain [51]. This indicates a path forward to better cross-modal training of foundation models, to increase the breadth of their capabilities.

### 8 Broader impact

Effective tools for iEEG analysis provide many practical use-cases [52], including supporting neurologists by processing large amounts of data collected through long-term implantations [53] and wearable devices [54]. Furthermore, models able to decode the information

|                 | EtCo | FaDe | HaWr | HaBe | JaVo | PEMS | SCP1 | SCP2 | $\operatorname{SpAr}$ | Uwave |
|-----------------|------|------|------|------|------|------|------|------|-----------------------|-------|
| MVPFormer       | 0.33 | 0.66 | 0.21 | 0.70 | 0.95 | 0.86 | 0.86 | 0.54 | 0.97                  | 0.80  |
| Transformer     | 0.29 | 0.64 | 0.20 | 0.70 | 0.91 | 0.84 | 0.83 | 0.54 | 0.95                  | 0.80  |
| PatchTST [3]    | 0.29 | 0.67 | 0.23 | 0.72 | 0.95 | 0.85 | 0.83 | 0.51 | 0.97                  | 0.82  |
| TimeMixer $[2]$ | N.A.                  | N.A.  |

Table 3: Accuracy on time-series classification tasks.

contained in iEEG may lead to a deeper understanding of the mechanisms underpinning neuronal activity [55, 4]. At the same time, the adaptability of MVPA extends beyond iEEG analysis, offering a promising solution for a wide range of applications in neuroscience and time-series forecasting across other domains.

Nonetheless, any tool developed for the clinical domain poses risks. Over-reliance on automated decision-making could lead to reduced oversight by human experts, diminishing the effectiveness of care. Moreover, when making large-scale iEEG datasets such as ours publicly available, great care must be placed to ensure full compliance with all regulatory and ethical guidelines to avoid putting patients in harm's way.

#### References

- M. Jin, S. Wang, L. Ma, Z. Chu, J. Y. Zhang, X. Shi, P.-Y. Chen, Y. Liang, Y.-F. Li, S. Pan, and Q. Wen, "Time-LLM: Time series forecasting by reprogramming large language models," in *The Twelfth International Conference on Learning Representations*, 2024.
- [2] S. Wang, H. Wu, X. Shi, T. Hu, H. Luo, L. Ma, J. Y. Zhang, and J. ZHOU, "Timemixer: Decomposable multiscale mixing for time series forecasting," in *International Conference* on Learning Representations (ICLR), 2024.
- [3] Y. Nie, N. H. Nguyen, P. Sinthong, and J. Kalagnanam, "A time series is worth 64 words: Long-term forecasting with transformers," in *International Conference on Learning Representations (ICLR)*, 2023.
- [4] P. L. Nunez and R. Srinivasan, *Electric Fields of the Brain*. Oxford University Press, Jan. 2006.
- [5] L. Kuhlmann, K. Lehnertz, M. P. Richardson, B. Schelter, and H. P. Zaveri, "Seizure prediction — ready for a new era," *Nature Reviews Neurology*, vol. 14, no. 10, 2018.
- [6] K.-O. Cho and H.-J. Jang, "Comparison of different input modalities and network structures for deep learning-based seizure detection," *Scientific Reports*, 2020.
- [7] P. Thuwajit, P. Rangpong, P. Sawangjai, P. Autthasan, R. Chaisaen, N. Banluesombatkul, P. Boonchit, N. Tatsaringkansakul, T. Sudhawiyangkul, and T. Wilaiprasitporn, "EEGWaveNet: Multiscale CNN-based spatiotemporal feature extraction for EEG seizure detection," *IEEE Transactions on Industrial Informatics*, vol. 18, no. 8, 2022.
- [8] C. Wang, V. Subramaniam, A. U. Yaari, G. Kreiman, B. Katz, I. Cases, and A. Barbu, "BrainBERT: Self-supervised representation learning for intracranial recordings," in *International Conference on Learning Representations (ICLR)*, 2023.
- [9] K. Saab, S. Tang, M. Taha, C. Lee-Messer, C. Ré, and D. L. Rubin, "Towards trustworthy seizure onset detection using workflow notes," *npj Digital Medicine*, vol. 7, no. 1, 2024.
- [10] H. Zhou, S. Zhang, J. Peng, S. Zhang, J. Li, H. Xiong, and W. Zhang, "Informer: Beyond efficient transformer for long sequence time-series forecasting," in *Proceedings* of the AAAI Conference on Artificial Intelligence, vol. 35, no. 12. AAAI Press, 2021, pp. 11106–11115.
- [11] H. Wu, J. Xu, J. Wang, and M. Long, "Autoformer: Decomposition transformers with Auto-Correlation for long-term series forecasting," in Advances in Neural Information Processing Systems (NeurIPS), vol. 34, 2021.
- [12] X. Liu and W. Wang, "Deep time series forecasting models: A comprehensive survey," *Mathematics*, vol. 12, no. 10, 2024.
- [13] Z. Yuan, F. Shen, M. Li, Y. Yu, C. Tan, and Y. Yang, "Brant-2: Foundation model for brain signals," arXiv preprint arXiv:2402.1025, 2024.
- [14] A. Vaswani, N. Shazeer, N. Parmar, J. Uszkoreit, L. Jones, A. N. Gomez, L. u. Kaiser, and I. Polosukhin, "Attention is all you need," in *Advances in Neural Information Processing Systems (NeurIPS)*, vol. 30, 2017.
- [15] A. Dosovitskiy, L. Beyer, A. Kolesnikov, D. Weissenborn, X. Zhai, T. Unterthiner, M. Dehghani, M. Minderer, G. Heigold, S. Gelly, J. Uszkoreit, and N. Houlsby, "An image is worth 16x16 words: Transformers for image recognition at scale," in *International Conference on Learning Representations (ICLR)*, 2021.
- [16] Y. Zhang and J. Yan, "Crossformer: Transformer utilizing cross-dimension dependency for multivariate time series forecasting," in *International Conference on Learning Representations (ICLR)*, 2023.

- [17] Q. Wen, T. Zhou, C. Zhang, W. Chen, Z. Ma, J. Yan, and L. Sun, "Transformers in time series: A survey," arXiv preprint arXiv:2202.07125, 2022.
- [18] Z. Dai, Z. Yang, Y. Yang, J. Carbonell, Q. Le, and R. Salakhutdinov, "Transformer-XL: Attentive language models beyond a fixed-length context," in *Proceedings of the 57th Annual Meeting of the Association for Computational Linguistics*, 2019.
- [19] R. Child, S. Gray, A. Radford, and I. Sutskever, "Generating long sequences with sparse transformers," arXiv preprint arXiv:1904.10509, 2019.
- [20] J. Ainslie, J. Lee-Thorp, M. de Jong, Y. Zemlyanskiy, F. Lebron, and S. Sanghai, "Gqa: Training generalized multi-query transformer models from multi-head checkpoints," in *Proceedings of the 2023 Conference on Empirical Methods in Natural Language Processing*, 2023.
- [21] T. Dao, D. Y. Fu, S. Ermon, A. Rudra, and C. Ré, "FlashAttention: Fast and memoryefficient exact attention with IO-awareness," in Advances in Neural Information Processing Systems (NeurIPS), vol. 35, 2022.
- [22] T. Dao, "FlashAttention-2: Faster attention with better parallelism and work partitioning," in *International Conference on Learning Representations (ICLR)*, 2024.
- [23] Y. LeCun, "A path towards autonomous machine intelligence version 0.9. 2, 2022-06-27," Open Review, vol. 62, no. 1, pp. 1–62, 2022.
- [24] J. Tack, J. Lanchantin, J. Yu, A. Cohen, I. Kulikov, J. Lan, S. Hao, Y. Tian, J. Weston, and X. Li, "LLM pretraining with continuous concepts," arXiv preprint arXiv:2502.08524, 2025.
- [25] M. Tschannen, C. Eastwood, and F. Mentzer, "GIVT: Generative infinite-vocabulary transformers," in *Computer Vision – ECCV 2024*. Cham: Springer Nature Switzerland, 2025, pp. 292–309.
- [26] S. Hao, S. Sukhbaatar, D. Su, X. Li, Z. Hu, J. Weston, and Y. Tian, "Training large language models to reason in a continuous latent space," arXiv preprint arXiv:2412.06769, 2024.
- [27] J. Geiping, S. McLeish, N. Jain, J. Kirchenbauer, S. Singh, B. R. Bartoldson, B. Kailkhura, A. Bhatele, and T. Goldstein, "Scaling up test-time compute with latent reasoning: A recurrent depth approach," arXiv preprint arXiv:2502.05171, 2025.
- [28] E. J. Hu, yelong shen, P. Wallis, Z. Allen-Zhu, Y. Li, S. Wang, L. Wang, and W. Chen, "LoRA: Low-rank adaptation of large language models," in *International Conference on Learning Representations (ICLR)*, 2022.
- [29] S. Schneider, A. Baevski, R. Collobert, and M. Auli, "wav2vec: Unsupervised pretraining for speech recognition," in *Interspeech 2019*, 2019.
- [30] H. Touvron, L. Martin, K. Stone, P. Albert, A. Almahairi, Y. Babaei, N. Bashlykov, S. Batra, P. Bhargava, S. Bhosale *et al.*, "Llama 2: Open foundation and fine-tuned chat models," *arXiv preprint arXiv:2307.09288*, 2023.
- [31] M. Shoeybi, M. Patwary, R. Puri, P. LeGresley, J. Casper, and B. Catanzaro, "Megatronlm: Training multi-billion parameter language models using model parallelism," arXiv preprint arXiv:1909.08053, 2019.
- [32] M. Tangermann, K.-R. Müller, A. Aertsen, N. Birbaumer, C. Braun, C. Brunner, R. Leeb, C. Mehring, K. J. Miller, G. R. Müller-Putz, G. Nolte, G. Pfurtscheller, H. Preissl, G. Schalk, A. Schlögl, C. Vidaurre, S. Waldert, and B. Blankertz, "Review of the BCI competition IV," *Frontiers in Neuroscience*, 2012.
- [33] A. Shoeb, "CHB-MIT scalp EEG database," 2010.
- [34] I. Obeid and J. Picone, "The temple university hospital eeg data corpus," Frontiers in Neuroscience, vol. 10, 2016.
- [35] P. Nejedly, V. Kremen, V. Sladky, J. Cimbalnik, P. Klimes, F. Plesinger, F. Mivalt, V. Travnicek, I. Viscor, M. Pail, J. Halamek, B. H. Brinkmann, M. Brazdil, P. Jurak, and G. Worrell, "Multicenter intracranial eeg dataset for classification of graphoelements and artifactual signals," *Scientific Data*, vol. 7, no. 1, 2020.
- [36] C. Wang, A. U. Yaari, A. K. Singh, V. Subramaniam, D. Rosenfarb, J. DeWitt, P. Misra, J. R. Madsen, S. Stone, G. Kreiman, B. Katz, I. Cases, and A. Barbu, "Brain treebank: Large-scale intracranial recordings from naturalistic language stimuli," in *The Thirtyeight Conference on Neural Information Processing Systems Datasets and Benchmarks Track*, 2024.
- [37] H. Danker-Hopfe, D. Kunz, G. Gruber, G. Klösch, J. L. Lorenzo, S. L. Himanen, B. Kemp, T. Penzel, J. Röschke, H. Dorn, A. Schlögl, E. Trenker, and G. Dorffner,

"Interrater reliability between scorers from eight european sleep laboratories in subjects with different sleep disorders," *Journal of Sleep Research*, vol. 13, no. 1, pp. 63–69, Feb. 2004.

- [38] A. Schlögl, F. Lee, H. Bischof, and G. Pfurtscheller, "Characterization of four-class motor imagery EEG data for the BCI-competition 2005," *Journal of Neural Engineering*, vol. 2, no. 4, pp. L14–L22, Aug. 2005.
- [39] M. L. McHugh, "Interrater reliability: the kappa statistic," Biochemia Medica, 2012.
- [40] J. R. Landis and G. G. Koch, "The measurement of observer agreement for categorical data," *Biometrics*, vol. 33, no. 1, p. 159, Mar. 1977.
- [41] J. Halford, D. Shiau, J. Desrochers, B. Kolls, B. Dean, C. Waters, N. Azar, K. Haas, E. Kutluay, G. Martz, S. Sinha, R. Kern, K. Kelly, J. Sackellares, and S. LaRoche, "Inter-rater agreement on identification of electrographic seizures and periodic discharges in icu eeg recordings," *Clinical Neurophysiology*, vol. 126, no. 9, 2015.
- [42] A. C. Grant, S. G. Abdel-Baki, J. Weedon, V. Arnedo, G. Chari, E. Koziorynska, C. Lushbough, D. Maus, T. McSween, K. A. Mortati, A. Reznikov, and A. Omurtag, "EEG interpretation reliability and interpreter confidence: A large single-center study," *Epilepsy & Behavior*, vol. 32, 2014.
- [43] M. Quigg, F. Sun, N. B. Fountain, B. C. Jobst, V. S. S. Wong, E. Mirro, S. Brown, and D. C. Spencer, "Interrater reliability in interpretation of electrocorticographic seizure detections of the responsive neurostimulator," *Epilepsia*, vol. 56, no. 6, 2015.
- [44] Conformer: Convolution-augmented transformer for speech recognition, 2020.
- [45] T. Zhou, Z. Ma, Q. Wen, X. Wang, L. Sun, and R. Jin, "FEDformer: Frequency enhanced decomposed transformer for long-term series forecasting," in *Proc. 39th International Conference on Machine Learning (ICML 2022)*, 2022.
- [46] D. Zhang, Z. Yuan, Y. Yang, J. Chen, J. Wang, and Y. Li, "Brant: Foundation model for intracranial neural signal," in Advances in Neural Information Processing Systems (NeurIPS), vol. 36, 2023.
- [47] N. Muennighoff, A. M. Rush, B. Barak, T. L. Scao, N. Tazi, A. Piktus, S. Pyysalo, T. Wolf, and C. Raffel, "Scaling data-constrained language models," in Advances in Neural Information Processing Systems (NeurIPS), 2023.
- [48] J. Hoffmann, S. Borgeaud, A. Mensch, E. Buchatskaya, T. Cai, E. Rutherford, D. d. L. Casas, L. A. Hendricks, J. Welbl, A. Clark, T. Hennigan, E. Noland, K. Millican, G. v. d. Driessche, B. Damoc, A. Guy, S. Osindero, K. Simonyan, E. Elsen, J. W. Rae, O. Vinyals, and L. Sifre, "Training compute-optimal large language models," Advances in Neural Information Processing Systems (NeurIPS), vol. 36, 2022.
- [49] J. Kaplan, S. McCandlish, T. Henighan, T. B. Brown, B. Chess, R. Child, S. Gray, A. Radford, J. Wu, and D. Amodei, "Scaling laws for neural language models," arXiv preprint arXiv:2001.08361, 2020.
- [50] G. Penedo, H. Kydlíček, L. B. allal, A. Lozhkov, M. Mitchell, C. Raffel, L. V. Werra, and T. Wolf, "The fineweb datasets: Decanting the web for the finest text data at scale," in *The Thirty-eight Conference on Neural Information Processing Systems Datasets and Benchmarks Track*, 2024.
- [51] F. S. Carzaniga, G. T. Hoppeler, M. Hersche, K. Schindler, and A. Rahimi, "The case for cleaner biosignals: High-fidelity neural compressor enables transfer from cleaner iEEG to noisier EEG," in *International Conference on Learning Representations (ICLR)*, 2025.
- [52] U. K. Patel, A. Anwar, S. Saleem, P. Malik, B. Rasul, K. Patel, R. Yao, A. Seshadri, M. Yousufuddin, and K. Arumaithurai, "Artificial intelligence as an emerging technology in the current care of neurological disorders," *Journal of Neurology*, vol. 268, no. 5, 2019.
- [53] G. K. Bergey, M. J. Morrell, E. M. Mizrahi, A. Goldman, D. King-Stephens, D. Nair, S. Srinivasan, B. Jobst, R. E. Gross, D. C. Shields, G. Barkley, V. Salanova, P. Olejniczak, A. Cole, S. S. Cash, K. Noe, R. Wharen, G. Worrell, A. M. Murro, J. Edwards, M. Duchowny, D. Spencer, M. Smith, E. Geller, R. Gwinn, C. Skidmore, S. Eisenschenk, M. Berg, C. Heck, P. Van Ness, N. Fountain, P. Rutecki, A. Massey, C. O'Donovan, D. Labar, R. B. Duckrow, L. J. Hirsch, T. Courtney, F. T. Sun, and C. G. Seale, "Long-term treatment with responsive brain stimulation in adults with refractory partial seizures," *Neurology*, 2015.
- [54] E. Donner, O. Devinsky, and D. Friedman, "Wearable digital health technology for epilepsy," *New England Journal of Medicine*, vol. 390, no. 8, 2024.

- [55] J. O. McNamara, "Emerging insights into the genesis of epilepsy," Nature, vol. 399, no. 6738, pp. A15–A22, Jun. 1999.
- [56] X. Si, Z. Yang, X. Zhang, Y. Sun, W. Jin, L. Wang, S. Yin, and D. Ming, "Patientindependent seizure detection based on long-term ieeg and a novel lightweight cnn," *Journal of Neural Engineering*, vol. 20, no. 1, 2023.
- [57] W. Cui, W. Jeong, P. Thölke, T. Medani, K. Jerbi, A. A. Joshi, and R. M. Leahy, "Neuro-GPT: Developing a foundation model for EEG," arXiv preprint arXiv:2311.03764, 2023.
- [58] N. Kitaev, L. Kaiser, and A. Levskaya, "Reformer: The efficient transformer," in International Conference on Learning Representations (ICLR), 2020.
- [59] B. Wang, "Mesh-Transformer-JAX: Model-Parallel Implementation of Transformer Language Model with JAX," 2021.
- [60] A. v. d. Oord, Y. Li, and O. Vinyals, "Representation learning with contrastive predictive coding," arXiv preprint arXiv:1807.03748, 2018.
- [61] M. P. Beenhakker and J. R. Huguenard, "Neurons that fire together also conspire together: Is normal sleep circuitry hijacked to generate epilepsy?" *Neuron*, 2009.
- [62] H. P. Zaveri, B. Schelter, C. A. Schevon, P. Jiruska, J. G. Jefferys, G. Worrell, A. Schulze-Bonhage, R. B. Joshi, V. Jirsa, M. Goodfellow, C. Meisel, and K. Lehnertz, "Controversies on the network theory of epilepsy: Debates held during the ictals 2019 conference," *Seizure*, vol. 78, 2020.
- [63] M. L. Martini, A. A. Valliani, C. Sun, A. B. Costa, S. Zhao, F. Panov, S. Ghatan, K. Rajan, and E. K. Oermann, "Deep anomaly detection of seizures with paired stereoelectroencephalography and video recordings," *Scientific Reports*, vol. 11, no. 1, 2021.
- [64] S. Ziyabari, V. Shah, M. Golmohammadi, I. Obeid, and J. Picone, "Objective evaluation metrics for automatic classification of EEG events," arXiv preprint arXiv:1712.10107, 2017.
- [65] V. Shah, I. Obeid, J. Picone, G. Ekladious, R. Iskander, and Y. Roy, "Validation of temporal scoring metrics for automatic seizure detection," in 2020 IEEE Signal Processing in Medicine and Biology Symposium (SPMB), 2020.
- [66] A. Van de Vel, K. Verhaert, and B. Ceulemans, "Critical evaluation of four different seizure detection systems tested on one patient with focal and generalized tonic and clonic seizures," *Epilepsy & Behavior*, vol. 37, 2014.
- [67] E. Bruno, P. F. Viana, M. R. Sperling, and M. P. Richardson, "Seizure detection at home: Do devices on the market match the needs of people living with epilepsy and their caregivers?" *Epilepsia*, 2020.
- [68] L. A. Shokooh, D. H. Toffa, P. Pouliot, F. Lesage, and D. K. Nguyen, "Intracranial eeg seizure onset and termination patterns and their association," *Epilepsy Research*, vol. 176, 2021.
- [69] J. Gotman, "A few thoughts on "what is a seizure?"," *Epilepsy & Behavior*, vol. 22, pp. S2–S3, Dec. 2011.
- [70] R. F. Betzel, J. D. Medaglia, A. E. Kahn, J. Soffer, D. R. Schonhaut, and D. S. Bassett, "Structural, geometric and genetic factors predict interregional brain connectivity patterns probed by electrocorticography," *Nature Biomedical Engineering*, 2019.
- [71] J. C. Pang, K. M. Aquino, M. Oldehinkel, P. A. Robinson, B. D. Fulcher, M. Breakspear, and A. Fornito, "Geometric constraints on human brain function," *Nature*, vol. 618, no. 7965, 2023.
- [72] H. H. Jasper, "Report of the committee on methods of clinical examination in electroencephalography: 1957," *Electroencephalography and Clinical Neurophysiology*, vol. 10, no. 2, pp. 370–375, 1958.

### A Details on multi-variate parallel attention (MVPA)

Algorithm 1 illustrates the multi-variate parallel attention algorithm.

Algorithm 1: Computation of MVPA **Input:**  $\boldsymbol{x}_{c.t} \in \mathbb{R}^{n_{\text{embed}}}$  output token of Encoder;  $n_{\text{embed}} = 768$ **Output:**  $o_{c,t} \in \mathbb{R}^{n_{\text{embed}}}$  output attention **Data:** t time encoding; c channel encoding; u, y, w biases;  $h \in [1, \ldots, n_{\text{heads}}]$ ;  $h_{k,v} \in [1, ..., n_{gqa}]$ 1 def MVPAttention $(x_{c,t})$ : # Compute query separately from key and value due to GQA  $m{q}_{c,t}^{h} \leftarrow ext{LINEARNOBIAS}(m{x}_{c,t}) \\ m{k}_{c,t}^{h_{k,v}}, m{v}_{c,t}^{h_{k,v}} \leftarrow ext{LINEARNOBIAS}(m{x}_{c,t}) \\ # ext{ Compute the three components of MVPA}$  $\mathbf{2}$ 3  $oldsymbol{g}^h_{(c,t),(c',t')} \leftarrow (oldsymbol{q}^h_{c,t} + oldsymbol{u}^{h_{k,v}})^{\hat{T}}oldsymbol{q}^h_{c',t'}$ 4 # Time and channel components are independent of the key content, so they do not need to be recomputed  $oldsymbol{s}_t^h \leftarrow (oldsymbol{q}_{c,t}^h + oldsymbol{y}^{h_{k,v}})^T oldsymbol{t}^{h_{k,v}}$ 5  $oldsymbol{l}_{c}^{h} \leftarrow (oldsymbol{q}_{c,t}^{h} + oldsymbol{w}^{h_{k,v}})^{T}oldsymbol{c}^{h_{k,v}}$ 6 # Shift the time and channel components to avoid recomputation, from Transformer-XL  $m{s}^h_{t,t'} \leftarrow ext{ShiftTime}_{t'}(m{s}^h_t) \ m{l}^h_{c,c'} \leftarrow ext{ShiftChannel}_{c'}(m{l}^h_c)$ 7 8 # Apply window and causal mask  $m{m}^h_{(c,t),(c',t')} \leftarrow ext{CAUSALMASK}(m{g}^h_{(c,t),(c',t')} + m{s}^h_{t,t'} + m{l}^h_{c,c'})$  $m{n}^h_{(c,t),(c',t')} \leftarrow ext{WINDOWMASK}_{10}(m{m}^h_{(c,t),(c',t')})$ 9 10 # Apply structured dropout  $egin{aligned} & egin{aligned} & egi$ 11 12 $oldsymbol{o}_{c,t}^h \leftarrow \sum_{c',t'} oldsymbol{a}_{(c,t),(c',t')}^h \cdot oldsymbol{v}_{(c,t),(c',t')}^h$ 13 return  $o_c^h$ 14

#### A.1 Further motivation of MVPA

Single-channel data can be treated equivalently to sentences, by dividing the signal into 1D patches, which form the tokens. This modality has attracted considerable interest [29, 44], frequently for speech recognition tasks that are once again related to the natural language domain.

There is no straightforward extension of vanilla attention to the 2D case. The Vision Transformer [15] processes images by extending the notion of the patches to the 2D case. It carves images into a collection of patches, which it then flattens into a 1D sequence. Each patch has 2D coordinates (i, j) which get flattened by an arbitrary function  $f : \mathbb{N} \times \mathbb{N} \to \mathbb{N}$  into a 1D index (l). This is a simple way to recover the 1D case, but it has several drawbacks. First, by flattening the patches we lose any notion of spatial structure, as nearby patches in space are no longer necessarily close in the sequence. Any information about the structure of the patches is lost. However, if the size of the images, the number of patches, and the flattening direction are kept constant, then the Transformer might autonomously learn it. If it learns the structure, then it cannot be exposed to different images as it would completely misinterpret them; if it does not learn the structure, then it is missing critical information. This leads to an inflexible model which cannot easily generalize to different images. One possible solution is to choose a bijective f, such as the Cantor pairing function, to have

a one-to-one correspondence between the position of the patch in the image and in the sequence. This solution is, however, quite unintuitive. Second, the Vision Transformer does not distinguish between the two dimensions of height and width, i.e., it does not distinguish between up, down, left, and right. For images this limitation is not too impactful, as most of the information is conveyed in the closeness of two patches and not their relative position in any dimension.

The patching schema of the Vision Transformer is unsuitable to multi-variate time-series, as the two dimensions of time and channels require delicate handling. Transformers for time-series are a well-known problem in the field [17]. Channel-independent approaches [3] reuse vanilla self-attention and discard all information content in the time dimension, while channel-mixing promises to preserve it by either fusing the channels [45] or processing them sequentially [16]. The second family of solutions is more promising in addressing the issue but is still limited either with respect to computational expense or expressiveness.

EEG signals are multi-variate recordings of the brain. Transformer-based approaches to EEG are sparse [56, 57], due to the often unmanageable complexity of the data. In iEEG recordings, the subjects are implanted with electrodes directly in multiple areas of the brain for the purpose of clinical diagnosis. There is no standardized location, or even number of electrodes, for intracranial implants. This makes iEEG an extremely heterogeneous data modality, intractable for conventional attention approaches. The channels present a fundamental source of information, as electric fields spread in different areas of the brain on different time-scales and with different intensities depending on the strength of the connection between the areas. Moreover, the relationship between brain regions is not always proportional to their spatial closeness, as distant areas might be more strongly connected than close ones. There is a tremendously intricate interplay between space and time, which the Transformer must exploit.

#### A.2 Efficient computation of MVPA

Vanilla attention is quadratic in the number of input elements, and this often represents a significant computational roadblock [58]. The input becomes intractable as the number of channels increases, especially for multi-variate time-series. At the same time, more channels imply more sources of information, and we cannot simply discard them.

MVPA is also quadratic, but we employ a number of techniques to significantly reduce the computational complexity and make the processing of very large signals feasible. Letting T be the number of time segments and C be the number of channels, the context length of the Transformer becomes  $T \times C$  and number of terms necessary to compute for vanilla attention  $O(T^2 \times C^2)$ . Given a reasonable estimation of 100 segments and 50 channels the context length would be 5,000, until recently intractable even for language models.

By dividing MVPA into three components we gain considerable advantages (see Table S1 for the complexity of each term). Computing the time- and channel-based terms efficiently requires two main techniques. First, we recognize that it is not necessary to compute the full square matrix, which would be quadratic in the context length (i.e., both time and space). By design, all elements of the time-based attention are the same for each channel, and all elements of the channel-based attention are the same for each time point. Hence, complexity is quadratic in one dimension and constant in the other. We then simply repeat the elements along the right dimension at no additional cost. Second, we employ the shifting operation described in Supplementary Section B of Transformer-XL [18] to compute all relative embeddings in one pass.

Time-based Channel-based Content-based (w/ window) Content-based (vanilla)

Complexity  $O(T^2 \times C) = O(T \times C^2)$   $O(L^2 \times C^2)$   $O(T^2 \times C^2)$ 

Table S1: Complexity of each component of MVPA. T is the number of time segments in the signal, C is the number of channels, and L is the size of the local window. Contentbased attention without window has the same complexity as vanilla attention.

Content attention, though stripped of positional encoding, remains the most expensive component. To further reduce computational cost, with little impact to performance, we make use of a local attention window [19] which focuses on the most recent L time points discarding ones that have little information content. Since time-based attention is not limited, the lookup window still spans the entire context (though it is affected, see Figure S17b). Thus, for  $L \ll T$ , the total complexity of MVPA is  $O(T^2 \times C + T \times C^2)$ , quadratic in each dimension but subquadratic in the context length. Combining all techniques, MVPA pushes the effective total context length to over 10,000.

Given the three components are independent of each other, it is possible to exclude any one and reduce computations even more. As an additional cost-saving measure, we use grouped query attention [20] to reduce the number of heads without loss of performance. In summary, MVPA correctly treats time and space as unrelated dimensions, forcing the model to consider them separately, all with little computational overhead.

#### A.3 Triton implementation of MVPA

MVPFormer's training effectiveness is heavily affected by batch size, as its training routine draws the negative samples from the batch. The bigger the batch size, the more variety in the negative samples and the better the model generalizes. Given the large context size of MVPFormer, up to 10k, a pure Python implementation of scaled dot product attention would consume too much VRAM to be useful. FlashAttention [21] and FlashAttention-2 [22] provide the blueprint to solve this problem, though they only apply to vanilla attention. Using tiling, FlashAttention makes VRAM consumption linear instead of quadratic in the context length, enabling training on much longer context.

We develop FlashMVPA using the same technique in the OpenAI Triton language, which gives lower-level access to CUDA primitives. While a CUDA implementation could likely deliver better raw performance, the choice of Triton is dictated by the much lower coding time, though Triton is less robust and more prone to unexpected behaviors at this point. The time-based and channel-based components of MVPA are computed using PyTorch's own matrix-multiply, but are then shifted (Transformer-XL trick) and added using Triton, while the content-based component is fully implemented in Triton. This is due to limitations in Triton. FlashMVPA reaches 20 TFlops on an A100.

#### Algorithm 2: Computation of FlashMVPA

**Input:**  $\boldsymbol{x}_{c,t} \in \mathbb{R}^{n_{\text{embed}}}$  output token of Encoder;  $n_{\text{embed}} = 768$ **Output:**  $o_{c,t} \in \mathbb{R}^{n_{\text{embed}}}$  output attention **Data:** t time encoding; c channel encoding; u, y, w biases;  $h \in [1, \ldots, n_{\text{heads}}]$ ;  $h_{k,v} \in [1,\ldots,n_{\text{gqa}}]$ 1 def FlashMVPAttention $(x_{c,t})$ : # Compute query separately from key and value due to GQA  $q_{c,t}^{h} \leftarrow \text{LINEARNOBIAS}(\boldsymbol{x}_{c,t})$  $\boldsymbol{k}_{c,t}^{h_{k,v}}, \boldsymbol{v}_{c,t}^{h_{k,v}} \leftarrow \text{LINEARNOBIAS}(\boldsymbol{x}_{c,t})$ # Need to compute time and channel components outside Triton  $\mathbf{2}$ 3  $\boldsymbol{s}^h_t \leftarrow (\boldsymbol{q}^h_{c,t} + v^{\hat{h}_{k,v}})^T \boldsymbol{t}^{h_{k,v}}$ 4  $\boldsymbol{l}_{c}^{h} \leftarrow (\boldsymbol{q}_{c,t}^{h} + w^{h_{k,v}})^{T} \boldsymbol{c}^{h_{k,v}}$  $\mathbf{5}$ # Triton MVPA combines all computations into one kernel  $\boldsymbol{o}_{c,t}^h \leftarrow \text{TritonMVPA}(\boldsymbol{q}_{c,t}^h, \boldsymbol{s}_t^h, \boldsymbol{l}_t^h, \boldsymbol{v}, \boldsymbol{u}, \boldsymbol{y}, \boldsymbol{w})$ 6 return  $o_{c,t}^h$ 

#### A.4 Relative shifting

By design, MVPA requires the computation of relative time and channel encodings, which can notably slow down the overall operation. While this does not affect vanilla attention, other relative attentions provide us with an elegant solution to this problem. In particular, the shifting operation from Transformer-XL provides us with an efficient alternative to recomputing the time- and channel-based attention components. To keep notation simple, let  $\boldsymbol{q}_t = \boldsymbol{x}_{:t}^T W_q^T$ ,  $\boldsymbol{p}_c = \boldsymbol{x}_c^T W_q^T$ ,  $\boldsymbol{T}_i = W_{k_t} \mathcal{T}_{T-1-i}$ , and  $\boldsymbol{C}_i = W_{k_c} \mathcal{C}_{C-1-i}$ . The shift in time can be performed as in the original implementation

$$\begin{pmatrix} q_0 T_0 & q_0 T_1 & \dots & \dots & q_0 T_{T-1} \\ q_1 T_0 & q_1 T_1 & \dots & \dots & q_1 T_{T-1} \\ \vdots & \vdots & \vdots & \vdots & \vdots \\ q_{T-1} T_0 & q_{T-1} T_1 & \dots & \dots & q_{T-1} T_{T-1} \end{pmatrix} \xrightarrow{\text{SHIFTTIME}} \begin{pmatrix} q_0 T_{T-1} & 0 & \dots & \dots & 0 \\ q_1 T_{T-2} & q_1 T_{T-1} & 0 & \dots & 0 \\ \vdots & \vdots & \vdots & \vdots & \vdots \\ q_{T-1} T_0 & \dots & \dots & q_T T_{T-1} \end{pmatrix}$$
(5)

The right triangular matrix is zeroed out as a requisite of autoregressive training, i.e., we cannot attend to keys in the future. The entire time shifting operation can be performed efficiently and quickly using tensor manipulation in PyTorch.

Thus, the time attention component does not require recomputation for each time position, i.e. each row in the matrix of the time component.

The shift in channels is more involved

$$\begin{pmatrix} p_{0}C_{0} & p_{0}C_{1} & \dots & \dots & p_{0}C_{C-1} \\ p_{1}C_{0} & p_{1}C_{1} & \dots & \dots & p_{1}C_{C-1} \\ \vdots & \vdots & \vdots & \vdots & \vdots \\ p_{C-1}C_{0} & p_{C-1}C_{1} & \dots & \dots & p_{C-1}C_{C-1} \end{pmatrix} \xrightarrow{\text{SHIFTCHANNEL}} \begin{pmatrix} p_{0}C_{C-1} & p_{0}C_{C-2} & \dots & \dots & p_{0}C_{0} \\ p_{1}C_{C-2} & p_{1}C_{C-1} & p_{1}C_{C-2} & \dots & p_{1}C_{1} \\ \vdots & \vdots & \vdots & \vdots & \vdots \\ p_{C-1}C_{0} & \dots & \dots & \dots & p_{C}C_{C-1} \end{pmatrix}$$

$$(6)$$

Here, no element is zeroed out, as all channels can attend to all other channels. The channel shifting operation does not (to our knowledge) enjoy an implementation which as efficient as the time shifting one in PyTorch, but requires relatively complex index manipulation which cannot be streamlined.

As before, thanks to this shifting operation the channel attention component does not require recomputation for each channel position.

We provide a Triton implementation for both operations which is much more efficient and must be preferred when training a model.

#### A.5 Structured attention dropout

Dropout is a common technique to improve the generalization performance of neural networks. In Transformers, it is often applied inside the attention block to randomly zero-out some query-key attentions, to avoid over-reliance of the model on specific connections.



Figure S1: **Structured dropout.** (a) Our structured dropout blanks entire channels and time steps, to reduce the number of correlated segments. The dropout rate is computed to maintain the same number of dropped out segments as conventional dropout. (b) Conventional dropout blanks segments randomly. This is less effective with time-series because adjacent segments in time or space contain much of the same information.

Dropout usually applies to all elements with equal probability and creates uniform holes in the attention matrix. This is not efficient in the case of multi-variate time-series, as for each hole the neighboring segments are likely to carry very similar information, reducing dropout's effectiveness. We introduce a structured dropout technique which blanks entire channels and time points instead of individual segments. This technique is in principle much more effective by removing all segments which are more likely to be strongly correlated. We keep the same parameters as in conventional Dropout and compute the channel-specific and time-specific dropout rates as

$$t_{\rm drop} = c_{\rm drop} = 1 - \sqrt{1 - r_{\rm drop}} \tag{7}$$

This ensures that approximately the same overall number of elements are zeroed (see Figure S1).

For the specific dropout rates and the location of the structured dropout layers refer to the description of the architecture in Appendix B.

### **B** MVPFormer architecture

MVPFormer is part of a family of predictive deep learning models with 74 million (MVPFormer-S) to 1.2 billion (MVPFormer-M, or simply MVPFormer for brevity) parameters based on the Transformer [14] architecture, capable of generating iEEG signals. See Tables S2 and S3 for a breakdown of the models sizes and hyperparameters.

Table S2: **Breakdown of the parameters of MVPFormer-S.** The dimensions are indicated for each of the components of MVPFormer-S.

| Transformer                         | Encoder                             | Signal                                       |
|-------------------------------------|-------------------------------------|--|
| $n_{\text{layers}} \leftarrow 12$   | $n_{\text{input}} \leftarrow 2560$  | $w_{\text{length}} \leftarrow 500  \text{s}$ |
| $n_{\text{heads}} \leftarrow 12$    | $n_{\mathrm{embed}} \leftarrow 768$ | $n_{\text{segments}} \leftarrow 100$         |
| $n_{\rm gqa} \leftarrow 4$          |                                     | $w_{\text{segment}} \leftarrow 5\text{s}$    |
| $n_{\mathrm{embed}} \leftarrow 768$ |                                     | $n_{\text{negatives}} \leftarrow 30$         |
| $n_{\text{inner}} \leftarrow 1728$  |                                     | $n_{\text{local}} \leftarrow 10$             |
| $r_{\rm drop} \leftarrow 0.1$       |                                     |  |

Table S3: **Breakdown of the parameters of MVPFormer-M.** The dimensions are indicated for each of the components of MVPFormer-M.

| Transformer                        | Encoder                              | Signal                                       |
|------------------------------------|--------------------------------------|--|
| $n_{\text{layers}} \leftarrow 24$  | $n_{\text{input}} \leftarrow 2560$   | $w_{\text{length}} \leftarrow 500  \text{s}$ |
| $n_{\text{heads}} \leftarrow 16$   | $n_{\mathrm{embed}} \leftarrow 1024$ | $n_{\text{segments}} \leftarrow 100$         |
| $n_{\rm gqa} \gets 8$              |                                      | $w_{\text{segment}} \leftarrow 5\text{s}$    |
| $n_{\text{embed}} \leftarrow 2048$ |                                      | $n_{\text{negatives}} \leftarrow 30$         |
| $n_{\text{inner}} \leftarrow 5362$ |                                      | $n_{\text{local}} \leftarrow 10$             |
| $r_{\rm drop} \leftarrow 0.1$      |                                      |  |

#### **B.1** Inference

The full end-to-end inference procedure is reported in Algorithm 3.

Algorithm 3: Full inference with MVPFormer Input:  $x \in \mathbb{R}^{C \times T}$  raw input; C, T number of channels and length resp.;  $n_{\text{layers}} = 12$ 

#### B.2 Encoder

The Encoder block is detailed in Section 3.1. The algorithmic overview is presented in Algorithm 4.

Algorithm 4: Encoder block of MVPFormer

**Input:**  $\boldsymbol{x}_{c,t} \in \mathbb{R}^{n_{\text{input}}}$  raw input segment;  $n_{\text{input}} = 2560; c \in [1, \dots, C]; t \in [1, \dots, T]$ **Output:**  $o_{c,t} \in \mathbb{R}^{n_{\text{embed}}}$  output token;  $n_{\text{embed}} = 768$ **Data:** l = 8 maximum decomposition level given  $n_{input}$ 1 def Encoder( $x_{c,t}$ ):  $d_{c,t} \leftarrow \text{DISCRETEWAVELETDECOMPOSITION}_{db4}(\boldsymbol{x}_{c,t}, l)$ 2  $\boldsymbol{z}_{c,t} \leftarrow \text{RMSNORM}(\boldsymbol{d}_{c,t})$ 3  $o_{c,t} \leftarrow \text{LINEAR}(\boldsymbol{z}_{c,t})$ 4 return  $o_{c.t}$  $\mathbf{5}$ 

#### **B.3** Decoder

The collection of vectors resulting from the Encoder block is flattened into a 1D sequence to provide a unified input interface to the Transformer decoder blocks, consistent with conventional Transformers. All the encoded segments corresponding to a window form the input to the Transformer module, which computes the MVPA among all the segments. The segments are sequentially processed by multiple Transformer layers, composed of attention and MLP blocks in a deep network configuration. The attention blocks are masked to guarantee that MVPFormer only has access to past segments to generate the target. The model produces one output embedding for each input segment. The algorithmic overview is presented in Algorithm 5, while the MLP block in Algorithm 6.

Algorithm 5: Decoder block of MVPFormer

**Input:**  $o_{c,t} \in \mathbb{R}^{n_{\text{embed}}}$  input tokens;  $n_{\text{embed}} = 768$ Output:  $o_{c,t} \in \mathbb{R}^{n_{ ext{embed}}}$ 1 def Decoder $(o_{c,t})$ :  $z_{c,t} \leftarrow \text{RMSNORM}(o_{c,t})$ 2 # Compute attention  $\boldsymbol{a}_{c,t} \leftarrow \text{MVPATTENTION}(\boldsymbol{z}_{c,t})$ 3  $d_{c,t} \leftarrow \text{Dropout}(\text{LinearNoBias}((a_{c,t})))$ 4 # Compute feedforward residuals in parallel with attention [59]  $\boldsymbol{s}_{c,t} \leftarrow \mathrm{MLP}(\boldsymbol{z}_{c,t})$ 5 # Sum residuals and attention  $oldsymbol{o}_{c,t} \leftarrow oldsymbol{o}_{c,t} + oldsymbol{d}_{c,t} + oldsymbol{s}_{c,t}$ 6

7 return  $o_{c,t}$ 

### Algorithm 6: MLP block of MVPFormer

**Input:**  $\boldsymbol{z}_{c,t} \in \mathbb{R}^{n_{embed}}$  normalised Decoder output;  $n_{embed} = 768$ Data:  $\boldsymbol{u}_{c,t}, \boldsymbol{g}_{c,t} \in \mathbb{R}^{n_{\text{inner}}}; n_{\text{inner}} = 1728$ Output:  $s_{c,t} \in \mathbb{R}^{n_{\text{embed}}}$ 1 def MLP $(z_{c,t})$ :  $\boldsymbol{u}_{c,t} \leftarrow \text{LINEARNOBIAS}(\boldsymbol{z}_{c,t})$  $\boldsymbol{g}_{c,t} \leftarrow \mathrm{SILU}(\mathrm{LINEARNOBIAS}(\boldsymbol{z}_{c,t}))$ 3  $s_{c,t} \leftarrow \text{LINEARNOBIAS}(u_{c,t} + g_{c,t})$  $\mathbf{4}$ 

5 return  $s_{c,t}$ 

2

### C Details on training

#### C.1 Generative pre-training

MVPFormer is used to generate neuronal activity while in the base prediction task. During training, the target for each output is the successive input segment in time, not in space. We sample random input segments from the rest of the batched windows to create the confounding targets  $Z = \{z_1, ..., z_n\}$ . These segments still represent actual iEEG signals, so they are plausible, but they are expected to be very different from the true target.

This scheme strikes the correct balance between too much similarity and too little. The objective of MVPFormer is to generate future iEEG signals, so we choose a contrastive loss to increase the cosine similarity of its output with the true target, while decreasing it with the confounding targets. As training progresses, MVPFormer starts to produce outputs that look like encoded segments, i.e., its inputs. MVPFormer becomes more and more capable of choosing the right target and thus is able to predict the future signal.

**Extraction of positive and negative examples** Out of the entire dataset, B windows are chosen at random to form a batch. Each window  $W_{i \in [1..B]}$  has an arbitrary sample rate and  $C_i$  channels. First, the sampling rate is normalized to 512 Hz, then the windows are divided into T non-overlapping segments per-channel, leaving us with  $C_i \times T$  segments per window. Each segment is passed in parallel through the encoder. For the sake of simplicity, suppose one window  $W^*$  (with  $C^*$  channels) is selected at random as the positive window, and all the others as the confounding windows. The embeddings of  $W^*$  form the input context E with length  $C^* \times T$ .

For each segment, n embeddings are selected at random from the confounding windows to form the negative samples Z. Each  $Z_{c,t}$  has n elements, thus Z has size  $C^* \times T \times n$ . Z is excluded from backpropagation.

MVPFormer processes the entire E at once and produces an output O also of size  $C^* \times T$ . We then compute the losses and iteratively optimize to train the model.

**Contrastive loss** We train MVPFormer using a contrastive loss [60] and an auxiliary loss. To compute the contrastive loss, we rely on having other windows in the batch, so a larger batch size leads in general to a more stable training and better generalization performance. Let  $e^i$ ,  $i \in [1, \ldots, B]$  be the outputs of the signal Encoder and  $o^i$ ,  $i \in [1, \ldots, B]$  the outputs of the Decoder stack, for B the batch size. For each  $i^*$ , we select at random  $n_{\text{negatives}}$  elements from  $e^i$ ,  $i \neq i^*$  to act as our negative samples  $n^{i*}$ . Clearly, the bigger the batch the greater the entropy. We compute the contrastive loss for each  $i^*$  as follows:

$$\mathcal{L}_{c,t}^{i} = -\log \frac{\exp(\sin(\boldsymbol{o}_{c,t}, \boldsymbol{e}_{c(t+1)})/\tau)}{\sum_{\boldsymbol{z}_{k} \in Z} \exp(\sin(\boldsymbol{o}_{c,t}, \boldsymbol{z}_{k})/\tau)}$$
(8)

Finally summing over every i, c, t gives us the optimization target for the generative task.

The loss is invariant to the channel c, which encourages all the outputs to be the same regardless of channel. The temperature  $\tau$  is 0.1.

#### C.2 Generation of neuronal activity

The generation of brain signals during inference proceeds analogously as during training. However, we do not have access to the same source of entropy as in training since the batches are limited to one subject at a time. This limitation implies that the evaluation scores of MVPFormer must be more punishing than the training objective, since we cannot reliably estimate the accuracy with which MVPFormer chooses the right target. For this reason, we measure the cosine similarity directly in a three-way reference scheme. First, we consider the cosine similarity of the output with the true target. Second, we consider the similarity with the maximally correlated target. Third, we measure the cosine similarity with the highest form of entropy available, random segments in the batch that are still close by in time. This measurement ensures that the difference in similarity between the true and confounding targets remains significant.

### D SWEC iEEG dataset

The SWEC iEEG dataset is presented in Section 4. Here, Table S4 shows the full details of the dataset in a subject-by-subject breakdown. Finally, Figure S2 shows two annotated seizures in the dataset.

Table S4: **Per-subject details of our SWEC iEEG dataset.** Ch. is the number of electrodes,  $f_s$  is the sampling frequency in Hz, Rec. [h] is the length of the recording in hours, and Ev. is the number of seizures. The entire dataset contains 68 subjects, 9328 hours of recording and 704 ictal events.

| Subject | Ch. | $f_s$ [Hz] | Rec [h] | Ev. | Subject | Ch. | $f_s$ [Hz] | Rec [h] | Ev. | Subject | Ch. | $f_s$ [Hz] | Rec [h] | Ev. |
|---------|-----|------------|---------|-----|---------|-----|------------|---------|-----|---------|-----|------------|---------|-----|
| ID01    | 88  | 512        | 293.4   | 2   | ID24    | 32  | 1024       | 40.7    | 14  | ID47    | 32  | 1024       | 330.4   | 3   |
| ID02    | 66  | 512        | 235.2   | 2   | ID25    | 128 | 512        | 109.4   | 4   | ID48    | 57  | 1024       | 28.4    | 6   |
| ID03    | 64  | 512        | 158.4   | 4   | ID26    | 34  | 1024       | 87.6    | 1   | ID49    | 60  | 512        | 140.4   | 6   |
| ID04    | 32  | 1024       | 40.7    | 14  | ID27    | 32  | 1024       | 146     | 8   | ID50    | 64  | 1024       | 177.2   | 2   |
| ID05    | 128 | 512        | 109.4   | 4   | ID28    | 75  | 512        | 69      | 4   | ID51    | 89  | 512        | 161.5   | 1   |
| ID06    | 32  | 1024       | 146     | 8   | ID29    | 61  | 1024       | 143.8   | 70  | ID52    | 69  | 512        | 112.6   | 2   |
| ID07    | 75  | 512        | 69      | 4   | ID30    | 48  | 1024       | 40.9    | 27  | ID53    | 22  | 1024       | 134.9   | 1   |
| ID08    | 61  | 1024       | 143.8   | 70  | ID31    | 32  | 1024       | 42.4    | 17  | ID54    | 54  | 1024       | 202     | 3   |
| ID09    | 48  | 1024       | 40.9    | 27  | ID32    | 32  | 1024       | 212.2   | 2   | ID55    | 24  | 1024       | 152.1   | 2   |
| ID10    | 32  | 1024       | 42.4    | 17  | ID33    | 104 | 512        | 53.6    | 1   | ID56    | 62  | 1024       | 130.5   | 3   |
| ID11    | 32  | 1024       | 212.2   | 2   | ID34    | 56  | 1024       | 191.4   | 9   | ID57    | 40  | 1024       | 90.7    | 12  |
| ID12    | 56  | 1024       | 191.4   | 9   | ID35    | 64  | 1024       | 104     | 7   | ID58    | 92  | 512        | 138.2   | 7   |
| ID13    | 64  | 1024       | 104     | 7   | ID36    | 24  | 1024       | 161.4   | 60  | ID59    | 54  | 1024       | 107.3   | 15  |
| ID14    | 24  | 1024       | 161.4   | 60  | ID37    | 98  | 512        | 195.9   | 2   | ID60    | 74  | 512        | 50.7    | 8   |
| ID15    | 98  | 512        | 195.9   | 2   | ID38    | 34  | 1024       | 177.1   | 5   | ID61    | 76  | 512        | 89.6    | 6   |
| ID16    | 34  | 1024       | 177.1   | 5   | ID39    | 60  | 1024       | 129.6   | 2   | ID62    | 60  | 1024       | 235.1   | 7   |
| ID17    | 60  | 1024       | 129.6   | 2   | ID40    | 42  | 1024       | 205.1   | 5   | ID63    | 64  | 512        | 179.8   | 4   |
| ID18    | 42  | 1024       | 205.1   | 5   | ID41    | 33  | 1024       | 82.7    | 3   | ID64    | 56  | 1024       | 36.3    | 20  |
| ID19    | 29  | 1024       | 21.7    | 25  | ID42    | 63  | 1024       | 87.8    | 2   | ID65    | 49  | 1024       | 139.7   | 8   |
| ID20    | 88  | 512        | 293.4   | 2   | ID43    | 126 | 512        | 63.2    | 2   | ID66    | 39  | 1024       | 212.3   | 2   |
| ID21    | 66  | 512        | 235.2   | 2   | ID44    | 60  | 1024       | 150.3   | 2   | ID67    | 63  | 512        | 111.7   | 4   |
| ID22    | 64  | 512        | 158.4   | 4   | ID45    | 47  | 1024       | 157.3   | 1   | ID68    | 32  | 1024       | 167.8   | 3   |
| ID23    | 32  | 1024       | 42.4    | 33  | ID46    | 86  | 512        | 140.5   | 21  |         |     |            |         |     |

#### a Patient 17 with typical seizure presentation

#### The seizure can be easily distinguished from the surrounding signal

|  | ~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~   |
|--|--|
|  |  |
|  |  |
|  |  |
| Alther and the state of the sta |  |
|  | ·  |
| while we wanted with the sector many sector and the sec   | hanness and the second se   |
|  |  |
|  |  |
|  |  |
| www.andersteelenselenselenselenselenselenselensel  |  |
| his way to go and a fair of the second and the second a   |  |
|  |  |
|  |  |
| and the second   | -  |
|  | I manual and the second s   |
|  | fur a second sec |
|  |  |
|  | ~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~   |
|  |  |
| · · · · · · · · · · · · · · · · · · ·  |  |
|  |  |
| xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx   | When when the work of the second of the seco |
| where we have the bird of the physical and the design of the the term of the first of the term of term | m-man MMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMM   |
| the manufactories of the second s   | from the second  |
| i i i i i i i i i i i i i i i i i i i  | 1  |
|  |  |
|  |  |
| head the first of the advector in the section of th |  |

b Patient 14 with atypical seizure presentation

The seizure is not immediately recognisable

| ไม่ไปการการการที่สาวที่ไม่มาสารที่สาวที่สาวไปไม่ไม่ไม่ได้ได้ได้มีสารการการที่สารไปไม่ไม่มาสารได้การที่ได้มากรับ<br>ที่ไม่มีการการที่สาวที่ไม่มาสารการการที่มีผู้ในการที่ได้ได้มีสารการการที่สารได้ได้ได้การที่ได้การที่มีการการการ<br>ที่มีที่มีการการที่สาวที่สารการการการการที่มีการการก็สารการการการการการที่ได้ได้การที่ได้การการการการการการการ | Yhele Yhele all and and an ar an<br>Ar an ar a<br>Ar an ar   | 19   |
|--|--|--|
| ามรูกการสารสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสาร<br>การสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการ<br>การสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสารการสา           | in a second s<br>in a second s<br>in a second se | yn fel fan en fan ste fan ste<br>fan yn ste fan ste fan<br>fan yn ste fan ste fan<br>fan yn ste fan ste fan<br>fan yn ste fan ste fan<br>fan ste fan ste   |
|  |  |  |
|  | n an   | a de la companya de<br>el companya de la comp<br>per serva de la companya de la comp<br>per serva de la companya de la companya<br>Recompanya de la companya de la comp |

Figure S2: **iEEG activity of two patients with different ictal patterns.** (a) Patient 17 of the SWEC iEEG dataset presents typical ictal events. The seizure can be clearly distinguished even by a non-expert, and MVPFormer performs very well on this patient. The number of channels is reduced from the original recording to facilitate comparison with the more difficult presentation. (b) Patient 14 of the SWEC iEEG dataset does not have typical events. The neuronal activity during seizures for this patient cannot be clearly distinguished, and assessment by experts would diverge considerably. As expected, MVPFormer has a high level of disagreement on this patient. All the channels of the original recording are presented to exclude the chance of some channels carrying additional information.

### **E** Generation of results

#### E.1 Episodic seizure post-processing

For episodic evaluation we apply three post-processing steps to the model output:

- Merge events happening within 5 minutes of each other
- Remove events shorter than 20 seconds in length
- Remove events with less than 5 positive responses

Moreover, when the subject has multiple seizures in one minute we merge them into one.

#### E.2 Online seizure thresholding

In the clinical evaluation setup we apply a simple thresholding to decide whether to report a seizure or not. We set 3 positive seconds out of 10 to be the lower limit for detecting a seizure, to deter false positives; events shorter than 3 seconds are thus not reported, and an additional latency of 3 seconds is to be considered. We find this trade-off has limited drawbacks in practice, as there is often large disagreement even among neurologists about very short events.

### E.3 Kappa score estimation

To estimate the Kappa score, we choose 300 random segments per subject to compare their classification from MVPFormer and the labels. We perform multiple iterations to ensure no bias in this computation. Figure S3 indicates that our choice of 250 iterations is sufficient for stable results.



Figure S3: Mean absolute error in Cohen's kappa estimation. Our estimation scheme for Cohen's kappa converges after very few iterations. The error is computed per-subject as the absolute difference between the running averages at each two consecutive iterations; the running average is the average of all preceding steps. The average and standard deviation across all subjects is reported here. We compute up to 250 random iterations to ensure precise reporting.

#### E.4 Landis and Koch criteria

Table S5 reports the commonly used Landis and Koch criteria for qualitative inter-rater agreement estimation from Kappa scores.

| Kappa       | Agreement      |
|-------------|----------------|
| 0 - 0.20    | Slight         |
| 0.21 - 0.40 | Fair           |
| 0.41 - 0.60 | Moderate       |
| 0.61 - 0.80 | Substantial    |
| 0.81 - 1.00 | Almost perfect |

Table S5: Landis and Koch criteria. Landis and Koch criteria [40] for evaluating Cohen's kappa in the context of inter-rater agreement between human experts on seizure classification.

### F Prediction of iEEG signals

We evaluate the effects of the size of the model and the attention mechanism on the iEEG prediction task. Figure S4 shows that both MVPA and the vanilla attention are effective at predicting the next brain states. Scaling up the model size from MVPFormer-S to MVPFormer-M has the effect of shortening the tail of the true distribution, effectively increasing the concentration of the cosine similarity towards the maximum.



Figure S4: **Performance of MVPFormer-M, MVPFormer-S, and MV-Llama on the iEEG prediction task.** We report the average two-step, random, and true target cosine similarities for the three different models. All three are effective at predicting iEEG activity, while the larger model takes advantage of the increased embedding size by increasing the concentration of cosine similarities towards 1.

Complete details are available for the iEEG prediction performance of MVPFormer-M (see Appendix F.1, MVPFormer-S (see Appendix F.2), and MV-Llama (see Appendix F.3).

### F.1 Prediction of iEEG and ictal activity

MVPFormer is primarily a neuronal prediction model, trained to generate neuronal activity regardless of whether such activity is pathological or physiological. To understand the behavior of MVPFormer with anomalous brain states, we evaluate its performance in generating ictal neuronal activity. The precise relationship between ictal and interictal states is a point of ongoing discussion [61, 62], but many consider an approach to seizures as anomalies [63] the most appropriate. The SWEC iEEG dataset contains many ictal events, so we are able to evaluate the performance of MVPFormer in generating anomalous activity. In particular, in this dataset the ratio between non-ictal and ictal states is approximately 500:1.

Figure S5 shows that ictal states are not anomalous for MVPFormer. In particular, the prediction similarity of MVPFormer does not degrade when generating ictal activity. Moreover, the prediction similarity in the ictal state is neither significantly different from the average similarity nor from the non-ictal similarity. This indicates that MVPFormer's understanding of the mechanisms of generation of neuronal activity encompasses the pathological ictal state as well. Therefore, MVPFormer must model patterns found both in physiological and pathological brain states. Finally, MVPFormer incorporates a model of seizure generation as a by-product of its predictive task, which is particularly noteworthy.



Figure S5: Performance of MVPFormer on the prediction task. (a) A three-reference evaluation scheme is used to assess MVPFormer's performance. The true target is the immediate future, i.e. the next five seconds of iEEG signal. The two-step target is the signal twice removed in the future, i.e. the five seconds of iEEG signal coming after the true target. Finally, the random target is chosen from iEEG signals which are close by in time with the true target. The distribution of the average similarity across the entire recording is shown together with the similarity within three representative subjects (with maximum, median, and minimum average similarity). (b) The prediction similarity is computed again for all three targets, distinguishing between targets which lie within an ictal event, without, or at the boundary. There is no significant difference in the performance of MVPFormer in predicting ictal or non-ictal activity, indicating that MVPFormer can encompass anomalous brain states as well, together with the transitions between physiological and anomalous.

### F.2 Effects of the scale of the model

Figure S6 shows the full details on the performance of MVPFormer-S on the SWEC iEEG dataset in the iEEG prediction task.



Figure S6: **Performance of MVPFormer-S on the prediction task.** Performance of MVPFormer-S in generating neuronal activity of unseen test subjects. The prediction similarity is computed for all three targets, distinguishing between targets which lie within an ictal event or without. There is no significant difference in the performance of MVPFormer-S in predicting ictal or non-ictal activity.

#### F.3 Effects of the attention mechanism

Figure S7 shows that vanilla attention is also effective in predicting the development of iEEG signal.



Figure S7: **Performance of MV-Llama on the prediction task.** The prediction similarity is computed for all three targets, distinguishing between targets which lie within an ictal event or without. There is no significant difference in the performance of MV-Llama in predicting ictal or non-ictal activity. This indicates that vanilla attention with a proper positional encoding scheme can effectively generate neuronal activity. However, this does not translate to improved performance in the seizure classification task (see Table S8 vs. Table S7)

#### F.4 Per-subject cosine similarity

We provide a detailed per-subject breakdown of the maximum cosine similarity measure for MVPFormer. Figure S8 shows the per-patient global similarity. Figure S9 shows the per-patient similarity within an anomaly. Figure S10 shows the per-patient similarity at the boundary of an anomaly.



Figure S8: **Breakdown of total cosine similarity per-patient.** Maximum cosine similarity of MVPFormer's output with the true, random, and two-step targets over the entire SWEC iEEG dataset. The data is shown patient-by-patient.



Figure S9: **Breakdown of anomaly cosine similarity per-patient.** Maximum cosine similarity of MVPFormer's output with the true, random, and two-step targets while within an anomaly (seizure) in the SWEC iEEG dataset. The data is shown patient-by-patient.



Figure S10: **Breakdown of boundary cosine similarity per-patient.** Maximum cosine similarity of MVPFormer's output with the true, random, and two-step targets while in the onset and offset zones of seizures in the SWEC iEEG dataset. The data is shown patient-by-patient.

### G Additional results

#### G.1 Conventional evaluation

In addition to our clinically motivated evaluation (see Section 5.2), we assess all our models using conventional machine learning metrics for seizure detection: F1-score, sensitivity, and false positive rate. These metrics are commonly used in benchmarking seizure detection models [64, 65], and allow comparison with prior work. The full seizure detection results of MVPFormer are shown in Figure S11.



Figure S11: **Performance of MVPFormer on the classification task.** (a) The F1-score, sensitivity, and fp/h are reported. Raw results are computed without any post-processing of MVPFormer's output, while episodic results follow a common post-processing procedure which merge close ictal classifications. (b) Cohen's Kappa is used to measure the agreement between the artificial assistant and the human expert. The average kappa is 0.57, competitive with the values obtained between human experts. The distribution of kappa values clearly indicates that a minority of subjects are the source of most disagreement, consistent with the variability of inter-rater agreement among human experts.

We evaluate against two baselines: Brant-2 [13], a SOTA Transformer model for iEEG, and MV-Llama, an ablation of MVPFormer-S that uses standard attention instead of MVPA (see Appendix G.3). Brant-2 is fine-tuned with its published pre-trained weights and protocol. MV-Llama is trained identically to MVPFormer-S.

We report both raw and episodic metrics. Episodic metrics reflect clinically meaningful detections by grouping predictions into events [64]. The detailed results are provided in Tables S7 and S6.

The similarity between raw and episodic F1-scores suggests that MVPFormer naturally learns to detect seizure episodes of realistic length and frequency. On the 50-subject SWEC test set, the false positive rate is 0.17 fp/h (0.11 for MVPFormer-S), comparable to commercial

Table S6: Details of seizure detection results of MVPFormer with 18 subject pre-training. Kappa is the inter-rater agreement. The classification metrics report the raw and episodic metrics relevant for the seizure classification task. The similarity reports the breakdown of the cosine similarity in each of the considered scenarios.

|         |       | Raw Episodic True |          | True        | rue  |         | Random |           | Two-st  |       | ep        |         |       |           |
|---------|-------|-------------------|----------|-------------|------|---------|--------|-----------|---------|-------|-----------|---------|-------|-----------|
| Subject | Kappa | f1-score          | f1-score | sensitivity | fp/h | average | ictal  | non-ictal | average | ictal | non-ictal | average | ictal | non-ictal |
| ID19    | 0.01  | 0.01              | 0.14     | 0.44        | 3.55 | 0.30    | 0.27   | 0.25      | 0.14    | 0.15  | 0.10      | 0.16    | 0.14  | 0.10      |
| ID20    | 1     | 0.95              | 1        | 1           | 0    | 0.89    | 0.95   | 0.91      | 0.17    | 0.19  | 0.12      | 0.14    | 0.17  | 0.12      |
| ID21    | 0.99  | 0.96              | 1        | 1           | 0    | 0.92    | 0.92   | 0.93      | 0.16    | 0.15  | 0.16      | 0.13    | 0.22  | 0.18      |
| ID22    | 0.99  | 0.95              | 1        | 1           | 0    | 0.89    | 0.92   | 0.90      | 0.15    | 0.15  | 0.10      | 0.12    | 0.15  | 0.08      |
| ID23    | 0.1   | 0.14              | 0.11     | 0.06        | 0.07 | 0.88    | 0.95   | 0.94      | 0.22    | 0.22  | 0.20      | 0.20    | 0.21  | 0.20      |
| ID24    | 0.92  | 0.92              | 0.93     | 0.93        | 0.02 | 0.80    | 0.93   | 0.89      | 0.13    | 0.13  | 0.09      | 0.14    | 0.11  | 0.14      |
| ID25    | 0     | 0                 | 0        | 0           | 0    | 0.84    | 0.90   | 0.90      | 0.17    | 0.24  | 0.14      | 0.15    | 0.20  | 0.13      |
| ID26    | 0.01  | 0                 | 0        | 0           | 0.59 | 0.88    | 0.96   | 0.94      | 0.16    | 0.21  | 0.09      | 0.13    | 0.16  | 0.10      |
| ID27    | 1     | 0.95              | 1        | 1           | 0    | 0.88    | 0.95   | 0.95      | 0.20    | 0.23  | 0.20      | 0.17    | 0.23  | 0.21      |
| ID28    | 0.75  | 0.85              | 0.75     | 0.75        | 0.01 | 0.85    | 0.89   | 0.61      | 0.26    | 0.14  | 0.07      | 0.16    | 0.15  | 0.10      |
| ID29    | 0.29  | 0.33              | 0.19     | 0.13        | 0.05 | 0.86    | 0.92   | 0.92      | 0.14    | 0.14  | 0.11      | 0.11    | 0.12  | 0.13      |
| ID30    | 0.96  | 0.92              | 0.98     | 0.96        | 0    | 0.89    | 0.95   | 0.94      | 0.17    | 0.17  | 0.13      | 0.14    | 0.16  | 0.14      |
| ID31    | 0.98  | 0.96              | 1        | 1           | 0    | 0.88    | 0.95   | 0.95      | 0.18    | 0.21  | 0.13      | 0.15    | 0.22  | 0.17      |
| ID32    | 1     | 0.97              | 1        | 1           | 0    | 0.93    | 0.93   | 0.94      | 0.17    | 0.15  | 0.12      | 0.13    | 0.16  | 0.11      |
| ID33    | 0     | 0                 | 0        | 0           | 0.37 | 0.95    | 0.96   | 0.95      | 0.39    | 0.38  | 0.31      | 0.39    | 0.45  | 0.46      |
| ID34    | 0.98  | 0.97              | 1        | 1           | 0    | 0.89    | 0.94   | 0.92      | 0.14    | 0.14  | 0.09      | 0.11    | 0.14  | 0.12      |
| ID35    | 0.98  | 0.94              | 1        | 1           | 0    | 0.86    | 0.93   | 0.94      | 0.18    | 0.15  | 0.25      | 0.17    | 0.17  | 0.35      |
| ID36    | 0.58  | 0.73              | 0.54     | 0.37        | 0    | 0.91    | 0.95   | 0.94      | 0.19    | 0.20  | 0.13      | 0.16    | 0.18  | 0.14      |
| ID37    | 1     | 0.97              | 1        | 1           | 0    | 0.77    | 0.74   | 0.66      | 0.16    | 0.12  | 0.29      | 0.14    | 0.15  | 0.49      |
| ID38    | 0.98  | 0.98              | 1        | 1           | 0    | 0.82    | 0.85   | 0.77      | 0.15    | 0.13  | 0.13      | 0.13    | 0.14  | 0.09      |
| ID39    | 1     | 0.98              | 1        | 1           | 0    | 0.83    | 0.70   | 0.83      | 0.16    | 0.12  | 0.10      | 0.14    | 0.12  | 0.10      |
| ID40    | 1     | 0.99              | 1        | 1           | 0    | 0.87    | 0.88   | 0.90      | 0.15    | 0.13  | 0.09      | 0.13    | 0.13  | 0.11      |
| ID41    | 0.68  | 0.72              | 0.67     | 1           | 0.04 | 0.86    | 0.79   | 0.86      | 0.16    | 0.15  | 0.10      | 0.13    | 0.15  | 0.10      |
| ID42    | 0.99  | 0.89              | 1        | 1           | 0    | 0.82    | 0.80   | 0.79      | 0.15    | 0.14  | 0.07      | 0.13    | 0.15  | 0.07      |
| ID43    | 0.92  | 0.69              | 1        | 1           | 0    | 0.80    | 0.77   | 0.84      | 0.13    | 0.14  | 0.12      | 0.10    | 0.15  | 0.11      |
| ID44    | 0     | 0                 | 0        | 0           | 0    | 0.82    | 0.89   | 0.89      | 0.11    | 0.14  | 0.49      | 0.08    | 0.09  | 0.65      |
| ID45    | 0.02  | 0.02              | 0.04     | 1           | 0.34 | 0.82    | 0.87   | 0.77      | 0.14    | 0.16  | 0.09      | 0.12    | 0.14  | 0.16      |
| ID46    | 0.5   | 0.5               | 0.56     | 0.43        | 0.01 | 0.85    | 0.87   | 0.61      | 0.18    | 0.14  | 0.28      | 0.16    | 0.11  | 0.51      |
| ID47    | 0.72  | 0.54              | 0.46     | 1           | 0.02 | 0.84    | 0.81   | 0.83      | 0.16    | 0.15  | 0.12      | 0.14    | 0.16  | 0.09      |
| ID48    | 0.78  | 0.65              | 0.91     | 0.83        | 0    | 0.80    | 0.83   | 0.79      | 0.12    | 0.12  | 0.08      | 0.10    | 0.11  | 0.09      |
| ID49    | 0.52  | 0.54              | 0.6      | 1           | 0.06 | 0.89    | 0.79   | 0.86      | 0.23    | 0.14  | 0.35      | 0.23    | 0.17  | 0.53      |
| ID50    | 0.67  | 0.57              | 0.44     | 1           | 0.03 | 0.86    | 0.90   | 0.84      | 0.16    | 0.24  | 0.11      | 0.14    | 0.21  | 0.10      |
| ID51    | 0.74  | 0.41              | 0.33     | 1           | 0.02 | 0.87    | 0.78   | 0.89      | 0.26    | 0.14  | 0.28      | 0.28    | 0.24  | 0.39      |
| ID52    | 0.44  | 0.48              | 0.36     | 1           | 0.06 | 0.82    | 0.84   | 0.79      | 0.21    | 0.18  | 0.11      | 0.19    | 0.16  | 0.14      |
| ID53    | 0     | 0.01              | 0        | 1           | 3.56 | 0.82    | 0.77   | 0.79      | 0.17    | 0.17  | 0.09      | 0.14    | 0.15  | 0.11      |
| ID54    | 0.01  | 0.01              | 0.02     | 0.67        | 1.21 | 0.82    | 0.71   | 0.83      | 0.16    | 0.14  | 0.12      | 0.13    | 0.14  | 0.09      |
| ID55    | 0.63  | 0.49              | 0.44     | 1           | 0.03 | 0.84    | 0.87   | 0.80      | 0.18    | 0.18  | 0.11      | 0.16    | 0.19  | 0.14      |
| ID56    | 0.62  | 0.52              | 0.57     | 0.67        | 0.02 | 0.87    | 0.81   | 0.79      | 0.19    | 0.18  | 0.12      | 0.17    | 0.20  | 0.18      |
| ID57    | 0     | 0                 | 0        | 0           | 0    | 0.87    | 0.87   | 0.81      | 0.16    | 0.21  | 0.14      | 0.13    | 0.17  | 0.17      |
| ID58    | 0.18  | 0.08              | 0.12     | 0.14        | 0.06 | 0.86    | 0.84   | 0.77      | 0.21    | 0.22  | 0.33      | 0.20    | 0.24  | 0.30      |
| ID59    | 0     | 0                 | 0        | 0           | 0    | 0.82    | 0.85   | 0.80      | 0.17    | 0.14  | 0.09      | 0.16    | 0.12  | 0.11      |
| ID60    | 0.16  | 0.08              | 0.55     | 0.38        | 0    | 0.87    | 0.88   | 0.93      | 0.15    | 0.14  | 0.56      | 0.12    | 0.11  | 0.37      |
| ID61    | 0     | 0                 | 0        | 0           | 0    | 0.89    | 0.87   | 0.89      | 0.18    | 0.16  | 0.11      | 0.15    | 0.13  | 0.11      |
| ID62    | 0.7   | 0.41              | 0.83     | 0.71        | 0    | 0.83    | 0.86   | 0.92      | 0.14    | 0.13  | 0.52      | 0.13    | 0.12  | 0.72      |
| ID63    | 0.41  | 0.43              | 0.57     | 1           | 0.03 | 0.82    | 0.86   | 0.86      | 0.15    | 0.16  | 0.09      | 0.13    | 0.15  | 0.06      |
| ID64    | 0.52  | 0.56              | 0.52     | 0.35        | 0    | 0.86    | 0.87   | 0.85      | 0.15    | 0.15  | 0.16      | 0.12    | 0.14  | 0.10      |
| ID65    | 0.72  | 0.54              | 0.75     | 0.75        | 0.01 | 0.83    | 0.87   | 0.85      | 0.13    | 0.15  | 0.09      | 0.10    | 0.13  | 0.08      |
| ID66    | 0.85  | 0.42              | 0.8      | 1           | 0    | 0.75    | 0.80   | 0.76      | 0.12    | 0.13  | 0.06      | 0.10    | 0.11  | 0.08      |
| ID67    | 0.04  | 0.05              | 0.03     | 0.5         | 1.06 | 0.83    | 0.82   | 0.81      | 0.20    | 0.15  | 0.07      | 0.19    | 0.12  | 0.10      |
| ID68    | 0.66  | 0.00              | 0.00     | 0.00        | 0.01 | 0.84    | 0.80   | 0.82      | 0.17    | 0.22  | 0.15      | 0.14    | 0.19  | 0.16      |

EEG devices used in clinical practice [66, 67]. As expected in a real-world dataset, false positive rates vary considerably across subjects, with 78% having fewer than 0.05 fp/h.

These results confirm that MVPFormer performs competitively on conventional seizure detection benchmarks, while also offering robust generalization to clinically realistic evaluation settings.

#### G.2 Effects of the scale of the model

The performance improvements of LLMs as a function of their model sizes have also been widely reported [48, 49]. According to Chinchilla's scaling law the training dataset is already not large enough to fully train MVPFormer-S (75M parameters), so we investigate whether a larger model (MVPFormer-M, 1.2B parameters) can provide any improvement in performance.

Figure S12 shows the seizure detection performance of MVPFormer-S on the SWEC iEEG dataset (see Table S7). As noted in the main results, MVPFormer-M marginally improves seizure detection results over MVPFormer-S. In particular, it reaches higher F1-score but higher fp/h rate as well, with small net improvement. Therefore, we have shown that the amount of iEEG data currently available is not sufficient to fully take advantage of the increase in model size of Transformers. We hope that making the SWEC iEEG dataset publicly available will increase overall availability and unlock further model scaling potential.



Figure S12: **Performance of MVPFormer-S on the classification task.** (a) Seizure detection results of MVPFormer-S on unseen subjects: the F1-score, sensitivity, and fp/h are reported. Raw results are computed without any post-processing of MVPFormer's output, while episodic results follow a common post-processing procedure which merge close ictal classifications. (b) Cohen's Kappa is used to measure the agreement between MVPFormer-S and the human expert. The average kappa is 0.54, competitive with the values obtained between human experts. The distribution of kappa values clearly indicates that a minority of subjects are the source of most disagreement, consistent with the variability of inter-rater agreement among human experts.

|              | Classification metrics |          |          |             |      | Similarity |       |           |         |       |           |         |       |           |  |
|--------------|------------------------|----------|----------|-------------|------|------------|-------|-----------|---------|-------|-----------|---------|-------|-----------|--|
|              | Kaw                    |          | Episodic |             |      | True       |       |           | Rando   | m     |           | Two-ste | ep    |           |  |
| Subject      | Kappa                  | f1-score | f1-score | sensitivity | fp/h | average    | ictal | non-ictal | average | ictal | non-ictal | average | ictal | non-ictal |  |
| ID19         | -0.05                  | 0        | 0        | 0           | 0.97 | 0.28       | 0.27  | 0.27      | 0.21    | 0.2   | 0.25      | 0.21    | 0.2   | 0.23      |  |
| ID20         | 0.98                   | 0.87     | 0.8      | 1           | 0    | 0.95       | 0.96  | 0.95      | 0.18    | 0.21  | 0.12      | 0.15    | 0.17  | 0.12      |  |
| ID21         | 0.92                   | 0.86     | 0.8      | 1           | 0    | 0.94       | 0.93  | 0.95      | 0.16    | 0.18  | 0.18      | 0.13    | 0.24  | 0.2       |  |
| ID22         | 0.98                   | 0.89     | 1        | 1           | 0    | 0.94       | 0.89  | 0.95      | 0.15    | 0.16  | 0.09      | 0.12    | 0.17  | 0.07      |  |
| ID23         | 0.26                   | 0.17     | 0.14     | 0.09        | 0.16 | 0.95       | 0.95  | 0.95      | 0.25    | 0.28  | 0.23      | 0.22    | 0.24  | 0.24      |  |
| ID24         | 0.91                   | 0.9      | 0.93     | 0.93        | 0.02 | 0.93       | 0.94  | 0.93      | 0.19    | 0.19  | 0.12      | 0.16    | 0.17  | 0.17      |  |
| ID25         | 0.35                   | 0        | 0        | 0           | 0.03 | 0.94       | 0.96  | 0.96      | 0.17    | 0.19  | 0.15      | 0.14    | 0.18  | 0.15      |  |
| ID26         | 0.08                   | 0.03     | 0.11     | 1           | 0.18 | 0.95       | 0.97  | 0.96      | 0.17    | 0.22  | 0.12      | 0.14    | 0.18  | 0.11      |  |
| ID27         | 0.99                   | 0.93     | 1        | 1           | 0    | 0.95       | 0.96  | 0.96      | 0.22    | 0.24  | 0.21      | 0.19    | 0.25  | 0.24      |  |
| ID28         | 0.71                   | 0.68     | 0.75     | 0.75        | 0.01 | 0.95       | 0.96  | 0.95      | 0.17    | 0.19  | 0.09      | 0.13    | 0.17  | 0.12      |  |
| ID29         | 0.25                   | 0.29     | 0.2      | 0.13        | 0.03 | 0.95       | 0.95  | 0.95      | 0.16    | 0.15  | 0.14      | 0.12    | 0.13  | 0.13      |  |
| ID30         | 0.92                   | 0.88     | 0.93     | 0.93        | 0.05 | 0.95       | 0.96  | 0.96      | 0.19    | 0.19  | 0.15      | 0.16    | 0.18  | 0.16      |  |
| ID31         | 0.98                   | 0.95     | 1        | 1           | 0    | 0.95       | 0.96  | 0.96      | 0.19    | 0.24  | 0.11      | 0.16    | 0.24  | 0.18      |  |
| ID32         | 1                      | 0.92     | 1        | 1           | ŏ    | 0.95       | 0.91  | 0.96      | 0.18    | 0.13  | 0.16      | 0.15    | 0.18  | 0.12      |  |
| ID33         | 0                      | 0        | 0        | 0           | 0.26 | 0.93       | 0.95  | 0.93      | 0.23    | 0.28  | 0.19      | 0.21    | 0.26  | 0.21      |  |
| ID34         | 0.89                   | 0.91     | ñ 9      | 1           | 0.01 | 0.95       | 0.95  | 0.95      | 0.15    | 0.15  | 0.15      | 0.21    | 0.15  | 0.13      |  |
| ID35         | 0.85                   | 0.87     | 0.82     | 1           | 0.01 | 0.92       | 0.94  | 0.92      | 0.17    | 0.10  | 0.19      | 0.15    | 0.17  | 0.22      |  |
| ID36         | 0.00                   | 0.62     | 0.51     | 0.37        | 0.03 | 0.92       | 0.04  | 0.95      | 0.11    | 0.10  | 0.14      | 0.10    | 0.21  | 0.17      |  |
| 1D27         | 0.74                   | 0.62     | 0.51     | 1           | 0.00 | 0.00       | 0.00  | 0.01      | 0.22    | 0.15  | 0.22      | 0.13    | 0.17  | 0.60      |  |
| 1D37         | 0.74                   | 0.05     | 0.01     | 1           | 0.02 | 0.92       | 0.91  | 0.91      | 0.10    | 0.15  | 0.33      | 0.14    | 0.17  | 0.09      |  |
| 1D30         | 0.94                   | 0.95     | 0.91     | 1           | 0.01 | 0.95       | 0.95  | 0.91      | 0.14    | 0.15  | 0.11      | 0.11    | 0.10  | 0.08      |  |
| ID 39        | 0.92                   | 0.69     | 0.0      | 1           | 0.01 | 0.95       | 0.9   | 0.92      | 0.10    | 0.15  | 0.11      | 0.14    | 0.15  | 0.11      |  |
| ID40<br>ID41 | 0.95                   | 0.95     | 0.91     | 1           | 0.04 | 0.94       | 0.95  | 0.95      | 0.15    | 0.10  | 0.14      | 0.13    | 0.14  | 0.12      |  |
| ID41<br>ID40 | 0.01                   | 0.05     | 0.67     | 1           | 0.04 | 0.95       | 0.91  | 0.94      | 0.16    | 0.10  | 0.14      | 0.13    | 0.10  | 0.11      |  |
| ID42         | 0.80                   | 0.84     | 0.8      | 1           | 0.01 | 0.93       | 0.91  | 0.92      | 0.15    | 0.17  | 0.08      | 0.12    | 0.10  | 0.08      |  |
| ID43         | 0.9                    | 0.64     | 1        | 1           | 0    | 0.93       | 0.9   | 0.95      | 0.14    | 0.14  | 0.1       | 0.11    | 0.16  | 0.12      |  |
| ID44         | 0                      | 0        | 0        | 0           | 0    | 0.93       | 0.95  | 0.92      | 0.12    | 0.12  | 0.51      | 0.09    | 0.1   | 0.64      |  |
| ID45         | 0.01                   | 0        | 0.01     | 1           | 1.06 | 0.94       | 0.95  | 0.93      | 0.15    | 0.18  | 0.13      | 0.12    | 0.14  | 0.14      |  |
| ID46         | 0.6                    | 0.56     | 0.69     | 0.52        | 0    | 0.94       | 0.95  | 0.76      | 0.18    | 0.15  | 0.36      | 0.16    | 0.12  | 0.53      |  |
| ID47         | 0.26                   | 0.29     | 0.16     | 1           | 0.09 | 0.94       | 0.93  | 0.94      | 0.18    | 0.2   | 0.12      | 0.15    | 0.18  | 0.12      |  |
| ID48         | 0.68                   | 0.62     | 0.8      | 0.67        | 0    | 0.93       | 0.93  | 0.92      | 0.12    | 0.11  | 0.06      | 0.1     | 0.12  | 0.09      |  |
| ID49         | 0.85                   | 0.79     | 0.86     | 1           | 0.01 | 0.95       | 0.92  | 0.93      | 0.22    | 0.14  | 0.26      | 0.21    | 0.19  | 0.51      |  |
| ID50         | 0.69                   | 0.6      | 0.44     | 1           | 0.03 | 0.93       | 0.95  | 0.94      | 0.17    | 0.19  | 0.12      | 0.14    | 0.23  | 0.11      |  |
| ID51         | 0.35                   | 0.23     | 0.15     | 1           | 0.07 | 0.93       | 0.91  | 0.93      | 0.23    | 0.15  | 0.2       | 0.23    | 0.28  | 0.26      |  |
| ID52         | 0.23                   | 0.35     | 0.18     | 1           | 0.16 | 0.93       | 0.94  | 0.95      | 0.16    | 0.18  | 0.1       | 0.13    | 0.16  | 0.1       |  |
| ID53         | 0.12                   | 0.21     | 0.09     | 1           | 0.16 | 0.93       | 0.9   | 0.93      | 0.19    | 0.17  | 0.13      | 0.15    | 0.16  | 0.12      |  |
| ID54         | 0.01                   | 0.01     | 0.01     | 0.67        | 1.72 | 0.94       | 0.88  | 0.94      | 0.18    | 0.19  | 0.14      | 0.14    | 0.16  | 0.1       |  |
| ID55         | 0.37                   | 0.35     | 0.29     | 1           | 0.07 | 0.94       | 0.94  | 0.94      | 0.2     | 0.18  | 0.14      | 0.17    | 0.2   | 0.17      |  |
| ID56         | 0.72                   | 0.6      | 0.67     | 0.67        | 0.01 | 0.94       | 0.91  | 0.88      | 0.19    | 0.18  | 0.13      | 0.17    | 0.19  | 0.16      |  |
| ID57         | 0                      | 0        | 0        | 0           | 0    | 0.94       | 0.94  | 0.89      | 0.14    | 0.13  | 0.14      | 0.11    | 0.09  | 0.16      |  |
| ID58         | 0.22                   | 0.25     | 0.18     | 0.29        | 0.09 | 0.94       | 0.94  | 0.86      | 0.21    | 0.22  | 0.18      | 0.2     | 0.24  | 0.15      |  |
| ID59         | 0                      | 0        | 0        | 0           | 0    | 0.92       | 0.94  | 0.91      | 0.14    | 0.14  | 0.11      | 0.12    | 0.12  | 0.09      |  |
| ID60         | 0.33                   | 0.23     | 0.46     | 0.38        | 0.04 | 0.95       | 0.96  | 0.98      | 0.14    | 0.14  | 0.61      | 0.11    | 0.12  | 0.37      |  |
| ID61         | 0                      | 0        | 0        | 0           | 0    | 0.96       | 0.96  | 0.96      | 0.19    | 0.17  | 0.16      | 0.16    | 0.16  | 0.14      |  |
| ID62         | 0.69                   | 0.45     | 0.77     | 0.71        | 0    | 0.92       | 0.93  | 0.97      | 0.14    | 0.14  | 0.43      | 0.11    | 0.11  | 0.73      |  |
| ID63         | 0.94                   | 0.9      | 1        | 1           | 0    | 0.93       | 0.94  | 0.95      | 0.15    | 0.17  | 0.09      | 0.12    | 0.16  | 0.07      |  |
| ID64         | 0.6                    | 0.64     | 0.62     | 0.5         | 0.06 | 0.95       | 0.95  | 0.95      | 0.16    | 0.17  | 0.13      | 0.13    | 0.14  | 0.11      |  |
| ID65         | 0.5                    | 0.35     | 0.59     | 0.62        | 0.03 | 0.94       | 0.95  | 0.94      | 0.14    | 0.17  | 0.08      | 0.11    | 0.14  | 0.09      |  |
| ID66         | 0.7                    | 0.36     | 0.57     | 1           | 0.01 | 0.87       | 0.92  | 0.92      | 0.13    | 0.13  | 0.08      | 0.1     | 0.13  | 0.08      |  |
| ID67         | 0.26                   | 0.15     | 0.33     | 0.25        | 0.01 | 0.94       | 0.94  | 0.94      | 0.22    | 0.17  | 0.11      | 0.21    | 0.13  | 0.1       |  |
| ID68         | 0.65                   | 0.00     | 0.00     | 0.00        | 0.01 | 0.94       | 0.92  | 0.95      | 0.15    | 0.13  | 0.31      | 0.10    | 0.10  | 0.10      |  |
|              |                        |          |          |             |      |            |       |           |         |       |           |         |       |           |  |

Table S7: Details of seizure detection results of MVPFormer-S with 18 subject pre-training. Kappa is the inter-rater agreement. The classification metrics report the raw and episodic metrics relevant for the seizure classification task. The similarity reports the breakdown of the cosine similarity in each of the considered scenarios.

Table S8: **Details of seizure detection results of MV-Llama with 18 subject pretraining.** Kappa is the inter-rater agreement. The classification metrics report the raw and episodic metrics relevant for the seizure classification task. The similarity reports the breakdown of the cosine similarity in each of the considered scenarios.

|              |       |          | assincari | E in l'     |              |         | There |           | 5       | D      | ity       |         | T       |           |
|--------------|-------|----------|-----------|-------------|--------------|---------|-------|-----------|---------|--------|-----------|---------|---------|-----------|
|              |       | Raw      | 64        | Episodic    | <i>c (</i> ) |         | True  |           |         | Randoi | n         |         | Two-ste | p         |
| Subject      | Kappa | f1-score | f1-score  | sensitivity | tp/h         | average | ıctal | non-ictal | average | ictal  | non-ictal | average | ictal   | non-ictal |
| ID19         | 0.00  | 0.01     | 0.12      | 0.06        | 0.00         | 0.32    | 0.31  | 0.29      | 0.23    | 0.23   | 0.23      | 0.23    | 0.22    | 0.22      |
| ID20         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.96    | 0.97  | 0.95      | 0.18    | 0.21   | 0.14      | 0.15    | 0.18    | 0.13      |
| ID21         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.90  | 0.96      | 0.16    | 0.15   | 0.15      | 0.13    | 0.23    | 0.20      |
| ID22         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.78  | 0.95      | 0.15    | 0.13   | 0.08      | 0.12    | 0.16    | 0.07      |
| ID23         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.96    | 0.96  | 0.96      | 0.25    | 0.26   | 0.20      | 0.22    | 0.24    | 0.24      |
| ID24         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.94  | 0.93      | 0.19    | 0.20   | 0.19      | 0.16    | 0.18    | 0.17      |
| ID25         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.96  | 0.96      | 0.17    | 0.19   | 0.16      | 0.14    | 0.18    | 0.15      |
| ID26         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.96    | 0.97  | 0.96      | 0.18    | 0.22   | 0.10      | 0.14    | 0.18    | 0.12      |
| ID27         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.96    | 0.96  | 0.96      | 0.22    | 0.25   | 0.19      | 0.19    | 0.25    | 0.24      |
| ID28         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.96    | 0.96  | 0.96      | 0.17    | 0.19   | 0.08      | 0.13    | 0.17    | 0.13      |
| ID29         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.95  | 0.96      | 0.16    | 0.15   | 0.09      | 0.12    | 0.13    | 0.13      |
| ID30         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.96    | 0.96  | 0.96      | 0.19    | 0.20   | 0.15      | 0.16    | 0.19    | 0.17      |
| ID31         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.96    | 0.96  | 0.96      | 0.20    | 0.24   | 0.16      | 0.16    | 0.24    | 0.18      |
| ID32         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.88  | 0.96      | 0.18    | 0.16   | 0.16      | 0.15    | 0.18    | 0.12      |
| ID33         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.93    | 0.95  | 0.92      | 0.23    | 0.24   | 0.19      | 0.21    | 0.25    | 0.20      |
| ID34         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.95  | 0.96      | 0.15    | 0.16   | 0.12      | 0.11    | 0.15    | 0.12      |
| ID35         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.93    | 0.94  | 0.92      | 0.17    | 0.18   | 0.16      | 0.14    | 0.17    | 0.18      |
| ID36         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.96    | 0.96  | 0.96      | 0.22    | 0.24   | 0.16      | 0.19    | 0.21    | 0.17      |
| ID37         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.92    | 0.88  | 0.87      | 0.16    | 0.15   | 0.36      | 0.14    | 0.16    | 0.57      |
| ID38         | 0.30  | 0.00     | 0.00      | 0.00        | 0.01         | 0.93    | 0.95  | 0.91      | 0.14    | 0.16   | 0.10      | 0.11    | 0.16    | 0.07      |
| ID39         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.93    | 0.87  | 0.94      | 0.16    | 0.14   | 0.08      | 0.13    | 0.13    | 0.11      |
| ID40         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.95  | 0.95      | 0.15    | 0.17   | 0.10      | 0.12    | 0.15    | 0.12      |
| ID41         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.96    | 0.92  | 0.95      | 0.16    | 0.15   | 0.13      | 0.13    | 0.15    | 0.11      |
| ID42         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.92  | 0.94      | 0.15    | 0.18   | 0.09      | 0.12    | 0.17    | 0.08      |
| ID43         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.91  | 0.96      | 0.15    | 0.11   | 0.12      | 0.11    | 0.17    | 0.12      |
| ID44         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.96  | 0.86      | 0.12    | 0.12   | 0.37      | 0.09    | 0.10    | 0.50      |
| ID45         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.95  | 0.93      | 0.15    | 0.15   | 0.12      | 0.12    | 0.14    | 0.15      |
| ID46         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.96  | 0.80      | 0.18    | 0.14   | 0.28      | 0.16    | 0.12    | 0.53      |
| ID47         | 0.74  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.94  | 0.95      | 0.18    | 0.20   | 0.14      | 0.15    | 0.18    | 0.13      |
| ID41<br>ID48 | 0.14  | 0.00     | 0.00      | 0.00        | 0.00         | 0.93    | 0.94  | 0.93      | 0.10    | 0.12   | 0.07      | 0.19    | 0.12    | 0.09      |
| ID49         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.01  | 0.92      | 0.21    | 0.14   | 0.29      | 0.20    | 0.19    | 0.46      |
| ID50         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.96  | 0.94      | 0.17    | 0.21   | 0.12      | 0.14    | 0.23    | 0.11      |
| ID51         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.93    | 0.90  | 0.92      | 0.22    | 0.15   | 0.12      | 0.14    | 0.28    | 0.21      |
| ID52         | 0.00  | 0.00     | 0.67      | 0.50        | 0.00         | 0.94    | 0.95  | 0.92      | 0.16    | 0.17   | 0.10      | 0.21    | 0.17    | 0.10      |
| ID52         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.03  | 0.94      | 0.10    | 0.20   | 0.11      | 0.16    | 0.17    | 0.10      |
| ID54         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.90  | 0.05      | 0.19    | 0.16   | 0.12      | 0.10    | 0.16    | 0.12      |
| ID55         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.05  | 0.95      | 0.10    | 0.10   | 0.12      | 0.15    | 0.10    | 0.10      |
| ID56         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.95  | 0.35      | 0.20    | 0.21   | 0.15      | 0.17    | 0.18    | 0.17      |
| ID57         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.05    | 0.01  | 0.00      | 0.13    | 0.10   | 0.10      | 0.10    | 0.10    | 0.13      |
| 1D59         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.95  | 0.90      | 0.14    | 0.13   | 0.10      | 0.10    | 0.09    | 0.14      |
| ID50         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.94  | 0.00      | 0.21    | 0.22   | 0.15      | 0.19    | 0.23    | 0.17      |
| ID59<br>ID60 | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.94  | 0.95      | 0.14    | 0.15   | 0.00      | 0.11    | 0.12    | 0.10      |
| ID61         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.90  | 0.98      | 0.14    | 0.15   | 0.02      | 0.11    | 0.12    | 0.41      |
| ID01         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.90    | 0.90  | 0.97      | 0.19    | 0.20   | 0.10      | 0.10    | 0.10    | 0.15      |
| 1D62         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.92    | 0.93  | 0.97      | 0.13    | 0.13   | 0.52      | 0.10    | 0.10    | 0.74      |
| 1D03         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.94    | 0.95  | 0.95      | 0.10    | 0.18   | 0.11      | 0.12    | 0.15    | 0.07      |
| 1D64<br>1D65 | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.95  | 0.96      | 0.16    | 0.18   | 0.12      | 0.13    | 0.15    | 0.11      |
| ID65         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.95  | 0.95      | 0.14    | 0.15   | 0.11      | 0.11    | 0.14    | 0.09      |
| 1D66         | 0.74  | 0.00     | 0.00      | 0.00        | 0.00         | 0.88    | 0.93  | 0.92      | 0.12    | 0.14   | 0.08      | 0.09    | 0.13    | 0.08      |
| ID67         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.95    | 0.95  | 0.94      | 0.22    | 0.17   | 0.09      | 0.20    | 0.14    | 0.10      |
| ID68         | 0.00  | 0.00     | 0.00      | 0.00        | 0.00         | 0.86    | 0.91  | 0.89      | 0.16    | 0.17   | 0.10      | 0.20    | 0.16    | 0.08      |

#### G.3 Effects of the attention mechanism

To assess the validity of our MVPA scheme, we train MV-Llama, a model almost equivalent to MVPFormer-S that uses vanilla attention instead of MVPA. While MV-Llama uses vanilla attention, it is still based on the SOTA Llama2 architecture. We also re-use the vanilla positional encoding, with a simple adjustment to recover a one-to-one correspondence between the positional encoding and the position of the patch in the time-series (the Cantor pairing function, see App. A.1).

Table 1 in the main text indicates that vanilla attention does not perform seizure detection at a level comparable to MVPA. In particular, the performance of MV-Llama is poor, indicating that it cannot generalize to this task. We argue this is due to higher flexibility of the internal representations generated by MVPA, which better lend themselves to further tasks, such as seizure classification.

#### G.4 Seizure detection with Brant-2

Figure S13 shows the full distribution of the results of Brant-2 on the seizure detection task. Table S9 presents the detailed subject-by-subject breakdown of the performance of Brant-2.



Figure S13: **Performance of Brant-2 on the classification task.** (a) The F1-score, sensitivity, and fp/h are reported. The results are computed using the same schema as MVPFormer. Brant-2, however, is not capable of generalizing to this dataset. (b) Cohen's kappa is used to measure the agreement between the artificial assistant and the human expert. The average kappa is 0.08, not competitive with MVPFormer.

| Table S9:  | Details of    | f seizure   | detection     | results o    | of Brant-   | -2 with | 18 sı  | ubject  | pre- |
|------------|---------------|-------------|---------------|--------------|-------------|---------|--------|---------|------|
| training.  | Kappa is th   | e inter-rat | er agreemen   | t. The clas  | ssification | metrics | report | the raw | and  |
| episodic m | etrics releva | ant for the | seizure class | sification t | ask.        |         |        |         |      |

| Classification metrics |                    |       |          |          |             |      |
|------------------------|--------------------|-------|----------|----------|-------------|------|
|                        |                    |       | Raw      |          | Episodic    |      |
|                        | $\mathbf{Subject}$ | Kappa | f1-score | f1-score | sensitivity | fp/h |
|                        | ID19               | N.A.  | N.A.     | N.A.     | N.A.        | N.A. |
|                        | ID20               | 0.47  | 0        | 0        | 0           | 0.02 |
|                        | ID21               | 0     | 0        | 0        | 0           | 0    |
|                        | ID22               | 0.4   | 0        | 0        | 0           | 0.01 |
|                        | ID23               | 0     | 0        | 0        | 0           | 0    |
|                        | ID24               | -0.02 | 0        | 0        | 0           | 0.32 |
|                        | ID25               | 0     | 0        | 0        | 0           | 0    |
|                        | ID26               | 0     | 0        | 0        | 0           | 0    |
|                        | ID27               | 0     | 0        | 0        | 0           | 0.12 |
|                        | ID28               | -0.01 | 0        | 0        | 0           | 0.49 |
|                        | ID29               | 0.01  | 0        | 0        | 0           | 0.01 |
|                        | ID30               | 0     | 0        | 0        | 0           | 0    |
|                        | ID31               | 0     | 0        | 0        | 0           | 0    |
|                        | ID32               | 0     | 0        | 0        | 0           | 0    |
|                        | ID33               | 0     | 0        | 0        | 0           | 0    |
|                        | ID34               | 0.64  | 0        | 0        | 0           | 0.01 |
|                        | ID35               | 0     | 0        | 0        | 0           | 0    |
|                        | ID36               | N.A.  | N.A.     | N.A.     | N.A.        | N.A. |
|                        | ID37               | 0     | 0        | 0        | 0           | 0    |
|                        | ID38               | 0.01  | 0        | 0        | 0           | 0.17 |
|                        | ID39               | 0     | 0        | 0        | 0           | 0    |
|                        | ID40               | 0     | 0        | 0        | 0           | 0.29 |
|                        | ID41               | 0     | 0        | 0        | 0           | 0.11 |
|                        | ID42               | 0.07  | 0        | 0        | 0           | 0.16 |
|                        | ID43               | 0     | 0        | 0        | 0           | 0    |
|                        | ID44               | 0.74  | 0        | 0        | 0           | 0.01 |
|                        | ID45               | 0     | 0        | 0        | 0           | 0    |
|                        | ID46               | 0.01  | 0        | 0        | 0           | 0.06 |
|                        | ID47               | 0.02  | 0        | 0        | 0           | 0.15 |
|                        | ID48               | 0     | 0        | 0        | 0           | 0    |
|                        | ID49               | 0.14  | 0        | 0        | 0           | 0.03 |
|                        | ID50               | 0     | 0        | 0        | 0           | 0.11 |
|                        | ID51               | 0     | 0        | 0        | 0           | 0.04 |
|                        | ID52               | 0     | 0        | 0        | 0           | 0.08 |
|                        | ID53               | N.A.  | N.A.     | N.A.     | N.A.        | N.A. |
|                        | ID54               | 0     | 0        | 0        | 0           | 0    |
|                        | ID55               | N.A.  | N.A.     | N.A.     | N.A.        | N.A. |
|                        | ID56               | 0.56  | 0        | 0        | 0           | 0.02 |
|                        | ID57               | 0     | 0        | 0        | 0           | 0.01 |
|                        | ID58               | 0     | 0        | 0        | 0           | 0    |
|                        | ID59               | 0.14  | Ő        | 0        | 0           | 0.01 |
|                        | ID60               | 0     | Ő        | 0        | Ő           | 0    |
|                        | ID61               | 0     | 0        | 0        | 0           | 0    |
|                        | ID62               | 0.16  | ů<br>0   | Ŭ<br>0   | Ŭ.          | 0.03 |
|                        | ID63               | 0     | 0        | 0        | Ő           | 0    |
|                        | ID64               | 0.07  | 0.02     | 0.09     | 0.08        | 0.22 |
|                        | ID65               | 0     | 0        | 0        | 0           | 0    |
|                        | ID66               | 0     | 0        | 0        | 0           | 0.29 |
|                        | ID67               | 0     | 0        | 0        | 0           | 0.25 |
|                        | ID68               | 0     | 0        | 0        | 0           | 0    |
|                        | 1000               |       | 0        |          | 0           | v    |

#### G.5 Effects of the scale of the pre-training dataset

The performance of LLMs as the size of their training dataset increases has been investigated quite thoroughly [48, 49], giving rise to a variety of scaling laws. Following Chinchilla's scaling law, a model with 75 million parameters like MVPFormer-S should be trained with around 2 billion tokens, while we only have 400 million at our disposal.

The architecture of the model and the nature of the training data, however, make it unclear whether such laws can be adopted for MVPFormer as well. We investigate this behavior by continuing the training of MVPFormer-S on 40 more subjects, to bring the total to 58 pre-training subjects for almost 7,000 hours of iEEG recordings. In particular, MVPFormer is initially trained on 304 ictal events, and then further on 323 more. Therefore, we are left with 10 unseen subject to test the downstream seizure detection task.



Figure S14: Seizure detection with 58 patient pre-training. (a) Seizure detection results of MVPFormer-S on the 10 remaining unseen patients: the F1-score, sensitivity, and fp/h are reported. All performance metrics are improved with respect to the original MVPFormer model. The raw and episodic F1-scores are significantly different here, indicating that MVPFormer benefits of the episode merging effect of post-processing on these patients. The false positive rate has decreased further with the scale of the pre-training dataset. (b) Inter-rater agreement of MVPFormer with the human expert: Cohen's kappa is used to measure the agreement between the artificial assistant and the human expert. The average kappa is increased to 0.48. The distribution of kappa values again indicates that there is considerable variability in the agreement.

Figure S14 shows the performance of the 58-subject MVPFormer-S on the 10 unseen subjects (for a detailed breakdown see Table S10). On the other hand, Figure S15 shows the results of the original MVPFormer-S model on those same 10 subjects (for a detailed breakdown see Table S11). All performance metrics improve with a growing pre-training dataset size, although on a small test cohort, indicating that increasing the number of subjects in the pre-training dataset has a net positive effect on the downstream classification task.



Figure S15: Seizure detection with 18 subjects pre-training on a selection of 10 unseen subjects. (a) Seizure detection results of MVPFormer-S on the 10 subjects excluded from the 58-subjects model: the F1-score, sensitivity, and fp/h are reported. The raw and episodic F1-scores are significantly different here, indicating that MVPFormer benefits of the episode merging effect of post-processing on these patients. These results are a subset of those presented in the Results section. (b) Inter-rater agreement of MVPFormer with the human expert: Cohen's kappa is used to measure the agreement between the artificial assistant and the human expert. The average kappa is 0.46, reduced from the overall results indicating that these subjects are more difficult than average. The distribution of kappa values again indicates that there is considerable variability in the agreement.

|         |       | Classification metrics |          |             |      |  |
|---------|-------|------------------------|----------|-------------|------|--|
|         |       | Raw                    | Episodic |             |      |  |
| Subject | Kappa | f1-score               | f1-score | sensitivity | fp/h |  |
| ID59    | 0.40  | 0.41                   | 0.44     | 0.86        | 0.86 |  |
| ID60    | 0.07  | 0.05                   | 0.10     | 0.07        | 0.07 |  |
| ID61    | 0.54  | 0.39                   | 0.47     | 0.50        | 0.50 |  |
| ID62    | 0.00  | 0.00                   | 0.00     | 0.00        | 0.00 |  |
| ID63    | 0.65  | 0.50                   | 0.71     | 0.86        | 0.86 |  |
| ID64    | 0.76  | 0.84                   | 0.67     | 1.00        | 1.00 |  |
| ID65    | 0.60  | 0.67                   | 0.62     | 0.45        | 0.45 |  |
| ID66    | 0.50  | 0.41                   | 0.52     | 1.00        | 1.00 |  |
| ID67    | 0.79  | 0.50                   | 0.57     | 1.00        | 1.00 |  |
| ID68    | 0.22  | 0.22                   | 0.15     | 0.75        | 0.75 |  |
|         |       |                        |          |             |      |  |

Table S10: Details of seizure detection results of MVPFormer-S with 58 patient pre-training. Kappa is the inter-rater agreement. The classification metrics report the raw and episodic metrics relevant for the seizure classification task.

Table S11: Details of seizure detection results of MVPFormer-S with 18 patient pre-training on a selection of 10 subjects. Kappa is the inter-rater agreement. The classification metrics report the raw and episodic metrics relevant for the seizure classification task.

|                    |       | Classification metrics |          |             |      |  |  |
|--------------------|-------|------------------------|----------|-------------|------|--|--|
|                    |       | Raw                    |          | Episodic    |      |  |  |
| $\mathbf{Subject}$ | Kappa | f1-score               | f1-score | sensitivity | fp/h |  |  |
| ID59               | 0.22  | 0.25                   | 0.18     | 0.29        | 0.09 |  |  |
| ID60               | 0.00  | 0.00                   | 0.00     | 0.00        | 0.00 |  |  |
| ID61               | 0.33  | 0.23                   | 0.46     | 0.38        | 0.04 |  |  |
| ID62               | 0.00  | 0.00                   | 0.00     | 0.00        | 0.00 |  |  |
| ID63               | 0.69  | 0.45                   | 0.77     | 0.71        | 0.00 |  |  |
| ID64               | 0.94  | 0.90                   | 1.00     | 1.00        | 0.00 |  |  |
| ID65               | 0.60  | 0.64                   | 0.62     | 0.50        | 0.06 |  |  |
| ID66               | 0.50  | 0.35                   | 0.59     | 0.62        | 0.03 |  |  |
| ID67               | 0.70  | 0.36                   | 0.57     | 1.00        | 0.01 |  |  |
| ID68               | 0.26  | 0.15                   | 0.33     | 0.25        | 0.01 |  |  |

#### G.6 Effects of the selection of channels

As discussed in previous sections, the number of channels can vary considerably across subjects. To speed up the seizure detection task, we quickly select a subset of the channels (up to 50) which visually appear least noisy and most relevant. We also test the effect of including all channels, expecting it to decrease both the speed and performance due to the decrease of the overall signal-to-noise ratio.

Figure S16 shows that the performance decreases when we use all channels (for a detailed breakdown see Table S12). This is expected, as the noise contained in the entire recording increases together with the information content. MVPFormer's ability to generalize is not affected by the number of channels, but the noise affects the performance. Therefore, the optimal real-world operation of MVPFormer is obtained by selecting a subset of channels for detection.

![](_page_42_Figure_3.jpeg)

Figure S16: Seizure detection with 18 patient pre-training and evaluation on all channels. (a) Seizure detection results of MVPFormer-S on the 10 remaining unseen patients evaluated on all channels: the F1-score, sensitivity, and fp/h are reported. The performance metrics are reduced with respect to the results obtained when selecting a subset of the channels. MVPFormer is not affected by the number of channels, but the increase of noise emerging from all the channels contributes to a reduction in performance. (b) The average kappa is 0.36, reduced from the evaluation on a subset of channels.

|         |       | C        | lassincati | on metrics  |      |
|---------|-------|----------|------------|-------------|------|
| ~ • •   |       | Raw      |            | Episodic    |      |
| Subject | Kappa | f1-score | f1-score   | sensitivity | fp/h |
| ID19    | -0.05 | 0        | 0          | 0           | 0    |
| ID20    | 0.84  | 0.47     | 0.5        | 0.5         | 0.5  |
| ID21    | 0.86  | 0.4      | 0.8        | 1           | 1    |
| ID22    | 0.91  | 0.81     | 0.89       | 1           | 1    |
| ID23    | 0.25  | 0.17     | 0.14       | 0.09        | 0.09 |
| ID24    | 0.91  | 0.9      | 0.93       | 0.93        | 0.93 |
| ID25    | 0.45  | 0.17     | 0.25       | 0.25        | 0.25 |
| ID26    | 0.25  | 0.04     | 0.22       | 1           | 1    |
| ID27    | 0.98  | 0.93     | 1          | 1           | 1    |
| ID28    | 0.64  | 0.42     | 0.67       | 0.5         | 0.5  |
| ID29    | 0.05  | 0.04     | 0.07       | 0.29        | 0.29 |
| ID30    | 0.2   | 0.31     | 0.19       | 0.11        | 0.11 |
| ID31    | 0.98  | 0.95     | 1          | 1           | 1    |
| ID32    | 0.98  | 0.92     | 1          | 1           | 1    |
| ID33    | 0.01  | 0        | 0          | 0           | 0    |
| ID34    | 0.04  | 0.06     | 0.04       | 1           | 1    |
| ID35    | 0.61  | 0.63     | 0.55       | 0.43        | 0.43 |
| ID36    | 0.5   | 0.62     | 0.51       | 0.37        | 0.37 |
| ID37    | 0.03  | 0.02     | 0.04       | 1           | 1    |
| ID38    | 0.97  | 0.97     | 1          | 1           | 1    |
| ID39    | 0.77  | 0.52     | 0.67       | 1           | 1    |
| ID40    | 0.37  | 0.5      | 0.31       | 0.8         | 0.8  |
| ID41    | 0.52  | 0.47     | 0.55       | 1           | 1    |
| ID42    | 0.63  | 0.54     | 0.57       | 1           | 1    |
| ID43    | 0.94  | 0.71     | 1          | 1           | 1    |
| ID44    | 0     | 0        | 0          | 0           | 0    |
| ID45    | Ő     | Ő        | 0          | 0           | Ő    |
| ID46    | 0.06  | 0.03     | 0.06       | 0.19        | 0.19 |
| ID47    | 0.23  | 0.29     | 0.16       | 1           | 1    |
| ID48    | 0.25  | 0.18     | 0.31       | 0.33        | 0.33 |
| ID49    | 0.4   | 0.47     | 0.4        | 1           | 1    |
| ID50    | 0.05  | 0.05     | 0.06       | 1           | 1    |
| ID51    | 0.06  | 0.05     | 0.07       | 1           | 1    |
| ID52    | 0     | 0.01     | 0.01       | 1           | 1    |
| ID53    | 0.13  | 0.21     | 0.09       | 1           | 1    |
| ID54    | 0.01  | 0.01     | 0.01       | 0.67        | 0.67 |
| ID55    | 0.37  | 0.35     | 0.29       | 1           | 1    |
| ID56    | 0.33  | 0.34     | 0.31       | 0.67        | 0.67 |
| ID57    | 0     | 0        | 0          | 0           | 0    |
| ID58    | 0.29  |          | 0          | 0           | 0    |
| ID59    | 0.02  | 0        | 0          | 0           | 0    |
| ID60    | 0.02  |          | 0          | 0           | 0    |
| ID61    | 0     | 0        | 0          | Ő           | 0    |
| ID62    | 0.04  | 0.03     | 0.07       | 0.29        | 0.29 |
| ID63    | 0.21  | 0.39     | 0.18       | 0.75        | 0.75 |
| ID64    | 0.25  | 0.3      | 0.29       | 0.25        | 0.25 |
| ID65    | 0.04  | 0        | 0          | 0           | 0    |
| ID66    | 0.76  | 0.31     | 0.4        | 0.5         | 0.5  |
| ID67    | 0.17  | 0.14     | 0.22       | 0.25        | 0.25 |
| ID68    | 0.61  | 0.00     | 0.00       | 0.00        | 0.01 |
|         |       |          |            |             |      |

Table S12: Details of seizure detection results of MVPFormer-S with evaluation on all channels. Kappa is the inter-rater agreement. The classification metrics report the raw and episodic metrics relevant for the seizure classification task. Classification metrics

#### G.7 Patient classification difficulty

In the medical practice each patient has unique seizure presentations, though they might be broadly grouped into different categories [68]. As an effect, some patients have seizures which can be considered more typical (Figure S2a), and hence easier to detect, while others might have very atypical events (Figure S2b). There might be broad disagreement among neurologists over these atypical seizures, and at the same time no disagreement at all over the typical patients [69].

This phenomenon intuitively creates a difficulty scale among the patients, which also affects MVPFormer and contributes to the spread of performance between the model and the human expert. To better assess the impact of this latent patient classification difficulty we performed a multiple correlation analysis using the total recording length, the number of seizures, and the frequency of seizures to predict the kappa score, yielding an  $R^2$  of 0.054. The model performance is thus independent of the three variables, and we believe the difficulty might help explain most of the variance. The literature supports this hypothesis, as the subjects themselves can account for up to 65% of the variance [42] while the clinical setup itself has no impact.

#### G.8 Channel connectivity map

Generating future iEEG signal embedding implicitly places the greatest emphasis on the dimension of time, but to do so it is necessary to consider the interactions between channels as well. MVPFormer thus takes into consideration all electrodes concurrently, as electric potentials flow across different areas and circuits in the brain following their intrinsic connections and constraints [70, 71]. The number of channels depends on the number of electrodes decided for a specific patient and the clinical setup. As iEEG implantations are decided on a case-by-case basis by a physician, there is no uniform standard on where to place the electrodes, in contrast with the 10-20 system [72] for EEG. Therefore, we cannot give MVPFormer any a priori knowledge of how the channels will interact in space, and the model has to learn it on its own. MVPA enables our model to dynamically learn these connections to build an internal map of the flow, becoming independent from a specific electrode configuration.

It is well-known that two neighboring brain regions might not be as strongly connected as two faraway regions. The relationship between the electrodes (and hence the channels) mirrors this behavior. To truly understand the link between two channels the model must build a map of the connection strength between different brain areas and how these connections impact the diffusion of electric fields across channels. In MVPA, this understanding is the underpinning of the channel-based component. In particular, the complex interplay between the query, the channel codebook, and the channel attention, acts as the first level of processing. Further, the deep structure with multiple layers provides more representational power, as is typical of deep models.

Figure S17 shows that the channel-based MVPA component encodes a form of the brain connectivity map. Initially, the map is random as MVPA is randomly initialized. As training progresses, the attention magnitude among the channels starts to differentiate, building a map of the connection strength. The map is dominated by the diagonal component, which indicates that in general neighboring channels are more related than distant channels. However, it is possible to clearly distinguish clusters of strongly connected electrodes and also skipped connections, which possibly refer to strong connections between distant regions. Since the channel distance is relative, it can apply to arbitrary clinical setups and is not limited to already seen channels. Moreover, as the channel attention is a function of both the query content and the channel distance, the combination of the two can effectively modulate the attention even on unseen subjects.

![](_page_46_Figure_0.jpeg)

Figure S17: Attention components before and after training. (a) The channelbased attention component is randomly initialized. At the end of training, it shows the diagonal structure which indicates that the relationship between the channels is mainly one of proximity. This is expected, as nearby channels are expected to be more closely related, and showcases MVPA's learning outcome. (b) The time-based attention component is also randomly initialized. At the end of training, it shows that segments which are close in time are more related. Particularly, few closest segments are attended to more strongly, as the content-based attention's lookback windows is limited to a few segments.

#### G.9 Results on the MAYO and FNUSA datasets

We test MVPFormer additionally against Brant-2 on the iEEG MAYO and FNUSA datasets. Both datasets are single-channel iEEG datasets containing both physiological and pathological activity. We use these datasets to evaluate MVPFormer's performance on extremely noisy data, and to compare its resilience with the SOTA.

The MAYO dataset contains 24 patients for a total of 130 hours of data. In particular, 36% of the data is non-ictal, 9% is ictal, and 53% is noise. Moreover, 18 subjects contain no ictal activity and 9 patients are fully noise, with 13 being majority noise. The FNUSA dataset contains 14 patients for a total of 160 hours. 48% of the dataset is non-ictal, 27% is ictal, and 23% is noise. Moreover, the data of 2 subjects is completely noise, and the data of 3 other subjects is fully ictal. Both datasets are much smaller in scale than the SWEC iEEG dataset, and are heavily dominated by noise, both artifact and powerline. In particular, the SWEC iEEG dataset is almost 2,000 times larger and is also carefully evaluated by an expert neurologist to remove channels which contain too much noise or artifacts.

As we wish to assess MVPFormer in a realistic, real-world scenario, we do not remove any noise from the dataset but test them as-is. In particular, we consider noise and physiological activity as one category, and pathological activity as another. However, kappa scores are not meaningful with such small datasets, therefore we provide the aggregate F1-score, and the average sensitivity and specificity. Specifically, given the fact that many patients do not contain ictal activity, we do not compute the average F1-score across subjects, but pool together all subject's results and compute the aggregate F1-score.

We use the same MVPFormer models pre-trained on our SWEC iEEG dataset, and train a specific classification head for either MAYO or FNUSA by fine-tuning on the first four patients. Then, we test on the remaining patients. We also use the Brant-2 model whose pre-trained weights are publicly available, and fine-tune in the same manner as MVPFormer using the fine-tuning code provided by the authors.

The results can be found in Tables S13 and S14. Given the very low signal-to-noise ratio of both datasets, overall performance is affected. On the FNUSA dataset, where the amount of noise is more moderate, all models perform similarly, with MVPFormer-M showing a higher specificity. However, MVPFormer has a clear advantage on the MAYO dataset, with almost double the F1-score with respect to Brant-2. The difference between MVPFormer-S and MVPFormer-M is minimal, as the sizes of the datasets involved are too small to fully train a very large model such as MVPFormer-M (see Appendix G.2 for more information).

| Table S13: S | Summary 6     | of seizure   | detection     | results of all m   | nodels on  | the MA   | YO iE  | $\mathbf{EG}$ |
|--------------|---------------|--------------|---------------|--------------------|------------|----------|--------|---------------|
| dataset. K   | appa is the   | inter-rater  | agreement.    | The classification | on metrics | report t | he raw | and           |
| episodic met | trics relevan | nt for the s | eizure classi | fication task.     |            |          |        |               |

| Model       | F1-score | Sensitivity | Specificity |
|-------------|----------|-------------|-------------|
| MVPFormer-M | 0.36     | 0.38        | 0.91        |
| MVPFormer-S | 0.35     | 0.41        | 0.88        |
| Brant-2     | 0.19     | 1.00        | 0.18        |

Table S14: Summary of seizure detection results of all models on the FNUSA iEEG dataset. Kappa is the inter-rater agreement. The classification metrics report the raw and episodic metrics relevant for the seizure classification task.

| Model       | F1-score | Sensitivity | Specificity |
|-------------|----------|-------------|-------------|
| MVPFormer-M | 0.46     | 0.94        | 0.10        |
| MVPFormer-S | 0.46     | 0.99        | 0.03        |
| Brant-2     | 0.46     | 0.99        | 0.02        |

#### G.10 Evaluation of maximum performance

To better characterize MVPFormer's ability to generalize to unseen subjects, we perform the seizure detection task on 40 subjects in two different scenarios. First, we use a model that is pre-trained on those 40 subjects (see Figure S18 and Table S15) to determine MVPFormer's maximum performance. Second, we use a model for which those 40 subjects are unseen (see Figure S19 and Table S16). As expected, with a Kappa score of 0.73 the model trained on the testing subjects achieves superior agreement even to human experts, and can therefore be seen as having learned the training set. On the other hand, as seen with previous results as well, in the unseen subject scenario MVPFormer reaches a Kappa score of 0.56, indicating a high degree of generalization.

![](_page_48_Figure_2.jpeg)

Figure S18: Seizure detection with 58 subject pre-training and evaluation on 40 pre-trained subjects. (a) Seizure detection results of MVPFormer-S on 40 previously seen subjects: the F1-score, sensitivity, and fp/h are reported. The performance metrics are notably improved due to testing on previously seen subjects. (b) Cohen's kappa is used to measure the agreement between the artificial assistant and the human expert. The average kappa is 0.73, notably improved from the baseline.

Table S15: Details of seizure detection results of MVPFormer-S with 58 subject pre-training and evaluation on 40 pre-trained subjects. Kappa is the inter-rater agreement. The classification metrics report the raw and episodic metrics relevant for the seizure classification task.

|                    |       | Classification metrics |          |             |      |  |  |
|--------------------|-------|------------------------|----------|-------------|------|--|--|
|                    |       | Raw                    |          | Episodic    |      |  |  |
| $\mathbf{Subject}$ | Kappa | f1-score               | f1-score | sensitivity | fp/h |  |  |
| ID19               | 0.36  | 0.29                   | 0.29     | 0.19        | 0.09 |  |  |
| ID20               | 0.92  | 0.92                   | 0.67     | 1.00        | 0.01 |  |  |
| ID21               | 0.75  | 0.62                   | 0.50     | 1.00        | 0.02 |  |  |
| ID22               | 0.77  | 0.78                   | 0.73     | 1.00        | 0.02 |  |  |
| ID23               | 0.36  | 0.39                   | 0.39     | 0.27        | 0.09 |  |  |
| ID24               | 0.69  | 0.64                   | 0.76     | 0.79        | 0.10 |  |  |
| ID25               | 0.43  | 0.00                   | 0.00     | 0.00        | 0.01 |  |  |
| ID26               | 0.98  | 0.80                   | 1.00     | 1.00        | 0.00 |  |  |
| ID27               | 0.80  | 0.74                   | 0.80     | 0.75        | 0.01 |  |  |
| ID28               | 0.62  | 0.54                   | 0.67     | 0.75        | 0.03 |  |  |
| ID29               | 0.24  | 0.30                   | 0.22     | 0.13        | 0.01 |  |  |
| ID30               | 0.85  | 0.79                   | 0.87     | 0.89        | 0.10 |  |  |
| ID31               | 0.97  | 0.94                   | 1.00     | 1.00        | 0.00 |  |  |
| ID32               | 0.99  | 0.92                   | 1.00     | 1.00        | 0.00 |  |  |
| ID33               | 0.00  | 0.00                   | 0.00     | 0.00        | 0.00 |  |  |
| ID34               | 0.88  | 0.88                   | 0.90     | 1.00        | 0.01 |  |  |
| ID35               | 0.92  | 0.89                   | 0.93     | 1.00        | 0.01 |  |  |
| ID36               | 0.26  | 0.42                   | 0.27     | 0.17        | 0.02 |  |  |
| ID37               | 0.93  | 0.85                   | 0.80     | 1.00        | 0.01 |  |  |
| ID38               | 0.98  | 0.95                   | 1.00     | 1.00        | 0.00 |  |  |
| ID39               | 0.87  | 0.84                   | 0.80     | 1.00        | 0.01 |  |  |
| ID40               | 0.77  | 0.82                   | 0.62     | 0.80        | 0.02 |  |  |
| ID41               | 0.99  | 0.94                   | 1.00     | 1.00        | 0.00 |  |  |
| ID42               | 0.99  | 0.95                   | 1.00     | 1.00        | 0.00 |  |  |
| ID43               | 0.83  | 0.79                   | 0.80     | 1.00        | 0.02 |  |  |
| ID44               | 0.98  | 0.87                   | 1.00     | 1.00        | 0.00 |  |  |
| ID45               | 0.91  | 0.67                   | 0.67     | 1.00        | 0.01 |  |  |
| ID46               | 0.80  | 0.76                   | 0.86     | 0.76        | 0.00 |  |  |
| ID47               | 0.93  | 0.84                   | 0.86     | 1.00        | 0.00 |  |  |
| ID48               | 0.97  | 0.93                   | 1.00     | 1.00        | 0.00 |  |  |
| ID49               | 0.67  | 0.74                   | 0.67     | 1.00        | 0.04 |  |  |
| ID50               | 0.81  | 0.75                   | 0.67     | 1.00        | 0.01 |  |  |
| ID51               | 0.71  | 0.53                   | 0.40     | 1.00        | 0.02 |  |  |
| ID52               | 0.36  | 0.63                   | 0.24     | 1.00        | 0.12 |  |  |
| ID53               | 0.99  | 0.91                   | 1.00     | 1.00        | 0.00 |  |  |
| ID54               | 0.55  | 0.35                   | 0.36     | 0.67        | 0.03 |  |  |
| ID55               | 0.98  | 0.85                   | 1.00     | 1.00        | 0.00 |  |  |
| ID56               | 0.85  | 0.87                   | 0.75     | 1.00        | 0.02 |  |  |
| ID57               | 0.20  | 0.38                   | 0.11     | 0.83        | 1.74 |  |  |
| ID58               | 0.40  | 0.41                   | 0.46     | 0.86        | 0.09 |  |  |

Table S16: Details of seizure detection results of MVPFormer-S with 18 subjects pre-training on a selection of 40 unseen subjects. Kappa is the inter-rater agreement. The classification metrics report the raw and episodic metrics relevant for the seizure classification task.

|                    |       | Classification metrics |          |             |      |  |  |
|--------------------|-------|------------------------|----------|-------------|------|--|--|
|                    |       | Raw                    |          | Episodic    |      |  |  |
| $\mathbf{Subject}$ | Kappa | f1-score               | f1-score | sensitivity | fp/h |  |  |
| ID19               | -0.05 | 0                      | 0        | 0           | 0.97 |  |  |
| ID20               | 0.98  | 0.87                   | 0.8      | 1           | 0    |  |  |
| ID21               | 0.92  | 0.86                   | 0.8      | 1           | 0    |  |  |
| ID22               | 0.98  | 0.89                   | 1        | 1           | 0    |  |  |
| ID23               | 0.26  | 0.17                   | 0.14     | 0.09        | 0.16 |  |  |
| ID24               | 0.91  | 0.9                    | 0.93     | 0.93        | 0.02 |  |  |
| ID25               | 0.35  | 0                      | 0        | 0           | 0.03 |  |  |
| ID26               | 0.08  | 0.03                   | 0.11     | 1           | 0.18 |  |  |
| ID27               | 0.99  | 0.93                   | 1        | 1           | 0    |  |  |
| ID28               | 0.71  | 0.68                   | 0.75     | 0.75        | 0.01 |  |  |
| ID29               | 0.25  | 0.29                   | 0.2      | 0.13        | 0.03 |  |  |
| ID30               | 0.92  | 0.88                   | 0.93     | 0.93        | 0.05 |  |  |
| ID31               | 0.98  | 0.95                   | 1        | 1           | 0    |  |  |
| ID32               | 1     | 0.92                   | 1        | 1           | 0    |  |  |
| ID33               | 0     | 0                      | 0        | 0           | 0.26 |  |  |
| ID34               | 0.89  | 0.91                   | 0.9      | 1           | 0.01 |  |  |
| ID35               | 0.85  | 0.87                   | 0.82     | 1           | 0.03 |  |  |
| ID36               | 0.5   | 0.62                   | 0.51     | 0.37        | 0.03 |  |  |
| ID37               | 0.74  | 0.69                   | 0.5      | 1           | 0.02 |  |  |
| ID38               | 0.94  | 0.95                   | 0.91     | 1           | 0.01 |  |  |
| ID39               | 0.92  | 0.89                   | 0.8      | 1           | 0.01 |  |  |
| ID40               | 0.95  | 0.95                   | 0.91     | 1           | 0    |  |  |
| ID41               | 0.61  | 0.65                   | 0.67     | 1           | 0.04 |  |  |
| ID42               | 0.86  | 0.84                   | 0.8      | 1           | 0.01 |  |  |
| ID43               | 0.9   | 0.64                   | 1        | 1           | 0    |  |  |
| ID44               | 0     | 0                      | 0        | 0           | 0    |  |  |
| ID45               | 0.01  | 0                      | 0.01     | 1           | 1.06 |  |  |
| ID46               | 0.6   | 0.56                   | 0.69     | 0.52        | 0    |  |  |
| ID47               | 0.26  | 0.29                   | 0.16     | 1           | 0.09 |  |  |
| ID48               | 0.68  | 0.62                   | 0.8      | 0.67        | 0    |  |  |
| ID49               | 0.85  | 0.79                   | 0.86     | 1           | 0.01 |  |  |
| ID50               | 0.69  | 0.6                    | 0.44     | 1           | 0.03 |  |  |
| ID51               | 0.35  | 0.23                   | 0.15     | 1           | 0.07 |  |  |
| ID52               | 0.23  | 0.35                   | 0.18     | 1           | 0.16 |  |  |
| ID53               | 0.12  | 0.21                   | 0.09     | 1           | 0.16 |  |  |
| ID54               | 0.01  | 0.01                   | 0.01     | 0.67        | 1.72 |  |  |
| ID55               | 0.37  | 0.35                   | 0.29     | 1           | 0.07 |  |  |
| ID56               | 0.72  | 0.6                    | 0.67     | 0.67        | 0.01 |  |  |
| ID57               | 0     | 0                      | 0        | 0           | 0    |  |  |
| ID58               | 0.22  | 0.25                   | 0.18     | 0.29        | 0.09 |  |  |

![](_page_51_Figure_0.jpeg)

Figure S19: Seizure detection with 18 subjects pre-training on a selection of 40 unseen subjects. (a) Seizure detection results of MVPFormer-S on 40 unseen subjects which are part of the training set for the 58-subjects model: the F1-score, sensitivity, and fp/h are reported. The raw and episodic F1-scores are notably lower here with respect to the 58-subjects model. This is expected given the 58-subject model is pre-trained on these subjects. These results are a subset of those presented in the Results section. (b) Cohen's kappa is used to measure the agreement between the artificial assistant and the human expert. The average kappa is 0.56, competitive with expert agreement but, as expected, reduced from the 58-subjects pre-trained model.

#### G.11 Ablation of the prediction task

We design MVPFormer with a two-phase training regime. First, during the generative pre-training task MVPFormer learns to predict the neuronal activity. Second, during the classification task it needs to correctly classify ictal periods. To determine the significance of the generative task on the classification task, we train MVPFormer only on the classification task and compare its performance with the full architecture. Figure S20 and Table S17 clearly indicate that the generative task is of fundamental importance to the overall architecture, with a Kappa score decrease to 0.52. This is below the original result of 0.54 and below the human agreement threshold. Moreover, the distribution of agreement has flattened, with an overall decrease of performance across the board and an increase of subjects with no agreement. This suggests that without pre-training the generalization capability of MVPFormer suffers. Therefore, the generative task is necessary and is a significant contributor to learning.

![](_page_52_Figure_2.jpeg)

Figure S20: Seizure detection with no generative pre-training and 18 subjects classification training. (a) Seizure detection results of MVPFormer-S on 40 unseen subjects which are part of the training set for the 58-subjects model: the F1-score, sensitivity, and fp/h are reported. The raw and episodic F1-scores are notably lower here with respect to the 58-subjects model. This is expected given the 58-subject model is pre-trained on these subjects. These results are a subset of those presented in the Results section. (b) Cohen's kappa is used to measure the agreement between the artificial assistant and the human expert. The average kappa is 0.56, competitive with expert agreement but, as expected, reduced from the 58-subjects pre-trained model.

|              |       |          | lassincati | on metrics  |      |
|--------------|-------|----------|------------|-------------|------|
| ~ • • •      |       | Raw      | 0.1        | Episodic    |      |
| Subject      | Kappa | f1-score | f1-score   | sensitivity | tp/h |
| ID19         | 0.37  | 0.46     | 0.3        | 0.19        | 0.05 |
| ID20         | 0.99  | 0.98     | 1          | 1           | 0    |
| ID21         | 0.47  | 0.3      | 0.29       | 1           | 0.04 |
| ID22         | 0.74  | 0.7      | 0.73       | 1           | 0.02 |
| ID23         | 0.37  | 0.51     | 0.36       | 0.24        | 0.07 |
| ID24         | 0.88  | 0.84     | 0.9        | 0.93        | 0.05 |
| ID25         | 0     | 0        | 0          | 0           | 0    |
| ID26         | 0.43  | 0.19     | 0.29       | 1           | 0.06 |
| ID27         | 0.95  | 0.84     | 1          | 1           | 0    |
| ID28         | 0.73  | 0.8      | 0.75       | 0.75        | 0.01 |
| ID29         | 0.17  | 0.28     | 0.17       | 0.1         | 0.01 |
| ID30         | 0.9   | 0.85     | 0.93       | 0.93        | 0.05 |
| ID31         | 0.98  | 0.95     | 1          | 1           | 0    |
| ID32         | 0.98  | 0.93     | 1          | 1           | 0    |
| ID33         | 0     | 0        | 0          | 0           | 0    |
| ID34         | 0.78  | 0.84     | 0.82       | 1           | 0.02 |
| ID35         | 0.71  | 0.77     | 0.74       | 1           | 0.05 |
| ID36         | 0.59  | 0.68     | 0.53       | 0.37        | 0.01 |
| ID37         | 0.66  | 0.61     | 0.5        | 1           | 0.02 |
| ID38         | 0.59  | 0.64     | 0.59       | 1           | 0.04 |
| ID39         | 0.76  | 0.71     | 0.57       | 1           | 0.02 |
| ID40         | 0.67  | 0.77     | 0.57       | 0.8         | 0.02 |
| ID41         | 0.75  | 0.72     | 0.75       | 1           | 0.02 |
| ID41<br>ID42 | 0.15  | 0.72     | 1          | 1           | 0.02 |
| ID42<br>ID43 | 0.15  | 0.14     | 0          | 0           | 0.06 |
| ID45<br>ID44 | 0.15  | 0 76     | 1          | 1           | 0.00 |
| ID44<br>ID45 | 0.03  | 0.03     | 0.04       | 1           | 033  |
| ID45<br>ID46 | 0.05  | 0.05     | 0.04       | 0.62        | 0.00 |
| ID40<br>ID47 | 0.0   | 0.02     | 0.08       | 1           | 0.03 |
| ID47<br>ID48 | 0.15  | 0.72     | 1          | 1           | 0.02 |
| ID40         | 0.08  | 0.10     | 0.11       | 1           | 0.71 |
| 1D49<br>ID50 | 0.08  | 0.08     | 0.11       | 1           | 0.71 |
| ID50<br>ID51 | 0.50  | 0.30     | 0.24       | 1           | 0.07 |
| 1D51<br>1D59 | 0.52  | 0.32     | 0.22       | 1           | 0.04 |
| 1D52<br>1D52 | 0.5   | 0.29     | 0.5        | 1           | 0.04 |
| 1D55<br>1D54 | 0.2   | 0.33     | 0.1        | 1           | 0.13 |
| ID54<br>IDFF | 0.01  | 0.51     | 0.29       | 0.33        | 0.01 |
| 1D55<br>1D56 | 0.35  | 0.25     | 0.25       | 1           | 0.08 |
| ID50<br>ID57 | 0.82  | 0.01     | 0.8        | 0.07        | 0    |
| ID57         | 0     | 0        | 0          | 0           | 0    |
| ID58<br>ID50 | 0.17  | 0.2      | 0.21       | 0.71        | 0.26 |
| ID59         | 0     | 0        | 0          | 0           | 0    |
| ID60<br>ID61 | 0.16  | 0.08     | 0.36       | 0.25        | 0.02 |
| ID61         | 0     | 0        | 0          | 0           | 0    |
| ID62         | 0.42  | 0.44     | 0.41       | 0.86        | 0.07 |
| ID63         | 0.54  | 0.25     | 0.57       | 0.5         | 0.01 |
| ID64         | 0.62  | 0.59     | 0.65       | 0.5         | 0.03 |
| ID65         | 0.31  | 0.24     | 0.39       | 1           | 0.18 |
| ID66         | 0.69  | 0.52     | 0.5        | 1           | 0.02 |
| ID67         | 0     | 0        | 0          | 0           | 0.13 |
| 1D68         | 0.78  | 0.33     | 0.5        | 0.33        | 0    |

Table S17: Seizure detection with no generative pre-training and 18 subjects classification training. Kappa is the inter-rater agreement. The classification metrics report the raw and episodic metrics relevant for the seizure classification task. Classification metrics

#### G.12 Evaluation on traditional long-term forecasting task

To provide a comprehensive evaluation of MVPA, we compare MVPFormer with existing SOTA architectures on a classical long-term forecasting task. Table S18 reports the results of MVPFormer, the vanilla Transformer [14], PatchTST [3], and TimeMixer [2] on the ETTh1, ETTh2, and Weather datasets [10, 11]. These datasets represent a well-known benchmark that allows us to decouple MVPFormer from the specific clinical setting. The lookback window is fixed to 96, while the forecast is performed at lengths of 96, 192, 336, and 720. These settings are well established in the literature [2]. MVPFormer notably surpasses the vanilla Transformer and is competitive with established architectures designed specifically for this task, achieving the best or second best result in most cases. Moreover, MVPFormer is on average 2x faster to train than the vanilla Transformer and 1.4x faster than PatchTST — still slower than TimeMixer, which a fully MLP-based model —, making it an excellent choice in many scenarios. Therefore, we have shown that MVPFormer and MVPA have a wide applicability and transfer their performance from the clinical task — for which they were designed — to more general time-series tasks as well.

Table S18: Classical time-series forecasting benchmark. MVPFormer is compared with SOTA architectures PatchTST and TimeMixer on the time-series forecasting task using the ETTh1, ETTh2, and Weather datasets. The lookback window is fixed at 96 and the forecasting length varies between 96 to 720. The vanilla Transformer is also included as a point of comparison. In bold are the best MSE results, in italics are the second best. MVPFormer notably outperforms the vanilla Transformer and is competitive with both PatchTST and TimeMixer, having either the best or second best result in most cases.

| Model                                       |        | MVP  | Former (ours) | Transf | tormer | Patch | IST  | TimeN | /lıxer |
|---|--------|------|---------------|--------|--------|-------|------|-------|--------|
| Metric                                      | Length | MSE  | MAE           | MSE    | MAE    | MSE   | MAE  | MSE   | MAE    |
| <b>Ե</b> ՐՐՆ 1                              | 96     | 0.38 | 0.40          | 0.83   | 0.72   | 0.38  | 0.40 | 0.37  | 0.40   |
|   | 192    | 0.45 | 0.44          | 0.96   | 0.78   | 0.43  | 0.43 | 0.44  | 0.43   |
|   | 336    | 0.49 | 0.46          | 1.04   | 0.83   | 0.47  | 0.46 | 0.50  | 0.46   |
|   | 720    | 0.49 | 0.48          | 1.16   | 0.86   | 0.52  | 0.51 | 0.49  | 0.48   |
|   | 96     | 0.30 | 0.35          | 2.64   | 1.30   | 0.31  | 0.35 | 0.29  | 0.35   |
| $\mathbf{F}\mathbf{T}\mathbf{T}\mathbf{b}9$ | 192    | 0.37 | 0.40          | 3.48   | 1.47   | 0.38  | 0.40 | 0.37  | 0.39   |
| E11112                                      | 336    | 0.42 | 0.43          | 4.07   | 1.62   | 0.43  | 0.44 | 0.43  | 0.44   |
|   | 720    | 0.43 | 0.45          | 3.28   | 1.52   | 0.43  | 0.45 | 0.47  | 0.47   |
|   | 96     | 0.17 | 0.22          | 0.33   | 0.38   | 0.17  | 0.22 | 0.16  | 0.21   |
| Weather                                     | 192    | 0.21 | 0.26          | 0.51   | 0.50   | 0.22  | 0.26 | 0.21  | 0.25   |
| weather                                     | 336    | 0.28 | 0.30          | 0.62   | 0.56   | 0.28  | 0.30 | 0.26  | 0.29   |
|   | 720    | 0.35 | 0.35          | 0.91   | 0.70   | 0.36  | 0.35 | 0.35  | 0.35   |